Biolon

Biolon:

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Biolon

Biolon:

Transparent Ø 120 mm thickness: 0,5mm, 0,75mm, 1 mm, 1,5 mm, 2 mm, 3,0 mm
Special sizes in round or squares are available on request

Values of the material:

Density:	1,27 g/cm ³
Rockwell Hardener:	R104 - 107
E-Module:	2050 N/mm ²
Notch resistance:	16 kJ/m ²
Viscose value:	70 °C
Water absorption:	≤ 0,1 %

Field of application:

Biolon is suitable for the manufacture of splints, support of medicine, orthodontic aids and bandage plates.

Processing:

The material can excellently be thermoformed condition. The blue protection layer has to be removed before thermoforming. The material is plastified in a suitable thermoforming unit (Vacformat / Druformat) until remaining traces can be seen after having touched it with a stump, metal instrument. The temperature that will be achieved is 160 - 170 °C. Compared to similar materials Biolon distinguishes by a heating time that is reduced by 25 %. The individual heating times for the different thicknesses of Biolon in combination with Dreve thermoforming units can be taken from the separate time table.

With the help of vacuum or pressure the plate is formed onto the object and distinguishes by its very good edge strength and high transparency. In order to maintain the transparency plaster models should be watered. In addition, an insulating agent on alginate basis can be used (Isolat film). Biolon is treated with hard metal cutters and grinding paper with rough grains in order to avoid a smearing of the edges by the influence of heat. Thin material that can also be cut with the scissors. It compounds with the resins on acrylate basis and with the help of a bonding agent it compounds with corresponding light curing composite materials (Lightdon Splint).

Important: Biolon has to be stored under dry conditions! Due to the chemical composition of the blank it might suck humidity from the surrounding air during storage. This could cause air bubbles in the material while preheating the blank.

Information on health:

Biolon is marked with a CE label and thus fulfils the general requirements of the medical device guideline. Biolon is successfully tested concerning biocompatibility.

Dreve Dentamid GmbH,
Tel.: ++49 2303 8807-0
info@dreve.de

Max-Planck-Str. 31,
Fax: ++49 2303 8807-55
www.dreve.com

59423 Unna / Germany