# Essix<sup>®</sup> Laminated Sports Mouthguard Material

DIRECTIONS FOR USE  
INSTRUCCIONES DE USO  
MODE D'EMPLOI  
ISTRUZIONI PER L'USO  
GEBRAUCHSANWEISUNG

Manufacturer:  
**Dentsply International Raintree Essix**  
7290 26th Court East  
Sarasota, FL 34243 USA  
Tel: 800.883.8733  
[www.essix.com](http://www.essix.com)

Biostar and Druformat are not trademarks of  
Dentsply Sirona. ©2018 Dentsply Sirona

#RT-DFU-LAMINATEDMG RTE-046-18 Rev.5 2/18

**R<sub>x</sub>**  
ONLY

  
SINGLE  
USE ONLY

# Essix® Laminated Sports Mouthguard Material

ENG

(Instructions for use on Essix® Vacuum Machines only)\*

## Intended Use:

Essix® Laminated Sports Mouthguard Material is thermoformed to fabricate custom sports mouthguards.

## Indications for Use:

Essix® Laminated Sports Mouthguard Material is indicated for the fabrication of sports mouthguards.

## Contraindications:

This product is contraindicated for patients with a history of allergic reactions to plastics.

## ⚠ Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

## Precautions:

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in a cool, dry place.

Single Use Only. Risks of reuse include infection and cross-contamination to patient and healthcare provider as well as other wear-related irregularities.

**Adverse Reactions:** Allergic reactions.

## Step by Step Instructions:

1. Sports mouthguards should be trimmed high to the vestibular margins, leaving the frenums. Full arch models should have the palate trimmed away.
2. **IMPORTANT:** Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
3. Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
4. Place the plastic into the frame on the vacuum machine and the model on the base plate.
5. Swing the heating element over the plastic. Observe the material as it heats.
6. Material will sag approximately ½" – this is the optimal thermoforming time.
7. Move the heating element away from the plastic.
8. Turn on the vacuum and lower the plastic over the model.
9. Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it “shrinks” the material around the cast, realizing a tighter fit.
10. Let vacuum continue to run until plastic is cool.
11. Trim and finish appliance as desired.

\*Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic “folds” when thermoformed, reduce heating time until material does not fold.

Thickness	Vacuum Heat Time**	Essix® SelectVac®	Drufomat
4mm (.160")	75 seconds	2:30	2:00
5mm (.200")	110 seconds	3:30	2:25

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