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rapporto metallo vergine/revert

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
©2018 Dentsply Sirona Inc. All rights reserved. NPK1817 Rev. 12 (07/18)

Indications For Use

Vitallium 2000 Alloy is a Cobalt Chromium Molybdenum Based Casting Alloy designed for use when casting removable partial denture frameworks.

Contraindications

Vitallium 2000 Alloy is contraindicated for patients and users with a history of sensitivity to cobalt, chromium, and molybdenum.

Warnings

- Prolonged exposure to alloy dust and/or fumes may lead to lung irritation, and/or pulmonary complications. Use appropriate engineering controls to limit exposure. For excessive inhalation of dust or fumes, seek medical advice.
- This alloy contains chromium. Some compounds of these elements are potential carcinogens. Advise your health care provider of exposure to these elements.
- This alloy contains elements that are known to the State of California to be carcinogenic.

Precautions

- When melting, grinding, or polishing alloys, use adequate ventilation, vacuum systems, protective eyewear, and protective masks and clothing.

Adverse Reactions

- Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications.
- Certain components of alloys are potential carcinogens. See Warnings.

PROCESSING INSTRUCTIONS	
Investment Recommendations	
Use Dentsply Sirona’s Ethyl Silicate Bonded Investment System	REF#
Vitallium 2000 InvestmentN022665
Vitallium 2000 Investment Liquid Set (A&B)N022608

Duplicating Material Recommendations
Use **Perflex**® Reversible Hydrocolloid Duplicating Material (REF# N021055) with **Vitallium 2000** Investment.

Silflex® III Polyvinyl Siloxane Duplicating Material (REF# N021501) is recommended for use for model duplication for cases with precision attachments and for combination cases.

NOTE: Follow the instructions for the investment and duplicating material that you are using for refractory model and investing steps.

EQUIPMENT
BURNOUT FURNACE REQUIREMENTS
For Vitallium 2000 Investment
High temperature gas fired burnout furnace capable of attaining 2000°F (1095°C).

CASTING MACHINES
Dentsply Sirona ECM IV Casting Machine or a Torch Casting Machine - i.e. Oxygen acetylene multi-orifice torch & spring activated centrifugal casting arm.

BURNOUT INSTRUCTIONS
Preheat the investment molds and a clean casting crucible to a final temperature of 2000 °F (1095°C) following the recommended burnout cycle for the investment being used. *Make sure that crucibles are only used for this alloy and that used crucibles are clean of slag and any remaining metal.*

CASTING INSTRUCTIONS
Virgin Metal/Revert Ratio: To maintain the physical and mechanical properties of this alloy, it is recommended that at least 50% new metal be used for each casting. Revert (buttons) must be sandblasted and cleaned prior to reuse.

By adhering to this recommendation, the chemical composition of the alloy will be maintained and excess oxide or slag build up in the alloy (a major cause of pits in castings) will be avoided.

DENTSPLY SIRONA'S ECM IV Casting Machine Settings
Casting Temperature Setting:..... 2750°F Vitallium Scale
Soak Timer:..... 5 seconds
Melt Selector Switch:..... On High

ACCELERATION REDUCER SETTING
Small Castings - 50
Large Casting - 75 to 100
RPM: 380-400
Adjust centrifugal arm counter weight to the proper position to provide for a smooth and balanced casting operation.

Follow the instructions for the machine that you are using. Melt and cast the alloy into the preheated investment mold and allow the molds to bench cool before further processing.

DO NOT FORCE COOL OR QUENCH THE MOLDS AS IT MAY CAUSE WARPAGE OF THE CASTINGS & UNDESIRABLE CHANGES IN THE MECHANICAL PROPERTIES OF THE ALLOY.

TORCH CASTING INSTRUCTIONS
***Warning:** Use extreme caution during torch melting operations. Follow the instructions and warnings of the torch manufacturer and wear eye protection designed for use when operating an oxygen acetylene torch.*

Use 3-5 turns on the casting arm depending on the tension of the spring. Set the counter weights at the correct position according to the size of the case to be cast.

Use a multi-orifice oxygen acetylene torch pressure regulators and flow meters.

Place the preheated mold and crucible in the casting machine. Place the required amount of alloy in the crucible.

Light the torch and with a circular motion move the flame slowly over the metal. Continue this motion to uniformly heat and melt the alloy. When the alloy begins to melt, the ingots and buttons will collapse. Add flux and continue to uniformly heat the alloy until the molten alloy can be seen to slump and move under the pressure of the flame. DO NOT overheat the alloy.

NOTE: The oxide film should not break.

Release the casting arm and allow the arm to coast to a complete stop. Remove the mold and allow it to bench cool before further processing.

DO NOT FORCE COOL OR QUENCH THE MOLDS AS IT MAY CAUSE WARPAGE OF THE CASTINGS & UNDESIRABLE CHANGES IN THE MECHANICAL PROPERTIES OF THE ALLOY.

INSTRUCTIONS FOR CLEANING PARTIAL DENTURES
Vitallium 2000 Alloy dentures can be cleaned using commercial denture cleansers or denture pastes which indicate on their labels that they are suitable for cleaning "bridges, orthodontic appliances and/or partial dentures." The use of cleansers having high hypochlorite contents and the use of household bleach is not recommended as they may cause tarnish or pitting of the framework.

INSTRUCTIONS FOR THE REPAIR OF VITALLIUM ALLOY FRAMEWORK
Frameworks can be repaired using Vitallium Welding Rod (REF# N041200 & 01) or White Gold Solder (REF# N085100 or N085200).

PHYSICAL & MECHANICAL PROPERTIES
0.2% Yield Strength..... 87,000 psi (600 MPa)
Ultimate Tensile Strength..... 124,000 psi (855 MPa)
Modulus of Elasticity29 X 10⁶ psi (200 GPa)

Percent Elongation..... 6.0%
Vickers Hardness410 HV5
Density 8.3 g/cm³
Melting Range..... 2282-2500 °F (1250-1370 °C)

Composition (% by mass)
Cobalt (Co).....63.1%
Chromium (Cr)28.5%
Molybdenum (Mo).....6.0%
Mn, Si, C<1.0%
Nickel and Beryllium Free Alloy

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Vitallium® Trademark is licensed to Dentsply Sirona by Howmedica Osteonics Corporation.
Perflex®, Polyflex®, and Silflex® are registered trademarks of Dentsply Sirona.