digit Power® Dispenser

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U.S.A.: Rx only.

Table of Contents-English

1. PRODUCT DESCRIPTION	2
1.1 Indications.	2
1.2 Contraindications	2
1.3 Delivery Forms	2
1.4 Compatible Materials.	2
2. GENERAL SAFETY NOTES	2
2.1 Warnings	2
2.2 Precautions	3
2.3 Adverse reactions	3
2.4 Storage conditions	3
3. STEP-BY-STEP INSTRUCTIONS	4
3.1 Description of digit Power® Dispenser	4
3.2 Unpacking the System	5
3.3 digit Power® Dispenser Installation	5
3.4. Charging the Handpiece and Optional Wireless Foot Control	5
3.5 Synchronizing the Handpiece and Optional Wireless Foot Control	6
3.6. Preparation for Use	6
3.7 Loading the digit Power® Dispenser	9
3.8 Dispenser operation	9
3.9 Dispenser disassembly/plunger retraction.	11
4. HYGIENE & MAINTENANCE	11
5. TECHNICAL SPECIFICATIONS	14
6. CLASSIFICATIONS	14
7. SYMBOL IDENTIFICATION	15
8. DISPOSAL OF UNIT	16
9. ELECTRO-MAGNETIC COMPATIBILITY PRECAUTIONS ..	16
10. TROUBLESHOOTING	18
11. SERIAL NUMBER	19
12. WARRANTY	19
13. TECHNICAL SUPPORT	19

1. PRODUCT DESCRIPTION

digit Power® Dispenser is a battery powered dental material dispensing handpiece for precision intraoral placement of dental materials delivered via digit Power® cartridges. digit Power® Dispenser is designed to improve the ergonomics of material delivery when compared to conventional, mechanical systems while allowing for easy and precise intraoral placement of dental materials. Portability and versatility are built into the design of digit Power® Dispenser with long battery life, multiple speed options on the handpiece and an optional Bluetooth® wireless foot control for full variable speed.

digit Power® Dispenser has a built-in, automatically activated timer, for convenient monitoring of work time (WORK TIME) of material while dispensing, and dispensed material set time/mouth removal time (SET TIME/MRT). Timer preset modes correspond to Aquasil® Ultra+ Smart Wetting® Impression Material FAST SET and REG SET. A convenient user-specified timer option is also built into the handpiece.

1.1 INDICATIONS

The digit Power® Dispenser is indicated for intraoral delivery of materials provided in digit Power® cartridges.

The digit Power® Dispenser is readily suitable for:

- Delivery of low viscosity Aquasil® Ultra+ Smart Wetting® Impression Material in digit Power® cartridges for single or multiple preparation impression technique.

1.2 CONTRAINDICATIONS

None known.

1.3 DELIVERY FORMS*

**Some delivery forms may not be available in all countries*

digit Power® Dispenser is available as a complete unit, with battery charging stand.

Product components:

- Handpiece Assembly
- Removable Autoclavable Cartridge Cap
- Charging Base Assembly
- Wall Plug-in Power Supply
- International Power Adapters

Available accessories:

- Wireless foot control for dispensing speed control
- Disposable Polyethylene Barrier Sleeves

1.4 COMPATIBLE MATERIALS

The digit Power® Dispenser is compatible with materials provided in digit Power® cartridges.

Materials not provided in digit Power® cartridges are not intended to be used with the digit Power® Dispenser.

See the Complete Directions for Use of selected material for compatibility and use.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

- To prevent damage, charge the handpiece using only the digit Power® Dispenser charging base for handpiece and power supply and cord.
- Sterilizing the handpiece will cause component damage to the handpiece and the sterilizing equipment, and may cause personal bodily injury.
- The cartridge cap must be steam sterilized before first use and between patients to prevent patient cross contamination. See **Section 4** for the Hygiene and Maintenance.
- It is the responsibility of the Healthcare Professional to determine the

appropriate uses of this product and to understand:


- the health of each patient
- the dental procedures being undertaken
- applicable industry and governmental agency recommendations for infection control in dental healthcare settings
- requirements and regulations for safe practice of dentistry
- these Directions for Use in their entirety
- Per FCC Part 15.21, changes or modifications not expressly approved by the party responsible for compliance (i.e., the manufacturer) could void the user's authority to operate this equipment.
- Failure to follow recommendations for environmental operating conditions (see **Section 5 and 9** for Specifications) could result in injury to patients or users.
- Inspect the handpiece system before each use for worn, loose or damaged parts.
- To prevent bodily injury and damage to the device, do not sterilize the handpiece, charging base for handpiece, optional wireless foot control or power supply and cord. Disinfect the handpiece and charging base for handpiece using only the tested and approved disinfectants listed in **Section 4**, Hygiene and Maintenance.
- The handpiece, optional wireless foot control, charging base for handpiece and power supply and cord are not waterproof. To prevent damage to the equipment, contamination or bodily injury, do not immerse any of these components in water or a chemical solution.
- Use only components and accessories listed in **Section 3.2** of this manual. Failure to do so will void the warranty, may decrease system performance, may lead to unsafe operation, may negatively affect Electro-Magnetic Compliance performance and result in non-compliance.
- There are no user-serviceable items in the handpiece, power supply and cord, cartridge cap, optional wireless foot control or charging base for handpiece. Opening any of these units may result in unsafe operation and will void the warranty.
- According to IEC 60601-1/UL60601-1, this device must not be used in the presence of a flammable anesthetic gas mixed with air, oxygen, or nitrous oxide. (Note: nitrous oxide by itself is not a flammable anesthetic gas.).
- User should not touch the patient and the accessible charging base contacts or USB contact simultaneously.

2.2 Precautions

- This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Always insert the cartridge/mixtip assembly into the cartridge cap first and then connect the cartridge cap to the handpiece, otherwise a malfunction may occur.
- Ensure plungers are fully retracted by visual inspection prior to placing cartridge cap/cartridge/mixtip assembly onto handpiece.
- Improper alignment of the cartridge within the cartridge cap, and/or the cartridge cap and handpiece may result in a device malfunction. An audible click will be heard when the cartridge cap is secured to the handpiece.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- As an additional precautionary measure, digit Power® Dispenser may be protected from gross debris but not from all contamination by applying a protective barrier. Reprocess reusable components after each use according to instructions in **Section 4**, Hygiene and Maintenance section.
- Do not autoclave the digit Power® Dispenser handpiece.
- Do not spray disinfectant or other fluids directly onto the handpiece, optional wireless foot control, charging base for handpiece, or power supply and cord to avoid liquid from pooling on the component. The user should spray solution onto a cloth or use a wipe to disinfect items per the instructions in **Section 4**.
- Prevent liquids from entering openings on the handpiece and charging base.
- Ensure handpiece and charging base are fully dry before placing handpiece into charging base.
- If plungers do not retract automatically at defined limits, or when manually

commanded (see Operation **Section 3.8**), allow cartridge to empty away from patient field, turn unit off (see Operation **Section 3.6**) and return unit for service. Do not attempt to disassemble or use further.

- Use only material in designated digit Power® cartridges. Suitable cartridges will be marked PWR on tab. Do not use digit Power® Dispenser with materials in digit® Targeted Delivery System cartridges designed for manual dispensing. Device and cartridge malfunction may occur.
- The handpiece is designed to be **lube-free**. Lubrication may cause damage to the handpiece.
- Oil and/or dirt may damage the motor, electronics and battery located inside the handpiece.
- The batteries are **not** user replaceable. If not functional, the units should be returned to the listed repair center for replacement.
- Do not place the system on or next to a radiator or other heat source. Excessive heat may damage the system's electronics.
- Inadvertent system shutdown may occur in the presence of strong non-compliant radio frequency-generating components.

2.3 Adverse reactions

None known.

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store at temperatures between -20°C/50°C.
- Use the product at room temperature.
- Protect from moisture.
- Store at relative humidity range 5 - 95% (non-condensing).

3. STEP-BY-STEP INSTRUCTIONS

3.1 Description of digit Power® Dispenser

Removable, Autoclavable Cartridge Cap

Charging Base for Handpiece

(See Section 3.4)

Charge Indicator

(Charging - See Section 3.4)

(Synchronization - See Section 3.5)

(Operation - See Section 3.8)

Wireless Foot Control (optional)

(Charging - See Section 3.4)

(Synchronization - See Section 3.5)

(Operation - See Section 3.8)

On/Off/Mode Button

Mode Indicator

Handpiece

(See Section 3.6)

3.2 Unpacking the System

As you unpack your digit Power® Dispenser, verify that the following components and accessories are included:

System:

- 1- Handpiece with 1 preassembled cartridge cap
- 2- Cartridge caps
- 1- Charging Base for Handpiece
- 1- Power Supply and Cord
- 1- Package of 50 Disposable Polyethylene Barriers (Single Use)
- 1- Package International Power Adapters

Optional Wireless Foot Control

(Available Separately):

- 1- Foot Control
- 1- Power Supply and Cord
- 1 Package- International Power Adapters

Cross-contamination

Do not reuse single use products. Dispose of in accordance with local regulations.
Reprocess reusable products as described in **Section 4, Hygiene and Maintenance** section.

*Cartridge caps should be steam-autoclave sterilized prior to first use.
Reprocessing is mandatory following each use thereafter.*

Additional Cartridge caps are available for purchase.

3.3 digit Power® Dispenser Installation

System Setup

- Both the handpiece and optional wireless foot control must be charged for at least 90 minutes prior to first use.

3.4. Charging the Handpiece and Optional Wireless Foot Control

- When the battery icon is empty (no bars), the handpiece must be charged.
- When the charge indicator light on the optional wireless foot control is solid orange, the foot control must be charged.
- A fully discharged handpiece or foot control will need approximately 3 hours to fully recharge.

i When the handpiece is not in use, storage in the charging base is recommended. The handpiece has built-in safety circuitry that prevents overheating or overcharging. The battery does not have “memory” and will always accept a full recharge.

1. Plug the power supply into a wall outlet.

Electrical Malfunction

- Connect to single phase AC power 100-240V power only. Otherwise, malfunction will occur.
- Do not unplug the power supply by pulling on the cord.
- Ensure Handpiece, charging base and optional wireless foot control are clean and fully dry before connecting power and/or inserting handpiece into charging base.

2. Plug the power supply cord into the charging base (for handpiece charging) or directly into the optional wireless foot control (for foot control charging).

3. Insert the handpiece into the charging base. Ensure contacts are clean and aligned. Handpiece should be stable in the charging base.

i A successful connection is shown as follows:
– The battery icon flashes and scrolls through all states of charge for handpiece. Battery icon will be solid when handpiece battery is at full charge.
– A solid Green Charge Indicator on top of the optional, wireless foot control appears.

i Unplug the power supply from the wall outlet (power mains) to remove all power to the unit being charged.

- The handpiece should be recharged after each full day's use.
 - If necessary, the handpiece may be "quick charged" for a single use in 15 minutes.
 - When in the charging base, the handpiece motor is not allowed to operate.
 - When fully charged, the battery icon in the display will be solid (not scrolling).
- The optional wireless foot control should be recharged when the charge state indicator light is solid orange.
- The foot control may be operated while charging.


CHARGE STATUS FOR FOOT CONTROL	
Solid green light	95-100%
Blinking green light	30-95%
Solid orange light	Less than 30%

3.5 Synchronizing the Handpiece and Optional Wireless Foot Control

Follow this procedure when:

- The optional wireless foot control or handpiece is replaced for any reason.
- The handpiece and optional wireless foot control do not seem to be communicating properly.

i Handpiece must be placed in charging base and connected to power supply in order to pair with the foot control.

1. Place the optional wireless foot control within 10 feet of the charging base for handpiece.
2. Invert the optional wireless foot control.
3. Ensure foot control is on. To turn foot control on, depress the power button for 1 second.
4. Depress the "Sync" button on the selected optional wireless foot control for at least 2 seconds and release.
5. Immediately depress the On/Off/Mode button on the handpiece once for 3 seconds and release to complete pairing.
6. A successful sync is indicated by the appearance of the Bluetooth® icon on the handpiece changing from flashing to solid, and the sync indicator light on the foot control is blinking green.
7. The Bluetooth icon will change to a "foot on the pedal" when the paired foot control is pressed.

Synchronization Indicator Status (Handpiece and Optional Wireless Foot Control Sync Mode)

Device Indicator	Ready to Pair	Successful Pairing	Failed Pairing
Handpiece Display	Bluetooth Icon Blinking	Solid Bluetooth Icon	Bluetooth Icon Gone
Foot Control LEDs	Blinking Blue	Blinking Green	Blinking Orange

3.6. Preparation for Use

1. Position the optional wireless foot control for use, ensuring that the floor is level.
2. Ensure that the handpiece (including cartridge cap) has been reprocessed according to **Section 4**, Hygiene and Maintenance.

Cross-contamination


Ensure the disposable polyethylene barrier has been correctly applied over the entirety of the handpiece and cartridge cap, with mixtip and intraoral tip protruding through hole of barrier, before beginning a procedure.

This Polyethylene Sleeve does not replace cleaning and disinfection of the dental instruments. Please clean and reprocess dental instruments after each patient as described in Section 4, Hygiene and Maintenance.

3. Description of Display Graphics

	On/Off Power Button, Timer Mode selection button
FAST SET	Fast Set Timer Mode
REG. SET	Regular Set Timer Mode
WORK TIME	User Set Timer Mode, Work Time adjustment screen
SET TIME	User Set Timer Mode, Set Time adjustment screen
	Battery level status icon
	Icon reminds user to use the ► button to increase User Set times and the ►► button to decrease User Set times
	Bluetooth logo, indicates pairing status with optional wireless foot control
SHUT DOWN	Device is powering off
ERROR CODE 1: TEMPERATURE	Unit is in overheated state
ERROR CODE 2: DRIVE TRAIN	Unit drive train is in error state
ERROR CODE 3: POWER SUPPLY	Unit power supply is in error state
ERROR CODE 4: SELF TEST	Unit self-diagnostic check is in error state

4. Powering the Handpiece On/Off

- To power on, press and hold the power/mode selector button for approximately 2 seconds.
- Handpiece will confirm proper function with slight advance and retraction of plungers.
- Refer to the table for handpiece display of battery charge level.


HANDPIECE DISPLAY BATTERY STATES		
Indicator Light	Percent of icon	Degree of Charge
	100 %	80-100%
	75 %	60-79%
	50 %	40-59%
	25 %	20-39%

- To power off, press and hold On/Off/Mode selector button.

i digit Power® Dispenser performs a self-diagnostic check upon start-up to ensure all systems are functioning properly.

1. If the display screen shows an error code (see Graphics description), power the handpiece off by pressing and holding the power/mode selector button.
2. Ensure the plungers are not fouled or blocked. Remove debris if necessary. Do not attempt to remove debris if plungers are in motion.
3. Attempt to power the handpiece on. If the unit display is normal, the handpiece has self-corrected the error and is fully functional.
4. If ERROR CODE 1: TEMPERATURE is still displayed, unit may be in an overheating condition. Power the handpiece off, and remove handpiece to a cool, dry area. Allow handpiece to cool for 10 minutes, then power the handpiece on. If the unit display is normal, the handpiece has self-corrected the error and is fully functional. Avoid exposing handpiece to conditions that caused the overheating.
5. If any error code (including ERROR CODE 1: TEMPERATURE after 10 minute cool-down as described above) is still present, note the code, power the handpiece off, and contact Dentsply Sirona Technical Support. Do not attempt any further disassembly or remedial actions. Do not attempt to use the handpiece.


5. Adjusting the Integrated Timer

- The integrated timer assists in ensuring intraoral placement of impression material is completed within material intraoral work time.
- The integrated timer automatically counts down the time from first dispensing until the recommended setting time/mouth removal time (SET TIME /MRT).
- There are three timer modes for digit Power® Dispenser with material work and set time notifications. "FAST SET" and "REG. SET" are preprogrammed to coincide with the work/set times for Aquasil® Ultra+ Smart Wetting® Impression Material in digit Power® cartridges. "USER SET" allows the user to define work and set time for the material being used.

First Use or to Change Settings

- Display and handpiece default to the "FAST SET" timer and audible signals on first use.
- For subsequent uses, the handpiece will default to the previously used setting.
- If the handpiece is paired with the optional wireless foot control, the Bluetooth icon will display.
- To move between each mode, press and release the power/mode select button.
 - o FAST SET - (intraoral WORK TIME = 0:35, SET TIME/MRT = 2:30) Corresponding to Aquasil Ultra+ Material Fast Set.
 - o REG. SET- (intraoral WORK TIME = 1:10, SET TIME/MRT = 5:00) Corresponding to Aquasil Ultra+ Material Regular Set.
 - o USER SET - User defined WORK TIME and SET TIME/MRT.

• User Set - user specifies Work Time and Set Time

1. Scroll by pressing and releasing the power/mode button until the "WORK TIME" displays. "0:00" will be flashing and the "+ -" adjust icon will appear.
2. Press the forward-located Slow ► to add time or Fast ►► speed buttons to subtract time to set the Work Time.
3. Press and release the power/mode selection button. "SET TIME" will display. The specified work time from the step above will be flashing and the "+ -" adjust icon will appear.

Note - pressing and holding either Slow ► or Fast ►► speed buttons will increment or decrement time by 10 seconds respectively.

4. Press the forward-located Slow ► to add time or Fast ►► speed buttons to subtract time to set the Set Time.

[Notice] - If the User defined Set Time is adjusted to be less than the User defined Work Time, the timer will disregard the longer user defined Work Time setting and use the Set Time for both times for audible notifications and automatic plunger retraction .


5. Press and release the power/mode selection button. "USER SET" will display. The time displayed will be the user defined Set Time from step 4 above.

i In the User Defined mode, the settings may be left to 0:00.


- The handpiece will function normally, with the audible signals disabled.
- The timer counts up from 0:00 showing elapsed time from initial dispensing.

3.7 Loading the digit Power® Dispenser


i To test functionality pre-operative, please follow directions below, extruding one cartridge, prior to use on patients.


1. Prepare a digit Power® cartridge for use by snapping the circular disk off the cartridge.


2. Insert unit dose cartridge into mixing tip. Ensure both cartridge spurs fully lock into the square holes on the side of the mixing tip.

3. With the cartridge cap removed from the handpiece, insert the digit Power® cartridge mixing tip assembly into the cartridge cap.

4. Ensure the digit Power® cartridge is locked into the cartridge cap by the retaining tab. If cartridge is misaligned (retaining tab not parallel to cartridge) cartridge cap cannot be seated onto handpiece.

5. Ensure plungers are fully retracted. Visually inspect and confirm plunger retraction before attempting to seat the cartridge cap.

6. Place cartridge cap with fully seated cartridge in place on handpiece so that the locks slide into place. An audible click will be heard when the cartridge cap is properly secured to the handpiece. Tug on the seated cartridge to confirm secure attachment.

Improper Assembly – Device Malfunction

- Ensure plungers are fully retracted before placing cartridge cap onto handpiece.
- Ensure the cartridge is fully seated into the cartridge cap.
- Ensure cartridge/handpiece assembly is fully seated onto the handpiece.
- Audible clicks may be heard to denote proper assembly.

7. Place a protective barrier sleeve over Dispenser.
8. Verify that all parts of the Dispenser are securely attached before use.
9. Charge Status: Handpiece Charge Level must display at least 1 bar. Refer to Section 3.4, Charging the Handpiece.

3.8 Dispenser operation

- The handpiece has two forward speed dispensing buttons for finger operation, operating discrete Slow and Fast speeds.
- The optional wireless foot control allows continuously variable speed control throughout the entire handpiece speed range.

i The wireless foot control connects to the handpiece via Bluetooth. Like all wireless devices, the connection is subject to interruption by environmental interference. If the connection with the handpiece is lost, complete the procedure using the manual speed buttons. When convenient, re-pair the wireless foot control and handpiece as outlined in Section 3.5.

1. Ensure the integrated timer is properly adjusted for the selected material.
2. Activate the Dispenser by depressing the fast or slow handpiece buttons, or depress the optional wireless foot control and, away from the patient, bleed a small amount through the tip to assure proper flow and a good mix.
3. Follow Directions for Use for selected material placement.
4. Speed may be changed during operation by repositioning finger from one button to the other. Speed may be varied with the wireless foot control by stepping harder or softer on the control at any time during operation.

i Start out using the slowest speed setting and adjust upward as desired after becoming familiar with the device.

5. Upon initiating a dispensing cycle (when material enters the mixing tip), a timer counts down the Work Time and provides an audible notification at three time points:

- o 10 seconds before WORK TIME **OR** if plungers are nearing full extension, cartridge nearly empty (whichever comes first), **3 short beeps**.
- o At WORK TIME, **1 long beep** – Intraoral placement should have been completed, tray should be seated.
- o At SET TIME, **3 long beeps**– Impression may be removed from the mouth.

(See Complete Directions for Use of selected material for material work and set times)

Notification Signal	Procedure Time
3 short beeps	Nearing Intraoral Worktime or end of cartridge
1 long beep	At Intraoral Worktime
3 long beeps	At Set Time/MRT

- If device is operating in user set mode, with WORK TIME/SET TIME = 0:00, there will be **no timing notification beeps** and the clock will count up once material begins dispensing, up to 9:59 (unless restarted by reloading and extruding another cartridge).

6. Stop the flow of material by releasing the button or stepping off of the wireless foot control. The digit Power® Dispenser will slightly retract the plungers to reduce continued material flow. Material dispensing may be restarted by re-depressing button or foot control.

Impression Defect-exceeding work time

- **Do not restart dispenser or attempt to place material after expiration of selected material work time.**
- **Do not attempt to remove and reload cartridge during a single clinical procedure.**
- **Ensure the loaded tray has been seated within the intraoral WORK TIME.**
- **Exceeding intraoral WORK TIME of syringed material before seating tray may lead to impression defects, requiring repeat of procedure.**
- **See Complete Directions for Use of selected material.**

7. Integrated timer – automatic plunger retraction

- **FAST SET, REG. SET, User Defined WORK TIME/SET TIME ≠ 0:00**
 - o These settings are designed for dispensing one cartridge only during a procedure.
 - o Once material is extruded, starting timer, the mode selection is disabled, and the timer may only be interrupted by power cycling the handpiece (power off, power back on).
 - o Plungers will automatically retract to fully retracted position when:
 - Timer reaches 1:00 (Fast Set) or 2:00 (Regular Set) or 0:00 (User defined ≠ 0:00) OR
 - 10 seconds after cartridge is emptied (plungers fully extended) regardless of timer OR
 - User simultaneously presses and releases Slow ► and Fast ►► speed buttons OR
 - Handpiece is power cycled off/on OR
 - Handpiece is placed in charging base
- **User Defined WORK TIME/SET TIME = 0:00**
 - o This setting will allow dispensing of more than one cartridge during a procedure.
 - o [NOTICE] – User must ensure material(s) are used within allowable work times. See complete Directions for Use of selected material.
 - o Once material is extruded, starting timer up, the mode selection is disabled. However, user may manually retract (simultaneously press and release Slow ► and Fast ►► speed buttons) and restart the timer by loading and dispensing another cartridge.

i If manually retracted, reloaded and second cartridge dispensed, the timer restarts counting up from 0:00, timing the second dispensed cartridge. Timing of first dispensed cartridge is no longer tracked.

- o Plungers will automatically retract to fully retracted position when
 - Timer reaches 9:59 OR
 - 10 seconds after cartridge is emptied (plungers fully extended) OR
 - User simultaneously presses and releases Slow ► and Fast ►► speed buttons OR
 - Handpiece is power cycled off/on OR
 - Handpiece is placed in charging base


- **ALL MODES/TIMER SETTINGS**

- o Plungers will automatically retract to fully retracted position when
 - 10 seconds after cartridge is emptied (plungers fully extended) OR
 - User simultaneously presses and releases Slow ► and Fast ►► speed buttons OR
 - Handpiece is power cycled off/on OR
 - Handpiece is placed in charging base OR
 - Timer reaches 1:00 (Fast Set) or 2:00 (Regular Set) or 0:00 (User defined ≠ 0:00)

8. Dispenser has a built-in safety mechanism to prevent cartridge over pressurization. If material is partially or completely set and/or cartridge is blocked, plungers will slow or retract. Manually retract plungers by simultaneously pressing and releasing Slow ► and Fast ►► speed buttons, and replace cartridge/mix tip assembly.

3.9. Dispenser disassembly/plunger retraction

1. Prior to disassembly, ensure Dispenser plungers are fully retracted.
 - Simultaneously press and release both Slow ► and Fast ►► speed forward speed buttons on the handpiece to retract plungers.
 - Plungers will automatically retract as outlined above without additional user input.

2. To disassemble: Press both side cartridge cap locks to release. Withdraw cartridge cap/cartridge/mix tip straight off. Do not twist the cartridge cap while removing. Slightly depress the retaining tab and withdraw the used cartridge/mix tip assembly.
3. To prepare/reprocess contaminated components for reuse, see Hygiene & Maintenance **Section 4**.

4. HYGIENE & MAINTENANCE

Cross-contamination - Infection

- Do not reuse single use products. Dispose of in accordance with local regulations.
- The barrier is designed for single use and must be disposed of after each use in accordance with local regulations. The barrier is not a replacement for cleaning, disinfection, and sterilization.
- Reprocess reusable products as described below.

NOTICE

Wrong cleaning or disinfection method.

Damage to the device.

1. Always disassemble digit Power® Dispenser before final reprocessing.
2. Do not clean or disinfect with aggressive agents (e.g. solutions based on orange oil or acetone).
3. Start reprocessing within 1 hour after use.
4. Repeated liquid contact may damage labels. **Vigorous wiping can destroy the label.** Wipe device labels gently.

For Manual Cleaning and Disinfection

- Discard used gloves according to local regulations.
- Disinfect hands with an appropriate bactericidal, virucidal, and fungicidal hand disinfectant solution according to local regulations. Use according to disinfectant solution manufacturer's Instructions for Use.
- Use a new clean pair of examination gloves.

INSTRUCTIONS FOR CLEANING, DISINFECTING AND STERILIZING CARTRIDGE CAP

Warnings	<ul style="list-style-type: none"> • These instructions are for use ONLY on the Cartridge Cap. The handpiece and charging base should be disinfected according to the procedures in the "Instructions for Cleaning and Disinfecting Other Parts" section. • Steam autoclaving sterilization is appropriate for the Cartridge Cap. • Do not allow the device to exceed 137° C.
Limitations on Reprocessing	<ul style="list-style-type: none"> • Repeated processing has minimal effect on the cartridge cap. • The cartridge cap can be subjected to the manual cleaning and disinfection process up to 800 times. • The cartridge cap can be subjected to the automated washer/disinfector process up to 800 times. • The cartridge cap can be subjected to the steam autoclave sterilization process up to 800 times. • Cold liquid immersion disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.
Point of Use	<ul style="list-style-type: none"> • Remove protective barrier and discard according to local regulations. • Remove excess soil with disposable cloth / paper wipe in combination with an alcohol-based, tuberculocidal, quaternary ammonium solution with a label claim for cleaning (e.g., VoloWipes® Disinfecting/Cleaning/Deodorizing Wipes). • Reprocess dispenser handpiece as outlined in the "Instructions for Cleaning and Disinfecting Other Parts" section. <p>Prior to Cleaning and Disinfection of the Cartridge Cap, the digit Power® Dispenser should be disassembled:</p> <ol style="list-style-type: none"> 1. Remove cartridge cap: Press both side cartridge cap locks to release. Withdraw cartridge cap/cartridge/mix tip straight off. Do not twist the cartridge cap while removing. 2. Slightly depress the retaining tab and withdraw the used cartridge/mix tip assembly. Do not attempt to refill or consolidate used cartridges. Dispose of used cartridge/mix tip assembly in accordance with local regulations. <p>It is recommended that the Cartridge cap and handpiece be reprocessed as soon as is reasonably practical following use.</p>

Containment and Transportation	No particular requirements. Handle with care.
Cleaning and Disinfection: Automated	Thermal disinfectant: Manufacturer's specification according to DIN EN ISO 15883. Cleaning program (A0 value > 600 or, at least 1 min at 90°C/194°F) as indicated by the manufacturer in the operating instructions. Follow manufacturer's recommendation for use of detergent and neutralizer, observing concentrations and contact times.
Cleaning and Disinfection (Low-level): Manual	As an alternative to automated cleaning and disinfection, the cartridge cap has to be cleaned by scrubbing with hot water and enzymatic detergent solution prepared according to the manufacturer's instructions (e.g., Resurge® Instrument Cleaning Solution). <ul style="list-style-type: none"> • Immerse cartridge cap in the detergent solution. Clean with a soft brush. • Rinse under running potable water. • Dry with a lint-free single-use cloth. Follow by disinfection as outlined below: <ul style="list-style-type: none"> • Wipe all device surfaces with a single-use cloth in combination with an alcohol-based, tuberculocidal, quaternary ammonium solution (e.g., VoloWipes® Disinfecting/Cleaning/Deodorizing Wipes). • Thoroughly disinfect the surface with the soaked cloth, paying special attention to any cavities and holes. Ensure the entire surface remains wet for the entirety of the required contact time in accordance with the disinfectant solution manufacturer's instructions for use. • Rinse for 30 seconds under running potable water. • Dry with a lint-free single-use cloth.
Packaging	Paper / plastic steam sterilization pouches (e.g., AssurePlus® Sterilization Pouches) may be used, but are not required.
Sterilization	After Cleaning and Disinfecting, steam autoclaving is required. Gravity Steam Sterilization: <ul style="list-style-type: none"> • Full Cycle: 132°C for 15 minutes. Prevacuum Steam Sterilization: <ul style="list-style-type: none"> • Full Cycle: 132°C for 4 minutes or 134°C for 3 minutes. Steam Flush Pressure Pulse Sterilization: <ul style="list-style-type: none"> • Full Cycle: 132°C for 3.5 minutes. Follow manufacturer's instructions for loading and operation cycle.
Drying	Use the drying cycle of the autoclave, minimum 30 minutes. Allow the cartridge cap to fully air dry before storage.
Maintenance	Visually inspect to ensure that all contamination has been removed. If the cartridge cap is discolored, damaged, worn, or distorted it should be discarded. No additional maintenance or lubrication is recommended.
Storage	Store the cleaned, disinfected and sterilized cartridge cap at room temperature, away from moisture or excessive humidity. Instruments steam autoclaved bagged should remain bagged until ready for use. Prior to subsequent re-use, inspect the sterilization pouch and the cartridge cap. If the integrity of the sterilization pouch has been compromised, the cartridge cap must be reprocessed prior to use. Instruments steam autoclaved unwrapped should be used immediately. Dispenser handpiece and Charging Base should be cleaned, disinfected and dried as outlined in the section below prior to storage.
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

INSTRUCTIONS FOR CLEANING AND DISINFECTING OTHER PARTS (DISPENSER HANDPIECE , CHARGING BASE)

WARNINGS	<ul style="list-style-type: none"> • The digit Power® Dispenser handpiece and charging base for the handpiece are not sterilizable by autoclave. • Intermediate-level disinfection is appropriate for the digit Power® Dispenser handpiece and charging base for the handpiece. • Do not autoclave in steam autoclave. • Do not clean/disinfect in automated washer/disinfector. • Do not immerse in liquid.
Limitations on Reprocessing	<ul style="list-style-type: none"> • Repeated reprocessing has minimum effect on these instruments. End of life is normally determined by wear and damage due to use. • The digit Power® Dispenser handpiece and charging base can be subjected to the manual cleaning and disinfection process up to 3120 times. • Cold liquid immersion disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use. <p>The Cartridge cap must be steam autoclave sterilized. Refer to cleaning, disinfection and steam sterilization procedures for the autoclavable cartridge caps in the section above.</p>
Point of Use	<ul style="list-style-type: none"> • Remove protective barrier and discard according to local regulations. • Remove excess soil with disposable cloth / paper wipe in combination with an alcohol-based , tuberculocidal, quaternary ammonium solution with a label claim for cleaning (e.g., VoloWipes® Disinfecting/Cleaning/Deodorizing Wipes). • Following handpiece Cleaning outlined below, the digit Power® Dispenser cartridge cap should be disassembled and reprocessed as outlined above. • Do not attempt further disassembly of handpiece or charging base. <p>It is recommended that the Cartridge cap and handpiece be reprocessed as soon as is reasonably practical following use.</p>
Containment and Transportation	No particular requirements. Handle with care.
Cleaning and Disinfection: Automated	Do not use automated washer/disinfectors for reprocessing digit Power® Dispenser handpiece or charging base. Component damage will occur.
Cleaning and Disinfection (Intermediate-level): Manual	<ul style="list-style-type: none"> • The digit Power® Dispenser handpiece and charging base have to be cleaned by scrubbing with a pre-impregnated wipe or disposable towel soaked with an alcohol-based , tuberculocidal, quaternary ammonium solution with a label claim for cleaning (e.g., VoloWipes® Disinfecting/Cleaning/Deodorizing Wipes). • Remove all visible soil, ensuring fluid penetrates all crevices. Do not allow solution to penetrate the casing. <p>Remove Cartridge Cap and reprocess as outlined above.</p> <ul style="list-style-type: none"> • After cleaning and disassembling Cartridge Cap, thoroughly wipe all device surfaces with a new single-use cloth in combination with an alcohol-based, tuberculocidal, quaternary ammonium solution (e.g., VoloWipes® Disinfecting/Cleaning/Deodorizing Wipes) approved according to local regulations, and use according to disinfectant solution manufacturer's Instructions for Use. Pay special attention to seams, areas around buttons, window, handpiece dock and crevices where the cap attaches. • Allow the handpiece and charging base to air dry. • For the handpiece, rinse by wiping with a sterile, clean, lint-free cloth dampened with distilled water.
Packaging	No particular requirements.
Sterilization	<p>Sterilization is not allowed. No methods have been validated.</p> <p>Do not subject components to steam autoclaving or liquid chemical sterilant immersion. Component damage will occur.</p> <p>The cartridge cap must be steam autoclave sterilized. Refer to steam sterilization procedures for the autoclavable cartridge caps in the section above.</p>
Drying	Dry the components with a sterile, clean, lint-free cloth. Allow the components to fully air dry before storage.
Maintenance	<p>Visually inspect to ensure that all contamination has been removed.</p> <p>Visually inspect power supply and cord for damage.</p> <p>Components that are damaged, worn, or distorted should be discarded.</p> <p>No additional maintenance or lubrication is recommended.</p>
Storage	<p>Store the digit Power® Dispenser handpiece and charging base at room temperature, away from moisture or excessive humidity.</p> <p>Cartridge cap should be cleaned, disinfected, sterilized and dried as outlined above prior to storage.</p>
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

Incidental surface contact of the optional wireless foot control, power supply and cord with water, soap or a water-based hospital-level disinfection solution will not damage the material of construction. Do not allow any solution to penetrate the casing.

The instructions provided above have been validated by Dentsply Sirona as being capable of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing is actually performed using equipment, materials and trained personnel in the processing facility to achieve the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.

5. TECHNICAL SPECIFICATIONS

Power Supply manufacturer: Power Supply model number:	GLOBTEK WR9QA1200USBNMEDRVB	
AC Input	Continuous (100-240 VAC)	
AC Input Current	Less than 0.7A	
AC Input Phases	Single	
AC Input Frequency	50-60 Hz	
AC Voltage Range for Power supply	120 - 240 V	
DC Output Power	5W	
DC Output Voltage and Current	+5VDC at 1A	
Output Regulation	+/- 10%	
Weight	Handpiece with cartridge cap = 120 g Foot Control = 200 g	
Dimensions	Handpiece with cartridge cap: L = 190 mm, W = 30 mm Optional Wireless Foot Control: W = 118 mm H = 40 mm	
Remote Communication Transmit/Receive	Frequency:	2405-2480 MHz
	Power:	1mW
	Channels:	40 (BLE)
	Bandwidth	2 MHz
Operating Environment	Ambient temperature:	10°C - 40°C
	Relative Humidity:	5-95% (non-condensing)
	Atmospheric Pressure:	80-106 kPa
	Altitude:	< 2,000 meters
Transport and Storage Conditions	Ambient temperature:	-20°C to 50°C
	Relative Humidity:	5-95% (non-condensing)
	Atmospheric Pressure:	50-106 kPa
	Altitude	< 5,000 meters
Optional Wireless Foot Control	Protection Class IPX1.	
Applied Parts	Single-use disposable polyethylene barrier sleeve Single-use disposable cartridge mixtip/intraoral tip assembly	

6. CLASSIFICATIONS

Type of protection against electric shock:	Class II
Degree of protection against electric shock:	Type B Applied Part
Mode of operation for Handpiece:	Operating, Charging, OFF
Mode of operation for Foot Control:	OFF, Continuous
Settings for Handpiece	3 (Fast Set, Regular Set, User Set)
Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:	Equipment not suitable for use in the presence of flammable mixtures
According to medical device directive:	I (Rule 12) (IEC 60601) UL 60601-1
Pollution Degree Classification	Pollution Degree 2
Overvoltage Category	Category II (connected to wall outlet)

7. SYMBOL IDENTIFICATION

	Class II Equipment
	Type B applied part
	MEDICAL EQUIPMENT WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA-C22.2 No. 601.1, ANSI/AAMI ES60601-1 (2005, 3rd ed.), CAN/CSA-C22.2 No. 60601-1 (2008),13VA
	Consult instructions for use
	Sterilizable up to the temperature specified
	Do Not Autoclave (Handpiece)
	Do not re-use
	Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union
IPX0	Protection Class IPX0 IPX0 Classification of ingress of water for Charger - not protected
IPX1	Protection Class IPX1 IPX1 Classification of ingress of water
IPX3	Protection Class IPX3 IPX3 Classification of ingress of water for Handpiece - Protected against falling spray.
	Serial Number
	Batch Code/Lot Number
	This symbol is a mandatory marking for devices entering the European market to indicate conformity with the essential health and safety requirements set out in European Directives. The symbol may be accompanied by a four-digit identification number of the notified body.
	Direct Current (DC) supply 5 Volts, 1 Amp

	Interference may occur in the vicinity of equipment marked with the following symbol
	Date of Manufacture

8. DISPOSAL OF UNIT

U.S. - Dispose of the system components in accordance with state and local laws.

EU - Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union.

9. ELECTRO-MAGNETIC COMPATIBILITY PRECAUTION

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The digit Power® Dispenser is intended for use in the electromagnetic environment specified below. The Customer or the user of the digit Power® Dispenser should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The digit Power® Dispenser uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The digit Power® Dispenser is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the digit Power® Dispenser or shielding the location.
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The digit Power® Dispenser is intended for use in the electromagnetic environment specified below. The customer or the end user of the digit Power® Dispenser should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electromagnetic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, ±4, ±8, ±15kV air	±8kV contact ±2, ±4, ±8, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Enclosure Port IEC 61000-4-3	Tested as specified in Table 9 as per IEC 61000-4-3	Tested as specified in Table 9 as per IEC 61000-4-3	
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±0.5, ±1, ±2kV mains No I/O lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±0.5, ±1kV Line to Line ±0.5, ±1, ±2kV Line to Earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) f or 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the digit Power® Dispenser requires handpiece charging during power mains interruptions, it is recommended that the digit Power® Dispenser charging base be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The digit Power® Dispenser is intended for use in the electromagnetic environment specified below. The customer or the end user of the digit Power® Dispenser should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Medical Electric Equipment (digit Power® Dispenser), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = [3.5/V1] \sqrt{P}$ $d = [3.5/E1] \sqrt{P}$ 80 MHz to 800 MHz $d = [7/E1] \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the digit Power® Dispenser is used exceeds the applicable RF compliance level above, the digit Power® Dispenser should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the digit Power® Dispenser. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

RECOMMENDED SEPARATION DISTANCE BETWEEN PORTABLE & MOBILE RF COMMUNICATIONS EQUIPMENT & THE MODEL

The medical electric equipment (digit Power® Dispenser) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the digit Power® Dispenser can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical electric equipment (digit Power® Dispenser) as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

POWER INPUT AND FREQUENCIES

Input	100-240V~, 50-60 Hz, 0.3A
Output	Direct Current (DC) supply 5 Volts, 1 Amp

- This device complies with Part 15 of the FCC Rules and with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesirable operation.
- This Class B digital apparatus complies with Canadian ICES-003.
- Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.
- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - o Reorient or relocate the receiving antenna.
 - o Increase the separation between the equipment and receiver.
 - o Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

10. TROUBLESHOOTING

PROBLEM	SOLUTION
Optional Wireless Foot Control does not charge.	<ol style="list-style-type: none">1. Ensure that the power cord is securely attached to the wall outlet and the optional wireless foot control.2. Remove the optional wireless foot control access door and look in the compartment; ensure that the LED in the bottom compartment is illuminated. If not, call Technical Support.
System intermittently shuts down.	<ol style="list-style-type: none">1. Handpiece may shutdown in the presence of strong non-compliant Radio Frequency generating components (interference).2. Call Technical Support.
Handpiece does not have the correct speed control.	<ol style="list-style-type: none">1. Handpiece or optional wireless foot control may be malfunctioning.2. Call Technical Support.
Difficulty removing cartridge cap from handpiece.	<ol style="list-style-type: none">1. Check the components for gross debris.2. Hold the handpiece securely, press both locking tabs and pull straight off. Do not twist cartridge cap while removing.3. Ensure both locking tabs are depressed before removing cartridge cap.4. Inspect parts for wear.5. Call Technical Support.
Charging base for the handpiece does not drain liquids.	Clear any debris from the hole at the bottom of the charging base for handpiece.
Handpiece does not charge.	<ol style="list-style-type: none">1. Clean the charge contacts on the handpiece and charging base for handpiece using one of the approved cleaning solutions described in section 4. Do not spray cleaning solutions directly onto parts.2. Verify that the power supply is properly connected to the charging base for handpiece and that the green LED on back lights up.3. Verify that the handpiece is able to properly sit inside the charging base for handpiece, and that there are no foreign obstructions.4. If still not charging, call Technical Support.
Handpiece does not hold charge.	<ol style="list-style-type: none">1. Verify that the handpiece properly charges (Battery icon flashes and scrolls when in the charger).2. Call Technical Support.
Power supply overheating.	Immediately unplug the unit and call Technical Support.
Power cord is frayed or damaged in any way.	Do not use. Call Technical Support.

11. SERIAL NUMBER

The digit Power® Dispenser, charging base and foot control have different serial numbers located on their respective bodies. These serial numbers should be included in all correspondence.

12. WARRANTY

The Dentsply Sirona digit Power® Dispenser is designed exclusively for dental use and this warranty is not applicable to other uses. This warranty extends to digit Power® Dispenser system purchased from an authorized Dentsply Sirona distributor, and only to the original purchaser.

All components of the digit Power® Dispenser are warranted for (1) year from the date of purchase against defects arising from faulty materials and workmanship.

Parts will be repaired or replaced at the discretion of Dentsply Sirona provided that the system has been operated and maintained as prescribed in these instructions and has not been subjected to apparent misuse, abuse or accident. Claims covered by this warranty will be honored when presented through your Dentsply Sirona distributor within thirty (30) days from discovery of defect within the applicable warranty period.

THERE ARE NO WARRANTIES, EXPRESS OR IMPLIED, WHICH EXTEND BEYOND THE DESCRIPTION ON THE FACE HEREOF. Dentsply Sirona neither assumes, nor authorizes any person to assume for it, any other liability in connection with the sale or use of its products. DAMAGES ARE LIMITED STRICTLY TO REPAIR OR REPLACEMENT OF PARTS. DENTSPLY EXPRESSLY DISCLAIMS LIABILITY FOR INCIDENTAL AND CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OF THE PRODUCTS.

13. TECHNICAL SUPPORT

For technical support, please contact the numbers below or your local Dentsply Sirona Sales Office directly:

United States of America
Dentsply Sirona Preventive
Technical Service and Repair Department
1301 Smile Way
York, PA 17404-1785
Phone: (800) 989-8826 or (717) 767-8502

UK
Dentsply Sirona Service Centre
Brunel Way, Stroudwater Business Park
Stonehouse, Gloucestershire GL10 3GB
United Kingdom
Phone: +44 (0)1453 793 190

Germany
Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz
Germany
Phone: 0 8000 735000 (for Germany)
Phone: 00800-00735000 (for Austria)

Canada
Dentsply Sirona Canada
161 Vinyl Court
Woodbridge Ontario
L4L 4A3 Canada
Phone: (905) 851-6060

Australia
Dentsply Sirona Australia
11-21 Gilby Road
Mount Waverly, Victoria 3149
Australia
Phone: (61) 3-9538-8280

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©2017 Dentsply Sirona. All Rights Reserved.
Bluetooth® is not a registered Trademark of Dentsply Sirona.
Dentsply LLC