
Midwest Stylus® Plus / Plus SH Midwest Tradition® TC / Pro TCF

Air-Driven-High-Speed Handpiece for connection: /
Turbine pour raccord: / Turbina para conexión:
Midwest®, Dentsply Sirona, KaVo, W&H


English / Français / Español

Operating instructions
Notice d'utilisation
Manual del operador

Midwest Stylus[®] Plus / Plus SH Midwest Tradition[®] TC / Pro TCF

High-speed handpiece for Midwest[®], Dentsply Sirona, KaVo and W&H coupling

English

Operating instructions

Table of Contents

1	Before you begin	6
1.1	Structure of the document	8
1.2	Service life	11
1.3	Trade names / registered trade marks of the manufacturers.....	12
2	Safety instructions	13
3	Technical description	17
3.1	Task	17
3.2	Structure of Midwest Stylus Plus	18
3.3	Structure of Midwest Stylus Plus SH	19
3.4	Structure of Midwest Tradition Pro TCF	20
3.5	Structure of Midwest Tradition TC	21

	3.6	Technical Data	22
	3.7	Available high-speed handpiece connection types	27
	3.8	Product labeling	29
4		Preparation	31
	4.1	Initial start-up and longer breaks in use	31
	4.2	Prior to starting the work day	31
	4.3	Before each patient.....	32
5		Operation	33
	5.1	Connecting the Dentsply Sirona quick coupling to the supply hose	34
	5.2	Replacing the high-speed handpiece	35
	5.3	Inserting and removing burs	38
	5.4	Adjusting the cooling spray.....	39

6	Postprocessing	40
6.1	After each treatment session	40
6.2	At the end of the work day	41
7	Reprocessing	42
7.1	Automated cleaning	42
7.2	Manual cleaning	44
7.3	Manual maintenance	47
7.4	Sterilizing	54
8	Maintenance	57
8.1	Testing the FG clamping system	57
8.2	Servicing Dentsply Sirona quick coupling	59
8.3	Servicing Midwest quick coupling	64

9	Spare parts and consumables	66
10	Disposal	70

1 Before you begin ...

The Midwest Stylus Plus / Midwest Tradition high-speed handpiece complies with the provisions of state-of-the-art technology.

The Midwest Stylus Plus / Midwest Tradition high-speed handpiece complies with the ISO 14457, ISO 9168 standards.

The units, the product is connected to, must comply with the requirements of IEC 60601-1.

1. Read the operating instructions before using the Midwest Stylus Plus / Midwest Tradition high-speed handpiece.
2. Use the Midwest Stylus Plus / Midwest Tradition high-speed handpiece only for the purposes described in the operating instructions.
3. Follow the hygiene regulations, work safety regulations, and accident prevention measures that apply to the Midwest Stylus Plus / Midwest Tradition high-speed handpiece in your country.

Intended use

The high-speed handpieces of the Midwest Stylus Plus / Midwest Tradition high-speed handpiece family are intended for:

- Preparation of cavities and crowns
- Removal of carious material
- Removal of fillings
- Processing of tooth and restoration surfaces
- Reducing hard tooth substance

Contraindications

none

Target group

This product is intended only for use by trained dental personnel in dental practices and laboratories.

Scope of the operating instructions

These operating instructions are applicable for the following high-speed handpieces:

- Midwest Stylus[®] Plus
- Midwest Stylus[®] Plus SH

- Midwest Tradition[®] TC
- Midwest Tradition[®] Pro TCF

1.1 Structure of the document

1.1.1 Labeling of information

Warnings

- To prevent injuries, please observe the warnings.

Warnings are labeled as follows:

 **DANGER!** indicates a danger **leading** to death or serious injury if not avoided.

 **WARNING!** indicates a danger that **may lead** to death or serious injury if not avoided.

 **CAUTION!** indicates a danger that **may lead** to injury if not avoided.

Instructions for use

- To prevent material damage and additional expenses, please observe all instructions for use.

Instructions for use are labeled as follows:

NOTE! indicates measures for the prevention of material damage.

IMPORTANT: indicates information on the avoidance of additional expenses and other important information.

Tip: indicates information for facilitating work.

1.1.2 Formats and symbols

The formats and symbols used in this document have the following meaning:

<ul style="list-style-type: none"> ✓ Prerequisite 1. First action step 2. Second action step or ➤ Alternative action ↪ Result ➤ Individual action step 	Requests you to do something.
Use of formats and symbols [→ 10].	Identifies a reference to another text passage and specifies its page number.
<ul style="list-style-type: none"> • List 	Identifies a list.

1.2 Service life

When used as intended:

- Non-moving parts have a typical service life of 5 years
- Moving parts have a typical service life of 3 years

No warranty claim can be inferred here, as wear may occur earlier or later than indicated above depending on use, frequency of sterilization, and frequency of maintenance.

Elastomers, e.g. O-rings, must be replaced depending on their degree of wear.

Check the clamping system of friction grip burs on a monthly basis [→ 57].

1.3 Trade names / registered trade marks of the manufacturers

The designations used in this operating instruction with reference to the instruments are trade names or registered trade names of the manufacturers:

KaVo, Biberach, Germany:

- MULTIflex LUX

Obligations of the user

Preventing the spread of infections and cross contamination

Prevention of eye damage

2 Safety instructions

- Use only fault-free materials that do **not** deviate from the specified data [→ 22].
- Protect yourself, patients, and others from danger. To do this, observe the safety information.
- Comply with the Intended use of the equipment.
- You should always keep these operating instructions within reach for further reference.

Prevent the spread of infections and cross contamination between patients, users, and third parties. Sterilize equipment after each patient.

Take the appropriate hygiene measures, e.g. wear protective gloves.

The LED is in risk class 2 according to the IEC 62471:2006 standard. The LED emits optical radiation that is potentially

Malfunction or damage

hazardous and may be harmful to the eyes! Potential damage to the retina from the blue light emission. Do not stare at the LED for longer periods of time while in operation.

Discontinue use immediately in case of malfunction, unusual or different sounds or damage. Damaged high-speed handpieces may cause injury. Notify the dental supplies distributor or the manufacturer.

Potential explosion hazard

Do **not** use this product in areas subject to explosion hazards.

Repair

Do **not** attempt to repair the high-speed handpiece yourself.

The replacement of bearings can lead to sudden failures and injuries to patients due to the different levels of wear of the newer and older components. Safe operation is **not** guaranteed if proper repairs have not been carried out. For repairs, contact the dental depot or the manufacturer.

Spare and accessory parts

First aid measures in the case of lubricant accidents

Use only original Midwest parts or parts approved by Midwest. Safe operation is **not** guaranteed for parts that have not been approved by Midwest.

General information:	Immediately remove any clothing soiled by the product.
After inhalation:	Supply fresh air; consult doctor in case of complaints.
After skin contact:	If skin irritation continues, consult a doctor.
After eye contact:	Rinse opened eye for several minutes under running water.
After swallowing:	If symptoms persist consult doctor.

For details download the Midwest Plus Spray Material Safety Data Sheet from the Midwest homepage: www.dentsplysirona.com

If you have any questions, please contact your dental supplies distributor or the manufacturer.


3 Technical description

3.1 Task

The Midwest Stylus Plus / Midwest Tradition high-speed handpiece is used to power dental instruments with an FG shank.


The quick coupling is used for media supply and connection to the supply hose and features convenient swivel.

3.2 Structure of Midwest Stylus Plus


A	High-speed handpiece head
B	Push button
C	Cooling spray outlet (4 spray nozzles)
D	Opening of the chuck system
E	Light aperture
F	Handpiece sleeve

3.3 Structure of Midwest Stylus Plus SH


A	High-speed handpiece head, small
B	Push button
C	Opening of the chuck system
D	Cooling spray outlet (1 spray nozzle)
E	Light aperture
F	Handpiece sleeve

3.4 Structure of Midwest Tradition Pro TCF

A	High-speed handpiece head
B	Push button
C	Cooling spray outlet (4 spray nozzles)
D	Opening of the chuck system
E	Light aperture
F	Handpiece sleeve

3.5 Structure of Midwest Tradition TC

A	High-speed handpiece head
B	Push button
C	Cooling spray outlet (4 spray nozzles)
D	Opening of the chuck system
E	Handpiece sleeve

3.6 Technical Data

Midwest Stylus Plus / Midwest Tradition high-speed handpiece

	Midwest Stylus Plus/ Plus SH / Tradition Pro TCF	Midwest Tradition TC
Light function	x	-
Back-suction stop in the head housing	x	x
Anti-retraction valve in the spray water area	x	x
Push-button cover clamping system	FG	FG

High-speed handpiece head

	Midwest Stylus Plus / Midwest Tradition
No-load speed in rpm	410 000 ± 40 000
Maximum torque in Nmm	~ 1.6

Friction grip bur instrument

	Midwest Stylus Plus / Midwest Tradition
Shank standard	ISO 1797-1, Type 3
Shank diameter in mm	1.59 - 1.60

	Midwest Stylus Plus / Midwest Tradition
Maximum total length in mm	21
Maximum working diameter in mm (ISO 2157)	2.1

Dental unit

	Midwest Stylus Plus / Midwest Tradition
Driving air pressure in bar (psi) (flowing, dry, and clean)	3.0 (43.5)
Driving air consumption in NI/min	50 ± 5
Return air pressure in bar (psi)	< 0.3 (< 4.35)
Spray air pressure in bar (psi)	~ 2.5 (36.3)

	Midwest Stylus Plus / Midwest Tradition
Spray air flow in l/min	> 2.5
Spray water pressure in bar (psi)	~ 2.0 (29)
Water content in spray ml/min	> 50

The indicated pressures can be checked by a service technician with a star tester that is placed between the supply hose and the quick coupling.

Quick couplings (Dentsply Sirona)

	6-Pin LED Coupler S-Style	5-Hole Coupler S-Style	4-Hole Coupler S-Style
Connection	6-pin	5-hole	4-hole
Standard	ISO 9168	ISO 9168	ISO 9168
Driving air return	to the hose	to the hose	to the hose
Spray water control	on the coupling	-	-
Back suction stop	in the spray water area	in the spray water area	in the spray water area
Light function	x	x	-
Lamp voltage (lamp with dark brown socket)	3.6 V (optional)	- (Lamp in the hose)	-
Lamp voltage (LED)	3.6V	- (Lamp in the hose)	-

3.7 Available high-speed handpiece connection types

	Coupling type available for	Midwest Stylus Plus	Midwest Stylus Plus SH	Midwest Tradition Pro TCF	Midwest Tradition TC
	Midwest quick coupling	-	x	-	-
	6-pin LED S style quick coupling (Dentsply Sirona)	x	-	-	-
	5-pin S style quick coupling (Dentsply Sirona)	x	-	-	-

	Coupling type available for	Midwest Stylus Plus	Midwest Stylus Plus SH	Midwest Tradition Pro TCF	Midwest Tradition TC
	4-pin S style quick coupling (Dentsply Sirona)	x	-	-	x
	Multiflex LUX quick coupling (KaVo)	x	x	x	-
	Roto Quick quick coupling (W&H)	x	-	-	-

3.8 Product labeling

Sterilizable in a steam sterilizer (autoclave) at the temperature specified

Serial Number

Date of manufacture

Order number

Data matrix code (here: an example)

Content of the data matrix code:

+E2768752050/\$\$+5192/16D20170309E

A

B

C

D

A Manufacturer code (here: E276)

B REF number (here: 875205)

C Serial no. (here: 5192)

D Manufacture date (YYYYMMDD)

4 Preparation

4.1 Initial start-up and longer breaks in use

- ✓ The high-speed handpiece can reach full power only if specified operating pressures are set [→ 22]. If necessary, have the pressures checked by a service technician with a star tester that is placed between the supply hose and the quick coupling.
- Sterilize the high-speed handpiece and accessory parts prior to startup [→ 54].
NOTE! Never sterilize the quick coupling.
- Clean and lubricate the high-speed handpiece after longer breaks in use.

4.2 Prior to starting the work day

- Purge the water and air channels for 30 seconds.

4.3 Before each patient

1. Purge the water and air channels for 30 seconds.
2. Attach the high-speed handpiece.
3. Insert the bur.
4. Set a sufficient amount of cooling water (> 50 ml/min) [→ 39].
5. Use filtered water only (< 50 µm).
6. Check the nozzles for blockages and lime deposits, for example, and clean the nozzles if necessary.

⚠ CAUTION! Insufficient cooling leads to overheating of the preparation site and to damage to the tooth substance. Ensure that the water content is ≥ 50 ml/min.

5 Operation

⚠ CAUTION! Insufficient cooling leads to overheating of the preparation site and to damage to the tooth substance. Ensure that the water content is ≥ 50 ml/min.


⚠ CAUTION! Never pull the patient's cheek back with the handpiece! This would actuate the pushbutton, thus creating a risk of burning the patient's oral mucosa.

⚠ CAUTION! Potentially hazardous optical radiation may cause harm to the eyes. Do **not** stare at the LED for longer periods of time while in operation.

IMPORTANT: For application details and operating data of the different preparation tools, please refer to the information provided by the manufacturer.

5.1 Connecting the Dentsply Sirona quick coupling to the supply hose

1. Attach the quick coupling (B) to the supply hose, noting the various tube diameters.
2. Screw the cap nut (C) tight. Secure the quick coupling using the key (A).

5.2 Replacing the high-speed handpiece

⚠ CAUTION! Do **not** remove the high-speed handpiece while the turbine is running.


5.2.1 High-speed handpiece with Midwest connection

Attaching the high-speed handpiece

- Insert the high-speed handpiece until it audibly clicks into place.

Removing the high-speed handpiece

- ✓ The high-speed handpiece is at a standstill.
- Press the quick coupling release buttons and hold the hose nut firmly and pull the high-speed handpiece while turning it slightly. Do **not** pull on or hold the supply hose.


5.2.2 High-speed handpiece with Dentsply Sirona or KaVo connection

Attaching the high-speed handpiece

- Insert the high-speed handpiece until it audibly clicks into place.

Removing the high-speed handpiece

- ✓ The high-speed handpiece is at a standstill.
- Hold the quick coupling and the hose nut firmly and pull the high-speed handpiece while turning it slightly. Do **not** pull the supply hose or hold it firmly.


5.2.3 Turbine with W&H connection

Attaching the high-speed handpiece

- Insert the high-speed handpiece until it audibly clicks into place.

Removing the high-speed handpiece

- ✓ The high-speed handpiece is at a standstill.
- Press the markings (A) on the W&H quick coupling and pull the high-speed handpiece while turning it slightly. Do **not** pull or hold onto the supply hose.


5.3 Inserting and removing burs

IMPORTANT: Check the push button to make sure it moves freely!

NOTE! Use only burs and diamond polishers that are sharp and undamaged. Use clean burs and diamond polishers to avoid dirt in the clamping system.

⚠ CAUTION! A loose or partially removed bur can detach itself from the head or break off. This may cause injury! Therefore, use the high-speed handpiece only if the bur is inserted at least 10 mm and is clamped securely in place.

5.4 Adjusting the cooling spray

- Adjust the flow rate of the cooling water using the control ring (A) (> 50 ml/min).

Tip: You can measure the amount of cooling water with a measuring cup and watch.

Water flow rate

Quick coupling F (Dentsply Sirona): The maximum water flow is set when the control ring clicks into place.

Multiflex LUX quick coupling (KaVo): The maximum water flow is set when the two markings face each other.

Roto Quick (W&H) quick coupling: The minimum water flow is set when the red line and the red dot are facing each other.

Turn the control ring to the left or right for the maximum water flow.

6 Postprocessing

6.1 After each treatment session

NOTE! Condition immediately, or at the latest, one hour after treatment.

NOTE! Use only Midwest Plus Spray.

- ✓ The bur must not be moving.
- ✓ Wear appropriate protective clothing.
- 1. Purge the water and air channels on the treatment center for 30 seconds.
- 2. Remove the bur with tweezers.
- 3. Predisinfect the high-speed handpiece and the quick coupling directly at the treatment center.
- 4. Take the high-speed handpiece and unscrew the quick coupling.

5. Transport the high-speed handpiece and the quick coupling to the hygiene room in a suitable transport container.
6. Conduct automatic reprocessing. Manual reprocessing is possible in exceptional cases if the national/local regulations are followed.
NOTE! Do not conduct automatic reprocessing of the quick coupling.
7. Lubricate the high-speed handpiece with spray [→ 47].
8. Sterilize the high-speed handpiece and the accessories [→ 54].
NOTE! Do not sterilize the quick coupling.

6.2 At the end of the work day

- Lubricate the high-speed handpiece [→ 47].

7 Reprocessing

7.1 Automated cleaning

The cleaning equipment used must be approved by its manufacturer for the cleaning of dental instruments.

For further details, refer to the operating instructions supplied with the unit.

NOTE! Do **not** carry out automated cleaning of the quick couplings.

- ✓ The high-speed handpiece is conditioned with a cleaning device.
- 1. Check whether the high-speed handpiece is clean after conditioning under good lighting (min. 500 lux) and color rendering index (min. 80 Ra).
- 2. If they are still dirty, repeat the process.

☞ For further processing, the high-speed handpiece should be dry and free of residue.

3. Blow the high-speed handpiece out with max. 3 bar.
4. Maintain mechanical parts manually [→ 47].
5. Maintain the push button chuck manually [→ 50].
6. Pack the high-speed handpiece in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
7. Conduct sterilization [→ 54].

7.2 Manual cleaning

IMPORTANT: Manual reprocessing is possible in exceptional cases if the national/local regulations are followed. The national/local regulations are to be checked before.

NOTE! Condition immediately, or at the latest, one hour after treatment.

NOTE! Never clean in an ultrasonic bath!


NOTE! Never immerse in disinfectants!

NOTE! Use only Midwest Plus Spray.

IMPORTANT: Use a soft, clean, and disinfected brush for cleaning.

IMPORTANT: The Midwest Handpiece Cleaner is a surfactant.

✓ Wear appropriate protective clothing.

- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Use only disinfectants with **no** protein-fixing properties.
- 1. Spray a sufficient amount of Midwest Plus Handpiece Cleaner onto a clean, lint-free cloth and use it to clean the entire external surface of the high-speed handpiece.
- 2. Brush the high-speed handpiece/quick coupling under running water (< 38 °C, < 100 °F, at least drinking water quality) for at least 10 seconds until no more dirt can be seen in good lighting (min. 500 lux) and color rendering index (min. 80 Ra).
- 3. Flush the drive channels of the high-speed handpiece with spray.
- 4. Lubricate mechanical parts manually [→ 47].
- 5. Maintain the push button chuck manually [→ 50].
- 6. Pack the high-speed handpiece in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.


- 7.** Sterilize the high-speed handpiece [→ 54].
- 8.** Disinfect the quick coupling with disinfectant wipes. Observe the manufacturer's instructions for using the instrument disinfectant.
- 9.** Wipe the disinfectant off with a clean, lint-free cloth.
- 10.** Purge the individual connections of the quick coupling with compressed air for at least 10 seconds.

7.3 Manual maintenance

7.3.1 Maintenance of mechanical parts

Intervals


- After each patient
- Prior to each sterilization
- During a longer preparation (> 10 minutes)

Required accessories

- A Spray adapter for Dentsply Sirona connection
- B Spray adaptor for KaVo connection
- C Spray adapter for W&H connection
- D Spray adapter for Midwest connection
- E Midwest Plus Spray

NOTE! Use only Midwest Plus Spray.

Procedure

- ✓ The spray adapter is disinfected.
 - ✓ The O-rings on the spray adapter are intact.
1. Fit the spray adapter onto the nozzle of the spray can.
 2. Insert the high-speed handpiece as far as it will go and hold it in place.
 3. Spray into the high-speed handpiece for 1-2 seconds.
IMPORTANT: Hold the spray can upright.
 4. Wipe away any spray that runs out with a clean, lint-free cloth.
 5. Repeat the process until the spray that runs out is clear.

7.3.2 Care of the push button chuck


Use Midwest Plus Spray on the push-button chuck to remove deposits and ensure proper functioning of the clamping system.

Interval


- At least once a week

Required accessories

A	Black nozzle insert
B	Midwest Plus Spray

NOTE! Use only Midwest Plus Spray.

Procedure


- ✓ The high-speed handpiece is cleaned and disinfected.
- 1. Press the high-speed handpiece head with the chuck firmly against the spray can nozzle.
- 2. Spray the chuck for 1 - 2 seconds.
IMPORTANT: Hold the spray can upright.
- 3. Wipe away any spray that runs out with a clean, lint-free cloth.

7.3.3 Care of cooling spray nozzle openings

If your tap water is very hard, lime deposits may constrict or completely block the cooling spray nozzles.

Interval


- At least once a week
1. Carefully clean the nozzle openings by running a cleaning wire through them.
 2. Let the high-speed handpiece run briefly with cooling spray.

7.3.4 Cleaning the light guide surfaces

Check the fiber optic surface and clean it if needed.

1. In order to avoid scratching the surfaces (A), blow off any dirt particles with air using a dry syringe.
2. Wipe the surfaces with a Q-tip or a soft cloth and alcohol.



7.4 Sterilizing

Intervals:

- Prior to initial operation
- Prior to every other use

Procedure

- ✓ The high-speed handpiece is cleaned.
- ✓ High-speed handpiece maintenance is complete.
- ✓ If necessary, the high-speed handpiece can be sterilized in packaging suitable for sterilization and storage, e.g. paper/plastic composite packaging or container.
- Sterilize the high-speed handpiece in the steam sterilizer with saturated water vapor.

Temperature: 132 °C (270 °F) / 135 °C (275 °F)

Overpressure: 1.87 bar (27.12 psi) / 2.13 bar (30.89 psi)

Wrapping	Steam sterilization	Holding time	Drying time
Wrapped	Gravitation process	15 minutes at 132 °C	30 minutes
Wrapped	Gravitation process	10 minutes at 135 °C	30 minutes
Wrapped	Pre-vacuum process	4 minutes at 132 °C	30 minutes
Wrapped	Pre-vacuum process	3 minutes at 135 °C	16 minutes

NOTE! Do not exceed 140 °C (284 °F), even during the drying phase.

After sterilizing

1. Remove the high-speed handpiece from the steam sterilizer immediately.
⚠ CAUTION! The high-speed handpiece will be hot. Risk of burns!
NOTE! Do **not** attempt to accelerate the cooling process by immersing the high-speed handpiece in cold water. This will damage your high-speed handpiece!
2. Store all high-speed handpieces so that they are protected from contamination.
3. Sterilize again once the storage period has elapsed.

8 Maintenance

Elastomers, e.g. O-rings, must be replaced depending on their degree of wear.

Check the clamping system of friction grip burs on a monthly basis [→ 57].


8.1 Testing the FG clamping system

Interval

Test the FG clamping system at least once a month.

Required accessories

Chuck tester

Process

- ✓ The expiry date of the chuck tester (A) has **not** elapsed.
- 1. Insert the chuck tester into the FG clamping system.
- 2. Tighten the chuck tester until the marking ring appears (withdrawal force: 22 N)

Does the chuck tester slide out of the chuck before the marking ring appears?


⚠ CAUTION! The chuck is defective and the secure fit of the bur is not guaranteed. This may cause injury!

1. Do not use the product.
2. Have the clamping system replaced by Midwest.

Tip: Record the time and result of the check for your own information.

8.2 Servicing Dentsply Sirona quick coupling

8.2.1 Replacing O-rings


1. Remove the defective O-ring.
2. Slide the tool (A), with the new O-ring attached, to just in front of the corresponding groove (B).
3. Insert the O-ring. Do **not** use any sharp objects.
4. Lightly oil the O-ring with Midwest Plus Spray.
NOTE! Do **not** use Vaseline or silicone grease on the O-rings.

8.2.2 Replacing the compression ring


In order to guarantee proper functioning of the “anti-retraction valve in the spray water duct”, we recommend replacing the compression ring (C) every six months.

1. Purge the spray water duct at the maximum water flow.
2. Slide the tool, with the new compression ring attached, to just in front of the corresponding groove (C).
3. Insert the compression ring. Do **not** use any sharp objects.
4. Check that the ring fits snugly all the way around and is **not** twisted.

8.2.3 Replacing the halogen lamp/LED

⚠ CAUTION! The lamp can be hot. Risk of burns! Allow the lamp to cool down.

1. Remove the high-speed handpiece.
2. Unscrew the cap (B).
3. Remove the defective halogen lamp/LED (A).
4. Insert the new halogen lamp/LED. Ensure the proper position of the contact surfaces.

NOTE! Pressing on the lens can destroy the LED. Use the assembly tool (C) to insert the LED.

5. Wipe the bulb of the halogen lamp with a clean cloth.
6. Unscrew the cap from the quick coupling.
7. Check the light function of the quick coupling.

If the LED does not light up:

- Remove the LED and re-insert it after rotating it 180° around its own axis.


8.2.4 Replacing the spray water cartridge


Required accessories

- Spray water cartridge
- Midwest Plus Spray

If the spray water cartridge leaks or is clogged.

1. Remove the quick coupling from the supply hose.
2. Remove the sealing washer (B).
3. Click the control ring into the maximum water flow. The spray water cartridge (A) can be replaced **only** in this position.
4. Remove the spray water cartridge.
5. Spray oil lightly onto the new spray water cartridge.
6. Insert the new spray water cartridge.

7. Note the position of the bulging side (C) when fitting the sealing washer.

8.2.5 Replacing the felt strip (6-Pin LED Coupler S-Style)

Replace the felt strip if the oil-laden return air that escapes is soiled. We recommend replacing the felt strip at least every three months.

1. Click the control ring into the maximum water flow.
2. Place the coupling upright on a firm, flat surface supported by the screw-on cap and firmly push the control ring (C) downward until it disengages.
NOTE! Note that the balls (A) and springs (B) may drop out.
3. Replace the contaminated felt strip (D).
4. Click the control ring back into position, taking care to keep the balls and springs in place.

8.3 Servicing Midwest quick coupling

8.3.1 Replacing the halogen lamp of the 6-pin coupling


⚠ CAUTION! The lamp can be hot. Risk of burns! Allow the lamp to cool down.

1. Remove the high-speed handpiece.
2. Insert blunt pick into recessed area on lamp base and slide lamp out of socket.
3. Remove the defective lamp (A).

IMPORTANT: Do not touch replacement lamp with bare fingers. Oils from the skin decrease lamp life. Use gloves when handling lamp.

4. Insert carefully the new lamp with contact pin end first into lamp socket.

Be sure to orient recessed area on lamp base to face out slot on coupling. If it is not possible to insert lamp, rotate lamp 180 degrees and try again.

5. Check the light function of the quick coupling.

9 Spare parts and consumables

Use only original Midwest parts or parts approved by Midwest.

	REF		REF
Midwest Plus Spray (6 cans, 250 ml each)	38 00 80M	Halogen lamp (dark brown base)	59 40 291
Chuck tester	33 27 793	Cleaning wire for spray nozzles	24 00 232
Rotor replacement package Stylus Plus / Tradition	67 27 346	Midwest Plus Handpiece Cleaner	38 01 40

High-speed handpiece with Dentsply Sirona connection

	REF		REF
Spray adapter	59 41 802	O-ring for spray adapter	18 91 840
Replacement kit for quick coupling: <ul style="list-style-type: none"> • 10 O-rings • 2 compressed rings • 1 tool for O-rings and compressed rings 	41 75 803	Tool for O-rings and compressed rings	41 74 343
		Screw cap for 5-hole/6-pin quick coupling	66 78 531
		Wrench for quick coupling	59 41 794
		LED for 6-pin quick coupling	63 14 558
Spring for 6-pin quick coupling	89 17 262	Control ring for 6-pin quick coupling	41 76 603
Ball Ø 2.5 mm for 6-pin quick coupling	34 20 742	Felt strip for 6-pin quick coupling	41 76 389

	REF		REF
Spray water cartridge for 6-pin quick coupling	77 43 444	Sealing washer for 6-pin quick coupling	89 16 629
Sealing washer for 5-pin quick coupling	41 79 946	Sealing washer for 4-pin quick coupling	58 61 765

High-speed handpiece with Midwest connection

	REF		REF
Spray adapter (Midwest)	38 00 85	Halogen lamp for 6-pin coupling, 2 pieces	79 02 57
Gasket for 6-pin coupling	79 02 58	O-ring for 6-pin coupling, 2 sets	79 02 56
O-ring for maintenance coupling, 2 sets	79 02 59		

High-speed handpiece with KaVo connection

	REF		REF
Spray adapter (KaVo)	54 56 954	O-rings (small) for spray adapter	18 90 842
O-ring (large) for spray adapter	18 91 444		

High-speed handpiece with W&H connection

	REF		REF
Spray adapter (W&H)	54 56 970	O-ring (small) for spray adapter	41 83 849
O-ring (large) for spray adapter	70 23 542	Tool for replacing lamps	54 57 515
Lamp holder	54 56 871		

10 Disposal

- According to current information, the product does not contain any substances that are hazardous to the environment.
- Disinfect the product prior to disposal.
- Observe the applicable disposal regulations for your area.

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Sirona Dental Systems GmbH