

Marginal bone maintenance with Astra Tech Implant System®

The Astra Tech Implant System is designed and proven to maintain marginal bone. A meta-analysis¹ concluded that the Astra Tech Implant System shows superior bone maintenance compared to the standard norm*. Another meta-analysis² shows that implants with the OsseoSpeed surface maintain the marginal bone better than competitor surfaces.

The features of the EV connection with the Conical Seal Design (internal conical implant-abutment connection) and MicroThread (minute thread on the implant neck) reduce peak stresses and transfer the load deeper into the bone, leading to positive biomechanical bone stimulation that maintains marginal bone levels long term.

Transparent systematic literature search

A systematic literature search of Astra Tech Implant System publications was conducted, applying the following inclusion criteria:

- Prospective studies with OsseoSpeed surface
- Bone level changes read from periapical intraoral radiographs
- ≥ 10 patients, followed for ≥ 1 year


In total, 125 publications—including 43 randomized studies, 27 long-term studies and 74 studies with radiographic baseline at implant placement—fulfilled the inclusion criteria, resulting in data from more than 5,600 patients with 11,500 implants. Any surgical approach was allowed; however, block augmentation or larger grafting were excluded from the analysis.

Overall, mean marginal bone level changes are small for the Astra Tech Implant System, averaging -0.3 mm after up to 10 years of follow-up.

Conclusion

The scientific literature presenting bone level data on the Astra Tech Implant System is extensive and shows remarkably well-maintained marginal bone levels in both short and long-term perspective. Preserved bone levels is a key prerequisite for the predictability and longevity of the implant treatment and quality of life for the patient.

* Standard norm: 1 mm bone loss during the first year, and an annual bone loss of less than 0.2 mm thereafter, Albrektsson 1986³



First author of publication	Mean MBLC ^a mm	Radiographic baseline ^b (months from IP)	Follow-up period (years)	No. of patients ^c	Restorations ^d	Implant survival rate (%) ^e
Five years or more follow-up						
Windael et al., 2018 ⁴	-0.49	0	10	25	F	100
Windael et al., 2020 ⁵	-0.97	0	10	121	F, S	86
Raes et al., 2018 ⁶	-0.49 /+0.98	0	8-10	23 /16	S	100 /93.8
Trabovic et al., 2019 ⁷	-1.29	0	8-11	25	F	81.9
Berberi et al., 2014 ⁸	-0.19 /-0.21	0	5	36	S	100 /90.9
Boven et al., 2017 ⁹	-0.23	3	5	25	OD	97
Boven et al., 2015 ¹⁰	-0.37 [#]	3	5	26	OD	99.1
Cooper et al., 2014 ¹¹	+0.10 /+2.06	1 to 3	5	58 /55	S	98.3 /94.6
Cooper et al., 2014 ¹²	-0.18 ^s	0	5	21	S	95.6
Donati et al., 2015 ¹³	-0.28 [#]	3	5	149	S	97.5
Doornewaard et al., 2020 ¹⁴	-0.36 [#]	0	5	26	OD	100
Galindo-Moreno et al., 2017 ¹⁵	-0.15	0	5	69	S	95.9
Guljé et al., 2019 ¹⁶	-0.14	0.5	5	21	S	100
Guljé et al., 2021 ¹⁷	-0.06 [#]	1.5	5	95	F	97.5
Guljé et al., 2019 ¹⁸	-0.13 [#]	4	5	38	S	97.5
Hosseini et al., 2019 ¹⁹	-0.12 [#]	6	5	19	S	100
Ioannidis et al., 2019 ²⁰	-0.18	0.5	5	29	F, S	96.1
Laass et al., 2019 ²¹	-0.09 [#]	No info	5	20	S	100
Lops et al., 2013 ²²	-0.30 [#]	5	5	85	S	100
Mertens et al., 2011 ²³	-0.09 ^{#s}	0 to 2	5	17	S, F	97
Matthys et al., 2019 ²⁴	-1.25 [#]	0	5	75	OD	100
Noelken et al., 2018 ²⁵	-0.70 ^{#s}	0	5	21	F, S	100
Salman et al., 2019 ²⁶	-0.55 [#]	0	5	30	OD	100
Schliephake et al., 2012 ²⁷	-0.08 [#]	0	5	44	F	100
Slot et al., 2016 ²⁸	-0.51 [#]	3	5	50	OD	99.6
Temmerman et al., 2019 ²⁹	-0.79	0	5	48	OD, F	96.5
Thoma et al., 2018 ³⁰	-0.40	0	5	101	S	99.3
Toljanic et al., 2016 ³¹	-0.44	0	5	51	F	93.2

Table text

^a MBLC, marginal bone level change in mm, from radiographic baseline to the end of the follow-up period (value for healed sites /extraction sockets).

^b IP = implant placement (0 is defined as <48 hours since surgery).

^c Number of patients at study start (patients having implants in healed sites /extraction sockets).

^d Restoration, S = single tooth; F = fixed restoration; OD = overdenture.

^e Cumulative implant survival rate (value for healed sites /extraction sockets).

[#] Value has been presented for different or different time periods, and a *mean change* value has been calculated; or a mean change have been interpreted from tables or graphs.

^s Value at implants placed in a mix of extraction sockets and healed sites.

First author of publication	Mean MBLC ^a mm	Radiographic baseline ^b (months from IP)	Follow-up period (years)	No. of patients ^c	Restorations ^d	Implant survival rate (%) ^e
Two to four years follow-up						
Güncü et al., 2016 ³²	-0.70	3 to 6	4	24	S	100
Arora et al., 2017 ³³	+0.26 #s	0	3	30	S	100
Barewal et al., 2012 ³⁴	-0.22	0	3	40	S	97.5
Basler et al., 2018 ³⁵	-0.04	4	3	23	S	no info
Berberi et al., 2014 ³⁶	-0.27 s	0	3	20	S	100
Bernard et al., 2019 ³⁷	-0.60 #	3	3	60	F, OD	100
Borges et al., 2021 ³⁸	-0.55 #	0	3	29	F	98.5
Clelland et al., 2016 ³⁹	+0.20 #	5	3	18	F	98.8
Cooper et al., 2019 ⁴⁰	-0.12	3	3	45	S	100
De Bruyn et al., 2013 ⁴¹	-0.40 /+1.56	0	3	58 /55	S	98.3 /94.5
Doornewaard et al., 2019 ⁴²	-0.28 #	0	3	26	OD	100
Eisner et al., 2018 ⁴³	-0.24 #	No info	3	20	S	100
Galindo-Moreno et al., 2017 ⁴⁴	-0.33 s	0	3	59	S	100
Geckili et al., 2011 ⁴⁵	-0.88	0	3	52	OD	100
Guida et al., 2020 ⁴⁶	-0.04 #	3	3	30	F	100
Hadzik et al., 2018 ⁴⁷	-0.28 #	0	3	30	S	no info
Han et al., 2018 ⁴⁸	-0.04	0	3	45	S	95.8
Hosseini et al., 2013 ⁴⁹	-0.22	8	3	59	S	100
King et al., 2016 ⁵⁰	-0.23	0	3	38	S	96.8
Kraus et al., 2019 ⁵¹	0.0 #	No info	3	44	S	97.2
Kutan et al., 2015 ⁵²	-0.89 #	3	3	28	S	100
Liu et al., 2019 ⁵³	+0.23	2	3	40	F	100
Maiorana et al., 2015 ⁵⁴	-0.09	0	3	69	S	95.9
Palmer et al., 2012 ⁵⁵	-0.10 #	4	3	29	S	96.6
Pieri et al., 2014 ⁵⁶	-0.22	6	3	50	S	100
Pohl et al., 2017 ⁵⁷	-0.44	0	3	101	S	100
Rossi et al., 2020 ⁵⁸	-0.16 #	3	3	10	S	100
Sanz et al., 2014 ⁵⁹	+0.25 #s	0	3	93	S	98.9
Tabrizi et al., 2013 ⁶⁰	-0.75	3	3	33	S	100
Tabrizi et al., 2016 ⁶¹	-0.42 #	3	3	23	F	100
Temmerman et al., 2015 ⁶²	-0.18	4	3	28	S, F, OD	100
Thor et al., 2014 ⁶³	-0.57	0	3	51	F	95.7
Toija et al., 2021 ⁶⁴	-0.04 #	2	3	56	F	99.5
Zadeh et al., 2018 ⁶⁵	-0.19	0	3	95	F	97.6
Collaert et al., 2011 ⁶⁶	-0.11	0	2	25	F	100
Emami et al., 2016 ⁶⁷	-0.40 #	0	2	18	OD	91.7
Ettl et al., 2020 ⁶⁸	-0.81 #	No info	2	39	OD	92.3
Ferrari et al., 2015 ⁶⁹	-0.49 #	8	2	47	S	100
Göcken-Röhlig et al., 2010 ⁷⁰	-1.17 /-1.36	4	2	10	OD	100/100
Koutouzis et al., 2015 ⁷¹	-0.31 #	0	2	29	S	100
Lee et al., 2016 ⁷²	+0.06	0	2	14	S	100
Mumcu et al., 2012 ⁷³	-0.43 #	6	2	48	OD	100
Noelken et al., 2020 ⁷⁴	+7.75 #s	0	2.5	50	S	98
Noelken et al., 2014 ⁷⁵	+0.89 #s	0	2	20	S, F	100
Pieri et al., 2012 ⁷⁶	-0.60	6	2	25	F	96.8
Raes et al., 2015 ⁷⁷	+0.06 #	2	2	85	S	96.5
Schegnitz et al., 2017 ⁷⁸	-0.30	0	2	184	S	99.2
Toia et al., 2019 ⁷⁹	+0.02	3	2	40	OD	100
Vervaeke et al., 2013 ⁸⁰	-0.11	0	2	25	F	100
Vervaeke et al., 2018 ⁸¹	-0.37 #	3 to 4	2	25	OD	100
Younes et al., 2019 ⁸²	-0.34	0	2	22	S	100

First author of publication	Mean MBLC ^a mm	Radiographic baseline ^b (months from IP)	Follow-up period (years)	No. of patients ^c	Restorations ^d	Implant survival rate (%) ^e
One year follow-up						
Barbier et al., 2019 ⁸³	-0.32	No info	1	10	F	90
Barbier et al., 2012 ⁸⁴	-0.21 s	6	1	20	F	100
Bashutski et al., 2013 ⁸⁵	-0.50 "#s	0	1	24	S	92
Borges et al., 2018 ⁸⁶	-0.71	0	1	33	F	98.5
Cooper et al., 2010 ⁸⁷	-0.40 /+1.30	1 to 3	1	60 /55	S	98.3 /94.5
Cooper et al., 2015 ⁸⁸	-0.22	0	1	48	S	100
D'Haese et al., 2013 ⁸⁹	-0.47	0	1	26	F	88.6
Donati et al., 2008 ⁹⁰	-0.32 #	0	1	151	S	97.4
Ebler et al., 2016 ⁹¹	-0.58 #	0	1	33	F, S	100
Galindo-Moreno et al., 2012 ⁹²	-0.07	0	1	69	S	95.9
Ghoveizi et al., 2013 ⁹³	-0.24 #	no info	1	10	S	no info
Guljé et al., 2011 ⁹⁴	-0.10	0	1	12	OD	96
Guljé et al., 2013 ⁹⁵	-0.13 #	0	1	95	F	98.6
Guljé et al., 2015 ⁹⁶	-0.14	4	1	21	S	100
Guljé et al., 2014 ⁹⁷	-0.10	4	1	41	S	100
Guljé et al., 2016 ⁹⁸	-0.13	< 1	1	37	S	100
Han et al., 2016 ⁹⁹	-0.13	0	1	45	S	95.8
Hosseini et al., 2011 ¹⁰⁰	-0.09 #	7	1	36	S	100
Ivanoff et al., 2021 ¹⁰¹	-0.03 #	3	1	92	No info	99.4
Kim et al., 2010 ¹⁰²	-0.06 #	3 to 6	1	12	F	100
Koutouzis et al., 2011 ¹⁰³	-0.19	0	1	18	S	95
Lee et al., 2014 ¹⁰⁴	-0.09 #	3 to 6	1	76	S	100
Lyngstadaas et al., 2015 ¹⁰⁵	-0.65 #	0	1	37	F	98.6
Malmström et al., 2016 ¹⁰⁶	-0.13 #	3	1	30	S, F	98.8
Marcelis et al., 2012 ¹⁰⁷	-0.13	0	1	29	S	no info
Noelken et al., 2014 ¹⁰⁸	-0.54	0	1	65	S	100
Norton 2017 ¹⁰⁹	-0.13 s	3	1	22	S	100
Pieri et al., 2012 ¹¹⁰	-0.37 #	0	1	40	F	98.7
Piero et al., 2012 ¹¹¹	-0.36	4	1	15	S, F	100
Raes et al., 2013 ¹¹²	-0.18 /+0.66	0	1	23 /25	S	100 /96
Raes et al., 2011 ¹¹³	-0.25 /+0.80 #	0	1	23 /16	S	100 /93.8
Rismanchian et al., 2011 ¹¹⁴	-0.48	0	1	10	F	100
Schepke et al., 2017 ¹¹⁵	+0.09 #	3	1	50	S	100
Schickaglia et al., 2015 ¹¹⁶	-0.30 #	0	1	101	S	100
Schickaglia et al., 2016 ¹¹⁷	-0.40 #	0	1	31	OD	96.5
Simmons et al., 2017 ¹¹⁸	-0.21 #	0	1	27	S	93.3
Slot et al., 2013 ¹¹⁹	-0.25 #	3	1	50	OD	99.7
Slot et al., 2014 ¹²⁰	-0.22	3	1	25	OD	98
Stanford et al., 2016 ¹²¹	-0.07 #	0	1	120	S, F	96.9
Temmerman et al., 2017 ¹²²	-0.66	0	1	48	F	99.3
Toia et al., 2019 ¹²³	-0.17	0	1	50	F	100
Toljanic et al., 2009 ¹²⁴	-0.50	0	1	51	F	96
Tu et al., 2021 ¹²⁵	-0.50 #	0	1	19	No info	100
Tsuda et al., 2011 ¹²⁶	+0.10 s	0	1	10	S	100
Vercruyssen et al., 2014 ¹²⁷	-0.47 #	0	1	72	F, OD	100
Zhou et al., 2016 ¹²⁸	-0.14	0	1	45	F	100

References

1. Laurell, L, Lundgren, D, Marginal bone level changes at dental implants after 5 years in function: a meta-analysis. *Clin Implant Dent Relat Res* 2011;13(1):19-28. [Abstract](#)
2. Norton, MR, Astrom, M, The Influence of Implant Surface on Maintenance of Marginal Bone Levels for Three Premium Implant Brands: A Systematic Review and Meta-analysis. *Int J Oral Maxillofac Implants* 2020;35(6):1099-111. [Abstract](#)
3. Albrektsson, T, Zarb, G, Worthington, P, et al., The long-term efficacy of currently used dental implants: a review and proposed criteria of success. *Int J Oral Maxillofac Implants* 1986;1(1):11-25. [Abstract](#)
4. Windael, S, Vervaeke, S, Wijnen, L, et al., Ten-year follow-up of dental implants used for immediate loading in the edentulous mandible: A prospective clinical study. *Clin Implant Dent Relat Res* 2018;20(4):515-21. [Abstract](#)
5. Windael, S, Vervaeke, S, De Buyser, S, et al., The Long-Term Effect of Smoking on 10 Years' Survival and Success of Dental Implants: A Prospective Analysis of 453 Implants in a Non-University Setting. *J Clin Med* 2020;9(4). [Abstract](#)
6. Raes, S, Cosyn, J, Noyelle, A, et al., Clinical Outcome After 8 to 10 Years of Immediately Restored Single Implants Placed in Extraction Sockets and Healed Ridges. *Int J Periodontics Restorative Dent* 2018;38(3):337-45. [Abstract](#)
7. Trbakovic, A, Toljanic, JA, Kumar, VV, et al., Eight to eleven-year follow-up of immediately loaded implants placed in edentulous maxillae with compromised bone volume and poor bone quality: A prospective cohort study. *Clin Implant Dent Relat Res* 2020;22(1):69-76. [Abstract](#)
8. Berberi, AN, Sabbagh, JM, Aboushelib, MN, et al., A 5-year comparison of marginal bone level following immediate loading of single-tooth implants placed in healed alveolar ridges and extraction sockets in the maxilla. *Front Physiol* 2014;5:29. [Abstract](#)
9. Boven, GC, Slot, JWA, Raghoobar, GM, et al., Maxillary implant-supported overdentures opposed by (partial) natural dentitions: a 5-year prospective case series study. *J Oral Rehabil* 2017;44(12):988-95. [Abstract](#)
10. Boven, GC, Meijer, HJ, Slot, W, et al., Does a large dehiscent implant surface at placement affect the 5-year treatment outcome? An assessment of implants placed to support a maxillary overdenture. *J Craniomaxillofac Surg* 2015;43(9):1758-62. [Abstract](#)
11. Cooper, LF, Reside, GJ, Raes, F, et al., Immediate provisionalization of dental implants placed in healed alveolar ridges and extraction sockets: A 5-year prospective evaluation. *Int J Oral Maxillofac Implants* 2014;29(3):709-17. [Abstract](#)
12. Cooper, LF, Reside, G, Raes, F, et al., Immediate provisionalization of dental implants in grafted alveolar ridges in the esthetic zone: a 5-year evaluation. *Int J Periodontics Restorative Dent* 2014;34(4):477-86. [Abstract](#)
13. Donati, M, La Scala, V, Di Raimondo, R, et al., Marginal bone preservation in single-tooth replacement: a 5-year prospective clinical multicenter study. *Clin Implant Dent Relat Res* 2015;17(3):425-34. [Abstract](#)
14. Doornewaard, R, Bruyn, H, Matthys, C, et al., The Long-Term Effect of Adapting the Vertical Position of Implants on Peri-Implant Health: A 5-Year Intra-Subject Comparison in the Edentulous Mandible Including Oral Health-Related Quality of Life. *J Clin Med* 2020;9(10). [Abstract](#)
15. Galindo-Moreno, P, Nilsson, P, King, P, et al., Clinical and radiographic evaluation of early loaded narrow-diameter implants: 5-year follow-up of a multicenter prospective clinical study. *Clin Oral Implants Res* 2017;28(12):1584-91. [Abstract](#)
16. Gulje, FL, Raghoobar, GM, Vissink, A, et al., Single crown restorations supported by 6-mm implants in the resorbed posterior mandible: A five-year prospective case series. *Clin Implant Dent Relat Res* 2019;21(5):1017-22. [Abstract](#)
17. Gulje, FL, Meijer, HJA, Abrahamsson, I, et al., Comparison of 6-mm and 11-mm dental implants in the posterior region supporting fixed dental prostheses: 5-year results of an open multicenter randomized controlled trial. *Clin Oral Implants Res* 2021;32(1):15-22. [Abstract](#)
18. Gulje, FL, Raghoobar, GM, Vissink, A, et al., Single crowns in the resorbed posterior maxilla supported by either 11-mm implants combined with sinus floor elevation or 6-mm implants: A 5-year randomised controlled trial. *Int J Oral Implantol (Berl)* 2019;12(3):315-26. [Abstract](#)
19. Hosseini, M, Worsaae, N, Gotfredsen, K, Tissue changes at implant sites in the anterior maxilla with and without connective tissue grafting: A five-year prospective study. *Clin Oral Implants Res* 2020;31(1):18-28. [Abstract](#)
20. Ioannidis, A, Heierle, L, Hammerle, CHF, et al., Prospective randomized controlled clinical study comparing two types of two-piece dental implants supporting fixed reconstructions-Results at 5 years of loading. *Clin Oral Implants Res* 2019;30(11):1126-33. [Abstract](#)
21. Laass, A, Sailer, I, Husler, J, et al., Randomized controlled clinical trial of all-ceramic single-tooth implant reconstructions using modified zirconia abutments: Results at 5 years after loading. *Int J Periodontics Restorative Dent* 2019;39(1):17-27. [Abstract](#)
22. Lops, D, Bressan, E, Chiapasco, M, et al., Zirconia and titanium implant abutments for single-tooth implant prostheses after 5 years of function in posterior regions. *Int J Oral Maxillofac Implants* 2013;28(1):281-7. [Abstract](#)
23. Mertens, C, Steveling, HG, Early and immediate loading of titanium implants with fluoride-modified surfaces: results of 5-year prospective study. *Clin Oral Implants Res* 2011;22(12):1354-60. [Abstract](#)
24. Matthys, C, Vervaeke, S, Besseler, J, et al., Five-year study of mandibular overdentures on stud abutments: Clinical outcome, patient satisfaction and prosthetic maintenance-Influence of bone resorption and implant position. *Clin Oral Implants Res* 2019;30(9):940-51. [Abstract](#)
25. Noelken, R, Moergel, M, Kunkel, M, et al., Immediate and flapless implant insertion and provisionalization using autogenous bone grafts in the esthetic zone: 5-year results. *Clin Oral Implants Res* 2018;29(3):320-27. [Abstract](#)
26. Salman, A, Thacker, S, Rubin, S, et al., Immediate versus delayed loading of mandibular implant-retained overdentures: A 40-month follow-up of a randomized clinical trial. *J Clin Periodontol* 2019;46(8):863-71. [Abstract](#)
27. Schliephake, H, Rodiger, M, Phillips, K, et al., Early loading of surface modified implants in the posterior mandible - 5 year results of an open prospective non-controlled study. *J Clin Periodontol* 2012;39 (2):188-95. [Abstract](#)
28. Slot, W, Raghoobar, GM, Cune, MS, et al., Maxillary overdentures supported by four or six implants in the anterior region: 5-year results from a randomized controlled trial. *J Clin Periodontol* 2016;43(12):1180-87. [Abstract](#)
29. Temmerman, A, Rasmussen, L, Kubler, A, et al., A Prospective, Controlled, Multicenter Study to Evaluate the Clinical Outcome of Implant Treatment in Women with Osteoporosis/Osteopenia: 5-Year Results. *J Dent Res* 2019;98(1):84-90. [Abstract](#)
30. Thoma, DS, Haas, R, Sporniak-Tutak, K, et al., Randomized controlled multicentre study comparing short dental implants (6 mm) versus longer dental implants (11-15 mm) in combination with sinus floor elevation procedures: 5-Year data. *J Clin Periodontol* 2018;45(12):1465-74. [Abstract](#)
31. Toljanic, JA, Ekstrand, K, Baer, RA, et al., Immediate loading of implants in the edentulous maxilla with a fixed provisional restoration without bone augmentation: A report on 5-year outcomes data obtained from a prospective clinical trial. *Int J Oral Maxillofac Implants* 2016;31(5):1164-70. [Abstract](#)
32. Guncu, MB, Cakan, U, Aktas, G, et al., Comparison of implant versus tooth-supported zirconia-based single crowns in a split-mouth design: a 4-year clinical follow-up study. *Clin Oral Investig* 2016;20(9):2467-73. [Abstract](#)
33. Arora, H, Khzam, N, Roberts, D, et al., Immediate implant placement and restoration in the anterior maxilla: Tissue dimensional changes after 2-5 year follow up. *Clin Implant Dent Relat Res* 2017;19(4):694-702. [Abstract](#)
34. Barewal, RM, Stanford, C, Weesner, TC, A randomized controlled clinical trial comparing the effects of three loading protocols on dental implant stability. *Int J Oral Maxillofac Implants* 2012;27(4):945-56. [Abstract](#)
35. Basler, T, Naenni, N, Schneider, D, et al., Randomized controlled clinical study assessing two membranes for guided bone regeneration of peri-implant bone defects: 3-year results. *Clin Oral Implants Res* 2018;29(5):499-507. [Abstract](#)
36. Berberi, AN, Noujeim, ZN, Kanj, WH, et al., Immediate placement and loading of maxillary single-tooth implants: a 3-year prospective study of marginal bone level. *J Contemp Dent Pract* 2014;15(2):202-8. [Abstract](#)
37. Bernard, L, Vercruyssen, M, Duyck, J, et al., A randomized controlled clinical trial comparing guided with nonguided implant placement: A 3-year follow-up of implant-centered outcomes. *J Prosthet Dent* 2019;121(6):904-10. [Abstract](#)
38. Borges, T, Montero, J, Leitao, B, et al., Periimplant bone changes in different abutment heights and insertion timing in posterior mandibular areas: Three-year results from a randomized prospective clinical trial. *Clin Oral Implants Res* 2021;32(2):203-11. [Abstract](#)
39. Clelland, N, Chaudhry, J, Rashid, RG, et al., Split-mouth comparison of splinted and nonsplinted prostheses on short implants: 3-year results. *Int J Oral Maxillofac Implants* 2016;31(5):1135-41. [Abstract](#)
40. Cooper, LF, Reside, G, Stanford, C, et al., Three-year prospective randomized comparative assessment of anterior maxillary single implants with different abutment interfaces. *Int J Oral Maxillofac Implants* 2019;34(1):150-58. [Abstract](#)
41. De Bruyn, H, Raes, F, Cooper, LF, et al., Three-years clinical outcome of immediate provisionalization of single Osseospeed® implants in extraction sockets and healed ridges. *Clin Oral Implants Res* 2013;24(2):217-23. [Abstract](#)
42. Doornewaard, R, Glibert, M, Matthys, C, et al., Improvement of Quality of Life with Implant-Supported Mandibular Overdentures and the Effect of Implant Type and Surgical Procedure on Bone and Soft Tissue Stability: A Three-Year Prospective Split-Mouth Trial. *J Clin Med* 2019;8(6). [Abstract](#)
43. Eisner, B, Naenni, N, Husler, J, et al., Three-year results of a randomized controlled clinical trial using submucosally veneered and unveneered zirconia abutments supporting all-ceramic single-implant crowns. *Int J Periodontics Restorative Dent* 2018;38(5):645-52. [Abstract](#)

44. Galindo-Moreno, P, Padial-Molina, M, Nilsson, P, et al., The influence of the distance between narrow implants and the adjacent teeth on marginal bone levels. *Clin Oral Implants Res* 2017;28(6):704-12. [Abstract](#)
45. Geckili, O, Bilhan, H, Mumcu, E, et al., Three-year radiologic follow-up of marginal bone loss around titanium dioxide grit-blasted dental implants with and without fluoride treatment. *Int J Oral Maxillofac Implants* 2011;26(2):319-24. [Abstract](#)
46. Guida, L, Annunziata, M, Esposito, U, et al., 6-mm-short and 11-mm-long implants compared in the full-arch rehabilitation of the edentulous mandible: A 3-year multicenter randomized controlled trial. *Clin Oral Implants Res* 2020;31(1):64-73. [Abstract](#)
47. Hadzik, J, Krawiec, M, Slawiecki, K, et al., The Influence of the Crown-Implant Ratio on the Crestal Bone Level and Implant Secondary Stability: 36-Month Clinical Study. *Biomed Res Int* 2018;2018:4246874. [Abstract](#)
48. Han, J, Tang, Z, Zhang, X, et al., A prospective, multi-center study assessing early loading with short implants in posterior regions. A 3-year post-loading follow-up study. *Clin Implant Dent Relat Res* 2018;20(1):34-42. [Abstract](#)
49. Hosseini, M, Worsaae, N, Schiodt, M, et al., A 3-year prospective study of implant-supported, single-tooth restorations of all-ceramic and metal-ceramic materials in patients with tooth agenesis. *Clin Oral Implants Res* 2013;24(10):1078-87. [Abstract](#)
50. King, P, Maiorana, C, Luthardt, RG, et al., Clinical and radiographic evaluation of a small-diameter dental implant used for the restoration of patients with permanent tooth agenesis (hypodontia) in the maxillary lateral incisor and mandibular incisor regions: a 36-month follow-up. *Int J Prosthodont* 2016;29(2):147-53. [Abstract](#)
51. Kraus, RD, Epprecht, A, Hammerle, CHF, et al., Cemented vs screw-retained zirconia-based single implant reconstructions: A 3-year prospective randomized controlled clinical trial. *Clin Implant Dent Relat Res* 2019;21(4):578-85. [Abstract](#)
52. Kutan, E, Bolukbasi, N, Yildirim-Ondur, E, et al., Clinical and radiographic evaluation of marginal bone changes around platform-switching implants placed in crestal or subcrestal positions: A randomized controlled clinical trial. *Clin Implant Dent Relat Res* 2015;17 Suppl 2:e364-75. [Abstract](#)
53. Liu, J, Huang, Q, Wang, X, et al., Early loading of splinted implants in posterior mandible: Three-year results of a prospective multicenter study. *Clin Oral Implants Res* 2019;30(10):1049-58. [Abstract](#)
54. Maiorana, C, King, P, Quasas, S, et al., Clinical and radiographic evaluation of early loaded narrow-diameter implants: 3 years follow-up. *Clin Oral Implants Res* 2015;26(1):77-82. [Abstract](#)
55. Palmer, RM, Howe, LC, Palmer, PJ, et al., A prospective clinical trial of single Astra Tech 4.0 or 5.0 diameter implants used to support two-unit cantilever bridges: results after 3 years. *Clin Oral Implants Res* 2012;23(1):35-40. [Abstract](#)
56. Pieri, F, Siroli, L, Forlivesi, C, et al., Clinical, esthetic, and radiographic evaluation of small-diameter (3.0-mm) implants supporting single crowns in the anterior region: a 3-year prospective study. *Int J Periodontics Restorative Dent* 2014;34(6):825-32. [Abstract](#)
57. Pohl, V, Thoma, DS, Sporniak-Tutak, K, et al., Short dental implants (6 mm) versus long dental implants (11-15 mm) in combination with sinus floor elevation procedures: 3-year results from a multicentre, randomized, controlled clinical trial. *J Clin Periodontol* 2017;44(4):438-45. [Abstract](#)
58. Rossi, AL, Palombo, D, Clivio, A, et al., Three Year Randomized Split Mouth Clinical Trial Comparing Marginal Bone Loss Around Implants with a Conical or Straight Collar Design Global Journal of Oral Science 2020;6:16-24.
59. Sanz, M, Cecchinato, D, Ferrus, J, et al., Implants placed in fresh extraction sockets in the maxilla: clinical and radiographic outcomes from a 3-year follow-up examination. *Clin Oral Implants Res* 2014;25(3):321-7. [Abstract](#)
60. Tabrizi, R, Pourdanesh, F, Zare, S, et al., Do angulated implants increase the amount of bone loss around implants in the anterior maxilla? *J Oral Maxillofac Surg* 2013;71(2):272-7. [Abstract](#)
61. Tabrizi, R, Arabion, H, Aliabadi, E, et al., Does increasing the number of short implants reduce marginal bone loss in the posterior mandible? A prospective study. *Br J Oral Maxillofac Surg* 2016;54(7):731-5. [Abstract](#)
62. Temmerman, A, Keestra, JA, Coucke, W, et al., The outcome of oral implants placed in bone with limited bucco-oral dimensions: a 3-year follow-up study. *J Clin Periodontol* 2015;42(3):311-8. [Abstract](#)
63. Thor, A, Ekstrand, K, Baer, RA, et al., Three-year follow-up of immediately loaded implants in the edentulous atrophic maxilla: A study in patients with poor bone quantity and quality. *Int J Oral Maxillofac Implants* 2014;29(3):642-9. [Abstract](#)
64. Toia, M, Stocchero, M, Corra, E, et al., Fixed full-arch maxillary prostheses supported by four versus six implants with a titanium CAD/CAM milled framework: 3-year multicentre RCT. *Clin Oral Implants Res* 2021;32(1):44-59. [Abstract](#)
65. Zadeh, HH, Gulje, F, Palmer, PJ, et al., Marginal bone level and survival of short and standard-length implants after 3 years: An Open Multi-Center Randomized Controlled Clinical Trial. *Clin Oral Implants Res* 2018;29(8):894-906. [Abstract](#)
66. Collaert, B, Wijnen, L, De Bruyn, H, et al., A 2-year prospective study on immediate loading with fluoride-modified implants in the edentulous mandible. *Clin Oral Implants Res* 2011;22(10):1111-6. [Abstract](#)
67. Emami, E, Cerutti-Kopplin, D, Menassa, M, et al., Does immediate loading affect clinical and patient-centered outcomes of mandibular 2-unsplinted-implant overdenture? A 2-year within-case analysis. *J Dent* 2016;50:30-6. [Abstract](#)
68. Ettil, T, Junold, N, Zeman, F, et al., Implant survival or implant success? Evaluation of implant-based prosthetic rehabilitation in head and neck cancer patients-a prospective observational study. *Clin Oral Investig* 2020;24(9):3039-47. [Abstract](#)
69. Ferrari, M, Cagidiaco, MC, Garcia-Godoy, F, et al., Effect of different prosthetic abutments on peri-implant soft tissue. A randomized controlled clinical trial. *Am J Dent* 2015;28(2):85-9. [Abstract](#)
70. Gökçen-Röhlig, B, Meric, U, Keskin, H, et al., Clinical and radiographic outcomes of implants immediately placed in fresh extraction sockets. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2010;109(4):1-7. [Abstract](#)
71. Koutouzis, T, Neiva, R, Lipton, D, et al., The effect of interimplant distance on peri-implant bone and soft tissue dimensional changes: a nonrandomized, prospective, 2-year follow-up study. *Int J Oral Maxillofac Implants* 2015;30(4):900-8. [Abstract](#)
72. Lee, PK, Siu, AS, et al., A two-year evaluation of a sloped marginal contour implant system placed in healed sites. *Int J Oral Maxillofac Implants* 2016;31(6):1423-28. [Abstract](#)
73. Mumcu, E, Bilhan, H, Geckili, O, et al., The influence of healing type on marginal bone levels of implants supporting mandibular overdentures: A randomized clinical study. *Indian J Dent Res* 2012;23(4):514-8. [Abstract](#)
74. Noelken, R, Pausch, T, Wagner, W, et al., Peri-implant defect grafting with autogenous bone or bone graft material in immediate implant placement in molar extraction sites-1- to 3-year results of a prospective randomized study. *Clin Oral Implants Res* 2020;31(11):1138-48. [Abstract](#)
75. Noelken, R, Neffe, BA, Kunkel, M, et al., Maintenance of marginal bone support and soft tissue esthetics at immediately provisionalized OsseoSpeed implants placed into extraction sites: 2-year results. *Clin Oral Implants Res* 2014;25(2):214-20. [Abstract](#)
76. Pieri, F, Aldini, NN, Fini, M, et al., Preliminary 2-year report on treatment outcomes for 6-mm-long implants in posterior atrophic mandibles. *Int J Prosthodont* 2012;25(3):279-89. [Abstract](#)
77. Raes, S, Rocci, A, Raes, F, et al., A prospective cohort study on the impact of smoking on soft tissue alterations around single implants. *Clin Oral Implants Res* 2015;26(9):1086-90. [Abstract](#)
78. Schiegnitz, E, Noelken, R, Moergel, M, et al., Survival and tissue maintenance of an implant with a sloped configurated shoulder in the posterior mandible-a prospective multicenter study. *Clin Oral Implants Res* 2017;28(6):721-26. [Abstract](#)
79. Toia, M, Wennerberg, A, Torrisi, P, et al., Patient satisfaction and clinical outcomes in implant-supported overdentures retained by milled bars: Two-year follow-up. *J Oral Rehabil* 2019;46(7):624-33. [Abstract](#)
80. Vervaeke, S, Collaert, B, De Bruyn, H, et al., The effect of implant surface modifications on survival and bone loss of immediately loaded implants in the edentulous mandible. *Int J Oral Maxillofac Implants* 2013;28(5):1352-7. [Abstract](#)
81. Vervaeke, S, Matthys, C, Nassar, R, et al., Adapting the vertical position of implants with a conical connection in relation to soft tissue thickness prevents early implant surface exposure: A 2-year prospective intra-subject comparison. *J Clin Periodontol* 2018;45(5):605-12. [Abstract](#)
82. Younes, F, Cosyn, J, De Bruyckere, T, et al., A 2-year prospective case series on volumetric changes, PROMs, and clinical outcomes following sinus floor elevation using deproteinized bovine bone mineral as filling material. *Clin Implant Dent Relat Res* 2019;21(2):301-09. [Abstract](#)
83. Barbier, L, Pottel, L, De Ceulaer, J, et al., Evaluation of quality of life after mandibular reconstruction using a novel fixed implant-supported dental prosthesis concept: A pilot study. *Int J Prosthodont* 2019;32(2):162-73. [Abstract](#)
84. Barbier, L, Abeloos, J, De Clercq, C, et al., Peri-implant bone changes following tooth extraction, immediate placement and loading of implants in the edentulous maxilla. *Clin Oral Investig* 2012;16(4):1061-70. [Abstract](#)
85. Bashutski, JD, Wang, HL, Rudek, I, et al., Effect of flapless surgery on single-tooth implants in the esthetic zone: a randomized clinical trial. *J Periodontol* 2013;84(12):1747-54. [Abstract](#)
86. Borges, T, Leitao, B, Pereira, M, et al., Influence of the abutment height and connection timing in early peri-implant marginal bone changes: A prospective randomized clinical trial. *Clin Oral Implants Res* 2018;29(9):907-14. [Abstract](#)
87. Cooper, LF, Raes, F, Reside, GJ, et al., Comparison of radiographic and clinical outcomes following immediate provisionalization of single-tooth dental implants placed in healed alveolar ridges and extraction sockets. *Int J Oral Maxillofac Implants* 2010;25(6):1222-32. [Abstract](#)
88. Cooper, LF, Reside, G, Stanford, C, et al., A multicenter randomized comparative trial of implants with different abutment interfaces to replace anterior maxillary single teeth. *Int J Oral Maxillofac Implants* 2015;30(3):622-32. [Abstract](#)
89. D'haese, J, Vervaeke, S, Verbanck, N, et al., Clinical and radiographic outcome of implants placed using stereolithographic guided surgery: a prospective monocenter study. *Int J Oral Maxillofac Implants* 2013;28(1):205-15. [Abstract](#)
90. Donati, M, La Scala, V, Billi, M, et al., Immediate functional loading of implants in single tooth replacement: a prospective clinical multicenter study. *Clin Oral Implants Res* 2008;19(8):740-48. [Abstract](#)

91. Ebler, S, Ioannidis, A, Jung, RE, et al., Prospective randomized controlled clinical study comparing two types of two-piece dental implants supporting fixed reconstructions - results at 1 year of loading. *Clin Oral Implants Res* 2016;27(9):1169-77. [Abstract](#)
92. Galindo-Moreno, P, Nilsson, P, King, P, et al., Clinical and radiographic evaluation of early loaded narrow diameter implants - 1-year follow-up. *Clin Oral Implants Res* 2012;23(5):609-16. [Abstract](#)
93. Ghozei, R, Alikhasi, M, Siadat, M-R, et al., A radiographic comparison of progressive and conventional loading on crestal bone loss and dentistry in single dental implants: A randomized controlled trial study. *J Dent (Tehran)* 2013;10(2):155-63. [Abstract](#)
94. Guljé, F, Raghoobar, GM, Ter Meulen, JW, et al., Mandibular overdentures supported by 6-mm dental implants: A 1-year prospective cohort study. *Clin Implant Dent Relat Res* 2012;14(Supplement 1):e59-e66. [Abstract](#)
95. Guljé, F, Abrahamsson, I, Chen, S, et al., Implants of 6 mm vs. 11 mm lengths in the posterior maxilla and mandible: a 1-year multicenter randomized controlled trial. *Clin Oral Implants Res* 2013;24(12):1325-31. [Abstract](#)
96. Guljé, FL, Raghoobar, GM, Vissink, A, et al., Single restorations in the resorbed posterior mandible supported by 6-mm implants: a 1-year prospective case series study. *Clin Implant Dent Relat Res* 2015;17 Suppl 2:e465-71. [Abstract](#)
97. Guljé, FL, Raghoobar, GM, Vissink, A, et al., Single crowns in the resorbed posterior maxilla supported by either 6-mm implants or by 11-mm implants combined with sinus floor elevation surgery: A 1-year randomised controlled trial. *Eur J Oral Implantol* 2014;7(3):247-55. [Abstract](#)
98. Guljé, FL, Raghoobar, GM, Erkens, WA, et al., Impact of crown-implant ratio of single restorations supported by 6-mm implants: A short-term case series study. *Int J Oral Maxillofac Implants* 2016;31(3):672-5. [Abstract](#)
99. Han, J, Zhang, X, Tang, Z, et al., A prospective, multicenter study assessing the DENTSPLY Implants, OsseoSpeed TX, length 6 mm in the posterior maxilla and mandible: a 1-year follow-up study. *Clin Oral Implants Res* 2016;27(4):452-7. [Abstract](#)
100. Hosseini, M, Worsaae, N, Schiodt, M, et al., A 1-year randomised controlled trial comparing zirconia versus metal-ceramic implant supported single-tooth restorations. *Eur J Oral Implantol* 2011;4(4):347-61. [Abstract](#)
101. Ivanoff, CJ, Lindhe, J, Ellner, S, et al., An open, randomised, multi-centre study, comparing straight and tapered apex implants design, in partially and totally edentulous maxillae. *Acta Odontol Scand* 2021;E-pub Mar 5 doi:10.1080/00016357.2021.1894352:1-7. [Abstract](#)
102. Kim, JJ, Lee, DW, Kim, CK, et al., Effect of conical configuration of fixture on the maintenance of marginal bone level: preliminary results at 1 year of function. *Clin Oral Implants Res* 2010;21(4):439-44. [Abstract](#)
103. Koutouzis, T, Koutouzis, G, Tomasi, C, et al., Immediate loading of implants placed with the osteotome technique: One-year prospective case series. *J Periodontol* 2011;82(11):1556-62. [Abstract](#)
104. Lee, DW, Lee, DW, Park, KH, et al., The effects of off-axial loading on periimplant marginal bone loss in a single implant. *J Prosthet Dent* 2014;112(3):501-7. [Abstract](#)
105. Lyngstadaas, SP, Verket, A, Pinholt, EM, et al., Titanium granules for augmentation of the maxillary sinus - a multicenter study. *Clin Implant Dent Relat Res* 2015;17 Suppl 2:e594-600. [Abstract](#)
106. Malmstrom, H, Gupta, B, Ghanem, A, et al., Success rate of short dental implants supporting single crowns and fixed bridges. *Clin Oral Implants Res* 2016;27(9):1093-8. [Abstract](#)
107. Marcelis, K, Vercruyssen, M, Naert, I, et al., Model-based guided implant insertion for solitary tooth replacement: a pilot study. *Clin Oral Implants Res* 2012;23(8):999-1003. [Abstract](#)
108. Noelken, R, Donati, M, Fiorellini, J, et al., Soft and hard tissue alterations around implants placed in an alveolar ridge with a sloped configuration. *Clin Oral Implants Res* 2014;25(1):3-9. [Abstract](#)
109. Norton, MR, The Influence of Low Insertion Torque on Primary Stability, Implant Survival, and Maintenance of Marginal Bone Levels: A Closed-Cohort Prospective Study. *Int J Oral Maxillofac Implants* 2017;32(4):849-57. [Abstract](#)
110. Pieri, F, Aldini, NN, Fini, M, et al., Immediate fixed implant rehabilitation of the atrophic edentulous maxilla after bilateral sinus floor augmentation: a 12-month pilot study. *Clin Implant Dent Relat Res* 2012;14 (Suppl 1):e67-82. [Abstract](#)
111. Piero, B, Mario, V, Niccolo, N, et al., Implant placement in combination with sinus membrane elevation without biomaterials: A 1-year study on 15 patients. *Clin Implant Dent Relat Res* 2012;14(5):682-9. [Abstract](#)
112. Raes, F, Cosyn, J, De Bruyn, H, Clinical, aesthetic, and patient-related outcome of immediately loaded single implants in the anterior maxilla: A prospective study in extraction sockets, healed ridges, and grafted sites. *Clin Implant Dent Relat Res* 2013;15(6):819-35. [Abstract](#)
113. Raes, F, Cosyn, J, Crommelinck, E, et al., Immediate and conventional single implant treatment in the anterior maxilla: 1-year results of a case series on hard and soft tissue response and aesthetics. *J Clin Periodontol* 2011;38(4):385-94. [Abstract](#)
114. Rismanchian, M, Fazel, A, Rakhshan, V, et al., One-year clinical and radiographic assessment of fluoride-enhanced implants on immediate non-functional loading in posterior maxilla and mandible: a pilot prospective clinical series study. *Clin Oral Implants Res* 2011;22(12):1440-5. [Abstract](#)
115. Schepke, U, Meijer, HJ, Kerdijk, W, et al., Stock versus CAD/CAM customized zirconia implant abutments - clinical and patient-based outcomes in a randomized controlled clinical trial. *Clin Implant Dent Relat Res* 2017;19(1):74-84. [Abstract](#)
116. Schincaglia, GP, Thoma, DS, Haas, R, et al., Randomized controlled multicenter study comparing short dental implants (6 mm) versus longer dental implants (11-15 mm) in combination with sinus floor elevation procedures. Part 2: clinical and radiographic outcomes at 1 year of loading. *J Clin Periodontol* 2015;42(11):1042-51. [Abstract](#)
117. Schincaglia, GP, Rubin, S, Thacker, S, et al., Marginal bone response around immediate- and delayed-loading implants supporting a locator-retained mandibular overdenture: a randomized controlled study. *Int J Oral Maxillofac Implants* 2016;31(2):448-58. [Abstract](#)
118. Simmons, DE, Maney, P, Teitelbaum, AG, et al., Comparative evaluation of the stability of two different dental implant designs and surgical protocols-a pilot study. *Int J Implant Dent* 2017;3(1):16. [Abstract](#)
119. Slot, W, Raghoobar, GM, Vissink, A, et al., Maxillary overdentures supported by four or six implants in the anterior region; 1-year results from a randomized controlled trial. *J Clin Periodontol* 2013;40(3):303-10. [Abstract](#)
120. Slot, W, Raghoobar, GM, Vissink, A, et al., Maxillary overdentures supported by anteriorly or posteriorly placed implants opposed by a natural dentition in the mandible: a 1-year prospective case series study. *Clin Implant Dent Relat Res* 2014;16(1):51-61. [Abstract](#)
121. Stanford, CM, Barwacz, C, Raes, S, et al., Multicenter clinical randomized controlled trial evaluation of an implant system designed for enhanced primary stability. *Int J Oral Maxillofac Implants* 2016;31(4):906-15. [Abstract](#)
122. Temmerman, A, Rasmussen, L, Kubler, A, et al., An open, prospective, non-randomized, controlled, multicentre study to evaluate the clinical outcome of implant treatment in women over 60 years of age with osteoporosis/osteopenia: 1-year results. *Clin Oral Implants Res* 2017;28(1):95-102. [Abstract](#)
123. Toia, M, Stocchero, M, Becktor, JP, et al., Implant vs abutment level connection in implant supported screw-retained fixed partial dentures with cobalt-chrome framework: 1-year interim results of a randomized clinical study. *Clin Implant Dent Relat Res* 2019;21(2):238-46. [Abstract](#)
124. Toljanic, JA, Baer, RA, Ekstrand, K, et al., Implant rehabilitation of the atrophic edentulous maxilla including immediate fixed provisional restoration without the use of bone grafting: a review of 1-year outcome data from a long-term prospective clinical trial. *Int J Oral Maxillofac Implants* 2009;24(3):518-26. [Abstract](#)
125. Tu, CC, Wang, SH, Hu, CY, et al., Comparison of regular- and slope-configured dental implants placed in the edentulous ridge with height discrepancy: A pilot randomized controlled trial. *J Formos Med Assoc* 2021;120(6):1386-93. [Abstract](#)
126. Tsuda, H, Rungcharassaeng, K, Kan, JY, et al., Peri-implant tissue response following connective tissue and bone grafting in conjunction with immediate single-tooth replacement in the esthetic zone: A case series. *Int J Oral Maxillofac Implants* 2011;26(2):427-36. [Abstract](#)
127. Vercruyssen, M, van de Wiele, G, Teughels, W, et al., Implant- and patient-centred outcomes of guided surgery, a 1-year follow-up: An RCT comparing guided surgery with conventional implant placement. *J Clin Periodontol* 2014;41(12):1154-60. [Abstract](#)
128. Zhou, J, Huang, Q, Wang, X, et al., Early loading of splinted implants in the posterior mandible: a prospective multicentre case series. *J Clin Periodontol* 2016;43(3):298-304. [Abstract](#)

THE DENTAL
SOLUTIONS
COMPANY™

