Lucitone[®] FRS[®]

Flexible Dental Resin

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Resina flessibile per uso dentale

Flexibler Dentalkunststoff

Resina dental flexible

Directions For Use With The SUCCESS[®] INJECTION SYSTEM

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Small (S), Medium (M) and Large (L)

Indications for use:

Lucitone FRS is intended for fabrication of partial or full removable dentures as well as occlusal splints and nightguards.

Product Information:

Lucitone FRS is a flexible and monomerfree thermoplastic dental resin. This product is intended for the fabrication of temporary removable partial dentures or small to medium size removable full dentures, as well as clasps, occlusal splints and nightguards. Available Sizes: Small (S), Medium (M) and Large (L).

Contraindications:

1. This product is not intended for long-term restorations and is intended only for provisional or temporary applications. Metal frame partial dentures remain the "standard" for long-term restorations.
2. This product is not intended to be used as a unilateral restoration. Restorations of the unilateral edentulous space must utilize a bilateral framework design.

Warnings:

1. When using Lucitone FRS, extreme caution must be used to avoid skin contact with the heated sleeve, cartridge furnace heating bay, hot cartridge, injection insert, piston head adapter, hot flasks, and heat lamps or alternate heat source. Always use heat resistant gloves with a rating of at least 525°F (267°C).

Precautions:

1. When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.
2. Store at 60°F - 90°F (16°C - 32°C),

away from moisture and direct sunlight. Excessive moisture content may cause voids (porosity) in processed dentures. Always store the Lucitone FRS cartridges in their foil bags until ready for use.

3. Accurate denture base shade is achieved only when processed as described below, and heating times are followed.
4. Dispose of contents/container in accordance with local and national regulations.

Adverse Reactions:

1. Particulates will be generated when grinding thermoplastic resins. Eye, skin and respiratory irritation may occur if appropriate safety precautions are not taken.

STEP BY STEP INSTRUCTIONS: **Properties and Recommended Procedures**

The following table describes typical properties of Lucitone FRS resin plus recommended processing temperatures and other important processing information. Detailed information is found later in the text.

Type	ISO 20795-1 Type 3
Material used to invest wax-up	Dental Stone
Storage temperature for powder and liquid	60°F - 90°F (16°C - 32°C)
Dosage	1 cartridge
Processing Time / Temp. (To reach injectable stage)	17 minutes @ 575°F

Injection Pressure	75 psi minimum (100 psi recommended) with Narrow Piston Head
Injection Time	1 minute
Cool time (prior to deflasking)	5 minutes

FABRICATION OF FULL AND PARTIAL DENTURES, CLASPS, OCCLUSAL SPLINTS AND NIGHTGUARDS

Note: Lucitone FRS resin should not be used on large full dentures. Lucitone FRS resin is indicated for partial dentures, small to medium full dentures, clasps, occlusal splints and nightguards.

For a firmer, resilient flexibility, finish the Lucitone FRS flexible appliance to no less than 2mm thickness. For the most flexibility, finish to no less than 1.5mm thickness.

I. INVESTMENT AND SPRUING

Soak the cast and wax-up in water for 25 minutes prior to investing. Mix investment and follow normal Success Instructions to embed the cast and wax-up in the lower half of the Success Flask.

IMPORTANT: REMOVE ALL EXCESS INVESTMENT FROM THE FLAT LAND AREA OF THE FLASK.

Use wax sticks, such as the D-shaped Success Sprue Wax Sticks (7mm X 180mm, Item #904584) to build the injection sprues. On full upper dentures, attach the sprue to the posterior border of the denture and make sure the sprue is

sufficiently wide. Extend the sprue over the palate area to allow adequate flow of the material throughout the palate area (Figure 1). For partial dentures, lower full dentures, occlusal splints and night guards, attach the sprue to both lingual extensions, as well as in the mid-line as shown in Figure 2. Following investment of the lower half of the case, embed the upper half of the flask following normal Success Instructions.


Figure 1:
Full Upper


Figure 2:
Partial Denture

II. WAX ELIMINATION AND SEPARATOR APPLICATION

Boil out: Place the flask in boiling water for 4 to 6 minutes to soften the wax. Loosen the bolts on the flask and remove the Metal Flask Brackets. Open the flask and remove the black Space Maintainer. Complete boil out procedure, discarding wax. Flush with clean boiling water. Bevel the stone with a lab knife around the mouth of the sprue. Check the flask margin and cavities to ensure that both flask halves fit together with intimate metal contact.

Separator: Apply a thin coat of Al-Cote® Separating Agent to model and allow the model to dry completely.

Retention: Add mechanical retention (diatorics) to the teeth. Use vent holes, if possible. Remove only enough tooth

material to enhance retention. Excessive removal of tooth material will weaken the tooth and may lead to fracture.

III. HEATING AND INJECTION

Important:

- Oven must be switched on with the Cartridge Sleeve in place and set to 575°F at least 15 minutes prior to processing.
- Keep the Cartridge Sleeve in the Cartridge Furnace at all times, except when injecting.
- If using Dentsply Sirona Heat Lamps (904860), remove protective plastic covering from base pan of heat lamps prior to use. Heat lamps must be used at their lowest height setting.
- The Success Narrow Piston Head Adaptor (reorder #904865) must be used for the injection of Lucitone FRS resin.

Heat Cartridge: Spray Dentsply Sirona Silicone Spray on a Lucitone FRS Cartridge. Using heat resistant gloves, insert the cartridge into the Cartridge Sleeve. The Cartridge Sleeve must remain in the furnace during this process. The nozzle of the cartridge must be facing inward. Insert the cartridge into the Cartridge Sleeve as far as possible. Use the angled end of the Knockout Rod to assist placement.

Heating Flask with Heat Lamps: Position Injection Insert on bolt side of the flask and place flask halves directly under heat lamps. Turn on heat lamps.

Heating Flask with Alternate Heat Source: Assemble the flask halves with brackets

and the injection insert. Heat the flask using a dry heat source at approximately 200°F.

Set Timer: Set timer for 17 minutes. Cartridges removed prior to 17 minutes may not properly inject. Heating for less than 17 minutes will result in a slightly darker shade, while more than 17 minutes will result in a slightly lighter shade. **DO NOT HEAT LONGER THAN 20 MINUTES.**

Inject: Caution flask will be hot. Use heat resistant gloves. After heating for 17 minutes, remove flask from heat source, and remove power. If heat lamps were used, assemble the warm flask halves with brackets and the injection insert. Place this flask assembly in front of the Success Injection System. Remove the Cartridge Sleeve and cartridge assembly from the furnace, keeping the cartridge assembly horizontal while transporting it to the flask assembly. Position the cartridge assembly on top of the flask so that the nozzle of the cartridge fits into the opening of the Injection Insert. Then grasp the flask assembly with one hand, the cartridge assembly with the other hand, and slide the combined assemblies into the Success Injection System. Make sure that the Narrow Piston Head is properly aligned with the Cartridge Sleeve. Engage the piston by depressing the activation switch.

Note: The time between removing the cartridge assembly from the oven and injection should be less than 45 seconds. If longer, the cartridge will begin to cool and may result in partial or no injection.


The injection piston should remain engaged for one (1) minute. After injection, immediately remove the flask assembly from the system. Hold the cartridge assembly with one hand, then

rotate the flask assembly 360° with the other hand. Then immediately disengage the cartridge assembly from the flask assembly. Finally, expel the used cartridge using the Knock Out Base and Knock Out Rod and return the Cartridge Sleeve to the furnace. The Cartridge Sleeve must remain in the furnace when not in use.

Bench cool the flask assembly for 5 minutes before deflasking. Caution should be used when handling the hot Cartridge Sleeve, Injection Insert, and used cartridge. Deflask carefully and remove the sprues with a cut off disk.

IV. CLEAR CLASPS

Complete investment and spruing as per normal Success Instructions. Eliminate wax as stated above. Inject the clasps using Lucitone FRS Clear resin with the heating and injection technique described above. Finish the clasps to final contour with carbide burs, leaving a tail section that will be incorporated into the denture base. Refit the clasps onto the master model. Relieve the tail area and grind retentive holes through it so that the denture base resin will completely encase the clasps. See illustration:

Note: At least 1mm of thickness is necessary on each side of the tail for the clear clasp to be properly incorporated. Complete denture wax-up and tooth setup. Repeat steps I through IV above, and

inject case with tissue shaded Lucitone FRS.

V. FINISHING AND POLISHING

Lucitone FRS resin should be finished and polished using normal procedures for acrylics. A sharp carbide bur can be used for bulk reduction. Use a bur followed by coarse pumice for finishing. For a firmer, resilient flexibility, finish the Lucitone FRS flexible partial to no less than 2mm thickness. For the most flexibility, finish no less than 1.5mm thickness.

VI. REPAIR TECHNIQUES

Voids or defects due to processing in **non-flexing** areas can be filled using traditional thermoplastic repair techniques. The area to be repaired should be ground out using a rough cut bur. Mechanical retention should be cut into the denture base. Use Lucitone® Repair Material or Fas-Por™+ pourable dental resin to fill the area. Cure in a pressure pot according to the manufacturer's recommended instructions. Finish repair using normal acrylic finishing techniques.

If repair of **flexing** areas is required, or tooth additions or relines are necessary, the existing case should be rebased with new Lucitone FRS dental resin. To rebase the case, mount the denture on the model, trimming away all denture base except the flange areas holding the teeth. Grind mechanical retention into these areas so that the newly injected material will bond securely. Perform tooth addition (if required) and eliminate wax as described above. Heat and inject as stated above.

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York, PA 17401 U.S.A.
EC REP DeguDent GmbH