

Sirona Connect

DIGITAL IMPRESSIONS TRAINING



THE DIGITALIZATION TREND IS PREVAILING.

THE CONVENTIONAL IMPRESSION



PRACTICE



LABORATORY

The production of an impression and the resulting working model is still one of the biggest sources of error between the practice and laboratory!

PROCESSES IN COMPARISON

PRACTICE

Tray
selection

Tray
preparation

Silicone
impression

Total
impression

Disinfect

Order form

Clean

Cut
impressions

Pour with
plaster

Trim the
dental arches

Set pins

Trim the
base

Saw cuts

Preparation
margins

LABORATORY

Intraoral
imaging

Feedback

Transfer data

Laboratory
receives data

Feedback

Production

CONNECT

SIRONA CONNECT WORKFLOW

PRACTICE

PORTAL

LABORATORY

HOW SIRONA CONNECT WORKS:



PROCESS
WITH inLab



Central
production
PROCESS
E.G. WITH InfiniDent



STL
PROCESS
WITH THIRD-PARTY SYSTEMS



GENERAL BENEFITS OF DIGITAL IMPRESSIONS

PRECISION

HIGH PRECISION

High quality standards thanks to modern technology

PROCESS

SIMPLE PROCESS

Less work steps – greater efficiency and time savings

COMFORT

NO MORE PHYSICAL IMPRESSIONS

Less stressed patient, without the feared regurgitation effect



BENEFITS FOR THE DENTIST: BECAUSE YOUR PATIENT DESERVES IT

SIMPLE

SIMPLER PRACTICE ROUTINE

Scan, check and send. The direct link to the laboratory

PRECISE

PRECISE IMAGES

Digital processes ensure consistent precision


PATIENT FOCUS

BECAUSE YOUR PATIENT DESERVES IT

High level of patient comfort
Keeping existing patients and gaining new ones (practice's prestige)



SIRONA CONNECT PORTAL



SIRONA CONNECT PORTAL

- ☐ Portal login for dentists and dental technicians
- ☐ Manage my accounts
- ☐ www.sirona-connect.com: what is on the homepage under Service/Downloads?
- ☐ Connect to external portals like infiniDent, Invisalign etc
- ☐ Connect Direct

REGISTRATION

WWW.SIRONA-CONNECT.COM

REGISTRATION

The screenshot displays the Sirona Corporate Website. At the top, the navigation bar includes 'SIRONA CORPORATE WEBSITE' on the left and 'LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN' on the right. The Sirona logo, 'The Dental Company', is positioned in the upper right. Below the navigation bar, a secondary menu lists 'ABOUT SIRONA CONNECT', 'FOR DENTAL PRACTICES', 'FOR DENTAL LABORATORIES', 'SUPPORT', and 'CONTACT'. The main content area features a large heading 'Digital impression with Sirona Connect' and the subtext 'New perspectives for dental prac-'. Below this is a 'Login' section with input fields for 'Username' and 'Password', and a 'Login' button. To the left of the login fields, a large orange arrow points to a list of links: 'Dentist Registration', 'Laboratory Registration', and 'Forgot password?'. To the right of the login section, there is a large image of two men in a dental office setting, one on the phone and one at a computer. A yellow box with the text 'FREE CEREC CONNECT SW 4.X DOWNLOAD' is overlaid on the image.

SIRONA CORPORATE WEBSITE

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

sirona.
The Dental Company

ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SUPPORT CONTACT

Digital impression with Sirona Connect

New perspectives for dental prac-

Login

Username

Password

[Dentist Registration](#)

[Laboratory Registration](#)

[Forgot password?](#)

Login

FREE CEREC CONNECT SW 4.X DOWNLOAD



**Set username
and
password**

STEP 1

STEP 2

STEP 3

DONE

SIRONA Connect
MY SIRONA CONNECT > DENTIST REGISTRATION

- Laboratory Registration
- Dentist Registration**
- Login

Welcome to SIRONA Connect



Dentist Registration

Step 1Step 2Step 3Done

Your access data

At first please choose your personal and unique user name, which you can easily remember. If you have questions, give us a call: +49 (0)6251 16-1670.

Your favorite username *

Your chosen password *

*marked fields are mandatory fields

STEP 1

STEP 2

STEP 3

DONE

SIRONA Connect

- Laboratory Registration
- **Dentist Registration**
- Login

TV SIRONA CONNECT > DENTIST REGISTRATION

Welcome to SIRONA
Connect

Dentist Registration

Step 1

Step 2

Step 3

Done

My user information

Salutation *

☒ Mr. ☐ Mrs.

Title

Dr.

First Name *

Mark

Last Name *

Homburger

Phone Number *

0177 69 60 311

E-Mail Address *

rob@F-Dental.com

Confirm e-Mail address *

rob@F-Dental.com

Skype username

DentFFM

My practice data

Company Name *

Praxis Dr Homburger

E-Mail Address

DrH@PraxisHomburger.de

Phone Number *

061 72 69 60 311

Website

www.Dental.de

Fax Number

My preferences

Preferred Language

english (US)

Timezone

(UTC) Coordinated Universal Time

*marked fields are mandatory fields

BACK

NEXT

Enter complete
user data

STEP 1

STEP 2

STEP 3

DONE

Enter billing and
delivery address

MY SIRONA.COM > DENTIST REGISTRATION

Welcome to SIRONA
Connect



Dentist Registration

Step 1 ✓ Step 2 ✓ Step 3 Done

My billing address

Address 1 *
Gutleutstr. 33

Address 2

ZIP Code *
60311

City *
Frankfurt

Country *
Germany

Region
Hessen

☒ different shipping address

My shipping address

Recipient
Praxis Dr. Homburger

Address 1 *
Gutleutstr. 19


Address 2

ZIP Code *
60311

City *
Frankfurt

Country *
Germany

Region
Hessen



* required fields are mandatory fields

BACK

DONE

STEP 1

STEP 2

STEP 3

DONE

MY SIRONA CONNECT > DENTIST REGISTRATION

Welcome to SIRONA Connect



Dentist Registration

Step 1 ▾

Step 2 ▾

Step 3 ▾

Done

Congratulations! You nearly did it.

Your new Sirona account was created. Prior to your first login, we would like to check whether your e-mail address is valid und working.

We have therefore sent you an e-mail with the subject **"Please confirm your e-mail address"** which you can use to activate your account. Therefore please follow the instructions inside the e-mail.

You will find a link in the e-mail

If you didn't receive the e-mail, please check your spam folder. If you have problems, call our support (phone number + 49 (0) 6251 / 16-1670) or send us an e-mail (support@sirona.com).

We wish you good luck with Sirona Connect!

**Login complete.
Confirmation sent
separately by
email.**

EMAIL CONFIRMATION FOR REGISTRATION AS A DENTIST

SIRONA Connect – Please activate your account

Von: cerec.connect@sirona.de

An: alex@drauz.net

Datum: 25.03.2014 19:16

Dear Mr. Homburger

Welcome at SIRONA Connect! You have been registered on 03/25/2014 6:16:57 PM.

Your details are:

Username: Praxis Dr Homburger

Dr. Mark Homburger

Phone: 061 72 69 60 311

Email: DrH@PraxisHomburger.de

Dentist data:

Praxis Dr Homburger

Gutleutstr 33

60311 Frankfurt

Germany

Fax:


Website: www.Dental.de

You can activate your account with the following [Activation Link](#)

Sincerely,
Your SIRONA Connect Team



CONFIRMATION OF ACTIVATION



LOGIN AFTER REGISTRATION

Enter username and
password in the login
field

SIRONA CORPORATE WEBSITE

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SUPPORT CONTACT

sirona
The Dental Company

Digital impression with Sirona Connect

New perspectives for dental practices and laboratories.

Login

Praxis Dr. Homburger

.....

[Dentist Registration](#)
[Laboratory Registration](#)
[Forgot password?](#)

Login

FREE CEREC CONNECT SW 4.X
DOWNLOAD

SIRONA Connect

- My Account
- User Administration
- Order List
- Log Off

MY SIRONA CONNECT > WELCOME

Welcome to SIRONA Connect



Welcome

Welcome to Sirona Connect!

Your last visit was on 03/26/2014 3:31 PM

> You can edit your personal account on the following pages.

**Activation was
successful!**

ADD YOUR SKYPE USER NAME UNDER MY ACCOUNT

Vorname *

Ronny

E-Mail-Adresse
bestätigen *

ronny.kucharczyk@sirona.com

Nachname *

Kucharczyk

Skype-Benutzername

ronny.kucharczyk

Chat-Avatar



Also a Chat-Avatar can
be uploaded

AS A LAB PLEASE INDICATE YOUR SERVICE OFFERING

SIRONA CORPORATE WEBSITE | SPRACHE DEUTSCH ▼

Änderungen erfolgreich gespeichert!

LABORREGISTRIERUNG | ZAHNARZTREGISTRIERUNG | ANMELDUNG

ÜBER SIRONA CONNECT | FÜR DIE PRAXIS | FÜR DAS DENTALLABOR | SERVICE | KONTAKT



sirona
The Dental Company

SIRONA Connect

- Mein Konto
 - Mein Benutzer
 - Mein Labor
 - Meine Rechnungsadresse
 - Meine Lieferadresse
 - Angebote/Services
 - Connect Direct
 - Externe Laborportale
 - Passwort ändern
- Benutzerverwaltung
- Auftragsliste
- Abmelden

MEIN SIRONA CONNECT > MEIN KONTO > ANGEBOTE/SERVICES

Willkommen bei
SIRONA Connect



Angebote/Services

Geben Sie die zusätzlichen Informationen zu Ihren Angeboten und Dienstleistungen ein.

Angebote	Services
<input checked="" type="checkbox"/> Ästhetische Frontzahnrestorationen	<input checked="" type="checkbox"/> Lieferung durch Fahrer
<input checked="" type="checkbox"/> Sirona Materialien	<input checked="" type="checkbox"/> Versand
<input checked="" type="checkbox"/> Veneers	<input checked="" type="checkbox"/> Farbnahme
<input checked="" type="checkbox"/> Ivoclar Vivadent Materialien	
<input checked="" type="checkbox"/> VITA Materialien	
<input checked="" type="checkbox"/> 3M COPE Materialien	
<input checked="" type="checkbox"/> Kieferorthopädische Apparaturen	

SPEICHERN

ABBRECHEN

CONNECT TO EXTERNAL PORTALS LIKE INVISALIGN

SIRONA CORPORATE WEBSITE | SPRACHE DEUTSCH ▾

LABÖRREGISTRIERUNG | ZAHNÄRZTREGISTRIERUNG | ANMELDUNG

ÜBER SIRONA CONNECT | FÜR DIE PRAXIS | FÜR DAS DENTALLABOR | SERVICE | KONTAKT

sirona
The Dental Company

SIRONA Connect

- Mein Konto
 - Mein Benutzer
 - Meine Praxis
 - Meine Rechnungsadresse
 - Meine Lieferadresse
 - Meine bevorzugten Labore
 - Externes Portal
 - Invisalign® Clear Aligner
 - Passwort ändern
- Benutzerverwaltung
- Auftragsliste
- Abmelden

MEIN SIRONA CONNECT

Willkommen bei
SIRONA Connect





Sie können die Schaltfläche „Verknüpfen“ verwenden, um eine Verbindung zum Invisalign®-Portal herzustellen. Mit der Schaltfläche „Verknüpfung aufheben“ können Sie die Verbindung trennen.

VERKNÜPFEN

MANAGE MY ACCOUNTS AS A USER

- The administrator can create more user accounts under Manage accounts
- If administrators are granted rights, then the role owner can change and save all data
- If a "user role" is granted, then the role owner can edit Sirona Connect cases, but he cannot change the user data (address, etc.)
- It is also possible to define whether the role owner receives confirmation emails



SIGNATURE

- Dentists can „sign“ the order sheet they send to the dental laboratory
- During the registration they can upload a certificate or an image of their signature and add a text that will be shown on the work ticket
- They can also add the information under My Account / My user

The screenshot shows the SIRONA Connect web interface. At the top, there's a navigation bar with 'SIRONA Connect' and links for 'Labordienstleistungen', 'Zahnarztregistrierung', and 'Anmeldung'. The main heading is 'Willkommen bei SIRONA Connect'. Below this, there's a 'Zahnarztregistrierung' section with buttons for 'Schritt 1', 'Schritt 2', 'Schritt 3', and 'Fertig'. The 'Meine Benutzerinformationen' section contains a form with fields for 'Anzahl' (with checkboxes for 'Hilf' and 'Frei'), 'Titel', 'Vorname', 'Nachname', 'E-Mail-Adresse', and 'Passwort'. A 'Chat-Avatar' field shows a small profile picture. Below the form, there's a 'My Workticket signature preferences' section with a checked box for 'Use Workticket signature'. To the right, there's a 'Signature File' field showing a sample signature 'John Doe' and a 'Workticket Signature Test' button.

IF YOU'RE A **LAB OWNER**

REGISTRATION AS A LABORATORY

The screenshot displays the Sirona Corporate Website. At the top, navigation links include "LABORATORY REGISTRATION", "DENTIST REGISTRATION", and "LOGIN". The main header features the Sirona logo and the tagline "The Dental Company". Below this, a secondary navigation bar lists "ABOUT SIRONA CONNECT", "FOR DENTAL PRACTICES", "FOR DENTAL LABORATORIES", "SUPPORT", and "CONTACT".

The central content area is titled "Digital impression with Sirona Connect" and "New perspectives for dental prac-". It includes a "Login" section with input fields for "Username" and "Password", and a "Login" button. To the left of the login fields, a large orange arrow points to a list of links: "Dentist Registration", "Laboratory Registration" (which is highlighted with a black border), and "Forgot password?".

Below the login section, there are two images of a man in a white shirt talking on a mobile phone while sitting at a desk with a computer. A yellow banner overlaid on the image reads "FREE CEREC CONNECT SW 4.X DOWNLOAD".

STEP 1

STEP 2

STEP 3

DONE

SIRONA Connect

MY SIRONA CONNECT » LABORATORY REGISTRATION

- Laboratory Registration
- Dentist Registration
- Login

Welcome to SIRONA
Connect



Laboratory Registration

Step 1

Step 2

Step 3

Done

Your access data

At first please choose your personal and unique user name, which you can easily remember. If you have questions, give us a call: +49 (0)6251 16-1670.

Your favorite username *

Your chosen password *

* marked fields are mandatory fields

CANCEL

NEXT

Set username
and
password

STEP 1

STEP 2

STEP 3

DONE

SIRONA Connect

- Laboratory Registration
- Dentist Registration
- Login

MY SIRONA CONNECT > LABORATORY REGISTRATION

Welcome to SIRONA
Connect

Laboratory Registration

Step 1 ▾

Step 2

Step 3

Done

My user information

Salutation *

☐ Mr.☐ Mrs.

Phone Number *

009 60311 09

Title

Dr

E-Mail Address *

rob@F-Dental.com

First Name *

Robin

Confirm e-Mail address *

rob@F-Dental.com

Last Name *

Frankfurter

Skype username

RobFFM

My laboratory data

Company Name *

Dentalab Frankfurter

E-Mail Address

RF@Frankfurter-Dental.com

Phone Number *

009 60311 71

Website

www.Dental-FFM.de

Fax Number

009 60311 77

My preferences

Preferred Language

english (GB) ▾

Timezone

(UTC) Casablanca ▾

* marked fields are mandatory fields

BACK

NEXT

Enter complete
user data

STEP 1

STEP 2

STEP 3

DONE

Enter billing and
delivery address

SIRONA Connect

- Laboratory Registration
- Dentist Registration
- Login

MY SIRONA CONNECT > LABORATORY REGISTRATION

Welcome to SIRONA
Connect



Laboratory Registration

Step 1 ▾ Step 2 ▾ Step 3 ▾ Done



My billing address

Address 1 * Goethestrasse 23b
Address 2
ZIP Code * 60311
City * Frankfurt
Country * Germany ▾
Region * Hessen ▾

☒ different shipping address

My shipping address

Recipient Dentalab Frankfurt
Address 1 * Goethestrasse 23b
Address 2
ZIP Code * 60311
City * Frankfurt
Country * Germany ▾
Region * Hessen ▾



* marked fields are mandatory fields

BACK

DONE

STEP 1

STEP 2

STEP 3

DONE

SIRONA Connect

- Laboratory Registration
- Dentist Registration
- Login

MY SIRONA CONNECT > LABORATORY REGISTRATION

Welcome to SIRONA Connect



Laboratory Registration

Step 1 ▾

Step 2 ▾

Step 3 ▾

Done

Congratulations! You nearly did it.

Your new Sirona account was created. Prior to your first login, we would like to check whether your e-mail address is valid und working.

We have therefore sent you an e-mail with the subject "**Please confirm your e-mail address**" to rob@F-Dental.com. You will find a link in the e-mail which you can use to activate your account. Therefore please follow the instructions inside the e-mail.

If you didn't receive the e-mail, please check your spam folder. If you have problems, call our support (phone number + 49 (0) 6251 / 16-1670) or send us an e-mail (support@sirona.com).

We wish you good luck with Sirona Connect!

**Login complete.
Confirmation sent
separately by
email.**

EMAIL CONFIRMATION FOR REGISTRATION AS A LABORATORY

☐ **SIRONA Connect – Please activate your account**

Von: cerec.connect@sirona.de

An: alex@drazu.net

Datum: 18:57

Dear Mr. Homburger

Welcome at SIRONA Connect! You have been registered on 25/03/2014 6:57:15 PM.

Your details are:

Username: Dentallabor Frankfurter

Dr Mark Homburger

Phone: 061 72 69 60 311

Email: DrH@PraxisHomburger.de

|

Laboratory data:

Dentallabor Frankfurter

Goethestrasse 3

60711 Frankfurt

Germany

Fax:

Website:

You can activate your account with the following [Activation Link](#)

Sincerely,

Your SIRONA Connect Team



CONFIRMATION OF ACTIVATION

The screenshot displays the SIRONA Connect website interface. At the top left, it says "SIRONA CORPORATE WEBSITE". A yellow notification box at the top center states "Your account has been activated successfully." The top right navigation bar includes links for "LABORATORY REGISTRATION", "DENTIST REGISTRATION", and "LOGIN". Below this is a main navigation bar with "ABOUT SIRONA CONNECT", "FOR DENTAL PRACTICES", "FOR DENTAL LABORATORIES", "SERVICE", and "CONTACT". The SIRONA logo, "The Dental Company", is on the right. The main content area features a "SIRONA Connect" section with a list: "Laboratory Registration", "Dentist Registration", and "Login". A breadcrumb trail shows "MY SIRONA CONNECT > LABORATORY REGISTRATION". A large "Welcome to SIRONA Connect" message is followed by two images of a man talking on a phone at a computer. A large orange arrow points to a yellow box labeled "Activation" with the text "Please log in and provide the additional data." The footer contains links for "SITEMAP", "TERMS & CONDITIONS", and "PRIVACY STATEMENT".

SIRONA CORPORATE WEBSITE

Your account has been activated successfully.

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SERVICE CONTACT

sirona
The Dental Company

SIRONA Connect

MY SIRONA CONNECT > LABORATORY REGISTRATION

- Laboratory Registration
- Dentist Registration
- Login

Welcome to SIRONA Connect

Activation
Please log in and provide the additional data.

SITEMAP | TERMS & CONDITIONS | PRIVACY STATEMENT

LOGIN AFTER REGISTRATION

SIRONA CORPORATE WEBSITE

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

sirona
The Dental Company

ABOUT SIRONA CONNECT | FOR DENTAL PRACTICES | FOR DENTAL LABORATORIES | SERVICE | CONTACT

Digital impression with Sirona Connect

New perspectives for dental practices and laboratories

Login



[Dentist Registration](#)

[Laboratory Registration](#)

[Forgot password?](#)

Login

FREE CEREC CONNECT SW 4.X
DOWNLOAD



Enter username and password in the login field

SIRONA Connect

- My Account
- User Administration
- Order List
- Log Off

MY SIRONA CONNECT > WELCOME

Welcome to SIRONA Connect



Welcome

Welcome **Ronny Kucharczyk** to Sirona Connect!

Your last visit was on 03/26/2014 3:31 PM


> You can edit your personal account on the following pages.

**Activation was
successful!**

ADD MORE ACCOUNT DATA

ADD MORE ACCOUNTS

- Laboratories can also add their **infiniDent** login data so that they can order SLA models directly from the Sirona Connect portal
- Just in Time via **Skype**: User data can be added so that dentists and dental technicians can consult with each other directly via Skype



The screenshot displays the 'My Account' page in the Sirona Connect portal. The left sidebar contains a menu with options: 'My user', 'My laboratory' (highlighted), 'My billing address', 'My shipping address', 'Change Password', 'Offerings / Services', 'Connect Direct', 'User Administration', 'Order List', and 'Log Off'. The main content area is titled 'Welcome to SIRONA Connect' and features two images of a man talking on a phone. Below the images, the 'My laboratory data' section contains input fields for 'Company Name *' (Dental Lab Frankfurter), 'Phone Number *' (0049 69 240524352), 'Fax Number', 'E-Mail Address' (info@DentalLabFrankfurter.com), and 'Website' (www.DentalLabFrankfurter.com). The 'My preferences' section includes dropdown menus for 'Preferred Language' (english (GB)), 'Timezone' (UTC+01:00 Amsterdam, Berlin...), 'List range' (30 Days), and 'infiniDent' login details (User Name: max.mustermann@sirona.com, Password: infinident). At the bottom, there are 'SAVE' and 'CANCEL' buttons, and a note stating 'marked fields are mandatory fields'.

My Account

- My user
- My laboratory
- My billing address
- My shipping address
- Change Password
- Offerings / Services
- Connect Direct
- User Administration
- Order List
- Log Off

Welcome to SIRONA Connect

My laboratory data

Company Name * Dental Lab Frankfurter E-Mail Address info@DentalLabFrankfurter.com

Phone Number * 0049 69 240524352 Website www.DentalLabFrankfurter.com

Fax Number

My preferences

Preferred Language english (GB) List range 30 Days

Timezone (UTC+01:00) Amsterdam, Berlin...

infiniDent

infiniDent User Name max.mustermann@sirona.com infinidEnt Password infinident

SAVE CANCEL

* marked fields are mandatory fields

HOW DOES CONNECT DIRECT WORK?

THE CONNECT DIRECT PRINCIPLE

- Exclusive customer loyalty: The intraoral scanner is exclusively connected to one laboratory, which means it is "locked" to this lab
- This means that the digital data from the practice can only be sent to one laboratory
- This takes place in the laboratory, with some manual processes carried out by Sirona
- Attractive business model offers new opportunities for laboratories and practices



LOGIN AS A LABORATORY

Access via login:
User ID +
password

My Sirona Connect [Home](#) > [My Sirona Connect](#) > [Login](#)

[▶ Laboratory Registration](#)
[▶ Dentist Registration](#)
[▶ Login](#)

Login

Please enter your **User ID** and **Password**.

User ID *

Password *

[Forgot your User ID or Password?](#)

ADD WORKSTATION INFORMATION

My Sirona Connect [Home](#) > [My Sirona Connect](#) > [Connect Direct](#)

Logged in as
lab


- ▶ [MyAccount](#)
 - ▶ Password
 - ▶ User
 - ▶ Laboratory
 - ▶ Billing Address
 - ▶ Shipping Address
 - ▶ Offerings / Services
 - ▶ **Connect Direct**
- ▶ [Order List](#)
- ▶ [Log Off](#)

Connect Direct

Please use the form below to add dentists to your Connect Direct customers.

Add Connect Direct workstation

Inventory Number *	<input type="text" value="000111222"/>
Workstation ID *	<input type="text" value="012345"/>
Dentist's contact person username *	<input type="text" value="dentist"/>



The image shows a Sirona product label. At the top is the Sirona logo. Below it is a barcode with the number 116212166512166. To the right of the barcode is a crossed-out symbol and a factory icon with the year 2011. Below the barcode, the text reads: Model-No.: 6212166 D3492, Serial.-No.: 12166. Further down are several certification logos: a triangle with 'DE', a 'C' in a circle, a 'CE' mark, a 'SP' in a circle with 'C' and 'US' below it, and a checkmark in a circle with 'N15657' below it. At the bottom, it says: Sirona Dental Systems GmbH, Fabrikstr. 31 D-64625 Bensheim, MADE IN GERMANY. A large yellow arrow points from the 'Inventory Number' field in the form to the barcode area on the label.

ADD WORKSTATION INFORMATION




Workstation was added.
More workstations can be added.

My Sirona Connect Home > My Sirona Connect > Connect Direct

Logged in as lab


- ▶ MyAccount
 - ▶ Password
 - ▶ User
 - ▶ Laboratory
 - ▶ Billing Address
 - ▶ Shipping Address
 - ▶ Offerings / Services
 - ▶ **Connect Direct**
- ▶ Order List
- ▶ Log Off

 Connect Direct successfully added.

Connect Direct

Please use the form below to add dentists to your Connect Direct customers.

Existing Connect Direct workstations

Inventory Number	Company Name	City	Activation	
000111222	Company Dentist	zip City	pending since 2/23/2012	Send activation e-Mail 

Add Connect Direct workstation


Inventory Number *

Workstation ID *

Dentist's contact person username *


[CANCEL](#) [ADD](#)

REQUEST TO ACCEPT CONNECT DIRECT



The dentist is requested to accept this direct connection


CONFIRMATION OF ACCEPTANCE




The dentist has accepted the connection and is sent an email containing the username and password

CONFIRMATION OF ACCEPTANCE


The dentist has accepted the connection and is sent an email containing the username and password




CONNECT DIRECT LOGIN



CONNECT DIRECT




THE DENTIST ACCEPTS THE CONNECTION



The laboratory receives confirmation that the dentist has accepted the connection.

THE DENTIST ACTIVATES THE CONNECTION



The laboratory receives confirmation that the dentist has accepted the connection.

CONNECT MORE WORKSTATIONS

My Sirona ConnectHome > My Sirona Connect > Connect Direct


Logged in as
lab

- ▶ MyAccount
 - ▶ Password
 - ▶ User
 - ▶ Laboratory
 - ▶ Billing Address
 - ▶ Shipping Address
 - ▶ Offerings / Services
 - ▶ **Connect Direct**
- ▶ Order List
- ▶ Log Off

Connect Direct

Please use the form below to add dentists to your Connect Direct customers.

Existing Connect Direct workstations

Inventory Number	Company Name	City	Activation	
000111222	Company Dentist	zip City	activated on 2/23/2012	

Add Connect Direct workstation

Inventory Number *

Workstation ID *

Dentist's contact person username *

CANCEL

ADD

The Dental Company

sirona.

INFORMATION ON THE MANUAL PROCESS WITH APOLLO DI

- The laboratory sends an email with all the information to Sirona at ronny.kucharczyk@sirona.com:
 - Usernames, addresses, email, etc. of the dentist and dental technician
 - Info on APOLLO DI (small number under the barcode)
- The dentist also has to log into the portal once with APOLLO DI.
- The process for APOLLO will be automated soon.




MORE INFORMATION ONLINE

THE SERVICE AREA


All information can be accessed online in the "Service" area:

- Free software updates/upgrades
- Marketing materials
- FAQ
- Videos
- Training sessions
- Tutorials




SOFTWARE UPDATES/UPGRADES

- All software updates and upgrades are available online free of charge




MARKETING MATERIALS

- Process video on digital impressions with Sirona Connect



Useful tutorials are available:

- APOLLO DI scanning strategy
- Clinical tutorials
 - How do you use Expasyl?
 - How does the double cord technique work?
 - How do I keep the mouth dry without using cotton rolls?
 - How can the assistant provide support during intraoral scanning?
 - How do I use the laser?



Sirona DI Hardware


APOLLO DI



GENERAL INFORMATION

- Use only original software or software which has been approved by Sirona
- Check restorations with trained personnel
- Check models with trained personnel
- Information on the location of the device or on the image brightness:
 - The image brightness is automatically controlled during imaging
 - The surroundings of the tooth to be imaged should be as weakly illuminated as possible: Switch off any operating lights and avoid sunlight due to reflections from the window

GENERAL INFORMATION

- Switch on the device using the main switch (right side at the top). This starts the operating system.
- Wait until the CodeMeter (symbol in the bottom right) turns green 
 - If the user attempts to start the software before the stick is activated, he will receive an error message that the license is not valid.
- Warm-up stage: Wait for approx. 15 min after starting up the system before beginning to scan. This allows the camera head to reach an optimum "operating temperature".
- Switching the device off at the ON button during operation may cause data loss and PC malfunctions.

WAIT UNTIL LICENSE STICK IS ACTIVATED



APOLLO DI



SET UP AN INTERNET CONNECTION

- APOLLO DI has a WLAN port
- Wireless network requirements
 - Practice network with internet connection
 - Recommended internet connection: Broadband speed, e.g., DSL, T1 or cellular phone network
- Touch the network symbol in the Windows bar
 - A list of available wireless networks is shown
 - Select the desired wireless network
 - If required, enter the password and click on OK.
Use the Windows keyboard to do this
- APOLLO DI is now connected to the wireless network!
- Microsoft Security Essentials is pre-installed.
Do not switch on firewalls or your own anti-virus programs



PREPARE THE CAMERA

- In two steps, the camera is cleaned (1) and disinfected (2).
- The camera can be removed from the device to facilitate cleaning and disinfection.
 - Remove the camera sleeve and wipe/disinfect the surfaces with the cloth
 - Dry the camera with a soft, clean cloth
 - Do not spray the camera with or immerse it in cleaning agents or disinfectants!
- Cannot be sterilized: do not sterilize the camera or video cable!
- Approved care and cleaning agents e.g., FD 312 from Dürr. Only use care and cleaning agents which have been approved by Sirona!

PREPARE THE CAMERA


- Camera sheets are used to protect the camera for each scan process
- Hold the sleeve firmly and insert the camera between the white flap and the paper. The optical side must face the right side
- Slowly insert the camera into the sleeve until the camera head is at the tip of the sleeve. Make sure that you do not exert too much pressure otherwise the camera sleeve may become damaged.
- Make sure that the optical path through the front window is not covered by the sleeve.
- Pull off the paper backing and carefully push the camera down so that the sleeve fits firmly.



APOLLO DI HYGIENIC PROTECTIVE SLEEVE

The hygienic protection sleeve must be pulled as far down as possible so that the optical path through the front window remains completely free.

Manufacturer	Product	ART. NO.	PRICE excl. VAT
TiDi®	Disposable sleeve 20 x 25 pcs	64 41 344	EUR 199



THE FIRST STEPS

○ START THE PROGRAM

- To start the Sirona Connect software, double click on the the Sirona Connect symbol

○ MULTI-TOUCH SCREEN

- You can use your finger to navigate and enter information
- You can open symbols by tapping them with your finger.



THE FIRST STEPS

○ NAVIGATE IN THE SOFTWARE

Example	Definition
Tap	Tap once to execute functions in the software
Call up context menus	Tap and hold the corresponding areas
Drag&Drop	Drag and drop

VIRTUAL KEYBOARD

Display the
keyboard by clicking
on the left edge

AP



ApolloDI

Recycle Bin Jennifer

www.greenfing.com

2:28 PM
1/21/2014

VIRTUAL KEYBOARD

APOLLO DI



APOLLO DI



Add Patient



Search Patient

Next

THE VIRTUAL KEYBOARD

- Activate it by touching the top left corner of the screen
- It also appears automatically when you click in a text field



- Useful key combinations

- Alt + Tab = Switches between various applications.
- Ctrl + Shift + Esc = Returns to the Desktop when APOLLO Connect software is open
- Windows logo + F = Opens Windows Task Manager
- Windows logo + F = Searches for a file on the hard disk
- Windows logo + D = Switch to desktop

Sirona DI Hardware

AC Connect





GENERAL INFORMATION

- Use only original software or software which has been approved by Sirona
- Check restorations with trained personnel
- Check models with trained personnel
- Information on the location of the device or on the image brightness:

HARDWARE BASIC KNOWLEDGE

Trackball

- To be used like a computer mouse
- Left button: tilt object
- middle: zoom
- Right button: move object



CAMERA HANDLING

- When switching on the system, the camera needs to warm up for 15 – 20 minutes. If the coated sapphire glass of the Omnicam is not sufficiently warm, it may steam up during the acquisition.
- As soon as the camera is pointed over a tooth or the gums, or the foot control is pressed, data acquisition begins.
- The distance must be between 0-15 mm (ideally: 5 mm). The camera does not rest on the teeth or the gums. If the distance is too great, no data will be obtained.

PRODUCT OPTIONS

CEREC AC Connect – Acquisition Center

- All user interface components are integrated centrally in a mobile unit.
- Scan and design made centrally on the AC.
- Optimum solution for most of the dental offices.



CEREC AF Connect – Acquisition Flexible

- Concept is flexible approach with camera mobility
- Scan and design can be flexible and separate from one another.
- This solution is suitable for group practices on different sites or on different floor levels and or practices with limited floor space.



PRODUCT OPTIONS

CEREC AI Connect – Acquisition Integrated

- All user interface components are integrated in TENEO dental chair
- CEREC Omnicam is put on a camera tray
- Tablet with Sirona remote SW



SIRONA REMOTE FOR CEREC

- The AC screen contents are transmitted to another monitor – on desktop, treatment center, a wall or a tablet
- Users work more ergonomically and have more freedom when placing the AC
- Customers in CEREC Club receive the software as part of their membership
- The AC screen contents can also be transmitted to a tablet (iOS and Android) -> CEREC AF – CEREC AI



GINGIVA MANAGEMENT

WHY GINGIVA AND LIQUID MANAGEMENT

- Making the preparation margins clearly visible
- Keep them free of liquids during scanning
- For adhesive technology dry preparations are required

USE OF OPTRAGATE

- Easier access
 - No additional instrument needed to keep away cheek
- Better visibility
- Higher patient comfort
 - Less lip irritations
- Three sizes: junior, small, regular
- Latex-free



Foto: Ivoclar Vivadent

LIQUID MANAGEMENT

During scanning with APOLLO DI it is recommended not to use cotton rolls or other movable parts which could create faulty scans.

Instead we recommend:

- DryTips® from Molnlycke
- DryTips® are glued to the Parotis and can absorb large amounts of liquid

Visit www.sirona-connect.de to see how the patient's mouth is kept dry




Photo: Molnlycke

LIQUID MANAGEMENT – DRYTIPS

- DryTips are available in 2 sizes
- They are applied with the rough surface to the cheek
- If an Optragate™ (Ivoclar) is used, the DryTips® are applied first
- Before scanning the scan area has to be blown dry with air. The assistant uses the syringe and mirror to keep the area clear of the tongue and saliva
- In case of anesthesia the gingival edge is supplied with less blood
- Remove the DryTips® by rinsing them with water

GINGIVA MANAGEMENT

There are mechanical, chemical and chirurgical methods to make the preparations margins visible.

- Retraction cords
- Expasyl®
- Gingivectomy by means of laser or electrotome

RETRACTION CORDS

- Large variety of retraction cords which can be used to open the sulcus
- Our recommendation: braided cords, as this does reduce defibration. When trimming, please consider sufficient length
- Optionally dip the retraction cord into retraction liquid
- Do not apply excess force when placing the cord to avoid damage to the sulcus
- 2-cord technique
- First apply the thin cord and then the thicker
- After placing the second cord wait some minutes until the hemostasis takes place
- After residence time the upper, thicker cord is removed cautiously. The lower, thinner remains in situ.
- Dry with air, powder and scan

EXPASYL®

- Expasyl® (Acteon, France) uses a combination of mechanical and chemical retraction
- It is achieved by the expansion of the material as well as due to the hemostatic properties of aluminium-chloride
- The paste is applied sub- and paragingival to the preparation with a blunt canulla into the sulcus
- After 2-3 minuten residence time the remains are thoroughly rinsed with water




Foto: [dental users](#)

GINGIVECTOMY WITH LASER

Visit www.sirona-connect.com
to see how to use the
SIROLaser

- In case of surgical methods the gingiva is not only ousted temporarily, but removed permanently
- It is common to use the Elektrotom / laser in addition to the before mentioned retraction methods
- The advantage of electro surgery is that cutting and hemostasis are one process
- We recommend for instance the SIROLaser by Sirona
<http://www.sirona.com/en/products/instruments/sirolaser-advance-xtend/?tab=309>

LASER APPLICATIONS



The mesial papilla and a small overlap in the distolabial area made it difficult to get an optimal digital impression.



Situation after removing a part of the mesial papilla, the circular chamfer is readily visible

Digital impression with CEREC



Source: Laser – international magazine of laser dentistry – SIROLaser Factbook; “Use of a Diode Laser for Gingival Troughing in Conservative and Prosthetic Dentistry”, Dr Almut Marsch

CEREC / PROSTHETICS

- Cosmetic Applications
- Desensitization
- Pain Therapy

LASER APPLICATIONS



Hemostasis: before treatment



Hemostasis: directly after treatment with the diode laser

Source: Dr. med. dent. Helmut Götte

CEREC / PROSTHETICS

-  Cosmetic Applications
-  Desensitization
-  Pain Therapy

Scan with Apollo DI



ACQUISITION TECHNIQUE

- Spray-based, flowing acquisition technique
- The technology that forms the basis of APOLLO DI measures the structure in the spray. The camera does not work without the spray.
- This means that fine black and white particles in the spray ensure high contrast and thus very precise images
- Do not use any cotton rolls or other movable parts when scanning
- Turn the surgery lamp away when scanning

APOLLO DI SPEEDSPRAY

- APOLLO DI SpeedSpray is used as the contrast agent
 - Place the spray head with cannula/nozzle onto the spray can
 - Check that the cannula/nozzle is seated correctly before each use by pulling it gently
 - Shake the container before use
 - Hold the can upright while spraying
 - The spray nozzle should be held approx. 10-15 mm away from the object



Product

APOLLO DI SpeedSpray

Reference

64 14 572

SPRAYING DOES NOT MEAN POWDERING!


- Even if the high-contrast spray is insensitive to moisture, there must not be any blood or saliva in the imaging area (see chapter on soft tissue management)
- After suctioning, the spray is applied to the teeth in a very thin layer. The rule is to apply as much as necessary but as little as possible.
- The aim is not to achieve 100% coverage but rather to apply the pattern of black and white particles evenly.
- If not enough spray has been applied, "the scanning process comes to a halt" and no data is collected.



Source: Dr. Martin Butz, Munich

SPRAYING DOES NOT MEAN POWDERING!

- The tooth color still shines through.
- Shake the can well and apply the spray to the teeth at a distance of 10-15 mm
- If the distance is less, the spray will be unpleasantly cold for the patient. Spray more often and short.
- The assistant should apply suction while the spray is being applied inter orally and after application. The assistant can quickly spray air over the surface so that no unnecessary coatings adhere.
- ☒ Enough spray: good scanning performance and results
- ☐ Not Enough spray: bad scanning performance



Source: Dr. Martin Butz, Munich


HOW DO I USE THE CAMERA?

- Press the black button to activate the LED and the live window in the Acquisition phase
- Enter the mouth, bring the camera into position of an occlusal surface and press the white button to start the scan
- You can pause the scan and switch off the LEDs by pressing the black button
- Pressing the white button will finalize this scan. A new scan (also known as “cluster”) can be added at any time
- The camera is calibrated before exiting the factory (ex works). Calibration on site is not necessary



DETAILED SCAN STRATEGY FULL JAW

- Scan the occlusal surface from distal to mesial, e.g. from 48 to 33.
- When going over the front, please make sure to scan in a way that not only the incisal is captured but also parts from labial and lingual.
- Next tilt the camera about 45 degrees and scan the lingual (or palatine) surfaces from 33 to 48. Scan in a way that your occlusal scan connects to the scan of the lingual side
- At 48 move the camera to buccal and scan the buccal/ vestibular side in way that you scan connects to your occlusal scan
- Stitch your scan by scanning 43 and 33 once again from vestibular to lingual
- Now check the model and fill the missing relevant areas
- All this is done in one image catalogue (cluster). Use a new cluster only if you need to re-spray during the scan
- Proceed the same way from 38 to 43 in a new scan (cluster)




Cluster

- Each cluster is optimized in itself, so try to scan as much as possible into a single cluster
- For the highest possible accuracy do not scan the prep in several clusters because there might be a loss in quality due to the stitching of the images
- In order to help the system please wait until the cluster is finalized (red or green state)

THE "DIFFICULT AREAS"

- Move the camera head up and down to make the most of the camera's scanning area
- With sufficient space, the camera can also move diagonally in the free area around the preparation. This movement helps to scan the proximal contacts
- Swivel the camera up to 90 degrees to view all areas




Scan with Omnica




CAMERA HANDLING

- Divide the acquisition into four consecutive sequences:
- Occlusal:
 - ✓ Move the CEREC Omnicam to the starting position. For this purpose, the CEREC Omnicam is in the occlusal view of the tooth, which is next to the prepared tooth in the distal direction.
 - ✓ Scan in the mesial direction. To do so, slowly move the CEREC Omnicam in the occlusal direction from the distal-positioned tooth over the prepared tooth to the mesial-positioned tooth.




CAMERA HANDLING

- Divide the acquisition into four consecutive sequences:
- Buccal:
 - ✓ Rotate the CEREC Omnicam between 45° to maximum 90° toward the buccal and guide it over the entire buccal distance in the distal direction over the prepared tooth.
 - ✓ Ensure that the CEREC Omnicam is held like a flute during buccal scans. Do not tilt it vertically to the direction of motion..



CAMERA HANDLING

- Divide the acquisition into four consecutive sequences:
- Lingual:
 - ✓ Rotate the CEREC Omnicam from 90° in the buccal direction to around 45° to maximum 90° in the lingual direction on the other side.
 - ✓ Guide the CEREC Omnicam over the entire lingual distance in the mesial direction over the prepared tooth.




CAMERA HANDLING

○ Divide the acquisition into four consecutive sequences:

○ Proximal:


- ✓ Move the CEREC Omnicam in the occlusal direction to the prepared tooth.
- ✓ Acquire the approximal surfaces in the distal and mesial direction by using a wave motion in the occlusal, buccal, and lingual direction over the prepared tooth.
- ✓ To do so, tilt the surface by 15° in the distal and mesial direction to gain a better view of the approximal contacts.



CAMERA HANDLING – QUADRANT OR FULL JAW

The first (fourth) quadrant is scanned up to the opposite second front tooth by moving the camera in parallel along the jaw arch.

Start the scanning process




➤ Position the camera occlusally above the last tooth on the right, to start the scanning process.

QUADRANT OR FULL JAW

The first (fourth) quadrant is scanned up to the opposite second front tooth by moving the camera in parallel along the jaw arch.


Start the scanning process



➤ Position the camera occlusally above the last tooth on the right, to start the scanning process.

QUADRANT OR FULL JAW

Completing the scanning process



1. Start as indicated above, on the occlusal surface of the right terminal tooth, and scan it occlusally.
Tilt the camera by 45° in a palatal direction (oral) and guide it from the distal to the mesial.
2. Tilt the camera another 45° in a palatal direction (oral) and move it in a distal direction.
3. Tilt the camera by 90° on to the occlusal surface and move it in a mesial direction.
4. Tilt the camera in a 45° buccal direction and move it back towards the distal.
5. Then tilt the camera a further 45° in a buccal direction to a total of 90° and move it in a mesial direction again.

TIPS & TRICKS


Scanning



PRACTICE TIPS & TRICKS

Where should the intraoral scanner ideally be positioned?

The dentist works with an assistant:




- The dentist sits in the 8-9 o'clock position
- The scanner is in the 11 o'clock position, level with the headrest
- The assistant sits on the assistant side and can operate the suction instruments efficiently

PRACTICE TIPS & TRICKS

Where should the intraoral scanner ideally be positioned?

The dentist works without the assistant:



- The dentist sits in the 12 o'clock position (behind the head)
- The scanner is on the assistant side
- The dentist can also operate the suction instruments himself

PRACTICE TIPS & TRICKS

- What actually can be delegated to the assistant?
 - Preparation: Add patient data, prepare the spray and camera, switch on the scanner, prepare the instruments
 - During spraying: Suction and blow air on the area after applying the spray
 - During scanning: Use the mirror to keep the tongue, cheeks and lips away, suction while scanning, clean the camera if necessary, make sure that the mirror or suction device are not in the scanning field

PRACTICE TIPS & TRICKS

- Do not use cotton rolls in the scan area as a rule.
- What happens with patients with extreme saliva flow? There are extreme cases in which cotton roll is used, such as in the lower jaw for example.
- The cotton roll is placed under the tongue and kept away with the mirror so that it does not appear in the image. It must not be visible!



PRACTICE TIPS & TRICKS

The patient can help during the scanning process by opening and closing his mouth. He has to open his mouth as wide as possible in the distal area on the occlusion. If the vestibular area is scanned, it can help if the mouth is slightly closed as this makes the cheeks less tense.

THE FIRST STEPS

- Like everything else, digital impressions require some learning
- Before working on patient cases, scanning should be tested on colleagues/assistants first.
- For the first three cases, a conventional impression should be taken in addition to the digital impression and both should be sent to the laboratory.




THE FIRST STEPS

- The dentist and laboratory consult with each other.
- The laboratory can connect to intra-oral scanner directly (using commercially available programs, such as TeamViewer), as long as the patient is still in the practice.
- The patient will not be sent away until the laboratory has validated the data.




BUCCAL REGISTRATION

- Always include some gingiva in the scan to achieve good automatic buccal registration
- Swivel the camera during scanning in order to receive additional 3D information.
- 2-3 teeth are sufficient
- If tongues and cheeks that have been scanned interfere with the image, you can remove them using the "correct" tool
- Try to make the buccal close to the preparation



Poor



Good

WHOLE JAW VERSUS QUADRANT

- It's up to the clinician and the dental laboratory how much is scanned for each particular case, below please find a rather general guideline
 - A whole jaw scan is always necessary for
 - wide spanning bridges
 - difficult occlusal conditions
 - canine teeth if there is canine guidance
 - A 3/4 scan that includes the two canine teeth is indicated for posterior teeth restorations and canine guidance
 - A quadrant or anterior tooth scan is usually sufficient for
 - incisors if there are adjacent teeth
 - single tooth restorations in the posterior tooth region if there is intact occlusal group guidance and if there are adjacent teeth

HOW MANY HOLES ARE ALLOWED IN THE OPTICAL IMPRESSION?

- This means that the preparation, adjacent teeth and antagonists must be completely scanned
- It is particularly important that all functional areas related to the restoration to be made are scanned in with the greatest care
- The software fills holes automatically, but make sure to capture the relevant areas

HOW MANY HOLES ARE ALLOWED IN THE OPTICAL IMPRESSION?

- These are in particular:
 - preparation margins
 - proximal areas with a contact to the restoration to be made
 - occlusal surfaces of the adjacent teeth and on the prepared tooth (for inlays and onlays)
 - labial surfaces of the adjacent teeth for anterior tooth restorations
- It must also be ensured that functional areas are not covered by artifacts (gingiva, tongue, suction device, etc.)

SIRONA Connect SW STEP-BY-STEP




THE 4-STAGE WORKFLOW



- Enter patient information
- Enter case information
- Generate 3D Data
- Evaluate 3D model
- Mark margin
- Complete order information
- Send order to the laboratory


Sirona Connect SW 4.4

- Sirona Connect SW 4.4 onwards the user will be informed about a new version, as soon as the application started and is connected to the internet
- The user can decide if he wants to download right now or later




Sirona Connect Center

- The new version will be provided by the Sirona Connect Center
- By clicking on the symbol the user can check the progress of the download
- Also CEREC SW and inLab SW will use this Sirona Connect Center




Installation

- After a successful download the user can install the Software
- All applications should be closed prior running the installation



Installation

- Installation runs
- Please read all disclaimers and confirm them



Start Screen




Search for patients



Display patient data base



Add new patient












Patient Data Base

- Patient data base is built up from left to right starting with „Patient“, „Cases“ in the center and „Preview“ on the right hand side
- The single columns can be moved (for example name and restoration type) and the size of the single windows can be adapted




User Interface

- The same symbols are used throughout the Sirona Connect workflow


	Edit
	Delete
	Save
	Cancel
	Open case
	Export case (with log files)
	Add new patient
	Add new case
	Add new 2D image with Omnicam




ADMINISTRATION



ACQUISITION



MODEL



CONNECT

Sirona
Connect

Sirona Connect



Search Patient Database



Show All Patients




Add New Patient



Patient Database Sirona Connect Portal

Add New Patient

- Create a new patient via entering the name and date of birth
- Due to data security we recommend to enter a patient ID only




The screenshot displays the Sirona software interface with a blue background. At the top, there is a navigation bar with icons for 'Patienten', 'Behandlung', 'Zahlung', and 'Einstellungen', along with the 'Sirona' logo. The main content area features a 'Neuer Patient' (New Patient) form. The form has a title bar and a sidebar with icons for 'Patienten', 'Behandlung', and 'Zahlung'. The form fields are as follows:

- Nachname:** A text field with a yellow highlight.
- Vorname:** A text field.
- Geburtsdatum:** A date field with a calendar icon.
- Patienten-ID:** A text field.
- Zahnarzt:** A text field.

At the bottom of the form, there are three buttons: 'Zurück' (Back), 'Neuer Patient' (New Patient), and 'Neuer Fall' (New Case). Below the buttons, there is a small icon of a person with a green checkmark and the text 'Neuer Patient'.




ADMINISTRATION



ACQUISITION



MODEL



CONNECT

Sirona
Connect

Sirona Connect



Search Patient Database



Show All Patients




Add New Patient




Patient Database Sirona Connect Portal

Add New Case

- First step is to choose a bridge or single tooth restoration and the design mode (full contour or veneered)
- Restorations can be added in both jaws





Add New Case

- Maryland bridges as well as a mix of veneered and full contour are possible
- All windows can be positioned free on the screen via drag and drop
- The active indication is highlighted in orange



Implant

- A Sirona scan body is needed for the Sirona implant system
- ScanPost was developed the intraoral scan (tiBase are smaller)
- Choose your implant system



IMPLANT SYSTEMS IN THE SIRONA CAD/CAM WORLD

Implantatsysteme		Sirona Komponente	
Implantathersteller	Implantatsystem	TiBasis Set	inCoris Zi meso Größe
Nobel Biocare	Replace® NP 3,5	NB RS 3,5 L	L
	Replace® RP 4,3	NB RS 4,3 L	
	Replace® WP 5,0	NB RS 5,0 L	
	Replace® 6,0	NB RS 6,0 L	
Nobel Biocare	Brånemark® 3,4	NB B 3,4 L	L
	Brånemark® 4,1	NB B 4,1 L	
Nobel Biocare	Nobel Active® 4,5	NB NA 4,5	L
	Nobel Active® 5,0	NB NA 5,0	
Straumann	Tissulevel® NH 3,5	S SO 3,5 L	L
	Tissulevel® RH 4,8	S SO 4,8 L	
	Tissulevel® WH 6,5	S SO 6,5 L	
Straumann	Bone Level® NH 3,3	S BL 3,3	L
	Bone Level® NH 4,1	S BL 4,1	
Astra Tech	Osseospec™ 3,5 S / 4,0 S	AT OS 3,5/4,0 L	L
	Osseospec™ 4,5 / 5,0	AT OS 4,5/5,0 L	
Friadent	FX480 / X480® 3,4	FX 3,4 S	S
	FX480 / X480® 3,6	FX 3,6 S	
	FX480 / X480® 4,5	FX 4,5 L	L
	FX480 / X480® 5,5	FX 5,5 L	
Biomet 3i	Ex. Hex® 3,4	BO 3,4 L	L
	Ex. Hex® 4,1	BO 4,1 L	
	Ex. Hex® 5,0	BO 5,0 L	
Biomet 3i	Cerian® 3,4	BC 3,4	
	Cerian® 4,1	BC 4,1	
	Cerian® 5,0	BC 5,0	
Zimmer	Tapered Screw-Vent® 3,6	ZTSV 3,6 L	L
	Tapered Screw-Vent® 4,5	ZTSV 4,5 L	
	Tapered Screw-Vent® 5,7	ZTSV 5,7 L	
Mendeliv® Implant	M Implant Tapered	MB 3,5/5,0	L
CAMLOG	CAMLOG® 3,5 / 3,6 / 4,3	—	E
	CAMLOG® 5,0 / 6,0	—	L
olytech®	Tube Line 3,4 / 3,6 / 5,0 (A3 P-Switch)	—	S
	Tube Line 4,3 / 5,0	—	L

* Für das CAMLOG Implantatsystem sind die TiBasisen ausschließlich bei ihrer CAMLOG Landesvertretung erhältlich. Die Scanbody sind bei Ihrem Fachhändler bestellbar.

** TiBasis ist! Scanbody ist ausschließlich beim Originalhersteller erhältlich.


Material and Color

- All mandatory steps are underlined in red, green means „done“ and no color means an optional step
- If a special material is requested choose “see additional notes” and specify it in the free text field in the Sirona Connect Portal




ACQUISITION PHASE

- Additional image catalogues are available for example to scan the pre-op situation or the gingiva or a reference tooth
- Size of the scan window can be adapted




ACQUISITION PHASE – APOLLO DI

- By clicking the green light in the middle of the screen all scanned catalogues become visible
- By drag and drop catalogues can be moved or deleted




ACQUISITION PHASE – APOLLO DI

- A blue light on the bottom of the screen means that the software tries to correlate image catalogues




ACQUISITION PHASE – APOLLO DI

- If a scan can't be connected it will be shown in red color. An exclamation mark is shown in the catalogue
- In such a case cut the disturbing areas away or scan a „bridge“ to have sufficient information to calculate the scans together



ACQUISITION PHASE – APOLLO DI


- To cut areas away, please activate the tool cut
- By clicking you set the start point of a line. Determine the curve of the line by more single clicks (no drawing as before)
- „Apply“ to end the cut process





ACQUISITION PHASE – APOLLO DI

- To cut areas away, please activate the tool cut
- By clicking you set the start point of a line. Determine the curve of the line by more single clicks (no drawing as before)
- „Apply“ to end the cut process



MODEL PHASE – BUCCAL BITE

- Buccal registration is automated for Apollo and Omnicam
- If it doesn't work the buccal scan can be correlated manually
- Click the buccal scan and move it by drag and drop
- Click the bottom and move it up, so the buccal scan turns upside down
- Click the top and move it up, so the buccal scan remains in the position



MODEL PHASE – CONTACT POINTS

- Show the contact points under tools
- If the bite is not correct, please go back and rescan the buccal




MODEL PHASE – MARGIN APOLLO DI

- Turn the model into position – margin marking tool is not activated
- Touch the screen on the model and slide your finger to the desired start point on the preparation
- Mark the margin by single clicks
- Corrections will be made in the manual mode




MODEL PHASE – EVALUATION TOOLS

- With a marked margin a new step is open
- In this step you can show undercuts of the preparation and the distance to the antagonist




MODEL PHASE – DISTANCE TO ANTAGONIST

- Wit this tool the user can show distances in different colors
- How the colors are shown can be set in the parameters



PARAMETER PREPARATION ANALYSIS

- Example: Set the ideal distance on 1500 μm and the tolerance on 500 μm
- Distances smaller than 1500 μm will be shown in red
- 1500 μm till 2000 μm in blue
- Larger distances will be shown in green color





Tools






Manual Trim



Analyzing Tools

Model
Contacts

Color Model




Full Model



Full Model


SIRONA CONNECT PORTAL

- Check the case once again
- Very fast upload with Sirona Connect SW 4.4




SIRONA CONNECT PORTAL

- Add additional information
- Mandatory fields are marked in red
- Additional text can be entered and files uploaded



SKYPE


- To use the new skype option the user needs to enter the skype name in the Sirona Connect portal



The screenshot shows a registration form for Skype within the Sirona Connect portal. The form includes the following fields and labels:

- Vorname ***: Input field containing "Ronny".
- Nachname ***: Input field containing "Kucharczyk".
- E-Mail-Adresse bestätigen ***: Input field containing "ronny.kucharczyk@sirona.com".
- Skype-Benutzername**: Input field containing "ronny.kucharczyk".
- Chat-Avatar**: A placeholder icon for a profile picture.
- Durchsuchen...**: A button to search for the user.

The Skype logo is displayed below the form fields.

- Click the speech bubble in the upper right side of the screen of the inLab SW 15.0 or Sirona Connect SW 
- A chat window will open that shows all chats and files send between dentist and dental technician

CONFIRMATION OF SENT DATA BY EMAIL




SIRONA CONNECT APP

- Just in Time consultation between the practice and laboratory
- Access to order administration for dentists and dental technicians
- Assessment of the digital model data by the dental technician
- Available to download free of charge from the Apple App Store




SIRONA CONNECT APP

- Evaluation and control of the quality of the scan data
- Margins are visible if drawn
- Picture of the patient can be seen



ORDER CAN BE DOWNLOADED AS A PDF


SIRONA Connect


Dentist:
Haselbauer - Praxis
Fabrikstr. 31
64625 Bensheim
Tel.: 06251-16-2078

Patient:
Gender: F
Age: 34
Insurance: GKV

Order Number:
232958

Date Due:
26/02/2014 10:00 AM



CC = Crown veneered | PO = Pontic | CM = Missing | C = Crown | I/O = Inlay
V = Veneer | A = Abutment | * = Margin line defined

Type:	Tooth Number:	Indication:	Material Class:	Material:	Stump Color:	Final Color:
Crown veneered	38	Crown veneered	See Additional Notes		B1	A2
Bridge veneered	31	Crown veneered	Zirconium Oxide		B1	A2
Bridge veneered	32	Pontic	Zirconium Oxide		B1	A2
Bridge veneered	33	Crown veneered	Zirconium Oxide		B1	A2

AC: Arges

Dentist Comments:
inCoris TZI

Lab Comments:
-

DESIGN SERVICE FOR CEREC USERS

- Starting from inLab SW 15.0 inLab users can offer a design service for CEREC users
- The dentist scans in the Sirona Connect SW and sends the file to the lab via the Sirona Connect Portal
- The lab designs in the inLab SW and sends a .dxd file incl. restoration back to the dentist
- The dentist opens the case in the CEREC SW 4.4 and mills it on his CEREC milling machine

ENJOY EVERY DAY. WITH SIRONA.



Treatment
centers



Hygiene
systems



CAD/CAM
systems



Imaging
systems



Instruments