Symbios®

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<sup>1</sup> Schopper C, Ziya-Ghazvini F, Goriwoda W, et al. HA/TCP compounding of a porous CaP biomaterial improves bone formation and scaffold degradation – a long-term histological study. *J Biomed Mater Res B Appl Biomater* 2005;74(1):458–67.

<sup>2</sup> Spassova E, Gintenreiter S, Halwax E, et al. Chemistry, ultrastructure and porosity of monophasic and biphasic bone forming materials derived from marine algae. *Mat.-wiss. u Werkstofftech.* 2007;38(12):1–8.

<sup>3</sup> Zhou AJ, Clokie CM, Peel SA. Bone formation in algae-derived and synthetic calcium phosphates with or without poloxamer. *J Craniofac Surg* 2013;24(2):354–9.

<sup>4</sup> Zhou AJ, Peel SA, Clokie CM. An evaluation of hydroxyapatite and biphasic calcium phosphate in combination with Pluronic F127 and BMP on bone repair. *J Craniofac Surg* 2007;18(6):1264–75.

<sup>5</sup> Schopper C, Moser D, Spassova E, et al. Bone regeneration using a naturally grown HA/TCP carrier loaded with rhBMP-2 is independent of barrier-membrane effects. *J Biomed Mater Res A* 2008;85(4):954–63.

### Frios Algipore

<sup>6</sup> Simunek A, Cierny M, Kopecka D, et al. The sinus lift with phycogenic bone substitute. *Clin Oral Implants Res* 2005;16(3):342–48.

<sup>7</sup> Schopper C, Moser D, Sabbas A, et al. The fluorohydroxyapatite (fha) Frios algipore is a suitable biomaterial for the reconstruction of severely atrophic human maxillae. *Clin Oral Implants Res* 2003;14(6):743–9.

<sup>8</sup> Schopper C, Moser D, Wanschitz F, et al. Histomorphologic findings on human bone samples six months after bone augmentation of the maxillary sinus with algipore. *J Long Term Eff Med Implants* 1999;9(3):203–13.

<sup>9</sup> Schopper C, Ziya-Ghazvini F, Goriwoda W, et al. Ha/tcp compounding of a porous cap biomaterial improves bone formation and scaffold degradation--a long-term histological study. *J Biomed Mater Res B Appl Biomater* 2005;74(1):458–67.

<sup>10</sup> Spassova E, Gintenreiter S, Halwax E, et al. Chemistry, ultrastructure and porosity of monophasic and biphasic bone forming materials derived from marine algae. *Mat.-wiss. u Werkstofftech.* 2007;38(12):1–8.

<sup>11</sup> Ewers R. Maxilla sinus grafting with marine algae derived bone forming material: A clinical report of long-term results. *J Oral Maxillofac Surg* 2005;63(12):1712–23.

<sup>12</sup> Scarano A, Degidi M, Perrotti V, Piattelli A, Iezzi G. Sinus augmentation with phycogene hydroxyapatite: Histological and histomorphometrical results after 6 months in humans. A case series. *Oral Maxillofac Surg* 2012;16(1):41–5.

<sup>13</sup> Neugebauer J, Iezzi G, Perrotti V, et al. Experimental immediate loading of dental implants in conjunction with grafting procedures. *J Biomed Mater Res B Appl Biomater* 2009;91(2):604–12.

## Membranes

### Membrane en collagène Symbios Collagen Membrane SR

<sup>19</sup> Behfarnia P, Khorasani MM, Birang R, Abbas FM. Histological and histomorphometric analysis of animal experimental dehiscence defect treated with three bio absorbable GTR collagen membrane. *Dent Res J (Isfahan)* 2012;9(5):574–81.

<sup>20</sup> Guda T, Walker JA, Singleton BM, et al. Guided bone regeneration in long-bone defects with a structural hydroxyapatite graft and collagen membrane. *Tissue Eng Part A* 2013;19(17–18):1879–88.

<sup>21</sup> Tovar N, Jimbo R, Gangolli R, et al. Evaluation of bone response to various anorganic bovine bone xenografts: an experimental calvaria defect study. *Int J Oral Maxillofac Surg* 2014;43(2):251–60.

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### Frios BoneShield

<sup>22</sup> Butz F, Bachle M, Ofer M, Marquardt K, Kohal RJ. Sinus augmentation with bovine hydroxyapatite/synthetic peptide in a sodium hyaluronate carrier (pepgen p-15 putty): A clinical investigation of different healing times. *Int J Oral Maxillofac Implants* 2011;26(6):1317–23.

<sup>23</sup> Kohal RJ, Straub LM, Wolkewitz M, Bachle M, Patzelt SB. Evaluation of guided bone regeneration around oral implants over different healing times using two different bovine bone materials: A randomized, controlled clinical and histological investigation. *Clin Implant Dent Relat Res* 2014;E-pub May 13, doi: 10.1111/cid.12206.

## Instruments

### Frios MicroSaw

<sup>24</sup> De Stavola L, Tunkel J. Results of vertical bone augmentation with autogenous bone block grafts and the tunnel technique: A clinical prospective study of 10 consecutively treated patients. *Int J Periodontics Restorative Dent* 2013;33(5):651–9.

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<sup>14</sup> Li ST, Chen HC, Yuen D. Isolation and Characterization of a Porous Carbonate Apatite From Porcine Cancellous Bone. *Science, Technology, Innovation*, Aug. 2014:1–13. Deligianni DD, Katsala ND, Koutsoukos PG, Missirlis YF, Effect of Surface Roughness of Hydroxyapatite on Human Bone Marrow Cell Adhesion, Proliferation, Differentiation and Detachment Strength. *Elsevier Biomaterials* 22 (2001) 87–96.

<sup>15</sup> Spense G., Patel N., Brooks R., Rushton N. 2009. Carbonate Substituted Hydroxyapatite: Resorption by Osteoclasts Modifies the Osteoblastic Response. *Journal of Biomedical Materials Research Part A* 217–224.

<sup>16</sup> Ellies LG, Carter JM, Natiella JR, Featherstone JDB, Nelson DGA. 1988. Quantitative Analysis of Early In Vivo Tissue Response to Synthetic Apatite Implants. *J Biomed Mater Res* 22:137–148.

<sup>17</sup> Landi E., Celotti G., Logroscino G., Tampieri A. 2003. Carbonated Hydroxyapatite as Bone Substitute. *Journal of the European Ceramic Society* 23: 2931–2937.

<sup>18</sup> Frank M. Klenke, Yuelian Liu, Huipin Yuan, Ernst B. Hunziker, Klaus A. Siebenrock, Willy Hofstetter. Impact of pore size on the vascularization and osseointegration of ceramic bone substitutes in vivo. *Journal of Biomedical Materials Research Part A*, 2007, 777–786

<sup>25</sup> De Stavola L, Tunkel J. A new approach to maintenance of regenerated autogenous bone volume: Delayed relining with xenograft and resorbable membrane. *Int J Oral Maxillofac Implants* 2013;28(4):1062–7.

<sup>26</sup> Khoury F, Hanser T. Mandibular bone block harvesting from the retromolar region: A 10-year prospective clinical study. *Int J Oral Maxillofac Implants* 2015;30(3):688–397.

### Clous de fixation Frios

<sup>27</sup> Butz F, Bachle M, Ofer M, Marquardt K, Kohal RJ. Sinus augmentation with bovine hydroxyapatite/synthetic peptide in a sodium hyaluronate carrier (pepgen p-15 putty): A clinical investigation of different healing times. *Int J Oral Maxillofac Implants* 2011;26(6):1317–23.

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Données non publiées, Collagen Matrix, Inc.