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1. Description

- 1.1 INTENDED USE
- 1.2 INDICATION
- 1.3 CONTRAINDICATION
- 1.4 ADVERSE EFFECTS
- 1.5 PATIENT TARGET GROUP
- 1.6 INTENDED USERS
- 1.7 SUITABLE FOR DENTSPLY-SIRONA-PRIMEPRINT-SOLUTION
- 1.8 SPECIAL MANUFACTURING REQUIREMENTS
- 1.9 SAFETY INFORMATION
- 1.10 STORAGE
- 1.11 MANUFACTURING PROCESS
- 1.12 FINISHING
- 1.13 SURFACE PROCESSING

1.1 Intended use

Resin for dental 3D printing

1.2 Indication

Indicated for the fabrication of temporary dental restorations in conjunction with extra-oral curing light equipment. Duration is less than 30 days in oral environment.

1.3 Contraindication

Contains (meth)acrylics and phosphine oxides. Some ingredients of Primeprint Temp cause allergic reactions in predisposed persons. In such cases refrain from using the product. Primeprint Temp only insert intraorally in completely polymerised state.

1.4 Adverse effects

Product may cause allergic reactions.

1.5 Patient target group

Persons being treated in the context of a dental procedure.

1.6 Intended users

Dentist, dental technician

1.7 Suitable for Dentsply-Sirona-Primeprint-Solution

Primeprint Temp is suitable for the 385 nm DLP printer Primeprint and the Primeprint PPU finishing unit by Dentsply Sirona.

1.8 Special manufacturing requirements

Printing machine and post-processing unit should be set-up, validated and maintained according to their labelling and instructions for use.

1.8.1 File format

blc (based on stl)


1.8.2 DLP printer and its operation software

Printer Manufacturer	SDS Primeprint
Light source	385 nm
Pixel size x, y	±70 µm
Light intensity	6,2 mW/cm ²
Operation Software	inLab CAM SW 22

Technical Data:
ISO 10477-Type 2, Class 2

Processing:
at 23 °C ± 2 °C

Storage:

Consult Instructions for Use

Catalogue Number

Batch Code

Use by Date

Rx only Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

Temperature limitation

Keep away from sunlight

Manufacturer

Medical Device

Warning

Warning

1.8.3 Printing Parameters

Parameter data set*	Primeprint Temp
Curing time	1,9 s
Height base plate	0,5 mm
Support point size Max. support size Support density	0,9 mm (tip) 1,5 mm 50 %
Connection between supports	diagonal struts
Slicing	50 µm
Build speed	30 mm/hr
Build path	surface layer drawing path
Recommended orientation angle (degree)	Printing mode optimized quality

*The set of parameter include all relevant material, printer and post-processing specific information

1.8.4 Note on processing

- ❑ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is very important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the chapter "Manufacturing process [1.11]").
- ❑ After storage, intensively shake and homogenize the material in the cartridge before use.
- ❑ Before inserting the cartridge into the material unit, check the same for integrity.
- ❑ Remove the seal label before inserting the cartridge into the material unit.
- ❑ Maximum curing depth* at direct post-exposure: 3 mm *In case of large objects and exposure on both sides, the material thickness can be up to 6 mm (Example Primeprint Temp – with a curing depth of 3 mm). Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- ❑ Primeprint Temp is applied for the production of temporary c & b in the anterior and posterior area and can be used for single crowns and bridges up to 4 elements containing one pontic. The occlusal minimum thickness amounts to 1.5 mm, circularly up to 0.8 mm. The connector cross-section for pontics in the anterior area is 12 mm² minimum, in the posterior area 15 mm². Any non-compliance of these requirements can lead to an undesirable outcome.
- ❑ 25 and 50 µm are the layer thickness for the printing process.
- ❑ Polish the surface mechanically. Prepolish by means of rotating brushes and prepolishing paste, highgloss by means of buffing wheels and highshine polishing paste for resins.

1.8.5 Environmental Conditions

Processing temperature 23 °C ± 2 °C.

1.8.6 Cleaning Kit

Rinse bath tubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, scraper.

1.8.7 Curing light equipment

SDS Primeprint PPU.

1.9 SAFETY INFORMATION

- ❑ Only for the specified intended use by trained specialists.
- ❑ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- ❑ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ❑ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ❑ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ❑ After contact with skin wash immediately with water and soap.
- ❑ Biocompatibility is only guaranteed with complete polymerisation.
- ❑ Refer to the relevant safety data sheet for hazard and safety information.

1.10 STORAGE

Primeprint Temp is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

1.11 Manufacturing process

Refer to the Primeprint operating instructions (REF 6745173) by Dentsply Sirona for the insertion, setup and correct handling of the Primeprint cartridge. Data processing and creation of the support structure according to the information of the Dentsply Sirona software.


Mixing

Before inserting the cartridge into the material unit the material must be shaken. Inadequate mixing could cause deviations of colour and print failures.

Prepare the printer

Inserting the material unit and the Primeprint-Box into the printer. Different batches of the same 3D printing resin may not be mixed. The 3D printing resin can remain in the vat for a maximum of five months.

Construction process

Visual inspection of print resin material: in case print resin material is visually contaminated, it may not be used and must be discarded. Generate a print job complying with printer and material parameters. For setting up the printer follow the instructions for use of the printer. Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer (inLab CAM).

Post-processing

If possible, the finishing must take place immediately after the building process in the Primeprint PPU. To do this, remove the Primeprint box (transport container with build platform) from the Primeprint including build job and transfer it to the Primeprint PPU. In the Primeprint PPU, the pre-cleaning and main cleaning steps are performed in two separate vessels with isopropanol and fresh isopropanol (purity >98%). The build job is then automatically fed to the drying and post-exposure process. These steps are performed in the system according to the material. No user interaction is required.

CAUTION: The use of Primeprint is only approved when applied with the compatible devices mentioned in this instructions for use. Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specification.

1.12 Finishing

The Primeprint box is removed from the Primeprint PPU following the steps that have been completed automatically in the Primeprint PPU. The objects are now removed from the build platform, and the support structures are detached.

1.13 SURFACE PROCESSING

Polish surface mechanically.

1.14 Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

1.15 Notes

- ❑ Always keep container tightly sealed, immediately close the container carefully after each use.
- ❑ Please follow the instructions on the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.de and to the competent authority in which the user and/or patient is established.

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Sirona Dental Systems GmbH
DETAX GmbH