INSTRUCTIONS FOR USE

INTENDED USE:

Essix® Dual Laminate is thermoformed to fabricate intra-oral appliances such as TMJ splints, nightguards, and clenched & grinders. Essix® Dual Laminate is suitable for use with bruxers/grinders.

INDICATIONS FOR USE:

Mouthguard and aligner materials are indicated for the fabrication of orthodontic and dental appliances such as aligners, bite planes, mouthguards, nightguards, snoring appliances, splints, retainers, repositioners, and temporary bridges.

PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS:

Performance characteristics of orthodontic plastics include a good combination of strength, durability and clarity, resistant to cracking or breaking, easily removeable from models and easy to trim.

Clinical benefits of orthodontic treatment may include overall improvement of tooth health and correction of malocclusion.

CONTRAINDICATIONS:

This product is contraindicated for patients with a history of allergic reactions to plastics.

WARNINGS:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

Report serious incidents that have occurred in relation to the device to device manufacturer and the Competent Authority in Member State of residence, or other applicable regulatory agency

PRECAUTIONS:

Patients with a history of allergic reactions to plastics should not use this product.

Do not reuse plastic sheet. Once thermoformed it cannot be reused.

Thermoformed device is for Single Patient Use. Risks of reuse include infection and contamination (i.e. bacterial, viral).

Store plastic sheets in a cool, dry place.

ADVERSE REACTIONS:

Allergic reactions.

STEP BY STEP INSTRUCTIONS FOR USE WITH ESSIX VACUUM MACHINES:

- All models should be trimmed to 3/4" high. Full arch models should have the palate trimmed away.
- IMPORTANT:** Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
- Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
- Place the plastic into the frame of the vacuum machine soft side up and the model on the base plate.
- Swing the heating element over the plastic. Heat for approximately 1-1½ minutes or until glaze is visible on the soft side. Material will smoke (normal).
- Swing the heating element away and flip the material over so the hard side is now facing up. Swing the heating element back over the plastic and heat for approximately 15-25 seconds (until it is glossy).
- Move the heating element away from the plastic.
- Turn on the vacuum and lower the plastic over the model.
- Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it "shrinks" the material around the cast, realizing a tighter fit.
- Let vacuum continue to run for approximately 1 minute. Allow 10 minutes for cooling.
- Trim and finish appliance as desired.

COMMON PROBLEMS:

Not heating the soft side long enough.

Heating the hard side too long - it will become yellow and burn.

SPECIAL NOTES:

The edges of Dual Laminate sheets may appear rough. This is normal and will not affect the function of the plastic. Cold cure acrylic may be applied to the hard side only. Roughen the surface of the desired area prior to applying the monomer. Beware that the monomer does not get on the soft side of the material - it will burn (yellow).

Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic "folds" when thermoformed, reduce heating time until material does not fold.

For step-by-step instructions refer to owner's manual of thermo-forming machine used.

Thickness	Vacuum Heat Time**	Biostar™ Code**	Essix® SelectVac®
2.0mm (.080")	Soft - 45 seconds	Soft - 151 (Triple Flip)	Soft - 2:45
	Hard - 15 seconds	Hard - 105 (Triple Flip)	Hard - 0:50
3.0mm (.120")	Soft - 90 seconds	Soft - 231 (Triple Flip)	Soft - 3:00
	Hard - 25 seconds	Hard - 115 (Triple Flip)	Hard - 1:00

**after preheating element

Note: (Triple Flip) is for Biostar™ process only: Heat soft surface, then flip and heat hard surface, then flip back to the soft surface before pressure forming.

Not recommended for use in a Drufoomat.

SAFE DISPOSAL:

Unused plastic sheets can be recycled.

Used dental appliances may be considered a biohazard. Follow your country-specific laws, directives, standards and guidelines for the disposal of used devices.

 Manufacturer: **Dentsply Sirona Orthodontics Inc.**
7290 26th Court East Sarasota, FL 34243 USA

 www.dentsplysirona.com/orthodontics-e-labeling

Emergo Europe
Prinsessegracht 20
2514 AP, The Hague
The Netherlands

MKT-DFU0037-ML Rev. 01 Issued 03-2020

Biostar and Drufoomat are not trademarks of Dentsply Sirona.

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©2020 Dentsply Sirona