



Ankylos®

Excellence
without exception



 Dentsply
Sirona

Excellence without exception


Ankylos is for people who demand excellent results. With more than 30 years' proven clinical use and development, this innovative, world-class implant technology is thoroughly established and highly dependable. It allows your expertise and eye for detail to shine through. Because it is your skills and your experience that are the real key to improving your patients' quality of life. And with Ankylos you can provide something that every single patient will continue to treasure – excellent anterior aesthetics.



Ankylos® Implanting TissueCare

The true value of an implant system reveals itself with time. For over 30 years, the Ankylos implant system has stood for reliable, lasting esthetics. Numerous publications and long term clinical experience clearly demonstrate that Ankylos maintains hard and soft tissue stability, resulting in natural looking aesthetics that last.

This success can be attributed to a number of key features in the Ankylos system design which, together, form the unique Ankylos TissueCare concept.




SoftTissue
Chamber™




One fits all
TissueCare
connection



Progressive thread



Friudent® plus surface



Connected with you for more than 30 years

The Ankylos implant system was designed in 1985, based on the idea that a prosthetic dental implant could be just as stable as a natural tooth. After more than 30 years in clinical use, the Ankylos system still represents exceptional quality, lasting hard and soft tissue maintenance and superior aesthetics.

1993

Launch of the newly named Ankylos implant system.

2005

The Friadent implant surface and micro-rough implant shoulder are launched.

2008

Launch of Ankylos C/X implants, for use with both non-indexed (C) and indexed (X) prosthetic components.

2012

Introduction of the SmartFix concept, a prosthetic solution for full arch restorations at an angle.



2011

Introduction of titanium base, two-piece CAD/CAM abutments.

1987


First clinical use of a groundbreaking new implant system. Its innovative features include tapered connection, significant horizontal offset and progressive thread design.

2001

Introduction of the unique Ankylos SynCone concept for prosthetic overdentures.

2005

The Ankylos TissueCare concept is launched, offering a unique collection of tissue-stabilising features.



2013

Launch of the 6.6 mm Ankylos, designed to reduce the need for vertical bone augmentation procedures.




2016

Launch of new Titanium Nitride (TiN) coated Ankylos drills without internal irrigation, for better contrast of depth marking.



2018

Azento is a convenient start-to-finish custom implant solution, enabling dentists to provide consistently excellent results.



2012

Introduction of Atlantis patient-specific CAD/CAM abutments.

2014

Launch of the world-exclusive WeldOne concept for stable, lasting chairside solutions.

2015

Introduction of the high cover screw.

2018

Launch of the Atlantis® Custom-Base solution with angulated screw access and Ankylos interface.

2018

The Acuris concept – fixed prosthetics (single crown) for the patient and removable for the dentist. Thanks to highly accurate conometric retention there is no need for cement or prosthetic screws. Introduced for Ankylos, Xive and Astra Tech Implant System EV.

The Ankylos[®] TissueCare Concept


Scientifically and clinically proven to maintain hard and soft tissue over time, the Ankylos TissueCare concept combines a number of key features from the Ankylos implant system.

SoftTissue Chamber™

The TissueCare connection and micro-rough implant shoulder make subcrestal placement possible. Together with a horizontal offset and concave abutment design, they create a three dimensional chamber for soft tissue and bone, which in turn increases tissue stability.

Friident[®] plus surface

The microstructure of the grit-blasted and high temperature-etched Friident plus implant surface provides superior bone-to-implant contact, promotes early osseointegration and results in rapid bone formation.



Where Ankylos[®]
goes, hard and
soft tissue follow...

One fits all TissueCare connection

All abutments, regardless of size, share the same tapered connection, which simplifies the treatment process. A steep connection taper promotes friction lock and eliminates micromovement, thereby providing non-indexed prosthetics with an indexed option.

Progressive thread

The progressive thread design means the functional load is transferred to the apical part of the implant, creating primary stability and enabling immediate loading.

The Ankylos TissueCare concept provides lasting tissue stability and excellent aesthetic results.


Stable peri-implant
with hard and soft
tissue after healing



Final restoration
in place



Two years after
prosthetic
restoration



Four years after
prosthetic
restoration




SoftTissue Chamber™

The Ankylos SoftTissue Chamber provides a highly favorable three dimensional space for biological growth around the abutment. Thanks to a wide horizontal offset design and subcrestal implant placement, this biological space fills with a three dimensional network of connective tissue and collagen fibers. Which means a higher interdental papilla, lasting hard and soft tissue stability and superior aesthetic results.

Positive tissue response

The SoftTissue Chamber creates enough space for bone to form over the implant shoulder, as well as plenty of connective tissue fibers for healthy, natural-looking soft tissue.

- An optimum biological space, thanks to a combination of subcrestal placement, a wide horizontal offset and a narrow, concave abutment emergence profile.
- Together with the unique properties of the TissueCare connection, such as the friction-lock design that eliminates all micromovement, the SoftTissue Chamber enables and optimises one abutment-one-time treatment.



SoftTissue Chamber - the home of lasting tissue stability

Horizontal offset (platform switching)

Thanks to the geometry of the Ankylos TissueCare connection, the transition between implant and abutment shifts to a more central position. This horizontal offset creates the perfect conditions for hard and soft tissue stability at the implant shoulder. And when combined with other benefits like no micromovement and no bacterial growth, it also provides lasting tissue maintenance.

possible. This leads to a positive bone response, enabling the desired emergence profile and transgingival healing. In other words, Ankylos offers the versatility of a two-piece system while maintaining healthy hard and soft tissue and keeping it free from irritation.

Micro-rough implant shoulder


The patented microstructure implant shoulder promotes bone formation up to the abutment. This supports the overlying soft tissue and creates the perfect conditions for stability and lasting tissue health.

Subcrestal placement

The unique friction-lock and keyed TissueCare connection creates such a close fit between the two components that subcrestal implant placement is



Immediate placement and loading, standard abutment




Healing after 6 months




Final crown


Courtesy of Dr Marco Degidi, Bologna, Italy



1



2



3

- 1 The subcrestal placement of Ankylos implants enables appositional growth even on the implant shoulder. Sulcus former in Ankylos implant three months post-op. (Courtesy of Dr Dietmar Weng, Starnberg, Germany)
- 2 The concave abutment design provides space for the tissue to heal.
- 3 An ideal biological situation thanks to the subcrestal placement, a narrow abutment transition and a rigid connection. (Courtesy of Dr Marco Degidi, Bologna, Italy)

One fits all TissueCare connection



- Friction-locked and keyed
- No micromovement
- Tight

One fits all connection



The friction-locked and keyed stability of the Ankylos TissueCare connection is similar to that of a one-piece implant. Its precision engineered design ensures optimum load bearing and eliminates any micromovement between the implant and abutment or mechanical irritation that might lead to bone resorption.

Designed for absolute tightness, the peri-implant bone and gingival structures do not treat the connection as a gap. This reduces bacterial colonisation and the risk of inflammation.


One connection for all implants, regardless of implant diameter, means treatment is much simpler with the TissueCare connection.

Stability and lasting tissue health

- Provides exceptional tissue stability and eliminates all micromovement
- Ensures lasting tissue health and positive aesthetic outcomes as a result



The TissueCare connection – unique freedom of non-indexed prosthetics with an indexing option.



Caring for tissue over time



One fits all – full surgical and prosthetic flexibility with interchangeable prosthetic components

Thanks to identical tapered connections, all Ankylos abutments are fully compatible with all Ankylos implants. This versatility simplifies the treatment process by reducing the number of prosthetic components and ensures optimum abutment fit, regardless of which implant is necessary to meet the surgical requirements.

Tapered connection with indexing option

The tapered TissueCare connection enables 360° alignment of the prosthetic abutments in any position. Abutments with a positioning aid index are also available, for those who prefer indexed prosthetics. In both cases,

the keyed and friction-locked taper acts as an anti-rotation device and once the prosthetic restoration is fixed into place, the implant and abutment form an extremely stable and rotation-free unit.

Wide range of prosthetic options

The Ankylos system offers a wide variety of prosthetics, for various indications and preferred approaches. Each range contains abutments in several different sizes and designs, with or without a positioning aid (index), to ensure optimum functional and aesthetic results for every patient.



Progressive thread

The progressive thread has been designed to match the bone structure and simplify insertion.

Design features

- Cylindrical geometry reduces load transfer to the cortical bone
- The gradual increase in thread depth improves the transfer of load to the cancellous bone

This results in a controlled loading force in the areas significant for the maintenance of crestal bone. Once the thread is in position it creates primary stability and enables immediate loading.



Friadent[®] plus surface




Implant macro and microstructure create the optimum conditions for the surrounding bone

All Ankylos C/X implants have a grit-blasted and high temperature-etched Friadent plus surface. The excellent properties of this microstructure provide superior bone-to-implant contact and promote early osseointegration, as documented in numerous in vitro and in vivo studies.


Surface and macro design – the perfect combination

- Extremely receptive to new cell growth as the surface is hydrophilic, which stimulates primary cell apposition
- Unique, three dimensional micro design that promotes the apposition of bone-forming cells and ensures optimum osseointegration
- Rapid new bone formation and increased bone maturation in the early stages creates excellent interface stability


In addition to the benefits of the progressive thread and its unique design, the Friadent plus surface provides superior bone-to-implant contact and primary stability, which enables immediate loading.



1 SEM (3000 x) of the Friadent plus surface structure – bimodular morphology with micropores (0.5–1µm) in macrostructure



2 Initial contact and anchorage of an osteoblast using thread-like extensions (filopodia) on the Friadent plus surface



3 Extracellular matrix on the Friadent plus surface



4 Histology (10 x) of bone-to-implant contact on the Friadent plus surface between the implant threads

Clinical experience

The Ankylos implant system is scientifically and clinically proven to maintain hard and soft tissue over time, with excellent aesthetics and high patient satisfaction. Extensive literature evaluating over 55 000 implants shows safe and predictable treatment outcomes for the Ankylos implant system, with >97% implant survival rate after 5 years or more of follow-up.

> 55 000
implants

> 17 000
patients¹

97%

success rate after 5 years
or more of follow-up

Long-term evaluation of peri-implant bone loss

(after up to 204 months of follow-up)¹

Implants	12 737
Patients	4 206
Horizontal	< 1 mm 85.7%
Vertical	< 1 mm 85.2%

Conclusion: High implant survival rate and low rates of peri-implant bone loss observed in 125 implants followed for 17 years.

Tapered connection

Single implants replacing molars

Implants	233
Minimum years	5
Maximum years	7.37
Average years	6.3
Failures (1)	0.5%
Loosening of abutment (3)	1.3%
Abutment fracture (0)	0%

Conclusion: High prosthetic reliability observed in 233 implants followed for at least 5 years.

6-year multicentered clinical study

Study centers	34
Dentists	80
Implants placed during 2.5 years	1 500
Patients	450
Study follow-up	3–5 years
Implant survival	97.5%


Conclusion: Survival rate of 97.5% in 1 500 implants followed for 3 to 5 years.

¹Krebs M, Schmenger K, Neumann K, Weigl P, Moser W, Nentwig GH.: Long-term evaluation of Ankylos dental implants, part I: 20-year life table analysis of a longitudinal study of more than 12,500 implants. Clin Implant Dent Relat Res 2015;17 Suppl 1:e275-86.

Weigl P: New prosthetic restorative features of Ankylos implant system. J Oral Implantol 2004;30(3):178-88.

²Morris HF, Ochi S, Crum P, et al.: AICRG, PART I: A 6-Year multicentered, multi-disciplinary clinical study of a new and innovative implant design. J Oral Implantol 2004;30(3):125-33.

Restoration of maxillary lateral incisor – five-year follow-up




1 Control image after final restoration of tooth 12 shows excellent bone apposition
2 & 3 The situation after almost a year (Mar 2006)

4 Clinical situation one year later (Feb 2007)
5 Perfectly maintained peri-implant bone, five years after insertion

Courtesy of Dr Patric Renner,
Goldbach, Germany

Two patients from a study carried out between 1991 and 2011 by the University of Frankfurt, involving over 4 000 patients and 12 000 implants



Male patient, born 1955


1a Image taken in 1997, after implant placement
1b Image taken in 2010, demonstrating 13 years of tissue stability

Male patient, born 1976

2a Image taken in 1993, after implant placement
2b Image taken in 2012, demonstrating 19 years of tissue stability

Courtesy of Dr Mischa Krebs et al.,
Frankfurt am Main, Germany

Immediate implant procedure after removal of a fractured lateral incisor – restoration of the Ankylos® implant using a Cercon® abutment



X-rays showing stable crestal bone level above the implant shoulder over time:

1 During the provisionalisation stage
2 After completion of the final restoration

3 After two years
4 After four years
5 After six years
6 & 7 After eight years of use

Courtesy of Dr Nigel Saynor,
Stockport, UK

Mandibular restoration using the SynCone® concept – thirteen years post-op



1 & 2 Clinical situation 13 years after implant placement

3 X-ray taken 13 years after the procedure

Courtesy of Dr Dittmar May,
Lünen, Germany

Ankylos[®] at a glance

Indications

- The Ankylos C/X implant system is for single-stage or two-stage surgical procedures and cemented, removable or screw-retained restorations.
- The Ankylos C/X implant system can be used for immediate placement and function on single tooth and/or multiple tooth applications when adequate primary stability is achievable, with appropriate occlusal loading, in order to restore chewing function. Multiple tooth applications may be splinted.


Implants

Implant geometry

- Cylinder screws

Implant material

- Pure titanium grade 2 (ISO 5832-2)





Implant surface

- Friadent plus microstructure (grit-blasted and high-temperature etched)

Implant-abutment connection

- Keyed and friction-locked tapered connection (TissueCare connection) with optional index



Implant diameters and lengths

Diameters	3.5 mm (A)	4.5 mm (B)	5.5 mm (C)	7.0 mm (D)
Lengths	6.6 mm	6.6 mm	6.6 mm	-
	8 mm	8 mm	8 mm	8 mm
	9.5 mm	9.5 mm	9.5 mm	9.5 mm
	11 mm	11 mm	11 mm	11 mm
	14 mm	14 mm	14 mm	14 mm
	17 mm	17 mm	17 mm	-


Surgical instrument set

The Ankylos modular plastic surgical kits contain all of the necessary instruments for standard surgical procedures. There are three versions of the Ankylos surgical kit: one for machine-driven procedures (A and B implants only), one for manual procedures and one for guided surgery.



Drill preparation

Externally irrigated drills with diameter and length markers.



Implant 900

Induction drive unit with preset programs for augmentation and implant placement.



Prosthetic restorations

	Single-tooth crowns	Fixed bridges	Removable dentures
Balance anterior C/ or /X	x	x	-
Cercon balance C/	x'	-	-
Titanium base C/ or /X	x	x	-
Atlantis CAD/CAM abutments	x	x	-
Regular C/ or /X	x	x	-
Balance posterior C/	x	x	-
Standard C/ abutment	x	x	x

	Single-tooth crowns	Fixed bridges	Removable dentures
Balance base abutment C/ and SmartFix concept	-	x	x
SynCone C/	-	-	x
Atlantis Conus	-	-	x
Locator C/	-	-	x
Snap attachment C/	-	-	x
Atlantis	-	x	x
WeldOne	-	x	-

Single-tooth restoration in front-tooth region only (regions 13-23 and 33-43)

Esthetic solutions

Balance anterior C/ or /X abutment Cercon balance C/ abutment

- Anatomically shaped abutments for challenging aesthetic indications



Titanium base C/ or /X

- Titanium base for individual metal-reinforced ceramic abutments



Atlantis patient-specific abutments for Ankylos

- Patient-specific abutments
- Available in titanium, gold-coloured titanium and zirconia



Crown and bridge abutments

Regular C/ or /X abutment

- For various types of fixed prosthetics



Balance posterior C/

- For various types of fixed prosthetics



Standard C/ abutment

- One abutment-one-time concept
- Transfer on abutment level



Full arch restorations

Balance base abutment C/ narrow

- Foundation for bar and bridge construction
- SmartFix concept with 15° and 30° angled abutments
- By using the SmartFix concept, augmentation and critical anatomical areas can be avoided



Snap attachment C/

- For fixing mandible dentures in edentulous patients



Locator C/ abutment

- For fixing mandible dentures in edentulous patients



SynCone C/ abutment

- For fast and economical construction of removable dentures based on prefabricated conical crowns



Atlantis conus abutment

- Consists of patient-specific Atlantis conus abutments together with prefabricated SynCone caps



Atlantis implant suprastructures

- CAD/CAM structures for retaining bar and bridge structures



Ankylos WeldOne

- Innovative intraoral welding concept
- Facilitates sturdy, lasting chairside solutions



Key references:

The Ankylos TissueCare concept

Degidi M, Nardi D, Piattelli A. One abutment at one time: non-removal of an immediate abutment and its effect on bone healing around subcrestal tapered implants. *Clin Oral Implants Res* 2011;22(11):1303-07.

Degidi M, Perrotti V, Shibli JA, et al. Equicrestal and subcrestal dental implants: a histologic and histomorphometric evaluation of nine retrieved human implants. *J Periodontol* 2011;82(5):708-15.

Degidi M, Daprile G, Nardi D, et al. Immediate provisionalization of implants placed in fresh extraction sockets using a definitive abutment: the chamber concept. *Int J Periodontics Restorative Dent* 2013;33(5):559-65.

Degidi M, Nardi D, Daprile G, Piattelli A. Nonremoval of immediate abutments in cases involving subcrestally placed post-extractive tapered single implants: a randomized controlled clinical study. *Clin Implant Dent Relat Res* 2014;16(6):794-805.

Donovan R, Fetner A, Koutouzis T, et al. Crestal bone changes around implants with reduced abutment diameter placed non-submerged and at subcrestal positions: a 1-year radiographic evaluation. *J Periodontol* 2010;81(3):428-34.

Schwartz Fo HO, Novaes AB, Jr., de Castro LM, et al. In vitro osteogenesis on a microstructured titanium surface with additional submicron-scale topography. *Clin Oral Implants Res* 2007;18(3):333-44.

Zipprich H, Weigl P, Lauer H-C, et al. Micromovements at the implantabutment interface: measurement, causes, and consequences. *Implantologie* 2007;15(1):31-46.

Soft and hard tissue maintenance

Chou CT, Morris HF, Ochi S, et al. AICRG, PART II: Crestal bone loss associated with the Ankylos implant: Loading to 36 month. *J Oral Implantol* 2004;30(3):134-43.

Degidi M, Nardi D, Piattelli A. Prospective study with a 2-year follow-up on immediate implant loading in the edentulous mandible with a definitive restoration using intra-oral welding. *Clin Oral Implants Res* 2010;21(4):379-85.

Degidi M, Nardi D, Piattelli A. The conometric concept: Coupling connection for immediately loaded titanium-reinforced provisional fixed partial dentures-a case series. *Int J Periodontics Restorative Dent* 2016;36(3):347-54.

Koutouzis T, Neiva R, Nonhoff J, et al. Placement of implants with platform-switched Morse taper connections with the implant-abutment interface at different levels in relation to the alveolar crest: a short-term (1-year randomized prospective controlled clinical trial). *Int J Oral Maxillofac Implants* 2013;28(6):1553-63.

Koutouzis T, Neiva R, Nair M, Nonhoff J, Lundgren T. Cone beam computed tomographic evaluation of implants with platform-switched Morse taper connection with the implant-abutment interface at different levels in relation to the alveolar crest. *Int J Oral Maxillofac Implants* 2014;29(5):1157-63.

Nentwig GH Ankylos implant system: concept and clinical application. *J Oral Implantol* 2004;30(3):171-7.

Romanos GE Nentwig GH Single molar replacement with a progressive thread design implant system: a retrospective clinical report. *Int J Oral Maxillofac Implants* 2000;15(6):831-6.

Romanos GE, Malmstrom H, Feng C, Ercoli C, Caton J. Immediately loaded platform-switched implants in the anterior mandible with fixed prostheses: a randomized, split-mouth, masked prospective trial. *Clin Implant Dent Relat Res* 2014;16(6):884-92.

Long-term follow-up

Frisch E, Ziebolz D, Ratka-Kruger P, Rinke S. Double crown-retained maxillary overdentures: 5-year follow-up. *Clin Implant Dent Relat Res* 2015;17(1):22-31.

Jesch P, Jesch W, Bruckmoser E, et al. An up to 17-year follow-up retrospective analysis of a minimally invasive, flapless approach: 18 945 implants in 7783 patients. *Clin Implant Dent Relat Res* 2018;E-pub Feb 15, doi: 10.1111/cid.12593

Krebs M, Schmenger K, Neumann K, et al. Long-term evaluation of ANKYLOS(R) dental implants, part i: 20-year life table analysis of a longitudinal study of more than 12,500 implants. *Clin Implant Dent Relat Res* 2015;17 Suppl 1:e275-86.

Romanos GE, Nentwig GH. Immediate functional loading in the maxilla using implants with platform switching: five-year results. *Int J Oral Maxillofac Implants* 2009;24(6):1106-12.

Romanos GE, May S, May D. Treatment concept of the edentulous mandible with prefabricated telescopic abutments and immediate functional loading. *Int J Oral Maxillofac Implants* 2011;26(3):593-7.

Romanos GE, Aydin E, Gaertner K, Nentwig GH. Long-term results after subcrestal or crestal placement of delayed loaded implants. *Clin Implant Dent Relat Res* 2015;17(1):133-41.

Sethi A, Kaus T, Sochor P, et al. Evolution of the concept of angulated abutments in implant dentistry: 14-year clinical data. *Implant Dent* 2002;11(1): 41-51.

Esthetic evaluation and patient satisfaction

Aboud M, Koeck B, Stark H, Wahl G, Paillon R. Immediate loading of single-tooth implants in the posterior region. *Int J Oral Maxillofac Implants* 2005;20(1):61-8.

Döring K, Eisenmann E, Stiller M. Functional and esthetic considerations for single-tooth Ankylos implant-crowns: 8 years of clinical performance. *J Oral Implantol* 2004;30(3):198-209.

Esposito M, Bressan E, Grusovin MG, et al. Do repeated changes of abutments have any influence on the stability of peri-implant tissues? One-year post-loading results from a multicentre randomised controlled trial. *Eur J Oral Implantol* 2017;10(1):57-72.

Morris HF, Ochi S, Rodriguez A, Lambert PM. AICRG, Part IV: Patient satisfaction reported for Ankylos implant prostheses. *J Oral Implantol* 2004;30(3):152-61.

About Dentsply Sirona Implants

Dentsply Sirona Implants offers comprehensive solutions for all phases of implant therapy, including Ankylos®, Astra Tech Implant System® and Xive® implant lines, digital technologies, such as Atlantis® patient-specific solutions and Simplant® guided surgery, Symbios® regenerative solutions, and professional and business development programs, such as STEPPS™. Dentsply Sirona Implants creates value for dental professionals and allows for predictable and lasting implant treatment outcomes, resulting in enhanced quality of life for patients.

About Dentsply Sirona

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with a 130-year history of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. As The Dental Solutions Company™, Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better, safer and faster dentistry. Dentsply Sirona's global headquarters is located in York, Pennsylvania, and the international headquarters is based in Salzburg, Austria. The company's shares are listed in the United States on NASDAQ under the symbol XRAY.

Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

ANZ [✉ clientservices@dentsplysirona.com](mailto:clientservices@dentsplysirona.com) [📍 11-21 Gilby Road, Mount Waverley, VIC 3149](https://www.dentsplysirona.com) [🌐 www.dentsplysirona.com](https://www.dentsplysirona.com)

Australia [📞 1300 552 929](tel:+6171300552929) | **New Zealand** [📞 0800 336 877](tel:+640800336877)

YOUR
CONFIDENCE
IS OUR
INSPIRATION.

