Document Detail

Type: 0060-GR DT
Document No.: 0060-90764[7]
Title: DFU, NUPRO Freedom
Comment
Status: CURRENT
Effective Date: 12-Dec-2018

Approvals

<u>Role</u>	<u>Sign-off By</u>	<u>Sign-off Date</u>	
DPD Document Control	Lucy Brenneman	05-Dec-2018 1:29 pm	GMT
0060-Preventive UDI Approver	Jean Gardner	06-Dec-2018 9:17 pm	GMT
DPD Marketing Management	Angela Donohue	07-Dec-2018 6:03 pm	GMT
DPD Research and Development	Michael Stephenson	07-Dec-2018 1:17 pm	GMT
DPD Operations-Manufacturing	Carl Burkhart	07-Dec-2018 1:16 pm	GMT
DPD Quality Engineer	Kalpesh Patel	07-Dec-2018 1:08 pm	GMT
DPD Quality Assurance/Regulatory Management	Donata Kline	10-Dec-2018 9:56 pm	GMT
DPD BPCS Coordinator	Natalie Weiss	11-Dec-2018 3:47 pm	GMT
0060-Preventive UDI Approver	Donata Kline	12-Dec-2018 11:13 am	GMT
DPD Document Control	Lucy Brenneman	12-Dec-2018 11:39 am	GMT

Artwork Specification

Rev. 7

Part Number:	90764
Title	DFU, Cordless Prophy System with SmartMode® Technology
Material:	Inside panels: 70# smooth white offset Cover: 100# Text
Dimensions:	8-1/2 x 11
Binding:	Saddle stitch or perfect binding (depending on recommendation by printer)
Print Colors:	4 Color Process, Full bleed on cover only
Revision	7 (11/18)

Notes:

1. The printer is responsible for adjusting color as needed to match the artwork on final print.

Revision History:

Rev	Description of Change	Author	Dept	Change No.
1	First release to production.	J. Cusick	Mktg	CHO-073651
2	Added French translation	J. Cusick	Mktg	CHO-076824
3	Insert grasp and adaptation technique; Change indicator light information; Update translations	J. Cusick	Mktg	CHO-081648
4	Correct Document Type from Literature to Labeling	J. Cusick	Mktg	CHO-089087
5	Updated the package contents list to include the foot pedal (only with 90741L05).	J. Cusick	Mktg	CHO-090452
6	Format change for page numbers	J. Cusick	Mktg	CHO-093078
7	Update sterilization instructions	J. Cusick	Mktg	CHO-114682

Nupro Freedom[®]

Cordless Prophy System with SmartMode[®] Technology

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Directions For Use / Mode d'emploi

Please read carefully and completely before operating unit.

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
Table of Contents

Grasp & Adaption Technique	4
Overview	5
1. Indications for Use	5
2. Contraindications	5
3. Warnings	5-6
4. Precautions	6
5. Adverse Reactions	6
6. Nupro Freedom® Cordless Prophy System with SmartMode® Technology Description	7
7. Unpacking the System	8
8. Infection Control Procedures	8-10
9. Charging the Cordless Handpiece and Sleep Mode	10-11
10. Powering and Charging the Optional Wireless Foot Pedal	11
11. Synchronizing the Handpiece and Foot Pedal	12-13
12. Preparation for Use	13-14
13. Using SmartMode® Technology and/or Optional Wireless Foot Pedal	14-15
14. Specifications	16
15. Classifications	16
16. Symbol Identification	17
17. Disposal of Unit	18
18. Electro Magnetic Compactibility Precautions	18-20
19. Troubleshooting	20-21
20. Accessories	22
21. Limited Warranty	22

Nupro Freedom® Cordless Prophy System with SmartMode® Technology

Grasp & Adaptation Technique

Using a modified pen grasp, rest the handpiece in the V-shaped area of your hand, between the index fingers and thumb.

To initiate cup rotation, gently press the cup against the tooth surface so there is slight flaring of the cup.


Stabilize the handpiece with a fulcrum or finger rest.

Polish using light, intermittent pressure, gliding the cup across the tooth without lifting the cup from the tooth surface.

Periodically evaluate your grasp to ensure correct balance and a light grasp to support optimal use of the intuitive Smartmode technology.

Remove excess paste or saliva buildup on the cup with the saliva ejector.

Position the handpiece so the cup is adapted to the facial surface of the most posterior tooth you are polishing.

This quick reference guide is not intended as a replacement for a Directions for Use Manual (DFU) for any products shown. Always check your manufacturer's guidelines for recommended use of equipment.

To learn more, visit dentsplysirona.com/nuprofreedomvideos

Overview

The Nupro Freedom® Cordless Prophy System with SmartMode® Technology offers a cordless design that eliminates cord drag and allows clinicians easy access and comfort during prophylaxis procedures. SmartMode® Technology controls the speed of the handpiece using the DPA's pressure against the tooth, eliminating the need for a foot pedal. Quieter than the original Midwest® RDH® handpiece, the system offers all-day battery life while delivering the same performance you would expect from a corded handpiece. The system includes a steam autoclavable outer sheath for infection control.

1. Indications for Use

The Nupro Freedom® Cordless Prophy System with SmartMode® Technology is a high-performance cordless prophy handpiece with an optional wireless foot pedal for use with Nupro Freedom® Disposable Prophy Angles in a hygiene operatory to perform cleaning and polishing procedures on teeth.

2. Contraindications

None Known.

Safety Conventions in This Document

WARNING	Use care to prevent personal and / or patient injury
PRECAUTION	Use care to prevent product damage and ensure safe and effective product use

3. Warnings

- To prevent damage, charge the cordless handpiece using only the Nupro Freedom® with SmartMode® Technology charging base for cordless handpiece and power supply and cord.
- Sterilizing the cordless handpiece inner module will cause component damage to the cordless handpiece and the sterilizing equipment, and may cause personal bodily injury.
- The outer sheath must be steam sterilized before first use and between patients to prevent patient cross-contamination. See Section 8 for the Infection Control Procedures.
- The Disposable Prophy Angles are designed for single-patient use only and should never be used more than once. Disposable Prophy Angles are not autoclavable or designed to withstand disinfection solutions. The risks of reuse of a Disposable Prophy Angle are damage to equipment and patient cross-contamination. Install a new Disposable Prophy Angle before each use.
- It is the responsibility of the Healthcare Professional to determine the appropriate uses of this product and to understand:
 - the health of each patient
 - the dental procedures being undertaken
 - applicable industry and governmental agency recommendations for infection control in dental healthcare settings
 - requirements and regulations for safe practice of dentistry
 - these Directions for Use in their entirety
- Per FCC Part 15.21, changes or modifications not expressly approved by the party responsible for compliance (i.e., the manufacturer) could void the user's authority to operate this equipment.
- Failure to follow recommendations for environmental operating conditions (see Section 14 for Specifications) could result in injury to patients or users.
- Inspect the cordless handpiece system before each use for worn, loose or damaged parts. Do not attempt to operate unless the Disposable Prophy Angle (DPA) is properly installed. A loose DPA could separate from the cordless handpiece causing bodily injury. Reinstall the DPA or replace any damaged parts as necessary.
- To prevent bodily injury and damage to the device, do not sterilize the Disposable Prophy Angle, cordless handpiece inner module, charging base for cordless handpiece, optional wireless foot pedal or power supply and cord. Disinfect the cordless handpiece inner module, charging base for cordless handpiece, optional wireless foot pedal or power supply and cord using only the tested and approved disinfectants listed in Section 8, Infection Control Procedures.
- The cordless handpiece inner module, optional wireless foot pedal, charging base for cordless handpiece and power supply and cord are not waterproof. To prevent damage to the equipment, contamination or bodily injury, do not immerse any of these components in water or a chemical solution.
- Do not spray disinfectant or other fluids directly onto the cordless handpiece inner module, optional wireless foot pedal, charging base for cordless handpiece, or power supply and cord to avoid liquid from pooling on the component. The user should spray solution onto a cloth or use a wipe to disinfect these items per the instructions in Section 8.

- Use only components and accessories listed in Section 7 of this manual. Failure to do so will void the warranty, may decrease system performance, may lead to unsafe operation, may negatively affect EMC performance and result in non-compliance.
- Never mount a Disposable Prophy Angle to the cordless handpiece body while it is operating.
- Dispose of Nupro Freedom® Disposable Prophy Angles after each patient according to CDC Guidelines for Infectious Waste and Federal, State and Local regulations.
- There are no user-serviceable items in the cordless handpiece inner module, power supply and cord, outer sheath, optional wireless foot pedal or charging base for cordless handpiece. Opening any of these units may result in unsafe operation and will void the warranty.
- According to IEC 60601-1/UL60601-1, this device must not be used in the presence of a flammable anesthetic gas mixed with air, oxygen, or nitrous oxide. (**Note:** nitrous oxide by itself is not a flammable anesthetic gas.)

4. Precautions

- Before using this product, carefully read and follow all instructions and save them for future reference. Observe all precautions and warnings.
- The cordless handpiece system can **only** be used with Nupro Freedom® Disposable Prophy Angles.
- As with all dental procedures, use universal precautions (i.e., wear face mask, eyewear, or face shield, gloves and protective gown).
- The cordless handpiece inner module motor is designed to be **lube-free**. Lubrication may cause damage to the cordless handpiece inner module.
- Oil and/or dirt may damage the motor, electronics and battery located inside the cordless handpiece inner module.
- The batteries are **not** user replaceable. When needed, the units should be returned to the listed repair center for replacement.
- Do not place the system on or next to a radiator or other heat source. Excessive heat may damage the system's electronics.
- Inadvertent system shutdown may occur in the presence of strong non-compliant radio frequency-generating components.
- This device complies with Part 15 of the FCC Rules and with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) This device must accept any interference received, including interference that may cause undesirable operation. Changes or modifications not expressly approved by the party responsible for compliance (i.e. the manufacturer) could void the user's authority to operate the equipment.
This Class B digital apparatus complies with Canadian ICES-003.
- Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.
- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

5. Adverse Reactions

There are no known adverse reactions.

6. Description of Nupro Freedom® Cordless Prophy System with SmartMode® Technology

7. Unpacking the System

As you unpack your Nupro Freedom® Cordless Prophy System with SmartMode® Technology, verify that the following components and accessories are included:

System:

- (1) Cordless Handpiece Inner Module
- (1) Outer Sheath
- (1) Charging Base for Cordless Handpiece
- (1) Power Supply and Cord
- (20) Nupro Freedom® Disposable Prophy Angles (DPAs)
- (1) Disposa-Shield® Disposable Barrier Trial Pack (25)
- (1) Cradle for Cordless Handpiece

Package:

- (1) Cordless Handpiece Inner Module
- (3) Outer Sheaths
- (1) Charging Base for Cordless Handpiece
- (1) Power Supply and Cord
- (200) Nupro Freedom® Disposable Prophy Angles (DPAs)
- (1) Box of Disposa-Shield® Disposable Barriers (500 Barriers)
- (1) Cradle for Cordless Handpiece
- (1) Wireless Foot Pedal* (90741L05 only)

8. Infection Control Procedures

The objective of the information provided in this section is to reduce the potential for patient cross contamination when using a Nupro Freedom® Cordless Prophy System with SmartMode® Technology during routine dental care. In the event any regulatory agency disagrees with this information, the agency requirements take precedence.

NOTE: Outer sheaths must be steam-autoclave sterilized prior to each use. Additional outer sheaths are available for purchase.

Instructions for Steam Sterilizing the Outer Sheath and Handpiece Cradle	
WARNINGS 	<p>These instructions are for use ONLY on the outer sheath and handpiece cradle. All other parts of the system should be disinfected according to the procedures in the "Disinfecting All Other Parts" section.</p> <p>The outer sheath and handpiece cradle for the Nupro Freedom® Cordless Prophy System with SmartMode® Technology are not sterile upon receipt and must be steam-autoclave sterilized prior to use in accordance with the following instructions.</p>
LIMITATIONS ON REPROCESSING	<p>The handpiece cradle may only be sterilized following the pre-vacuum sterilization procedure. The gravity steam sterilization procedures have not been validated for the handpiece cradle.</p> <p>Do not use an automatic washer/disinfector for the outer sheath. Discoloration of material will occur.</p> <p>Repeated cleaning and steam sterilization cycles have minimum effect on these instruments. End of life is normally determined by wear and damage due to use.</p> <p>Do not use chemical disinfectants prior to sterilization or rapid deterioration of the material may occur.</p> <p>Cold liquid disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.</p> <p>Do not immerse the outer sheath or handpiece cradle in an ultrasonic bath.</p>
POINT OF USE	<p>Remove excess soil with disposable cloth or paper wipe.</p> <p>It is recommended that instruments are reprocessed as soon as is reasonably practical following use.</p>
CONTAINMENT AND TRANSPORTATION	<p>Protect the outer sheath and handpiece cradle from contact with other dental instruments that may cause damage.</p>
PREPARATION FOR DECONTAMINATION	<p>Remove the outer sheath from the cordless handpiece inner module of the cordless handpiece. Only the outer sheath and handpiece cradle may be steam-autoclave sterilized.</p>

Instructions for Steam Sterilizing the Outer Sheath and Handpiece Cradle (Cont.)

CLEANING: MANUAL	<ol style="list-style-type: none"> 1. Disassemble by removing the disposable prophy angle (DPA) from the Outer Sheath, and the Outer Sheath from the inner module. 2. Prepare enzymatic solution per manufacturer's recommendation. 3. Using a lint-free cloth dampened with the prepared detergent, thoroughly wet the Outer Sheath with the prepared detergent. <ol style="list-style-type: none"> 3.1 Using a soft bristled brush, brush the Outer Sheath until all visible soil is removed, paying attention to the disposable prophy angle area. 4. Rinse the Outer Sheath under potable running tap water for a minimum of 30 seconds to remove detergent residue. 5. Dry the Outer Sheath using clean lint-free cloth. 6. Visually inspect each Outer Sheath for visible soil. If soil is seen, repeat the cleaning procedure.
DISINFECTION	<p>Do not use disinfectant solution on sterilizable outer sheaths or handpiece cradle.</p>
PACKAGING	<p>Place each instrument in a separate paper/plastic steam-sterilizing pouch. If using a sterilizing cassette, ensure that the sterilizer's maximum load is not exceeded.</p>
STEAM STERILIZATION	<p>Place bagged instruments into a steam autoclave, per the autoclave manufacturer's instructions.</p> <p>Gravity Steam Sterilization - Outer Sheath Full Cycle: 135°C (275°F) for 10 minutes</p> <p>Pre-vacuum Steam Sterilization - Outer Sheath and Handpiece Cradle Full Cycle: 132°C (270°F) for 4 minutes</p> <p>Alternate Method: Place non-bagged instruments into the steam autoclave and run at the above listed cycles.</p> <p>NOTE: Instruments steam sterilized unbagged should be used immediately.</p>
DRYING	<p>To dry, use the drying cycle of the steam autoclave. Set cycle for 20 to 30 minutes. Do not exceed 137°C (278°F).</p>
MAINTENANCE	<p>Visually inspect to ensure that all contamination has been removed.</p> <p>Check for distortion, damage or wear. Discard damaged, worn or corroded instruments.</p>
STORAGE	<p>To maintain sterility, instruments should remain bagged until ready for use.</p>
MANUFACTURER CONTACT	<p>In the United States, contact DENTSPLY Professional Customer Service Technical Support at 800-989-8826. For areas outside the United States, contact your local DENTSPLY Division.</p>

Instructions For Disinfecting All Other Parts (Charging Base, Inner Module, Foot Pedal and Power Supply)

WARNINGS 	<p>The charging base for cordless handpiece, cordless handpiece inner module, optional wireless foot pedal, power supply and cord are not sterilizable by autoclave, but can be disinfected following the procedures listed below.</p> <p>Only use water based non-immersion type disinfectant solutions.</p> <p>Per the Centers for Disease Control and Prevention (CDC), chemical germicide registered with the EPA as a "hospital disinfectant" and labeled for "tuberculocidal" (i.e., mycobactericidal) activity is recommended for disinfecting surfaces that have been soiled with patient material.</p> <p>Recommended Disinfectant:</p> <p>The following tuberculocidal disinfectant has been found to be effective in achieving a desired level of disinfection</p> <ul style="list-style-type: none"> • VoloWipes® Surface Disinfectant (manufactured by Sultan Healthcare) or equivalent <p>Do not immerse in liquid.</p>
LIMITATIONS ON REPROCESSING	<p>Repeated cleaning has minimum effect on these instruments. End of life is normally determined by wear and damage due to use.</p> <p>Do not use disinfectant solution on sterilizable outer sheaths. Refer to steam sterilization procedures for sterilizable outer sheaths.</p>

POINT OF USE	Remove excess soil with disposable cloth or paper wipe. Discard wipe after use.
CONTAINMENT AND TRANSPORTATION	Handle with care.
CLEANING	Generously spray disinfectant solution on a clean cloth. Wipe the outer surfaces of the charging base for cordless handpiece, cordless handpiece inner module, optional wireless foot pedal, power supply and cord. Discard used cloth. Wipe dry with a clean cloth.
DISINFECTION	Generously spray disinfectant solution on a clean cloth or use a disinfectant wipe. Wipe the outer surfaces of the charging base for cordless handpiece, cordless handpiece inner module, optional wireless foot pedal, and the power supply and its cord. Discard used cloth. Allow disinfectant to air dry.
DRYING	When cleaning, wipe surfaces dry with a clean cloth. To achieve disinfection, allow surfaces to air dry.
MAINTENANCE	Visually inspect to ensure that all contamination has been removed. Visually inspect power supply and cord for damage.
STORAGE	Temperature range: -20°C to 50°C Relative humidity range: 45 - 95% (non-condensing)
MANUFACTURER CONTACT	In the United States, contact DENTSPLY Professional Customer Service Technical Support at 800-989-8826. For areas outside the United States, contact your local DENTSPLY Division.

The instructions provided above have been validated by DENTSPLY as being capable of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing is actually performed using equipment, materials and trained personnel in the processing facility to achieve the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

9. Charging the Cordless Handpiece and Sleep Mode

The cordless handpiece should be charged for at least 2.5 hours prior to first use. The cordless handpiece should then be recharged after each full day's use. If necessary, the cordless handpiece may be "quick charged" for a single use in 15 minutes. Please note that the cordless handpiece motor will not operate when charging.

Use only the Nupro Freedom® Cordless Prophy System Power Supply and cord. Failure to do so could cause your system to malfunction and void your warranty.

9.1 When the bottom charge indicator light on the cordless handpiece inner module is solid Orange, the cordless handpiece must be charged.

9.2 Plug the power supply into the wall outlet and the mini USB connector into the charging base. A successful connection is shown with a Solid Green power-on indicator on back of the charging base.

PRECAUTION: Connect to a single phase AC power 100-240V power only. Otherwise, malfunction will occur. Do not unplug the power supply by pulling on the cord.

9.3 Place the cordless handpiece inner module on the charging base, aligning the charging contacts on each.

9.4 Refer to the figure and table below to determine the cordless handpiece's degree of charge.

Charge Indicator Color (when handpiece is in charging base)

Indicator Lights	Degree of Charge
Three Solid Green	95-100%
Two Solid Green and one Blinking Green	75-94%
One Solid Green and two Scrolling Green	50-74%
Three Scrolling Green	25-49%
Scrolling Orange and Green	< 24%

9.5 Sleep Mode (if cordless handpiece is not in a powered charging base)

- After 1 minute of non-use, the cordless handpiece will enter a Sleep Mode and the charge indicators and mode indicator will turn off. Upon picking up the handpiece within 1 hour, the Mode Indicator and charge indicators will illuminate and the unit will enter Standby Mode.
- After 1 hour of inactivity, the cordless handpiece will turn off to conserve battery charge. After this, the On/OFF/Mode Button will need to be pressed to bring the handpiece into Standby Mode.

10. Powering and Charging the Optional Wireless Foot Pedal

The optional wireless foot pedal should be charged for at least 2.5 hours prior to first use. The optional wireless foot pedal should then be recharged every month.

Use only the Nupro Freedom® Cordless Prophy System Power Supply and cord. Failure to do so could cause your system to malfunction and void your warranty.

10.1 If plugged in using the Nupro Freedom Power Supply, the optional wireless foot pedal may be synchronized with the cordless handpiece and depressed to control the DPA speed. If receiving power from the Nupro Freedom Power Supply, the Charge indicator on the top of the optional wireless foot pedal will illuminate according to the charge level per the table below:


Indicator Lights	Degree of Charge
Solid Green	95-100%
Blinking Green	30-95%
Blinking Orange	< 30%

10.2 If the foot pedal is tapped without the power supply connected, a solid Orange indicator light that stays lit for 1 second indicates that <30% charge is remaining.

10.3 Note that the optional wireless foot pedal will not charge unless it is ON.

10.4 To determine if the optional wireless foot pedal is ON, press down the foot pedal without it connected to the Nupro Freedom Power Supply. If it is ON, the foot pedal's charge indicator light will turn Green for a second (if the indicator turns Orange the foot pedal's charge is less than 30%). If there is no light when the foot pedal is depressed, the foot pedal is OFF.

10.5 If the optional wireless foot pedal is not turned ON, invert the optional wireless foot pedal and Press the POWER button for 1 second. Then test as in 10.3 above to confirm that the foot pedal is ON.

11. Synchronizing the Cordless Handpiece and Optional Wireless Foot Pedal

Follow this procedure when:

- the optional wireless foot pedal or cordless handpiece inner module is replaced for any reason
- the cordless handpiece and optional wireless foot pedal do not seem to be communicating properly

11.1 Place the optional wireless foot pedal within 10 feet of the charging base for cordless handpiece.

11.2 Verify that the foot pedal is ON (when you press the foot pedal the charge light should turn either Green or Orange depending on its level of charge). If it is OFF press the power button approximately 1 second to turn it on.

11.3 Invert the optional wireless foot pedal.

11.4 Verify the charging base is plugged in and turned on (the power indicator on the back of the charging base should be green).

11.5 Place the cordless handpiece inner module on the charging base for cordless handpiece, ensuring that at least one of the three Charge Indicator lights on the handpiece are Green. If the bottom charge indicator is solid Orange, the inner module may not have enough charge to be synchronized.


11.6 Depress the On/Off/Mode button on the cordless handpiece inner module twice within 3 seconds and release.

11.7 When the Mode indicator and bottom charge Indicator begin blinking Blue alternately, the cordless handpiece inner module is ready to sync.

11.8 Depress the Sync button on the selected optional wireless foot pedal for at least 2 seconds and release when handpiece mode indicator and bottom charge indicator blink green.

11.9 A successful sync is indicated by the cordless handpiece inner module Mode Indicator and bottom Charge Indicator blinking Green simultaneously and the Sync Indicator on top of the optional wireless foot pedal turning solid Blue.

11.10 A failed sync is indicated by the cordless handpiece inner module Mode Indicator and bottom Charge Indicator blinking Orange simultaneously.

11.11 Synchronization Indicator Status (Cordless Handpiece and Optional Wireless Foot Pedal Sync Mode)

Mode Indicator Color	Functional Light Description
Blinking Blue	Ready to Sync: Mode Indicator and bottom Charge Indicator on cordless handpiece and Sync Indicator on top of selected foot pedal will blink Blue
Blinking Green	Sync was successful: Mode Indicator and bottom Charge Indicator on cordless handpiece will blink Green
Blinking Orange	Sync was not successful: Mode Indicator and bottom Charge Indicator on cordless handpiece will blink Orange

12. Preparation for Use

12.1 Position the optional wireless foot pedal for use, ensuring that the floor is level.

12.2 Ensure that the outer sheath has been sterilized according to the Infection Control Procedures (Section 8).


12.3 When it is sufficiently charged, insert the cordless handpiece inner module into a Disposa-Shield® disposable barrier. Ensure that the nose of the cordless handpiece inner module completely clears the opening of the Disposa-Shield® disposable barrier.

PRECAUTION: Do not install the covered cordless handpiece inner module into the outer sheath without ensuring the Disposa-Shield® disposable barrier is clear of the cordless handpiece inner module nose. Do NOT apply the Disposa-Shield® disposable barrier to the outer sheath, which must be sterilized before each use.

The Polyethylene Barrier does not replace cleaning and disinfection of the dental instruments. Please clean and disinfect dental instruments after each patient as described by the instrument manufacturer.

12.4 Install the cordless handpiece inner module into the outer sheath by aligning the angle on each and snapping in place.

12.5 Attach a Nupro Freedom® Disposable Prophy Angle to the cordless handpiece by aligning the swoops on each and snapping in place.

12.6 Verify that all parts of the cordless handpiece are securely attached before use.

Fully-Assembled Handpiece

12.7 Charge Status: Cordless Handpiece Charge Level

Charge Indicator Color (when outside charging base)	Light Functional Description
Three Solid Green	Cordless Handpiece Charge is between 75%-100%
Two Solid Green	Cordless Handpiece Charge is between 50%-74%
One Solid Green	Cordless Handpiece Charge is between 25%-49%
Solid Orange	Cordless Handpiece Charge is less than 24%
No Color	Cordless Handpiece is - turned off, or - in Sleep mode

13. Using SmartMode® Technology and/or Optional Wireless Foot Pedal

Follow standard prophylaxis procedures as you would with any corded prophy device.

Applying excessive force to the polishing surface may slow or stop the rotation of the DPA and result in ineffective polishing. You only need to apply a light pressure during polishing, enough to flare the lip of the rubber cup slightly subgingivally and into the interproximal spaces. Avoid direct contact on gingival tissues as speed increases with pressure.

The DPA (motor) will only rotate in one direction.

Sleep Mode (if cordless handpiece is not in a powered charging base):

After 1 minute of non-use, the cordless handpiece will enter a Sleep Mode and the charge indicators and mode indicator will turn off. Upon picking up the handpiece within 1 hour, the Mode Indicator and charge indicators will illuminate and the unit will enter Standby Mode.

After 1 hour of inactivity, the cordless handpiece will turn off to conserve battery charge. After this, the On/OFF/Mode Button will need to be pressed to bring the handpiece into Standby Mode.

13.1 Attach the sterilized Outer Sheath to the cordless handpiece inner module and then attach the DPA.

13.2 To operate the cordless handpiece in SmartMode:
With the cordless handpiece in Standby Mode (Mode Indicator solid Green), press the On/Off/Mode button quickly and verify that the Mode Indicator changes from solid Green to solid Purple to indicate that the unit is in SmartMode. Note that the DPA cup/tip/brush will begin to rotate at a low baseline speed of 500 rpm.

13.2.1 Press the cup of the DPA into the prophy paste applying a downward force and light pressure.

13.2.2 SmartMode will allow you to control the speed of the DPA cup based upon the pressure applied from the DPA cup to the tooth. By slightly increasing the pressure, the speed of the DPA cup will increase. By reducing pressure, the speed of the DPA cup will decrease. When the DPA cup is removed from the tooth surface, the DPA cup will return to a low baseline speed.

13.2.3 You may exit SmartMode (stop the DPA cup from spinning) by pressing the On/Off/Mode button or by quickly tapping the optional wireless foot pedal. A solid Green Mode Indicator shows that the cordless handpiece is in Standby Mode. Should the optional wireless foot pedal be synchronized to the handpiece, you may then control the rotation speed using the optional wireless foot pedal.

13.2.4 To enter SmartMode again, press the On/Off/Mode button quickly. The Mode Indicator will change from Green to Purple and the DPA cup will begin spinning at a low baseline speed.

13.3 To operate the cordless handpiece with the optional wireless foot pedal:
Set the handpiece in Standby Mode (with the Mode Indicator solid Green).


13.3.1 Apply pressure on the optional wireless foot pedal to control the rotation speed of the DPA throughout the procedure (mode Indicator turns solid Blue). By slightly increasing the pressure, the speed of the DPA cup will increase. By reducing pressure, the speed of the DPA cup will decrease.

13.3.2 To return the cordless handpiece to Standby Mode, remove pressure from the optional wireless foot pedal. A solid Green Mode Indicator will show on the cordless handpiece.

13.3.3 To move the cordless handpiece to SmartMode mode, remove pressure from the optional wireless foot pedal. A solid Green Mode Indicator will show on the cordless handpiece. Press the On/Off Mode button on the cordless handpiece, and the Mode Indicator will turn solid Purple to show that the unit has entered SmartMode mode.

13.4 Mode Status (Cordless Handpiece)

Whenever the cordless handpiece is in normal operational mode (i.e., not being charged or synchronized), the Mode Indicator on the cordless handpiece displays information about the system according to the table below:

NOTES:

- If the cordless handpiece seems to be losing power, recharge as soon as possible.
- If the cordless handpiece battery is completely discharged, the battery indicator will blink and the unit will shut down. The handpiece must be put in the cordless handpiece charging base before it's able to operate again.

Cordless Handpiece Mode	Mode Indicator Color	Functional Light Description
Power Off (when outside of charging base)	No Color	Unit is - turned off, or - in Sleep mode
Standby	Solid Green	Unit is ready for use and is in Standby Mode
SmartMode® Operation	Solid Purple	Unit is in SmartMode®
Foot Pedal Operation	Solid Blue	Foot Pedal is depressed and the prophy cup is spinning in response to pressing a synchronized optional wireless foot pedal
Temperature Indicator	Solid Orange	Internal temperature of the Handpiece is higher than normal. The light will turn off when the temperature returns to normal. Operator should refrain from using the handpiece for several minutes until the solid Orange Mode Indicator goes out. If the solid Orange Mode Indicator persists, call technical support

13.5 To turn the handpiece off:
With the handpiece in Standby Mode (solid Green Mode Indicator), press and hold the On/Off/Mode button for approximately one and one half second until the charge indicators turn off, and then release the button.

14. Specifications



Power Supply manufacturer: Power Supply model number:	GLOBTEK WR9QA1200USBNMEDRVB	
AC Input	Continuous (100-240 VAC)	
AC Input Current	Less than 0.7A	
AC Input Phases	Single	
AC Input Frequency	50-60 Hz	
DC Output Power	5W	
DC Output Voltage and Current	+5VDC at 1A	
Output Regulation	+/- 10%	
Weight	Handpiece with metal sheath = 120 g Foot Pedal = 200 g	
Dimensions	Handpiece with sheath & disposable angle L = 190 mm, W = 30 mm Foot Pedal W = 118 mm H = 40 mm	
Optional Wireless Foot Pedal	Protection Class IPX1. Not for operating theatres.	
Remote Communication	Frequency: Power: Channels:	2405-2480 MHz 1mW 40 BLE
Operating Environment	Ambient temperature: Relative Humidity: Atmospheric Pressure: Altitude:	10°C-40°C 45-95% (non-condensing) 80-106 kPa ≤ 2,000 meters
Transport and Storage Conditions	Temperature: Relative Humidity: Atmospheric Pressure: Altitude:	-20°C to 50°C 45-95% (non-condensing) 50-106 kPa ≤ 5,000 meters

15. Classification

Type of protection against electric shock:	Class II
Degree of protection against electric shock:	Type B Applied Part
Mode of operation for Handpiece:	Non-continuous: 5 minutes ON, 25 minutes OFF
Mode of operation for Foot Pedal:	Continuous
Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:	Equipment not suitable for use in the presence of flammable mixtures
According to medical device directive:	IIA (Rule 9) (IEC 60601) UL 60601-1
Pollution Degree Classification	Pollution Degree 2
Overvoltage Category	Category II (connected to wall outlet)

16. Symbol Identification

The following standard symbols appear on the device label.

	Class II Equipment
	Type B applied part
	MEDICAL EQUIPMENT WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA-C22.2 No. 601.1, ANSI/AAMI ES60601-1 (2005, 3rd ed.), CAN/CSA-C22.2 No. 60601-1 (2008),13VA
	Consult instructions for use
	Sterilizable up to the temperature specified
	Do not re-use (for DPAs)
	Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union
IPX0	Protection Class IPX0 IPX0 Classification of ingress of water for Charger - not protected
IPX1	Footswitch not for operating theatres Protection Class IPX1 IPX1 Classification of ingress of water
IPX3	Protection Class IPX3 IPX3 Classification of ingress of water for Inner Module - Protected against falling spray.
	Duty Cycle for handpiece: 5 minutes ON 25 minutes OFF
SN	Serial Number
LOT	Batch Code/Lot Number
CE	This symbol is a mandatory marking for devices entering the European market to indicate conformity with the essential health and safety requirements set out in European Directives. The symbol may be accompanied by a four-digit identification number of the notified body.
 5V-1A	Direct Current (DC) supply. 5 Volts, 1 Amp
	Do Not Autoclave (Inner Module)
	Interference may occur in the vicinity of equipment marked with the following symbol

17. Disposal of Unit

U.S. - Dispose of the system components in accordance with state and local laws.
 EU - Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union.

18. Electro-Magnetic Compatibility Precautions

Guidance and manufacturer's declaration - electromagnetic emissions		
The Nupro Freedom® is intended for use in the electromagnetic environment specified below. The Customer or the user of the Nupro Freedom® should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Nupro Freedom® uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Nupro Freedom® is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class B	Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Nupro Freedom® or shielding the location.
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Nupro Freedom® is intended for use in the electromagnetic environment specified below. The customer or the end user of the Nupro Freedom® should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electromagnetic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±2, ±4, ±6kV contact ±2, ±4, ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±0.5, ±1, ±2kV mains No I/O lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±0.5, ±1kV Line to Line ±0.5, ±1, ±2kV Line to Earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5 sec	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Nupro Freedom® requires continued operation during power mains interruptions, it is recommended that the Nupro Freedom® be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U _T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic emissions			
The Nupro Freedom® is intended for use in the electromagnetic environment specified below. The customer or the user of the Nupro Freedom® should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the MidWest RDH Freedom, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = [3.5/V1] fflP d = [3.5/E1] √P 80 MHz to 800 MHz d = [7/E1] √P 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Nupro Freedom® is used exceeds the applicable RF compliance level above, the Nupro Freedom® should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Nupro Freedom®. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

Recommended separation distance between portable and mobile RF communications equipment and the model _____			
The Nupro Freedom® is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Nupro Freedom® can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Nupro Freedom® as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

19. Troubleshooting

Problem	Solution
DPA is not revolving.	<ol style="list-style-type: none"> 1. Ensure the outer sheath and the DPA are snapped together securely. 2. Ensure that the cordless handpiece inner module and the outer sheath are snapped together securely. 3. Verify the cordless handpiece is powered up and properly charged. If the cordless handpiece lights do not illuminate, place the cordless handpiece in the charging base for cordless handpiece for a minimum of 5 seconds, and then remove it to use. 4. If using an optional wireless foot pedal and it is not plugged into the Nupro Freedom Power Supply, make sure that the foot pedal is not off or discharged (see section 10). Or, you may plug the optional wireless foot pedal into the Nupro Freedom Power Supply and the wireless foot pedal will operate via the power supply. Then (with the foot pedal powered) confirm that the foot pedal is synchronized to the cordless handpiece by synchronizing them (see section 11). The optional wireless foot pedal will operate while charging. 5. Ensure the DPA is not damaged by removing the DPA and spinning the cup between your fingers. The cup should spin freely. 6. Ensure that the outer sheath is not damaged by removing the cordless handpiece inner module from the outer sheath, leaving the DPA connected, and spinning the DPA cup between your fingers. The cup should spin freely. If the cup does not spin freely, place 1-2 drops of an FDA-cleared lubricant into the outer sheath nose and try to spin again. If the cup still does not spin freely, call Technical Support. If lubricant is applied, the outer sheath must be sterilized prior to use. Follow sterilization instructions in Section 8 (Infection Control Procedures). 7. If DPA still does not spin freely, remove the outer sheath from the cordless handpiece inner module and verify the cordless handpiece inner module motor spins when the optional wireless foot pedal is depressed. If the cordless handpiece inner module motor does not spin, resynchronize the optional wireless foot pedal with the cordless handpiece inner module as per instructions contained in this manual. 8. Resynchronize the units. See Section 11 Synchronizing the Cordless Handpiece and Optional Wireless Foot Pedal. If the synchronization is not successful and/or the cordless handpiece still does not spin, call Technical Support.
Excessive noise or vibration during operation.	<ol style="list-style-type: none"> 1. Ensure that the outer sheath is aligned correctly with the cordless handpiece inner module. 2. Try a new disposable prophylaxis angle (DPA). 3. Check the components for gross debris or contamination and adhere to all infection control procedures. 4. Check for damaged, worn or broken components. 5. Call Technical Support.
Optional Wireless Foot pedal does not charge.	<ol style="list-style-type: none"> 1. Ensure that the power cord is securely attached to the wall outlet and the optional wireless foot pedal. 2. Remove the optional wireless foot pedal access door and look in the compartment; ensure that the LED in the bottom compartment is illuminated. If not, call Technical Support.
System intermittently shuts down.	<ol style="list-style-type: none"> 1. Handpiece may shutdown in the presence of strong non-compliant Radio Frequency generating components (Interference). 2. Call Technical Support.

Cordless Handpiece does not have the correct speed control.	<ol style="list-style-type: none"> 1. Excessive force is being applied to the polishing surface, reducing the DPA rotation. Use a lighter force and follow standard prophylaxis procedures. 2. Cordless handpiece or optional wireless foot pedal may be malfunctioning. 3. Call Technical Support.
Difficulty removing outer sheath from cordless handpiece inner module.	<ol style="list-style-type: none"> 1. Check the components for gross debris. 2. Hold the cordless handpiece outer sheath securely and twist the cordless handpiece inner module. 3. Inspect parts for wear. 4. Call Technical Support.
Charging base for cordless handpiece does not drain liquids.	Clear any debris from the hole at the bottom of the charging base for cordless handpiece.
Cordless handpiece does not charge.	<ol style="list-style-type: none"> 1. Clean the charge contacts on the cordless handpiece and charging base for cordless handpiece using one of the approved cleaning solutions described in section 8. Do not spray cleaning solutions directly onto parts. 2. Verify that the power supply is properly connected to the charging base for cordless handpiece and that the green LED on back lights up. 3. Verify that the cordless handpiece is able to properly sit inside the charging base for cordless handpiece, and that there are no foreign obstructions. 4. If still not charging, call Technical Support.
Cordless handpiece does not hold charge.	<ol style="list-style-type: none"> 1. Verify that the cordless handpiece properly charges (LEDs scroll when in the charger). 2. Call Technical Support.
Orange Light on Cordless Handpiece Mode Indicator illuminates.	<ol style="list-style-type: none"> 1. Handpiece is heating up due to excess ON time or load. Do not use handpiece for several minutes and allow it to cool down. Avoid excessive load by applying a light pressure and following standard prophylaxis procedures. 2. If LED still illuminates even after allowing handpiece to cool down, set aside and call Technical Support.
Power supply overheating.	Immediately unplug the unit and call Technical Support.
Power cord is frayed or damaged in any way.	Do not use. Call Technical Support.

NOTE: For current Technical Support contact information, see "Manufacturer contact" in the control procedures (section 8).

20. Accessories

Nupro Freedom® Cordless Prophy System with SmartMode® Technology

Item Description	DENTSPLY Part Number
Outer Sheath	90796
Cordless Handpiece Inner Module, Lavender	90799
Charging Base for Cordless Handpiece, Lavender	90797
Wireless Foot Pedal	90742
Power Supply-Domestic	90798
Cradle for Cordless Handpiece	9070801
Disposa-Shield for Cordless Prophy System Inner Module, 500 ct. box	A88019

Nupro Freedom® Disposable Prophy Angles and Nupro Freedom® Prophy Packs

The Nupro Freedom® Cordless Prophy System with SmartMode® Technology is compatible with a wide range of Nupro Freedom® Disposable Prophy Angles and Nupro Freedom® Prophy Packs.

For a full listing, please see our web site at www.dentsply.com or contact your local DENTSPLY Customer Support Center.

21. Limited Warranty

The DENTSPLY Professional Nupro Freedom® Cordless Prophy System with SmartMode® Technology is designed exclusively for dental use and this warranty is not applicable to other uses. This warranty extends to Nupro Freedom® systems purchased from an authorized DENTSPLY distributor, and only to the original purchaser.

All components of the Nupro Freedom® Cordless Prophy System with SmartMode® Technology, except the Disposable Prophy Angle, which is a single use only item, and the Cordless Handpiece Inner Module are warranted for (1) year from the date of purchase against defects arising from faulty materials and workmanship. The Cordless Handpiece Inner Module has a (3) three year warranty.

Parts will be repaired or replaced at the discretion of DENTSPLY Professional provided that the system has been operated and maintained as prescribed in these instructions and has not been subjected to apparent misuse, abuse or accident. Claims covered by this warranty will be honored when presented through your DENTSPLY Professional distributor within thirty (30) days from discovery of defect within the applicable warranty period.

There are no warranties, express or implied, which extend beyond the description on the face hereof DENTSPLY neither assumes, nor authorizes any person to assume for it, any other liability in connection with the sale or use of its products. Damages are limited strictly to repair or replacement of parts. DENTSPLY expressly disclaims liability for incidental and consequential damages resulting from the use of the products.

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