

hurrimix™ Automatic Mixer
for Alginate and Stone

Mix with Confidence

The hurrimix™ increases office efficiencies
by standardizing the mixing process
to deliver a top quality mix every time.

1.800.645.5530 gacintl.com



The Ultimate in Style and Durability

The hurrimix™ features a stylish, stainless steel body, electronic microprocessor controller, and waterproof controls protected by a polycarbonate membrane.

Capacity and Convenience

The maximum container capacity is 80g powder and water for alginates and 100g for stones. The mixing containers have rounded internal surfaces to easily remove the mixed material and for cleaning.

Software with 20 Mixing Programs

The hurrimix™ is programmable with 10 options for alginates and 10 programs for stones. Users can modify any programs at any time. You can set programs for specific alginate and stone products resulting in a fast, reliable mixing process, with no time wasted on manually programming the mixer, and guaranteed correct mixing times.

Door Locking Device

A built in safety mechanism ensures the lid is properly closed to prevent accidental opening and any loss of material during mixing.

Large Digital Display

An extremely intuitive, large digital display for simple, rapid operation.



The hurrimix™ provides you with highly consistent quality mixes, while reducing mix time by as much as 80% vs hand mixing — a significant savings. The number of rotations achievable with an automatic mixer are not achievable by mixing manually.



Mix Alginate and Stone in 5 Easy Steps

hurrimix™ delivers consistent mixes with minimal effort from the operator.



1. Add Powder



2. Add Water



3. Insert Container



4. Start Machine



5. Remove Mix

Consistent Mix Quality — Manual Mixing vs. Automatic Mixing

Manual mixing can leave bubbles in the mixture and can cause distortion and other imperfections, compromising the quality of the finished work.



Manual Mixing:
Revolution max 180 rpm



Automatic Mixing with the hurrimix™
Revolution \approx 600 rpm
Rotation \approx 3000 rpm



Schedule Your FREE Lunch and Learn with GAC

We'll bring lunch to your office and provide a hands-on demonstration for you and your staff.*

To sign up for your FREE GAC Lunch and Learn, visit:

dentsplysirona.com/orthodontics/ortholearn

1-800-645-5530

The hurrimix™ Includes: 1 spatula, 1 pre-mix stick, 1 spatula holder, 4 mixing containers (2 Alginate and 2 Stone), 1 water bottle with a graduated measuring cup, Powercord and DFU.

*Some geographical restrictions may apply. RTE-020-17

Dentsply GAC International • One CA Plaza, Suite 100 • Islandia, NY 11749