In-Ovation[®] C In-Ovation[®] Czr

Orthodontic Self-Ligating Ceramic Brackets

DIRECTIONS FOR USE
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
GEBRUIKSAANWIJZING

Manufacturer:
DENTSPLY GAC International
One CA Plaza, Suite 100
Islandia, NY 11749 USA
Tel: 800-645-5530
Fax: 631-357-8793
www.dentsplysirona.com/orthodontics

SINGLE
USE ONLY

Rx
ONLY

Emergo Europe
Prinsessegracht 20
2514 AP, The Hague
The Netherlands

Made in USA
#DFU-0115-US Rev. 01
Issued 06/17

In-Ovation® C / In-Ovation® Czr

ENG

Orthodontic Self-Ligating Ceramic Brackets

INTENDED USE:

Intended to be used in orthodontic treatment of natural teeth.

INDICATIONS FOR USE:

Indicated for patients requiring orthodontic treatment by movement of natural teeth.

In-Ovation C brackets are indicated for use on natural teeth, excluding mandibular bicuspid teeth.

In-Ovation Czr brackets are indicated for use on natural, mandibular bicuspid teeth only.

CONTRAINDICATIONS:

These brackets should not be bonded to restorative or compromised enamel surfaces.

In-Ovation C brackets are contraindicated for patients with an allergy to Nickel.

⚠️ WARNINGS:

In-Ovation C/In-Ovation Czr ceramic brackets are for single use only. Risks of reuse include infection and cross-contamination to patient and healthcare provider as well as other wear-related irregularities.

Due to the hardness of ceramic brackets, bonding brackets in occlusion should be avoided to prevent wearing of enamel surfaces during all phases of treatment.

Deep bite cases should be opened prior to bonding with brackets to prevent wearing of tooth enamel by the ceramic bracket. Bite ramps and other intrusion mechanics may also be necessary to prevent bracket to tooth contact.

Bonding of ceramic brackets to compromised teeth (i.e., with large restorations, to peg laterals, or preexisting pathological conditions) can increase the risk of tooth damage.

Device contains nickel and may cause allergic reactions.

PRECAUTIONS:

For professional use only. U.S. Federal law restricts this device to sale by or on the order of a Dentist.

In-Ovation C/In-Ovation Czr ceramic brackets are capable of withstanding all normal torque requirements. However, care should be taken when making torquing activations. Large corrections with full size stainless steel wires may result in bracket failure and should be avoided.

Instruct patients not to chew or bite hard substances such as hard candy, ice or carrots. Careful and thorough patient instruction is a key to avoiding appliance or enamel damage.

Bonding to porcelain crowns or facings may cause chipping or breakage of the crown or facing during treatment or debonding.

Give tooth brushing instruction to patients, because food residue and tartar attached to the bracket may cause clip movement problems.

Open or close the clip after removing the food debris from the clip groove by water pick, or clip may be deformed.

Do not exert excessive force on brackets, it may cause deformation of the bracket or bond failure.

Please remove the color ID after placing the bracket.

Take care so as not to have the incisal edge or cusp of the counter tooth contact the ceramic bracket, or the enamel may get worn off.

Debonding of In-Ovation C/In-Ovation Czr brackets can be done using a direct bond bracket remover such as an ODG1001 (Pin & Ligature Cutter). No other existing ceramic bracket-debonding instrument should be used to debond brackets since it may result in fractured tie wings or tooth damage.

To ensure successful bonding, always use gloved hands and/or tweezers when handling brackets. The natural oils in skin may compromise the properties of the bonding base.

Used bracket is likely to cause deformation, deterioration, or reduction of strength and induce infection.

ADVERSE REACTIONS:

Allergic reaction.

ORTHODONTIC CERAMIC BRACKET IDENTIFICATION:

*In-Ovation C/In-Ovation Czr brackets are easily identifiable by a color-coded dot on the distolingual tie wing, as well as the site gauge found in the bracket slot.


UR & UL CENTRALS:	BLUE
UR & UL LATERALS:	PINK
UR & UL CUSPIDS:	LT. GREEN
UR & UL BICUSPIDS:	VIOLET
LOWER ANTERIORS:	ORANGE
LR & LL CUSPIDS:	LT. BLUE
LR & LL FIRST BICUSPIDS:	BLACK
LR & LL SECOND BICUSPIDS:	RED

BRACKET BONDING:

- To ensure successful bonding, always use tweezers when handling brackets. The natural oils in skin may compromise the properties of the bonding base.
- In-Ovation C/In-Ovation Czr brackets feature a mechanical lock base. No special primers or pretreatments are necessary.
- Excessive adhesive flash should be removed when the bracket is bonded to minimize staining of the tooth by the adhesive.
- In-Ovation C/In-Ovation Czr brackets can be used with traditional direct or indirect bonding methods. No change in bonding technique is necessary. For bonding procedures please follow the adhesive manufacturer's recommendation.
- When placing In-Ovation C/In-Ovation Czr brackets, hold the bracket mesio-distally. If held occluso-lingually the clips may close.
- When bonding, take precaution not to let the bonding material adhere to the clip.

SITE GAUGE REMOVAL:

- ONCE THE ADHESIVE HAS CURED, REMOVE THE SITE GAUGE PRIOR TO INSERTING THE ARCH WIRE.
- To remove the site gauge, use a scaler type instrument to lift up one end of the site gauge. (Fig. 1)
- Once one end has been lifted out from the slot, continue to peel the site gauge out of the slot. The site gauge will come out of the slot uniformly & completely. (Fig. 2)

(Site gauge removal Fig. 1)

(Site gauge removal Fig. 2)

BRACKET REBONDING PROCEDURE:

If a spontaneous bond failure occurs, it may be necessary to rebond a bracket. The following steps are recommended:

- Carefully inspect the bracket for any damage.
- Contour cured adhesive on bracket base with a bur.
- Leave a thin layer of adhesive on the bracket base when removing any excess adhesive. An adhesive-to-adhesive bond is necessary when rebonding In-Ovation C/In-Ovation Czr Czr brackets.
- If bracket has been contaminated (e.g. moisture), rinse the bracket in isopropyl alcohol or acetone and allow it to dry.
- Prepare the tooth surface and bond the bracket using the procedure as described by the adhesive manufacturer.


OPERATING INSTRUCTION:

- For a wire change, open the clip by inserting the tip of any pointed dental instrument such as a scaler into the facial hole in the clip of the bracket. Apply a gingival to occlusal sliding force.
- Alternatively, open the clip by applying gingival to occlusal sliding pressure to the rounded end of the clip with an In-Ovation C/In-Ovation Czr opening instrument.

SUGGESTED DE-BONDING PROCEDURE:

- NOTE: Failure to follow this debonding procedure may result in tooth damage.
- Remove adhesive flash from the occlusal and gingival debonding channels found on the In-Ovation C/In-Ovation Czr ceramic brackets.
- Visually check instrument tips for sharpness. To ensure durability and sharpness of instrument, it is recommended that the debonding instrument should be used exclusively with In-Ovation C/In-Ovation Czr brackets. We recommend using the ODG1001 (pin & ligature cutter) tool for debonding In-Ovation C/In-Ovation Czr Brackets.
- Place the debonding instrument occluso-lingually in the debonding channels and squeeze. Torquing of bracket is not indicated during debonding procedure.
- NOTE: Remove any debris from the pliers' jaws after debonding bracket to ensure even contact and force distribution on the next removal.

ODG1001

STORAGE:

- Keep away from non-professionals.
- Keep in location without risk of negative effect by air containing a large amount of dust, sulfur, salt, and moisture, and without external pressure against the products.
- Store away from direct sunlight, heat and moisture.

WARRANTY AND LIMITED REMEDY

Manufacturer warrants that this product is free from defects in materials and manufacture.

Manufacturer's sole obligation and customer's sole remedy in the event of a claimed defect shall be limited to replacement merchandise or refund of the purchase price.

All implied warranties, including implied warranties of merchantability or fitness for any particular purpose or use are disclaimed.

LIMITATION OF LIABILITY

Except where prohibited by law, Manufacturer will not be liable for any loss or damages arising from this product, whether direct, indirect, special, incidental, or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability. This limitation does not apply to third party personal injury claims.

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