52.5 mm

fold

EN BioBiteCorrector™ MS Directions for Use

The BioBiteCorrector™ MS (BBC) is a fixed rigid Class II appliance working with the functional principles of the Herbst-appliance.

INDICATIONS:

- Correction of Class II malocclusions and functional mandibular midline shifts
- Use as anchorage appliance for closing a gap in the posterior teeth of the lower jaw while simultaneously correcting a Class II malocclusion
- Avoidance of dysgnathic surgeries

CONTRAINDICATIONS:

- Severe protrusion of lower incisors

MATERIAL:

- The BioBiteCorrector is only manufactured of titanium grade 5
- The fixing screws and the distance spacers are manufactured of titanium
- The screwdriver is manufactured of stainless steel (sterilizable and re-usable)

The BioBiteCorrector is determined for single use only.

THE FIXATION OF THE BBC:

Bands on the upper molars are not required to fix the BioBiteCorrector.

1. Preconditions

- Precondition for the integration of the BioBiteCorrector is a multi-bracket appliance with a steel archwire with a minimum size of .017" x .025" in the upper and lower jaws.
- In the lower jaw the steel archwire must have a cinch back to avoid a proclination of the incisors. For the cinch back we recommend a Twister-Instrument. With the Twister you can bend the arch very simply without risking to remove the bracket from the tooth.

2. Choice of size

- Standard size is used in non-extraction cases.
- Only in extraction cases in the upper jaw small size is used.

3. Fixation of the BBC

- The BBC is already pre-assembled. The appliance is distinguished in left and right sides. First the BBC is fixed in the upper jaw and afterwards in the lower jaw.

4. Fixation in the upper jaw

- The upper jaw connecting element is fixed between the first molar and second premolar. The connecting element should have no contact with the bracket of the first molar. The screwdriver should be assured with dental floss to avoid aspiration or swallowing.
- Slide the body on the archwire from lateral and then fix the screw.

5. Fixation in the lower jaw

- After fixing the upper jaw connecting element, the lower jaw connecting element must be fixed between canine and first premolar. Slide the body on the archwire from lateral and then fix the screw. The connecting element should be situated as distal as possible mesially to the bracket of the first premolar. The connecting element must not have any contact with the canine bracket at all. Fix the other side in the same way.

6. Adjustment in neutral occlusion

- After both sides are fixed, guide the lower jaw back as far as possible. If there is still distal occlusion the lower jaw can be shifted forward by using distance spacers. These spacers can be crimped on the appliance, whereby the lower jaw can be shifted into neutral occlusion.
- The distance spacers are attached to the thinnest telescope pipe. Fine corrections can be carried out by shifting the connecting element.

The product range description is included on the product website.

The average duration of treatment with the BioBiteCorrector is approximately 4-7 months.

LIABILITY

- There is no liability for unsatisfying results. Do not use the appliance for other than dental or orthodontic purposes. The dentist or orthodontist is responsible for the use of the appliance.

SAFETY INSTRUCTIONS

- The firmness of the screws must be checked by the doctor in regular intervals to recognize loosened screws and to prevent a loss of the appliance with the danger of swallowing or aspirating the appliance or parts of it. The patient must be instructed for self-control.

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