Lucitone Digital Print Denture

Materials for Carbon® M-Series Printers

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Lucitone Digital Print™ 3D Denture Resin is a light-cured resin indicated for the fabrication of denture bases in dental laboratories, including full dentures, and implant overdentures.

Lucitone Digital Try-In™ 3D Trial Placement Resin is used as a try-in material for evaluation prior to fabrication of the final restoration.

- Fabrication of these prostheses require a computer-aided design and manufacturing (CAD/CAM) system using an additive printer.

Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent is indicated for use in enhancing the bond of denture teeth to denture base.

Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin is utilized as an aid in bonding denture base to denture teeth as well as repair using traditional techniques.

Lucitone Digital Fuse™ Step 3 – 3D Denture Sealer is a light-cured sealant that produces a smooth, glossy surface finish on the denture.

Requirements for Denture Fabrication

- > File in STL Format
- > Lucitone Digital Print™ 3D Denture Resin
- > Lucitone Digital Try-In™ 3D Trial Placement Resin
- > IPN 3D™ Digital Denture Teeth
- > Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent
- > Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin
- > Lucitone Digital Fuse™ Step 3 – 3D Denture Sealer
- > Ultrasonication Bath or Orbital Shaker
- > Isopropyl Alcohol (≥ 99%)
- > UV Tack-Cure Light
- > inLab Speedcure™ Processing Unit
- > Electric Hot Plate
- > Metal Container Kit
- > Flocked Applicator Tips
- > Biobond Brush, #4 Flat
- > 4x4 Gauze Pads
- > Pipette
- > Small Plastic Container (see Figure 5)
- > Carbon® 3D M-Series Printer Software
- > Carbon® 3D M-Series Printer

CAUTION: Any unauthorized changes to the process, workflow, materials, equipment, parameters, or software may result in a device that is out of specification and is not recommended. Contact Dentsply Sirona for a list of compatible components.

Contraindications

Lucitone Digital Print Denture components are contraindicated for patients with known hypersensitivity or severe allergic reactions to any of the components.

Warnings

1. Lucitone Digital Print Denture materials contain components which may cause skin dryness or irritation, sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wear protective gloves, eye protection and protective clothing when handling.
Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.
 2. Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent - liquid is highly flammable. Keep container tightly closed and use/store away from sources of ignition. Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization or other allergic reactions.
 3. Lucitone Digital Fuse Step 3 - 3D Denture Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. Suspected of damaging fertility or the unborn child. Harmful to aquatic life with long lasting effects. Avoid release into the environment. Dispose contents.
 4. Lucitone Digital Print Denture materials - Avoid inhalation or ingestion. High vapor concentration might induce headache, nausea and irritation of eyes and respiratory system. Work in a well-ventilated area. Avoid contact with eyes. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effect.
Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
Ingestion: Contact your regional Poison Control Center immediately.
 5. Proper ventilation and personal protective equipment should be used when grinding devices as the particulate generated during grinding may cause respiratory, skin and/or eye irritation.
 6. Isopropyl Alcohol is highly flammable. Keep containers tightly closed and use/store away from sources of ignition.
-

Precautions

1. These products are intended to be used only as specifically outlined in the Instructions for Use (IFU). Any use of these products inconsistent with the IFU is at the discretion and sole responsibility of the practitioner.
2. Store resins and system components at room temperature 60°F - 80°F (16°C - 27°C) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date. Dispose of contents and containers in accordance with local and national regulations.
3. Lucitone Digital Print™ 3D Denture Resin, Lucitone Digital Print™ 3D Trial Placement Resin, and Lucitone Digital Fuse™ Step 2 - 3D Denture Bonding Resin:
When washing a 3D printed appliance with a solvent or polishing the appliance, it should be done in a properly ventilated environment with proper protective masks and gloves. Expired or unused resins shall be disposed in accordance with local regulations.
4. Lucitone Digital Fuse™ Step 1 - 3D Tooth Conditioning Agent:
Use in a ventilated workspace. Handle teeth with solvent-resistant gloves or tweezers. Product is light sensitive. Minimize exposure to light and keep the lid on the container to maintain its effective use. Product is very volatile and flammable. Keep away from heat, sparks and flame. Do not add Lucitone Digital Fuse™ Step 1 to a warm metal container or

- while the container is on the electric hot plate.
- Lucitone Digital Fuse™ Step 3 – 3D Denture Sealer:
Use in a ventilated workspace. Keep bottle tightly closed when not in use. Lucitone Digital Fuse™ Step 3 is volatile and highly flammable. Keep away from heat, sparks and flames. Product is light sensitive. Minimize exposure to light and keep the lid on the container to maintain its effective use.
 - The following disinfectants have been used successfully with Lucitone Digital Print Denture materials: Glutaraldehyde-based, Quaternary ammonium chlorides, Quaternary ammonium chlorides/alcohol based and Phenolic sprays. Use of other disinfectant products may not be suitable. Common cleansers that have been used successfully include Efferdent Denture Cleanser and Polident Denture Cleanser.

Adverse Reactions

- Lucitone Digital Print Denture materials contain components which may cause skin dryness or irritation, sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wear protective gloves, eye protection and protective clothing when handling.
Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.
- Lucitone Digital Fuse Step 1 – 3D Tooth Conditioning Agent liquid is highly flammable. Keep containers tightly closed and use/store away from sources of ignition. Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization or other allergic reactions.
- Lucitone Digital Fuse Step 3 - 3D Denture Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. Suspected of damaging fertility or the unborn child. Harmful to aquatic life with long lasting effects. Avoid release into the environment. Dispose contents.
- Lucitone Digital Print Denture materials: Avoid inhalation or ingestion. High vapor concentration might induce headache, nausea and irritation of eyes and respiratory system. Work in a well-ventilated area. Avoid contact with eyes. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects mouth must be considered by the physician/dentist.
Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
Ingestion: Contact your regional Poison Control Center immediately.

Symbols on Product Labels

	Manufacturer		Keep Away from Sunlight
	Batch Code		CE European Conformity Marking
	Reorder Number		Medical Device
	Manufacture Date		U.S. Prescription Only
	Expiration Date		Store Between
	Consult Instructions for Use		Keep Dry
	Caution		

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

Step-by-Step Instructions

A. Procedure for Fabrication: PRINT

1. .stl file(s) should be created per scan and design manufacturer's instructions. Create Carbon M2 new project and upload .stl file(s) for the print job.

Select **"Lucitone Digital Print"** from the resin drop down menu when printing denture bases with pockets.

Select **"Lucitone Digital Try-In"** when printing denture try-ins.

Contact Dentsply Sirona Lab Technical Team for print orientation, support design and nesting guidelines at 800-243-1942 ext. 54212 or dsdigitaldentures@dentsplysirona.com

2. Based on software recommendation, pour resin volume slowly into the cassette.
3. Once cassette and platform are ready and door is closed, start 'Print'. The screen will change to display the print progress and estimate print time.

B. Procedure for Fabrication: WASH (TRY-INS AND BASES)

Note: Section B should be done under proper ventilation. Use of personal protective equipment is recommended.

1. Once print is complete, remove build platform from printer.
2. Remove appliances from build platform by snapping off each appliance by hand. Moderate pressure is required. (Figure 1).
3. Break off fence supports from the appliance (Figure 2). **Leave the bar supports in place.** (Figure 3 & 4). **Bar supports will be removed after post-cure.**

Figure 1

Figure 2

Figure 3

Figure 4

4. Add water to ultrasonication bath.

5. Place appliance pocket-down, if applicable, in a plastic container and wash with Isopropyl Alcohol ($\geq 99\%$ IPA) using a two-bath approach with fresh $\geq 99\%$ IPA each time. Using an ultrasonication bath, wash for two (2) minutes (1st cycle), followed by one (1) minute (2nd cycle). **Container should be closed when placed in the ultrasonication bath.** Wash each appliance individually. (Figure 5 & 6)

CAUTION: Prolonged exposure to $\geq 99\%$ IPA could cause changes to physical properties of the appliance. Follow wash times carefully.

When washing bases with pockets, remove from first bath cycle and use a wet brush (wet with $\geq 99\%$ IPA) to ensure tooth pockets are clean of residual resin. (Figure 7)

CAUTION: Residual resin in tooth pockets could affect fit and bond.


Figure 5

Figure 6

Figure 7

When washing a try-in, remove from first bath cycle and use a wet brush (wet with $\geq 99\%$ IPA) to wash occlusal surface and interproximal of teeth. (Figure 8)

CAUTION: Residual resin on occlusal surfaces could affect occlusion.

CAUTION: Do not wash dentures in “used” IPA.

Note: If using an Orbital Shaker, perform washing for three (3) minutes (1st cycle) followed by two (2) minutes (2nd cycle).

Figure 8

6. Blow dry the appliance using compressed air to ensure all parts are thoroughly dried before moving to tooth fusing/bonding.

C. Procedure for Fabrication: FUSE (BOND)

For post-curing a Try-In, skip to section E.

Note: Section C should be done under proper ventilation. Use of personal protective equipment is recommended.

1. Plug in the electric hot plate and set the temperature to 90°C.

Note: If the electric hot plate is not digital, set to 100°C.

Allow the electric hot plate to preheat for at least ten (10) minutes.

Note: The temperature setting may need to be adjusted slightly from the recommended setting due to several factors.

2. Place IPN 3D Digital Denture Teeth into respective pockets to check for a dry fit.
3. Shake the bottle of Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent to ensure contents are thoroughly mixed. Pour enough Fuse Step 1 into the round metal container to fill to top of the foam insert. (Figure 9) Refer to Precautions.

CAUTION: Using an inadequate amount of Fuse Step 1 will compromise tooth bond strength.

CAUTION: Keep the metal container and the foam insert free of any solid buildup. This can compromise the tooth bond strength.

4. Use tweezers to place up to fourteen (14) IPN 3D Digital Denture Teeth into the holes so that the teeth are positioned with the ridge lap down and in contact with the liquid. If the liquid is not level with the top of the foam insert, add additional Fuse Step 1 to metal container. Immediately replace the lid on the container. (Figure 10)
5. Place the tooth-filled container inside the metal ring of the electric hot plate until the 40°C rectangle on the temperature strip turns green (normally in four (4) to six (6) minutes). Refer to Precautions. (Figure 11)

CAUTION: Tooth bond strength will be compromised if 40°C is not reached or if higher warming temperatures are used to accelerate the process.

6. Carefully remove the metal container from the electric hot plate. Remove the lid from the metal container and use tweezers to retrieve the teeth. Allow teeth to air dry for minimum of two (2) minutes. Do not wipe teeth.

Figure 9

Figure 10

Figure 11

7. Place the lid back on the metal container immediately. Do not return the metal container to the electric hot plate until ready to treat the next set of teeth.
8. Apply a small amount of Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin in the tooth pockets (two (2) teeth at a time). (Figure 12) Be careful not to place resin in adjacent pockets. Insert teeth into tooth pockets (Figure 13) and apply pressure to occlusal surface. Hold in place while removing excess resin material from buccal side using 4x4 gauze.

Note: Fuse 2 is available in shade matched equivalents to Lucitone Digital Print 3D Denture Resin. Fuse 2 original shade may be used as a universal fusing option for all denture resin shades, with the exception of dark reddish pink.

9. While continuing to hold teeth in place, use the UV light to tack-cure buccal sides of the two (2) teeth for ten (10) seconds. (Figure 14)

CAUTION: Place finger over the adjacent tooth socket to keep excess resin from curing.

10. Remove excess resin from lingual side using 4x4 gauze. While continuing to hold teeth in place, tack-cure for ten (10) seconds. (Figure 15)

CAUTION: Place finger over the adjacent tooth socket to keep excess resin from curing.

11. Repeat this process to fuse the remaining teeth into the denture base (two (2) teeth at a time).
12. Use a Flocked Applicator Tip to apply a thin layer of the Lucitone Digital Fuse Step 3 - 3D Denture Sealer around the margins and the interproximal areas of the teeth where Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin was applied. (Figure 16)

The denture is now ready for post-cure.

Figure 12

Figure 13

Figure 14

Figure 15

Figure 16

D. Procedure for Fabrication: CURE FUSED DENTURE

For post-curing a Try-in, skip to section E.

CAUTION: Do not cure more than two (2) appliances at a time during post-cure step. This can cause poor tooth bonding and poor fit in an appliance.

1. Carefully place up to two (2) dentures arches, occlusal/tooth side facing up (Figure 17), on the turntable inside the inLab Speedcure Processing Unit.

Note: When curing two (2) appliances at one time, position anterior teeth of each appliance toward the center of the turntable. (Figure 17)

CAUTION: Do not stack appliances or position appliances in a way to allow for direct contact with post processing unit bulbs.


Figure 17

2. Select the menu option “Lucitone Digital” and then “Denture Base”. The program will run for ten (10) minutes followed by a 3-minute cool down. After one side is cured, flip the arch(es) and cure the intaglio/tissue side using the same menu option. The program will run for ten (10) minutes followed by a 3-minute cool down.

CAUTION: Handling the appliance before the 3-minute cool down may compromise fit accuracy.

E. Procedure for Fabrication: TRY-IN POST-CURE

CAUTION: Do not cure more than two (2) appliances at a time during post-cure step.

1. Carefully place up to two (2) try-in arches, occlusal/tooth side facing up, on the turn table inside the Dentsply Sirona inLab Speedcure.

Note: When curing two (2) appliances at one time, position anterior teeth of each appliance toward the center of the turntable. (Refer to Figure 17)

2. Select the menu option “Lucitone Digital” and then “Try-In.” The program will run for ten (10) minutes followed by a 3-minute cool down. After one side is cured, flip the arch(es) and cure the intaglio/tissue side using the same menu option. The program will run for ten (10) minutes followed by a 3-minute cool down.

CAUTION: Handling the appliance before the 3-minute cool down may compromise fit accuracy.

F. Procedure for Fabrication: FINISH AND POLISH

1. Remove bar supports using cut-off disc or similar tool. Minimal finishing of the arch(es) may be necessary.
2. For final dentures, follow your laboratory’s recommended technique for finishing, festooning, stippling and polishing devices.

3. The following disinfectants have been used successfully with Lucitone Digital Print Denture materials: Glutaraldehyde-based, Quaternary ammonium chlorides, Quaternary ammonium chlorides/alcohol based and Phenolic sprays. Use of other disinfectant products may not be suitable. Common cleansers that have been used successfully include Efferdent Denture Cleanser and Polident Denture Cleanser.

CAUTION: Do not steam clean Lucitone Digital Print appliances.

G. Recycle Excess Lucitone Digital Print Resin

After each print, left-over resin in the cassette can be poured back into the original container (Figure 18). The recycled resin can be used for up to five (5) months.

1. Place the original resin container, for the material used, by the pour-back unit.

Note: Different batches, shades or containers should not be mixed when recycling resin.

2. Position the cone funnel and disposable paper filter. The filter will help ensure no printed particles or foreign objectives are poured back into the container.
3. Place the cassette with excess resin on top of the stand and tilt so that the excess resin pours into the funnel.
4. Follow the Carbon Printer instructions to properly clean the cassette.

Note: Resin should not be left in the cassette for an extended period of time. Recycle immediately after the print job is finished.

Figure 18

H. Surface Repair

1. Prepare the repair area by grinding and roughening the surface with a carbide bur. (Figure 19)

Figure 19

2. Apply Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent to the repair area and air dry for 1 minute. (Figure 20)

Figure 20

3. Apply Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin to the repair area with a brush. Allow the material to flow into the repair area to avoid air pockets - do not brush. (Figure 21)

Figure 21

4. Using the UV Tack-Cure Light cure for 10 seconds. (Figure 22)


Figure 22

5. Apply Lucitone Digital Fuse™ Step 3 - 3D Denture Sealer to the repair area. (Figure 23)

Figure 23

6. Next, place the appliance repair side up in the inLab Speedcure Processing Unit. Select "Lucitone Digital" and then "Repair." The program will run for 10 minutes, followed by a 3-minute cool down. (Figure 24)

Note: Only cure one repair denture at a time. Only cure one side of the arch.


Figure 24

7. Finish and polish the denture using conventional techniques.

CAUTION: Do not steam clean Lucitone Digital Print appliances.

I. Tooth Repair

1. Prepare the repair area by grinding and roughening the surface with a carbide bur. (Figure 25)
2. Use Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent to treat/condition the replacement tooth/teeth. Refer to section C for step-by-step instructions.
3. Apply Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent to the repair area and air dry for 1 minute. (Figure 26)
4. Use sticky wax to position the tooth/teeth for repair. (Figure 27)

Figure 25


Figure 26

Figure 27

5. Apply Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin to the repair area with a brush. Allow the material to flow into the repair area to avoid air pockets - do not brush. (Figure 28)

Note: Lucitone Digital Fuse Step 2 is available in shade-matched equivalents to Lucitone Digital Print 3D Denture Resin. Fuse 2 original shade may be used as a universal fusing option for all denture resin shades, with the exception of dark reddish pink.

Figure 28

6. Using the UV Tack-Cure Light, cure each side of the repair for 10 seconds. (Figure 29)

Figure 29

7. Use a Flocked Applicator Tip to apply a thin layer of Lucitone Digital Fuse Step 3 - 3D Denture Sealer to the repair area including around the margins and the interproximal areas. (Figure 30)

Figure 30

8. Place the appliance repair side up in the inLab Speedcure Processing Unit and select "Lucitone Digital" and then "Repair." The program will run for 10 minutes followed by a 3-minute cool down. (Figure 31)

Note: Only cure one repair denture at a time. Only cure one side of the arch.

Figure 31

9. Finish and polish the denture using conventional techniques.

CAUTION: Do not steam clean Lucitone Digital Print appliances.

J. Fracture Repair

Lucitone Digital Fuse™ Step 2 - 3D Denture Bonding Resin (1 of 2 Techniques)

1. Tack the fractured denture together using an adhesive glue or sticky wax. (Figures 32A & B)

Figure 32A

Figure 32B

2. Fabricate a putty model using a silicone putty. (Figure 33)

Figure 33

3. After the model hardens remove the denture from the model. (Figure 34)

Figure 34

4. Prepare the repair area by grinding and roughening the surface with a carbide bur. (Figure 35)

Figure 35

5. Apply Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent to the repair surface and allow it to dry. (Figure 36)

Figure 36

6. Secure the denture on the putty model with sticky wax or similar type material. (Figure 37)

Figure 37

7. After the denture is secured, apply Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin to the repair area and tack cure for 10 seconds using the UV Tack-Cure Light. Repeat this process as needed to achieve the desired repair thickness (Figures 38A, B, C & D)

Figure 38A

Figure 38B

Figure 38C


Figure 38D

8. Apply Lucitone Digital Fuse™ Step 3 - 3D Denture Sealer to the repair area. (Figure 39)


Figure 39

9. Place the appliance repair side up in the inLab Speedcure Processing Unit and select "Lucitone Digital" and then "Repair." The program will run for 10 minutes followed by a 3-minute cool down. (Figure 40)


Figure 40

Note: Only cure one repair denture at a time. Only cure one side of the arch.

10. Finish and polish the denture using conventional techniques. (Figures 41A, B, C & D)


Figure 41A

Figure 41B

Figure 41C

Figure 41D

CAUTION: Do not steam clean Lucitone Digital Print appliances.

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Refer to the Lucitone HIPA - High Impact Pour Acrylic Instructions for Use “Repair Procedure” section for step-by-step fracture instructions.

CAUTION: Do not steam clean Lucitone Digital Print appliances.

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Refer to the Lucitone HIPA - High Impact Pour Acrylic Instructions for Use “Indirect Reline Procedure” section for step-by-step instructions.

CAUTION: Do not steam clean Lucitone Digital Print appliances.

L. Soft Reline

Dentsply Sirona recommends the following soft reline materials for use with Lucitone Digital Print Dentures:

- Tokuyama SOFRELINER TOUGH® Soft and SOFRELINER TOUGH® Medium†
- GC Corporation GC RELINE™ II Soft, RELINE™ II Extra Soft and GC RELINE™ II Extra Extra Soft††

CAUTION: Do not steam clean Lucitone Digital Print appliances.

Questions and Support

Please contact the Dentsply Sirona technical team:

- Email: [dsgdigitaldentures@dentsplysirona.com](mailto:digitaldentures@dentsplysirona.com)
- Phone: 800-243-1942

†The trademarks SOFRELINER TOUGH® Soft and SOFRELINER TOUGH® Medium are owned by Tokuyama Corporation.

††The trademarks GC RELINE™ II Soft, RELINE™ II Extra Soft and GC RELINE™ II Extra Extra Soft are owned by GC Corporation.

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Carbon is a registered trademark of Carbon, Inc.
**DeguDent GmbH**
Carbon er et registrert varemerke som tilhører Carbon, Inc.
†Varemerkene SOFRELINER TOUGH® Soft and SOFRELINER TOUGH® Medium eies av Tokuyama Corporation.
††Varemerkene GC RELINE™ II Soft, RELINE™ II Extra Soft and GC RELINE™ II Extra Extra Soft eies av GC Corporation.