Lucitone Digital Print™ Denture Materials For Carbon M Series Printer

Directions for Use
Mode d'emploi

Lucitone Digital Print™ 3D Denture Resin is a light-cured resin indicated for the fabrication of denture bases in dental laboratories, including full dentures.

Lucitone Digital Try-In™ Resin can be used as a try-in material for evaluation prior to fabrication of the final restoration.

Fabrication of these prostheses require a computer-aided design and manufacturing (CAD/CAM) system using an additive printer.

Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent is indicated for use in enhancing the bond of denture teeth to denture base.

Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin can be utilized as an aid in bonding denture base to denture teeth as well as repair using traditional techniques.

Lucitone Digital Fuse™ Step 3 – 3D Denture Sealer is a light-cured sealant. It is intended to produce a smooth, glossy surface finish on the denture.

B. Requirements for Denture Fabrication

- > File in STL Format
- > Lucitone Digital Print™ 3D Denture Resin
- > Lucitone Digital Try-In™ 3D Trial Placement Resin
- > IPN 3D™ Teeth
- > Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent
- > Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin
- > Lucitone Digital Fuse™ Step 3 – 3D Denture Sealer
- > Ultrasonication Bath or Orbital Shaker
- > Isopropyl Alcohol (≥ 99%)
- > UV Tack Cure Light
- > inLab® Speedcure™
- > Electric Hot Plate
- > Metal Container Kit
- > Flock Tip Applicators
- > Dentsply Sirona Brush
- > 4x4 Gauze Pads
- > Pipette
- > Small Plastic Container (see Figure 5)
- > Carbon** 3D Printer M Series Software
- > Carbon** 3D Printer M Series

CAUTION: Any unauthorized changes to the process, workflow, materials, equipment, parameters, or software may result in a device that is out of specification and is not recommended. Contact Dentsply Sirona for a list of compatible components.

*Carbon® is a registered mark owned by Carbon, Inc.

C. Contraindications

Lucitone Digital Print™ 3D System components are contraindicated for patients with known hypersensitivity or severe allergic reactions to any of the components.

D. Warnings

1. Lucitone Digital Print 3D System Materials contain components which may cause skin dryness or irritation, sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wear protective gloves, eye protection and protective clothing when handling.
Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.
2. Lucitone Digital Fuse Step 1 - Liquid is highly flammable. Keep containers tightly closed and use/store away from sources of ignition. Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization or other allergic reactions.
3. Lucitone Digital Fuse Step 3 - 3D Denture Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. Suspected of damaging fertility or the unborn child. Harmful to aquatic life with long lasting effects. Avoid release into the environment. Dispose contents.
4. Lucitone Digital Print 3D System Materials- Avoid inhalation or ingestion. High vapor concentration might induce headache, nausea and irritation of eyes and respiratory system. Work in a well-ventilated area. Avoid contact with eyes. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects mouth must be considered by the physician/dentist.
Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
Ingestion: Contact your regional Poison Control Center immediately.
5. Proper ventilation and personal protective equipment should be used when grinding devices as the particulate generated during grinding may cause respiratory, skin and/or eye irritation.
6. Isopropyl Alcohol is highly flammable. Keep containers tightly closed and use/store away from sources of ignition.

E. Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use (DFU). Any use of this product inconsistent with the DFU is at the discretion and sole responsibility of the practitioner.
2. Store denture resin and system components at room temperature 60°F - 80°F (16°C - 27°C) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date. Dispose contents and containers in accordance with local and national regulations.

3. Lucitone Digital Print™ 3D Denture Resin & Try-In Resin and Lucitone Digital Fuse™ Step 2 - 3D Denture Bonding Resin
When washing 3D printed appliance with a solvent or polishing the appliance, it should be done in a properly ventilated environment with proper protective masks and gloves. Expired or unused resins shall be disposed in accordance with local regulations.
4. Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent
Use in a ventilated workspace. Handle teeth with solvent-resistant gloves or tweezers. Product is light sensitive. Minimize exposure to light and keep the lid on the container to maintain its effective use. Product is very volatile and flammable. Keep away from heat, sparks and flame. Do not add the Lucitone Digital Fuse™ Step 1 Agent to a warm metal container or while the container is on the electric hot plate.
5. Lucitone Digital Fuse™ Step 3 – 3D Denture Sealer
Use in a ventilated workspace. Keep bottle tightly closed when not in use. Lucitone Digital Fuse™ Step 3 is volatile and highly flammable. Keep away from heat, sparks and flames. Product is light sensitive. Minimize exposure to light and keep the lid on the container to maintain its effective use.
6. The disinfectants that have been used successfully with Lucitone Digital Print systems are : Glutaraldehyde-based, Quaternary ammonium chlorides, Quaternary ammonium chlorides/ alcohol based and Phenolic sprays. Use of other disinfectant products may not be suitable. Common cleansers that have been used successfully include Efferdent Denture Cleanser and Polident Denture Cleanser.

F. Adverse Reactions

1. Lucitone Digital Print 3D System Materials contain components which may cause skin dryness or irritation, sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wear protective gloves, eye protection and protective clothing when handling.
Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.
2. Lucitone Digital Fuse Step 1 - Liquid is highly flammable. Keep containers tightly closed and use/store away from sources of ignition. Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization or other allergic reactions.
3. Lucitone Digital Fuse Step 3 - 3D Denture Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. Suspected of damaging fertility or the unborn child. Harmful to aquatic life with long lasting effects. Avoid release into the environment. Dispose contents.
4. Lucitone Digital Print 3D System Materials- Avoid inhalation or ingestion. High vapor concentration might induce headache, nausea and irritation of eyes and respiratory system. Work in a well-ventilated area. Avoid contact with eyes. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects mouth must be considered by the physician/dentist.
Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
Ingestion: Contact your regional Poison Control Center immediately.

G. Procedure for Fabrication: PRINT

1. .stl file(s) should be created per your scan and design manufacturer's instructions. Create Carbon M2 new project and upload .stl file(s) for your print job.

Select **"Lucitone Digital Print"** from the resin drop down menu when printing denture bases with pockets.

Select **"Lucitone Digital Try-In"** when printing denture try-ins.

Contact Dentsply Sirona Lab Technical Team for print orientation, support design and nesting guidelines at 800-243-1942 ext. 54212 or dsdigitaldentures@dentsplysirona.com

2. Based on software recommendation, pour resin volume slowly into the cassette.
3. Once cassette and platform are ready and door is closed, start 'Print'. The screen will change to display the print progress and estimate print time.

H. Procedure for Fabrication: WASH (TRY-INS AND BASES)

Note: Section H should be done under proper ventilation.

1. Once print is complete, remove build platform from printer.
2. Remove appliances from build platform by snapping off each using your hand. Moderate pressure is required. (Figure 1).
3. Break off fence supports from the appliance (Figure 2). Leave the bar supports in place. (Figure 3 & 4). They will be removed after post-cure.

Figure 1

Figure 2

Figure 3

Figure 4

4. Place appliance pocket-down in a plastic container and wash with Isopropyl Alcohol ($\geq 99\%$ IPA) using a two-bath approach with fresh $\geq 99\%$ IPA each time. Using an ultrasonication bath, wash for 2 minutes (1st cycle), followed by 1 minute (2nd cycle). Container should be closed when placed in the ultrasonication bath. Wash each appliance individually. (Figure 5 & 6)

CAUTION: Prolonged exposure to $\geq 99\%$ IPA could cause changes to physical properties of the appliance. Follow wash times carefully.


Figure 5

When washing bases with pockets, remove from first bath cycle and use a wet brush (wet with \geq 99% IPA) to ensure tooth pockets are clean of residual resin. (Figure 7)

CAUTION: Residual resin in tooth pockets could affect fit and bond.

When washing a try-in, remove from first bath cycle and use a wet brush (wet with \geq 99% IPA) to wash occlusal surface and interproximal of teeth. (Figure 8)

CAUTION: Residual resin on occlusal surfaces could affect occlusion.

Note: If using an Orbital Shaker, perform washing for 3 minutes (1st cycle) followed by 2 minutes (2nd cycle).

5. Blow dry the appliance using compressed air to ensure all parts are thoroughly dried before moving to tooth fusing/bonding.

Figure 6


Figure 7

Figure 8

I. Procedure for Fabrication: FUSE (BOND)

For post-curing a Try-In, skip to section K.

Note: Section H should be done under proper ventilation.

1. Plug in the electric hot plate and set the temperature to 90°C.

Note: If the electric hot plate is not digital, set to 100.

Allow the electric hot plate to preheat for at least ten (10) minutes.

Note: The temperature setting may need to be adjusted slightly from the recommended setting due to several factors.

2. Place IPN 3D teeth into respective pockets to check for a dry fit. When necessary use dental burs or tools to remove 'sprues' from teeth.

3. Shake the bottle of Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent to ensure contents are thoroughly mixed. Pour enough Fuse Step 1 into the round metal container to fill to top of foam insert. (Figure 9) Refer to Precautions.

CAUTION: Using inadequate amount of Fuse Step 1 will compromise tooth bond strength.

CAUTION: Keep the metal container and the foam insert free of any solid buildup. This can compromise the tooth bond strength.


Figure 9

4. Use tweezers to place up to 14 IPN 3D denture teeth into the holes so that the teeth are positioned with the ridge lap down and in contact with the liquid. If liquid is not level with the top of the foam insert, add additional Fuse Step 1 to metal container. Immediately replace the lid on the container. (Figure 10)

Figure 10

5. Place the tooth-filled container inside the metal ring of the electric hot plate until the 40°C rectangle on the temperature strip turns green (normally in four (4) to six (6) minutes). Refer to Precautions. (Figure 11)

CAUTION: Tooth bond strength will be compromised if 40°C is not reached.

6. Carefully remove the metal container from the electric hot plate. Remove the lid from the metal container and use tweezers to retrieve the teeth. Allow teeth to air dry for minimum of 2 minutes. Do not wipe teeth.

Figure 11

7. Place the lid back on the metal container immediately. Do not return the metal container to the electric hot plate until ready to treat the next set of teeth.
8. Apply a small amount of Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin in the tooth pockets (two (2) teeth at a time). (Figure 12) Be careful not to place resin in adjacent pockets. Insert teeth into tooth pockets (Figure 13) and apply pressure to occlusal surface. Hold in place while removing excess resin material from buccal side using 4x4 gauze.

Figure 12

Figure 13

9. While continuing to hold teeth in place, use the UV light to tack-cure buccal sides of the two (2) teeth for ten (10) seconds. (Figure 14)

CAUTION: Place finger over the adjacent tooth socket to keep excess resin from curing.

Figure 14

10. Remove excess resin from lingual side using 4x4 gauze. While continuing to hold teeth in place, tack-cure for 10 seconds. (Figure 15)

CAUTION: Place finger over the adjacent tooth socket to keep excess resin from curing.

Figure 15

11. Repeat this process to fuse the remaining teeth into the denture base (two (2) teeth at a time).
12. Use a flock tip applicator to apply thin layer of the Lucitone Digital Fuse Step 3 - 3D Denture Sealer around the margins and the interproximal areas of the teeth where Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin was applied (Figure 16).

Denture is now ready for post-cure.

Figure 16

J. Procedure for Fabrication: CURE FUSED DENTURE

For post-curing a Try-in, skip to section K.

CAUTION: Do not cure more than two (2) appliances at a time during post-cure step. This can cause poor tooth bonding and poor fit in an appliance.

1. Carefully place up to two (2) dentures arches occlusal/tooth side facing up (Figure 17) on the turntable inside the Dentsply Sirona inLab Speedcure.

Note: When curing two (2) appliances at one time, position anterior teeth of each appliance toward the center of the turntable. (Refer to Figure 17)

Figure 17

2. Select the menu option “Lucitone Digital” and then “Denture Base”. The program will run for 10 minutes followed by a 3 minute cool down. After one side is cured, flip the arches and cure the intaglio/tissue side using the same menu option. The program will run for 10 minutes followed by a 3 minute cool down.

CAUTION: Handling the appliance before the 3 minute cool down may compromise fit accuracy.

K. Procedure for Fabrication: TRY-IN POST-CURE

CAUTION: Do not cure more than two (2) appliances at a time during post-cure step.

1. Carefully place up to two (2) try-in arches occlusal/tooth side facing up on the turn table inside the Dentsply Sirona inLab Speedcure.

Note: When curing two (2) appliances at one time, position anterior teeth of each appliance toward the center of the turntable. (Refer to Figure 17)

2. Select the menu option “Lucitone Digital” and then “Try-In”. The program will run for 10 minutes followed by a 3 minute cool down. After one side is cured, flip the arch(es) and cure the intaglio/tissue side using the same menu option. The program will run for 10 minutes followed by a 3 minute cool down.

CAUTION: Handling the appliance before the 3 minute cool down may compromise fit accuracy.

L. Procedure for Fabrication: FINISH AND POLISH

1. Remove bar supports using cut-off disc or similar tool. Minimal finishing of the arches may be necessary.
2. For final dentures, follow your laboratory’s recommended technique for finishing, festooning, stippling and polishing devices.
3. The disinfectants that have been used successfully with Lucitone Digital Print systems are : Glutaraldehyde-based, Quaternary ammonium chlorides, Quaternary ammonium chlorides/ alcohol based and Phenolic sprays. Use of other disinfectant products may not be suitable. Common cleansers that have been used successfully include Efferdent Denture Cleanser and Polident Denture Cleanser.

M. Procedure for Fabrication: REPAIR

Note: Repairs will be performed using traditional method, not 3D printing. Please refer to our website, dentsplysirona.com, for detailed description on our repair procedure.

N. Procedure for Recycling Excess Lucitone Digital Resin

After each print, left over resin in the cassette can be poured back into the original container (Figure 18). The recycled resin can be used for up to 5 months.

1. Place the original resin container for the material used by the pour back unit.
Note: Different batches, shades or containers should not be used.
2. Position the cone funnel and disposable paper filter. These will ensure no printed particles or foreign objectives are poured back into the container.
3. Place the cassette with excess resin on top of the stand and tilt so that the excess resin pours into the funnel.

Figure 18

Questions and Support

Please contact the Dentsply Sirona technical team:

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York, PA 17401 U.S.A.