

Celtra® Press

FactFile

Properties of zirconia-reinforced lithium silicate ceramics (ZLS)

What is Celtra Press?

Celtra Press is a press ceramic, available as a pellet in two sizes - 3 g and 6 g. This ceramic material was specially developed for the production of indirect ceramic restorations at the dental laboratory. Celtra Press is characterized by its excellent pressing behavior with all restoration types.

The special power firing step before veneering or glazing will "healing" microscopic faults within the glass matrix

Characteristics

Celtra Press is a multiphase ceramic consisting of a glass matrix and lithium disilicate crystals having a crystal length of about 1.5 µm plus nano-scale lithium phosphate (see the SEM images). In addition to Li₂O and SiO₂, Celtra Press contains about 10% zirconia (ZrO₂), which is dissolved completely in the glass phase as in Celtra Duo rather than in crystalline form.

Celtra Press is characterized by a high strength of about 500 MPa (after power firing) and excellent flow properties during pressing. In conjunction with the newly developed investment, only a minimal reaction layer is formed,

caused by sandblasting or finishing. This ensures the consistently high strength of more than 500 MPa of the material.

The simplified sprueing of bridges, the low viscosity and good flowability of the ceramic material during pressing, and faster divesting/sandblasting steps will shorten the time needed to produce the restorations.

resulting in an excellent fit and shortening the time needed for finishing, because the surface acidification step is eliminated.

Celtra Press is easy to polish, should minor corrections in the patient's mouth turn out to be necessary. Though no post-processing in the laboratory is required and excellent esthetics is nevertheless ensured.

Technical specifications for Celtra® Press

CTE (25-500°C):	9,7 *10 ⁻⁶ /K ⁻¹
Softening Point Sp:	820°C
Glass-transition temperature Tg:	560°C
Flexural strength after power firing:	>500MPa
Chemical solubility:	<30 µg/cm ²

Celtra® Press FactFile

Cementing using the Celtra® cementation system

Depending on the indication, Celtra Press restorations may be cemented self-adhesively, fully adhesively, or (in the case of crowns) with Ketac Cem glass-ionomer cement.

Compatible time-proven adhesive cementing materials are available as part of the Dentsply Sirona range of products. Cements are available separately.

	Self-adhesive	Fully adhesive	Glass-ionomer ¹⁾
Inlays	R	HR	-
Onlays	R	HR	-
Crowns	HR	HR	R
Veneers	-	HR	-
Bridges	HR	HR	R

R = recommended
HR = highly recommended
¹⁾ Ketac*

* Ketac is a trademark of 3M ESPE

Esthetic properties

Celtra Press is available in selected VITA shades that can be used to produce almost any restoration. The individual pellets are assigned either to the HT (high translucency), the MT (medium translucency), or the LT (low translucency) groups, intended for different indications:

- 1 HT ingots** are high translucency ingots and are available in values of I1, I2 and I3. HT ingot can be used for restorations in the incisal region like inlays, onlays and veneers.
- 2 MT ingots** are ingots with a medium translucency and are available in A-D shades. MT ingots can be used for crowns and bridges designed as a full contour restoration. MT ingots can also be used for restoration with a cut-back of the incisal region, pressing and veneering with Celtra Ceram to complete the esthetic restoration.
- 3 LT ingots** are shaded to VITA® A-D shades in dentin. These are used for veneers or crowns and bridges with a cut-back of the incisal region, pressing and veneering with Celtra Ceram to complete the esthetic restoration.

The shades of Celtra Press are adapted to the VITA dentin shades. Thanks to the opalescent properties of Celtra Press ceramics, the incisal-edge or cusp regions will appear very natural even in the case of restorations that have been pressed to full anatomic contour.

Celtra Press is fluorescent; the intensity of the fluorescence may be enhanced by using the system's fluorescent glaze.

Celtra stains can be used to individually characterize any restoration.

Celtra® Press FactFile

In-vitro studies

Chewing simulation studies and aging tests conducted by the universities of Heidelberg and Regensburg as well as company-internal studies have shown that the performance of Celtra Press is at least as good as that of existing lithium disilicate (LS₂) pressable ceramics, such as e.max Press.

In-vitro studies on the wear behavior of Celtra Press and enamel-like antagonist have shown that the abrasion of Celtra Press and existing LS₂ pressable ceramics has the same order of magnitude as that of natural tooth enamel.

Shade philosophy

SMART SHADE SYSTEM

Celtra Press ingot portfolio

Type of restoration	Translucency	Shade			Individualization technique				
Incisal (inlay, onlay, veneer)	HT		I1	I2	I3			Glaze	
Full Contour (posterior)	MT	BL2*	A1	A2	A3	B1	C1	D2	Stain & Glaze
Cut-back (anterior)	LT	BL2*	A1	A2	A3	B1	C1	D2	Build-up Stain & Glaze

* coming soon


The three different translucency levels of Celtra Press (HT, MT, and LT) have been defined such that a true translucency gradient is obtained throughout for all shades within a translucency class. The chart (next page) shows relative translucency as a function of translucency classes and the A series of shades. While in conventional LS₂ glass ceramics, translucencies are not uniform (for example, MT can be more opaque than LT and, for example, A3 HT less translucent than A3 LT), Celtra Press possesses a consistent translucency gradient that is independent of the selected shade within the VITA A group.

The second chart (next page) shows the classification of all levels of translucency and associated shades compared to conventional LS₂ glass ceramics and an enamel of the Duceram Kiss range (S03) and one dentin of the Duceram Kiss range (DA3). Within a shade (e.g., VITA A), translucency decreases (e.g., A1 MT to A3) but always keeps the same distance in terms of translucency from to the next, less translucent shade group (A1 LT to A3 LT). This makes it easier for the dental technician to benefit from experience gained within one shade group to one of the other shade groups (e.g., from the A group to the B, C, or D group). Decreasing translucency within a shade group (e.g., VITA A) is caused by the higher chroma, as in natural teeth in an age-related series.

LIGHT TRANSMISSION

L [%]


Relative translucency



This conventional Lithium Disilicate ceramic material (HT & LT) is not consistent in the nomenclature regarding translucency.

Also a too high translucency without enough opacity results in a greying effect in the mouth

Relative translucency



Microstructures

The first pair of images shows the pellets of Celtra Press and a conventional LS₂ glass-ceramic material.

The scanning electron microscope (SEM) images show the microstructural differences of zirconia-reinforced glass ceramics compared to conventional glass ceramics. This pair of images shows the crystal size of the lithium disilicate within the pellet, as delivered. It can be clearly seen that the crystals in the Celtra Press pellet are smaller, which results in better compressibility and flowability (lower viscosity)

during the pressing process. This means that it is easier to press thin-walled sections of the restoration without needing a lot of sprues.

The final pair of images shows the crystallite size of lithium disilicate in the final restoration. Here, too, smaller crystalline structure of Celtra Press is obvious, which facilitates simpler and faster polishing.



Celtra® Press FactFile

Celtra® Press



Conventional LS₂ glass ceramic



Ingot



Ingot microstructure



Final restoration microstructure

Celtra® Press FactFile

Celtra® Press Investment

The Celtra Press investment was specially developed for Celtra Press minimizes the formation of a hard reaction layer between the investment and the pressed restoration. Along with then lower pressing temperature of Celtra Press compared to conventional LS₂ glass ceramics (50-60°C), the Celtra Press Investment and Celtra Press pair eliminates the need to used cleaning fluid containing hydrofluoric acid. Any reaction layer residue that still exists can easily be removed by sandblasting, not least because the hardness of the reaction layer is greatly reduced by the by the lower

pressing temperature and thanks to the boron nitride contained in the investment.

The left sample illustrates the divesting of Celtra Press from the Celtra Press Investment by mere sandblasting; no reaction layer is discernible. The right sample shows the corresponding situation when conventional lithium disilicate pressable ceramics are used; here, too, only sandblasting was used. Residues of the reaction layer are clearly visible.



The left sample illustrates the divesting of Celtra Press from the Celtra Press Investment by mere sandblasting; no reaction layer is discernible. The right sample shows the corresponding situation when conventional lithium disilicate pressable ceramics are used; here, too, only sandblasting was used. Residues of the reaction layer are clearly visible.

Celtra® Press FactFile

Dentsply Sirona Celtra Ceram

Description

Celtra Ceram a ceramic veneering material for all-ceramic frameworks, a low-fusing, leucite-reinforced feldspathic ceramic material optimized for veneering and characterizing Celtra Press zirconia-reinforced lithium silicate (ZLS) frameworks produced using heat-pressing technology. With a CTE 25–500 °C of $9.0 \times 10^{-6} \text{ K}^{-1}$ and a firing temperature of 770°C (1st dentine firing), Celtra Ceram is also suitable for lithium disilicate frameworks.

Benefits

Primary

- Comprehensive veneering porcelain system optimized for bonding to Celtra Press and lithium disilicate substructures, eliminating the need for wash firing
- Exceptional performance derived from proprietary feldspathic base formulation that delivers an incredibly stable and robust result – every time, even after multiple firings
- Shade confidence with out-of-the-bottle shade accuracy

Secondary

- Color-coding (organic pigments) of different types of porcelain provides dental technicians the convenience of identifying/distinguishing each layer of porcelain during build-up
- Optimized particle size distribution of dentin/enamel porcelain facilitates the wetting with modeling liquid for smooth application and seamless build-up
- Correction/Add-on can be applied and fired simultaneously with glaze/stain porcelain, or after glaze/stain firing.

Indications

- Celtra Press zirconia-reinforced lithium silicate (ZLS) frameworks: CTE (25–500°C): $9,7 \times 10^{-6} \text{ K}^{-1}$
- Lithium disilicate frameworks: CTE (100–500°C): $10,5 \times 10^{-6} \text{ K}^{-1}$

Physical properties

CTE (25–500°C): $10,5 \times 10^{-6} \text{ K}^{-1}$. Firing temperature: 770°C (1st dentin firing)

Flexural strength (three-point bending test): 108 N/mm². Vickers hardness: 5000 N/mm². Chemical solubility: 28 µg/cm².

celtra-dentsplysirona.com