Sirona Connect

TRAINING: DIGITALE ABFORMUNG

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ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SUPPORT CONTACT

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FERTIG

SIRONA Connect

MY SIRONA CONNECT > DENTIST REGISTRATION

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MY SIRONA CONNECT > DENTIST REGISTRATION

Welcome to SIRONA
Connect

Dentist Registration

Step 1

Step 2

Step 3

Done

My user information

Salutation *	<input checked="" type="radio"/> Mr. <input type="radio"/> Mrs.	Phone Number *	<input type="text" value="0177 69 60 311"/>
Title	<input type="text" value="Dr."/>	E-Mail Address *	<input type="text" value="rob@F-Dental.com"/>
First Name *	<input type="text" value="Mark"/>	Confirm e-Mail address *	<input type="text" value="rob@F-Dental.com"/>
Last Name *	<input type="text" value="Homburger"/>	Skype username	<input type="text" value="DentFFM"/>

My practice data

Company Name *	<input type="text" value="Praxis Dr Homburger"/>	E-Mail Address	<input type="text" value="DrH@PraxisHomburger.de"/>
Phone Number *	<input type="text" value="061 72 69 60 311"/>	Website	<input type="text" value="www.Dental.de"/>
Fax Number	<input type="text"/>		

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MY SIRONA CONNECT > DENTIST REGISTRATION

Welcome to SIRONA
Connect

Dentist Registration

Step 1 ▾ Step 2 ▾ **Step 3** ▾ Done ▾

My billing address

Address 1 *

Address 2

ZIP Code *

City *

Country *

Region

different shipping address

My shipping address

Recipient

Address 1 *

Address 2

ZIP Code *

City *

Country *

Region

*marked fields are mandatory fields

BACK

DONE

SCHRIIT 1

SCHRIIT 2

SCHRIIT 3

FERTIG

MY SIRONA CONNECT > DENTIST REGISTRATION

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We have therefore sent you an e-mail with the subject "**Please confirm your e-mail address**" which you can use to activate your account. Therefore please follow the instructions inside the e-mail.

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SIRONA Connect – Please activate your account

Von: cerec.connect@sirona.de

An: alex@drauz.net

Datum: 25.03.2014 19:16

Dear Mr. Homburger

Welcome at SIRONA Connect! You have been registered on 03/25/2014 6:16:57 PM.

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Sincerely,

Your SIRONA Connect Team

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The screenshot shows the SIRONA Connect website interface. At the top left, it says "SIRONA CORPORATE WEBSITE". In the top center, a yellow box contains the message "Your account has been activated successfully.". At the top right, there are links for "LABORATORY REGISTRATION", "DENTIST REGISTRATION", and "LOGIN". Below the navigation bar, the "sirona" logo is displayed with the tagline "The Dental Company". The main content area is divided into two columns. The left column, titled "SIRONA Connect", lists "Laboratory Registration", "Dentist Registration", and "Login". The right column, titled "MY SIRONA CONNECT", features a "Welcome to SIRONA Connect" message and two images of a man talking on a phone while working at a computer. A large yellow arrow points to a yellow box labeled "Activation" with the text "Please log in and provide the additional data.". At the bottom right, there are links for "SITEMAP", "TERMS & CONDITIONS", and "PRIVACY STATEMENT".

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ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SUPPORT CONTACT

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ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SUPPORT CONTACT

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New perspectives for dental practices and laboratories.

Login

Praxis Dr Homburger

•••••

[> Dentist Registration](#)
[> Laboratory Registration](#)
[> Forgot password?](#)

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SIRONA Connect

- My Account
- User Administration
- Order List
- Log Off

MY SIRONA CONNECT > WELCOME

Welcome to SIRONA Connect

Welcome

Welcome to Sirona Connect!

Your last visit was on 03/26/2014 3:31 PM

> You can edit your personal account on the following pages.

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Telefon-Nr.:
 E-Mail-Adresse:
 E-Mail-Adresse bestätigen:
 Skype-Dienstname:

My Workticket signature preferences

Use Workticket signature

Signature File: Workticket Signature Text:

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REGISTRIERUNG ALS LABOR

The screenshot displays the Sirona Corporate Website with the following elements:

- Top Navigation:** SIRONA CORPORATE WEBSITE | LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN
- Secondary Navigation:** ABOUT SIRONA CONNECT | FOR DENTAL PRACTICES | FOR DENTAL LABORATORIES | SUPPORT | CONTACT
- Logo:** sirona. The Dental Company
- Main Content:**
 - Text:** Digital impression with Sirona Connect. New perspectives for dental prac-
 - Form:** Login section with fields for Username and Password, and a Login button.
 - Links:** > Dentist Registration, > Laboratory Registration (highlighted with a yellow arrow), > Forgot password?
 - Image:** A man in a white polo shirt talking on a mobile phone, with a computer monitor displaying a 3D dental model in the background.
 - Call to Action:** FREE CEREC CONNECT SW 4.X DOWNLOAD
 - Image:** A man in a light blue shirt talking on a mobile phone in a dental laboratory setting, with a computer monitor displaying a 3D dental model.

SCHRITT 1

SCHRITT 2

SCHRITT 3

FERTIG

SIRONA Connect

MY SIRONA CONNECT > LABORATORY REGISTRATION

- Laboratory Registration
- Dentist Registration
- Login

Welcome to SIRONA
Connect

Laboratory Registration

Step 1

Step 2

Step 3

Done

Your access data

At first please choose your personal and unique user name, which you can easily remember. If you have questions, give us a call: +49 (0)6251 16-1670.

Your favorite username *

Your chosen password *

* marked fields are mandatory fields

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Step 1 ▾ Step 2 Step 3 Done

My user information

Salutation *	<input checked="" type="radio"/> Mr. <input type="radio"/> Mrs.	Phone Number *	<input type="text" value="069 60311 69"/>
Title	<input type="text" value="Dr"/>	E-Mail Address *	<input type="text" value="rob@F-Dental.com"/>
First Name *	<input type="text" value="Robin"/>	Confirm e-Mail address *	<input type="text" value="rob@F-Dental.com"/>
Last Name *	<input type="text" value="Frankfurter"/>	Skype username	<input type="text" value="RobFFM"/>

My laboratory data

Company Name *	<input type="text" value="Dentallab Frankfurter"/>	E-Mail Address	<input type="text" value="RF@Frankfurter-Dental.com"/>
Phone Number *	<input type="text" value="069 60311 71"/>	Website	<input type="text" value="www.Dental-FFM.de"/>
Fax Number	<input type="text" value="069 60311 77"/>		

My preferences

Preferred Language	<input type="text" value="english (GB)"/>	Timezone	<input type="text" value="(UTC) Casablanca"/>
--------------------	---	----------	---

* marked fields are mandatory fields

BACK

NEXT

SCHRITT 1

SCHRITT 2

SCHRITT 3

FERTIG

SIRONA Connect

- Laboratory Registration
- Dentist Registration
- Login

MY SIRONA CONNECT > LABORATORY REGISTRATION

Welcome to SIRONA Connect

Laboratory Registration

Step 1 ▾ Step 2 ▾ Step 3 ▾ Done ▾

My billing address

Address 1 *

Address 2


ZIP Code *

City *

Country *

Region

different shipping address

My shipping address

Recipient

Address 1 *


Address 2

ZIP Code *

City *

Country *

Region

* marked fields are mandatory fields

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FERTIG

SIRONA Connect

- Laboratory Registration
- Dentist Registration
- Login

MY SIRONA CONNECT > LABORATORY REGISTRATION

Welcome to SIRONA
Connect

Laboratory Registration

Step 1 ▾

Step 2 ▾

Step 3 ▾

Done

Congratulations! You nearly did it.

Your new Sirona account was created. Prior to your first login, we would like to check whether your e-mail address is valid und working.

We have therefore sent you an e-mail with the subject "Please confirm your e-mail address" to rob@F-Dental.com. You will find a link in the e-mail which you can use to activate your account. Therefore please follow the instructions inside the e-mail.

If you didn't receive the e-mail, please check your spam folder. If you have problems, call our support (phone number + 49 (0) 6251 / 16-1670) or send us an e-mail (support@sirona.com).

We wish you good luck with Sirona Connect!

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☒ **SIRONA Connect – Please activate your account**

Von: cerec.connect@sirona.de

An: alex@drauz.net

Datum: 18:57

Dear Mr. Homburger

Welcome at SIRONA Connect! You have been registered on 25/03/2014 6:57:15 PM.

Your details are:

Username: Dentallabor Frankfurter

Dr Mark Homburger

Phone: 061 72 69 60 311

Email: DrH@PraxisHomburger.de

|

Laboratory data:

Dentallabor Frankfurter

Goethestrasse 3

60711 Frankfurt

Germany

Fax:

Website:

You can activate your account with the following [Activation Link](#)

Sincerely,

Your SIRONA Connect Team

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The screenshot displays the SIRONA Connect website interface. At the top left, it says "SIRONA CORPORATE WEBSITE". A yellow notification box in the top center contains the text "Your account has been activated successfully.". On the top right, there are links for "LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN". Below this is a navigation menu with "ABOUT SIRONA CONNECT", "FOR DENTAL PRACTICES", "FOR DENTAL LABORATORIES", "SERVICE", and "CONTACT". The SIRONA logo "sirona. The Dental Company" is on the right. The main content area features a "SIRONA Connect" section with a list: "Laboratory Registration", "Dentist Registration", and "Login". A breadcrumb trail reads "MY SIRONA CONNECT > LABORATORY REGISTRATION". A large banner area contains the text "Welcome to SIRONA Connect" and two images of a man talking on a phone at a computer workstation. A large yellow arrow points to a yellow box labeled "Activation" with the text "Please log in and provide the additional data.". At the bottom right, there are links for "SITEMAP | TERMS & CONDITIONS | PRIVACY STATEMENT".

SIRONA CORPORATE WEBSITE

Your account has been activated successfully.

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SERVICE CONTACT

sirona.
The Dental Company

SIRONA Connect

- Laboratory Registration
- Dentist Registration
- Login

MY SIRONA CONNECT > LABORATORY REGISTRATION

Welcome to SIRONA Connect

Activation
Please log in and provide the additional data.

SITEMAP | TERMS & CONDITIONS | PRIVACY STATEMENT

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SIRONA CORPORATE WEBSITE

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES SERVICE CONTACT

sirona
The Dental Company

Digital impression with Sirona Connect

New perspectives for dental practices and laboratories

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- My Account
 - My user
 - My laboratory
 - My billing address
 - My shipping address
 - Change Password
 - Offerings / Services
 - Connect Direct
- User Administration
 - Order List
 - Log Off

Welcome to SIRONA Connect

My laboratory data

Company Name* E-Mail Address

Phone Number* Website

Fax Number

My preferences

Preferred Language List range

Timezone

infiniDent

infiniDent User Name infiniDent Password

* marked fields are mandatory fields

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My Sirona Connect Home > My Sirona Connect > Login

- Laboratory Registration
- Dentist Registration
- Login**

Login

Please enter your **User ID** and **Password**.

User ID *

Password *

[Forgot your User ID or Password?](#)

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My Sirona Connect Home > My Sirona Connect > Connect Direct

Logged in as **lab**

- ▶ MyAccount
 - ▶ Password
 - ▶ User
 - ▶ Laboratory
 - ▶ Billing Address
 - ▶ Shipping Address
 - ▶ Offerings / Services
 - ▶ Connect Direct
- ▶ Order List
- ▶ Log Off

Connect Direct

Please use the form below to add dentists to your Connect Direct customers.

Add Connect Direct workstation

Inventory Number *	<input type="text" value="000111222"/>
Workstation ID *	<input type="text" value="012345"/>
Dentist's contact person username *	<input type="text" value="dentist"/>

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My Sirona Connect Home > My Sirona Connect > Connect Direct

Logged in as
lab

- ▶ MyAccount
 - ▶ Password
 - ▶ User
 - ▶ Laboratory
 - ▶ Billing Address
 - ▶ Shipping Address
 - ▶ Offerings / Services
 - ▶ **Connect Direct**
- ▶ Order List
- ▶ Log Off

 Connect Direct successfully added.

Connect Direct

Please use the form below to add dentists to your Connect Direct customers.

Existing Connect Direct workstations

Inventory Number	Company Name	City	Activation	
000111222	Company Dentist	zip City	pending since 2/23/2012	Send activation e-Mail 

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

My Sirona Connect Home > My Sirona Connect > Connect Direct

- ▶ Laboratory Registration
- ▶ Dentist Registration
- ▶ Login

 Connect Direct workstation accepted successfully!

Connect Direct

Hello Name1 Dentist Name2 Dentist,

You have accepted Connect Direct for your AC with the serial number 012345.
Please open the portal to the AC and enter the following login data:

Username: 806E24764B43BD5CDC725CA5167056715A040788
Password: 1234


With this step Connect Direct is activated on your AC. For questions, please contact your laboratory.

Company Lab
adress
city, FL zip
United States
Phone: Phone Lab

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



My Sirona Connect Home > My Sirona Connect > Connect Direct

Logged in as **lab**

- ▶ MyAccount
 - ▶ Password
 - ▶ User
 - ▶ Laboratory
 - ▶ Billing Address
 - ▶ Shipping Address
 - ▶ Offerings / Services
 - ▶ **Connect Direct**
- ▶ Order List
- ▶ Log Off

Connect Direct

Please use the form below to add dentists to your Connect Direct customers.

Existing Connect Direct workstations

Inventory Number	Company Name	City	Activation	
000111222	Company Dentist	zip City	activated on 2/23/2012	

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

The screenshot shows a web page for Sirona software updates. The page has a navigation menu on the left with 'Support' and sub-items: 'Software', 'FAQ', 'Sirona Connect Video', and 'Video tutorials'. The main content area is titled 'Software' and includes a breadcrumb trail 'WELCOME > SUPPORT > SOFTWARE'. Below the title, there are two images of a man talking on a phone while working at a computer. The page is divided into two sections: 'APOLLO DI' and 'CEREC Connect Software'. Each section provides instructions on how to find or download the current software and its manual, with download links for various file types and sizes.

Support

- Software
- FAQ
- Sirona Connect Video
- Video tutorials

WELCOME > SUPPORT > SOFTWARE

Software

APOLLO DI

Here you find the current APOLLO Connect Software and its operator's manual.

- APOLLO Connect Software
exe file | 1,4 GB
- Update description SW 4.5.5
pdf file | 88 KB
- Operator's manual APOLLO Connect Software
pdf file | 1,3 MB

CEREC Connect Software

Here you can download the current CEREC Connect Software.

Please note the installation instructions from the following link www.sirona.com/firmwareupdate if your previous version is smaller than CEREC Connect SW 4.2.1.

- CEREC Connect SW 4.2.4:
exe file | 393 MB
- Update description CEREC Connect SW 4.2.4:
pdf file | 88 KB

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SIRONA CORPORATE WEBSITE

LABORATORY REGISTRATION | DENTIST REGISTRATION | LOGIN

ABOUT SIRONA CONNECT FOR DENTAL PRACTICES FOR DENTAL LABORATORIES **SUPPORT** CONTACT

sirona
The Dental Company

Support

- Software
- FAQ
- Sirona Connect Video**
- Video tutorials

WELCOME > SUPPORT > **SIRONA CONNECT VIDEO**

Sirona Connect Video

JOIN US

DIGITAL IMPRESSION WITH SIRONA CONNECT.

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The first (fourth) quadrant is scanned up to the opposite second front tooth by moving the camera in parallel along the jaw arch.

Start the scanning process

- Position the camera occlusally above the last tooth on the right, to start the scanning process.

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The first (fourth) quadrant is scanned up to the opposite second front tooth by moving the camera in parallel along the jaw arch.


Start the scanning process

- Position the camera occlusally above the last tooth on the right, to start the scanning process.

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Completing the scanning process

1. Start as indicated above, on the occlusal surface of the right terminal tooth, and scan it occlusally.
Tilt the camera by 45° in a palatal direction (oral) and guide it from the distal to the mesial.
2. Tilt the camera another 45° in a palatal direction (oral) and move it in a distal direction.
3. Tilt the camera by 90° on to the occlusal surface and move it in a mesial direction.
4. Tilt the camera in a 45° buccal direction and move it back towards the distal.
5. Then tilt the camera a further 45° in a buccal direction to a total of 90° and move it in a mesial direction again.

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The screenshot shows the Sirona Connect software interface. At the top, there are navigation icons for ADMINISTRATIVEN, AUFNAHME, MODELL, and CONNECT. The main area is divided into three sections: Patient, Fall, and Vorschau. The Patient section contains a table with the following data:

Zahnarzt	Nachname	Vorname	Patienten-ID	Geburtsdatum
			4444	
FH			05444	
Multiple			087098787	
Bar			76544	
Abutments			977655 (1)	
Telescopes			977655	
Connect	Angulo	Lucy		11.05.2015
Service	beta		1	
Full Arch Scan	Omniscan	Omniscan-Full Arc...		02.02.2015
Whisker	Doug			03.03.1966
yau	mak88			09.09.1999

The Fall section contains a table with the following data:

Name

The Vorschau section contains a table with the following data:

Zahnum...	Restaurationstyp	Designmodus	Material

At the bottom right of the Vorschau section, there is a summary of system information:

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The screenshot displays the Sirona Connect software interface. At the top, there is a navigation bar with icons for ADMINISTRATION, AUFNAHME, MODEL, and CONNECT, along with the Sirona CONNECT logo. The main content area features a 'Neuer Patient' form with the following fields:

- Nachname: [Empty text input field]
- Vorname: [Empty text input field]
- Geburtsdatum: dd.MM.yyyy [Date selection field]
- Patienten-ID: [Empty text input field]
- Zahnarzt: [Empty text input field]

At the bottom of the form, there are three buttons: 'Zurück' (Back), 'Neuer Patient' (New Patient), and 'Neuen Fall hinzufügen' (Add New Case). A 'Neuer Patient' icon is also visible at the bottom center of the interface.

New Patient

	Last Name:	De
	First Name:	
	Date of Birth:	
	Patience:	
	Dental:	

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

Date Due:

26/02/2014 10:00 AM

CC = Crown veneered | PO = Pontic | CM = Missing | C = Crown | I/O = Inlay |
V = Veneer | A = Abutment | * = Margin line defined

Type:	Tooth Number:	Indication:	Material Class:	Material:	Stump Color:	Final Color:
Crown veneered	38	Crown veneered	See Additional Notes		B1	A2
Bridge veneered	31	Crown veneered	Zirconium Oxide		B1	A2
Bridge veneered	32	Pontic	Zirconium Oxide		B1	A2
Bridge veneered	33	Crown veneered	Zirconium Oxide		B1	A2

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© SIRONA Dental Systems GmbH - SIRONA Connect Server Version 1.0.0.0