# Protraction Facemask

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Masque facial de protraction

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DIRECTIONS FOR USE

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## INTENDED USE

Intended to be used in orthodontic treatment of natural teeth.

## INDICATIONS FOR USE

This protraction facemask is used as a treatment option for the skeletal class III maxilla forward while stabilizing or retracting the mandible. It is also effective in closing space and correcting negative overjet

## CONTRAINDICATIONS

Do not wear the device while playing sports or other demanding physical activities.

Contraindicated for patients with a known Nickel allergy.

## WARNINGS

Metal components of facemask contain Nickel that may induce an allergic reaction in predisposed subjects

U.S. Federal law restricts this device to sale by or on the order of a health care professional.

## PRECAUTIONS

Recommended for home use only. Extreme care should be taken to avoid hitting or pulling the appliance while being worn.

While wearing or removing the device, strictly follow your dental professional's instructions.

## ADVERSE REACTIONS

Allergic Reactions

## STEP BY STEP INSTRUCTIONS

1. Position the forehead rest near the center of the patient's forehead and tighten the hexagon screw. The metallic support will oscillate for better adaptation and comfort.
2. Place the chin cup on the patient's chin and tighten the hexagon screw to the desired fit.
3. Adjust the crossbar by loosening the hexagon screw and sliding the bar along the vertical main frame. Once the desired position is achieved, tighten the crossbar with the hexagon wrench. The crossbar can also be reversed and placed on the outside of the metal frame for additional elastic force.

4. Elastics can be placed on both sides of the crossbar and attached intraorally to either bracket or band hooks.

#### WARRANTY AND LIMITED REMEDY

Manufacturer warrants that this product is free from defects in materials and manufacture.

Manufacturer's sole obligation and customer's sole remedy in the event of a claimed defect shall be limited to replacement merchandise or refund of the purchase price.

All implied warranties, including implied warranties of merchantability or fitness for any particular purpose or use are disclaimed.

#### LIMITATION OF LIABILITY

Except where prohibited by law, Manufacturer will not be liable for any loss or damages arising from this product, whether direct, indirect, special, incidental, or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability. This limitation does not apply to third party personal injury claims.

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