

ATLANTIS™ suprastructures compatibility chart



| Manufacturer/ Company | Implant system | ATLANTIS™ suprastructures prosthetic screw/ ASA screw availability | Screw Order REF number from Manufacturer/ Company* | Torque Ncm |
|------------------------------------|--|--|--|------------|
| BioComp | IMPLANT LEVEL | | | |
| | D4 Implant Level | | BC-1819-225220 | |
| | D5.2 Implant Level | | BC-1819-225220 | |
| | ABUTMENT LEVEL | | | |
| BioHorizons | D4/D5.2 Abutment Level | | BC-1824-282015 | |
| | IMPLANT LEVEL | | | |
| | Internal/ Tapered 3.5, 4.5, 5.7 | x | PXAS | 30 |
| | Internal/ Tapered/ Single 3.5, 4.5, 5.7 | x | PXAS | 30 |
| | External Implant 3.5 | | 130-300 | |
| | External Implant 4.0, 5.0, 6.0 | | 140-300 | |
| | ABUTMENT LEVEL | | | |
| | Internal D3.5/D4.5/D5.7 Abutment for screw | | PXABS | |
| | D3.5/D4.5/D5.7 Angled Abutment for screw | | PXABS | |
| | External D3.5/D4/D5 Abutment for screw | | 222-110 | |
| Biomet 3i | IMPLANT LEVEL | | | |
| | Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 | x | IUNIHT | 20 |
| | Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 | x | IUNIHT | 20 |
| | Certain™ XP 4/5, 5/6 | x | IUNIHT | 20 |
| | MicroMiniplant™ 3.25 | x | UNIHT | 20 |
| | Miniplant™ 3.25 | x | UNIHT | 20 |
| | Osseotite™ 3.75, 4.0, 5.0, 6.0 | x | UNIHT | 20 |
| | Osseotite™ XP 4/5, 5/6 | x | UNIHT | 20 |
| | ABUTMENT LEVEL | | | |
| | Standard Abutment | | GSH30 | |
| | Conical Abutment D3.4 | | GSH30 | |
| | Conical Abutment D4.1/D5 | | GSH30 | |
| | TG Hex Abutment | | GSH30 | |
| | IOL Abutment | | GSH30 | |
| Low Profile Abutment | x | LPCTSH | 10 | |
| Camlog | IMPLANT LEVEL | | | |
| | Screw-line Implant 3.3 | x | J4005.1601 | 20 |
| | Screw-line/Root-line Implant 3.8, 4.3 | x | J4005.1601 | 20 |
| | Screw-line/Root-line Implant 5.0, 6.0 | | J4005.2001 | |
| | CONOLOG D3.3, D3.8, D4.3 Implant Level | | C4005.1601 | |
| | CONOLOG D5 Implant Level | | C4005.2001 | |
| | ABUTMENT LEVEL | | | |
| | D3.3/D3.8/D4.3 Abutment level | | J4005.1602 | |
| D5/D6 Abutment level | | J4005.2002 | | |
| Vario SR Abutment D3.8/D4.3, D5/D6 | | J4005.2004 | | |
| DENTSPLY Implants | IMPLANT LEVEL | | | |
| | ASTRA TECH Implant System™ | | | |
| | OsseoSpeed™ EV 3.0 Green | x | 25203 | 25 |
| | OsseoSpeed™ EV 3.6 Purple | x | 25204 | 25 |
| | OsseoSpeed™ EV 4.2 Yellow | x | 25205 | 25 |
| | OsseoSpeed™ EV 4.8 Blue | x | 25206 | 25 |
| | OsseoSpeed™ EV 5.4 Brown | x | 25207 | 25 |
| | OsseoSpeed™ EV Profile 4.2 | x | 25205 | 25 |
| | OsseoSpeed™ EV Profile 4.8 | x | 25206 | 25 |
| | OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow | x | 24189 | 15 |
| | OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua | x | 24449 | 20 |
| | OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac | x | 24209 | 25 |
| | OsseoSpeed™ TX Profile 4.5/5.0 | x | 24209 | 25 |
| | OsseoSpeed™ 3.5 – Yellow (Previous) | | 24208 | |
| | OsseoSpeed™ 4.0 – Blue (Previous) | | 22486 | |
| | XiVE® | | | |
| | XiVE® S 3.0 | x | 46-4304 | 24 |
| | XiVE® S 3.4, 3.8, 4.5, 5.5 | x | 46-4305 | 24 |
| | Frialit™ | | | |
| | Frialit™ 3.4, 3.8, 4.5, 5.5 | x | 46-4305 | 24 |
| | IMZ | | | |
| | IMZ TwinPlus D3.3 Implant Level (After 1995) | | 51-1209 | |
| | IMZ TwinPlus D4 Implant Level (After 1995) | | 51-1208 | |

* Note: We recommend to check screw REF number and torque Ncm with the manufacturers, excluded DENTSPLY Implants.

Note: ATLANTIS suprastructures assortment can vary from country to country. Please check ATLANTIS WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.

| Manufacturer/ Company | Implant system | ATLANTIS™ suprastructures prosthetic screw/ ASA screw availability | Screw Order REF number/ from Manufacturer/ Company * | Torque Ncm |
|--------------------------|-----------------------------------|---|---|---------------|
| | ABUTMENT LEVEL | | | |
| | ASTRA TECH Implant System™ | | | |
| | UniAbutment EV | x | 25481 | 15 |
| | Angled Abutment EV | x (ASA not available) | 25481 | 15 |
| | 20° UniAbutment | x | 22435 | 15 |
| | 45° UniAbutment | x | 22435 | 15 |
| | Angled Abutment 20° | x (ASA not available) | 22435 | 15 |
| | ANKYLOS® | | | |
| | ANKYLOS® Balance Base Narrow D4.2 | x | 3105 6022 | 10 |
| | ANKYLOS® Balance Base D5.5 | x | 3106 6022 | 10 |
| | XiVE® | | | |
| | XiVE® TG 3.4, 3.8, 4.5 | x | 46-4300 | 24 |
| | XiVE® MP 3.4, 3.8, 4.5 | x | 46-4300 | 24 |
| | XiVE® MP 5.5 | x | 46-4305 | 24 |
| | IMZ | | | |
| | IMZ TwinPlus D3.3 Abutment Level | | 51-1209 | |
| | IMZ TwinPlus D4 Abutment Level | | 51-1208 | |
| Dyna | IMPLANT LEVEL | | | |
| | D3 | | 82E50 | |
| | D3.6 | | 31AS1 | |
| | D4 | | 31AS2 | |
| | Octa Lock D4 | x | 81AS0 | 30 |
| | Helix TM Octa | x | 81AS0 | 30 |
| | Octa Extension Level | | 82E50 | |
| Keystone Dental | IMPLANT LEVEL | | | |
| | PrimaConnex SD 3.3,3.5 | x | 45060K | 30 |
| | PrimaConnex RD 4.0,4.1 | x | 45060K | 30 |
| | PrimaConnex WD 5.0 | x | 45060K | 30 |
| | Genesis 3.8, 4.5, 5.5, 6.5 | x | G71125 | 30 |
| | Renova SDI | | 43116 | |
| | Renova RDI | | 43116 | |
| | Restore SD 3.3, 3.5 | | R9203-34-48K | |
| | Restore RD 3.75, 4.0 | | R9202-40-48K | |
| | Restore WD 5.0, 5.5, 6.0 | | R9203-50-48K | |
| | MultiUnit RD | | 45164K | |
| | Stage1 RDS Abutment Level | | S2408 | |
| | Detachable Abutment | | R9460-48K | |
| | Standard Abutment | | R9200-40-48K | |
| Nobel Biocare | IMPLANT LEVEL | | | |
| | NobelActive NP 3.5 | x | 36917 | 35 |
| | NobelActive RP 4.3, 5.0 | x | 36918 | 35 |
| | NobelReplace NP 3.5 | x | 29474 | 35 |
| | NobelReplace RP 4.3 | x | 29475 | 35 |
| | NobelReplace WP 5.0 | x | 29475 | 35 |
| | NobelReplace 6.0 | x | 29475 | 35 |
| | Branemark System NP 3.3 | | 29282 | |
| | Branemark System RP 3.75, 4.0 | x | 29283 | 35 |
| | Branemark System WP 5.0, 5.5 | | 29284 | |
| | ABUTMENT LEVEL | | | |
| | Multi-Unit Abutment RP | x | 29285 | 15 |
| | Multi-Unit Abutment WP | | 29286 | |
| | Standard Abutment | | 32983 | |
| | Estheticone Abutment | | 32981 | |
| | Steri-Oss PME Abutment | | 2346 | |

* **Note:** We recommend to check screw REF number and torque Ncm with the manufacturers, excluded DENTSPLY Implants.

Note: ATLANTIS suprastructures assortment can vary from country to country. Please check ATLANTIS WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.

Please turn over →

| Manufacturer/ Company | Implant system | ATLANTIS™ suprastructures prosthetic screw/ ASA screw availability | Screw Order REF number from Manufacturer/ Company * | Torque Ncm |
|--------------------------------------|---|---|--|---------------|
| Straumann | IMPLANT LEVEL | | | |
| | Bone Level 3.3 NC | x | 0254.900 | 35 |
| | Bone Level 4.1, 4.8 RC | x | 0254.900 | 35 |
| | NNC | | 048.314 | |
| | Standard Plus 3.3 NN (Narrow Neck) | | 049.177 | |
| | Standard/Standard Plus 4.8 RN (Regular Neck) | x | 048.356 | 35 |
| | Standard/Standard Plus 4.8 WN (Wide Neck) | x | 048.356 | 35 |
| | ABUTMENT LEVEL | | | |
| | Monotype 8° Conus | x | 048.350 | 15 |
| | RN Abutment Level | x | 048.350 | 15 |
| | RN Abutment Level (Old) | x | 048.350 | 15 |
| | WN Abutment Level | x | 048.350 | 15 |
| | Bone Level Multi-Base Abutment D3.5, 4.5, 6.5 | x | 025.0900-04 | 15 |
| | Bone Level Multi-Base Angled Abutment | x | 025.0900-04 | 15 |
| | Screw-Retained Abutment 3.5 | x | 023.4763 | 15 |
| | Screw-Retained Abutment 4.6 | x | 023.4763 | 15 |
| Zimmer Dental | IMPLANT LEVEL | | | |
| | Tapered Screw-Vent™ 3.5, 4.5, 5.7 | x | MHLAS | 30 |
| | Screw Vent™ 3.3, 3.7, 4.5 | x | MHLAS | 30 |
| | SwissPlus™ 3.8 Narrow | x | MHLAS | 30 |
| | SwissPlus™ 4.8 | | GPCAS | |
| | Advent D4.5 | | AVHLS | |
| | Spline 3.25, 4, 5mm | | 1537 | |
| | ABUTMENT LEVEL | | | |
| | Screw Vent™ TSI Abutment | | TSF | |
| | Tapered Abutment | x | SCTS | 20 |
| | Spectra Cone Abutment | x | SCTS | 20 |
| Spline Shouldered Abutment D4.5, 6.5 | | 1621 | | |

* **Note:** We recommend to check screw REF number and torque Ncm with the manufacturers, excluded DENTSPLY Implants.

Note: ATLANTIS suprastructures assortment can vary from country to country. Please check ATLANTIS WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.