

Enhance® Finishing System

Coarse when you need it, fine when you want it.

The Enhance® Finishing System accomplishes the entire process of excess removal, contouring and finishing simply by varying the pressure. Push harder to remove excess composite, ease it off for final results. In just one step the Enhance® Finishing System delivers a smooth, contoured surface with a natural looking finish. It's all about pressure. **Finish in one step.**



NEW






Enhance® mini

Improved access and versatility

Fewer Steps. Less Time. Impressive Results.


More than 70% of dentists say that simplicity is one of the most important factors for a finishing and polishing system. That's why we designed the Enhance® Finishing System. Achieve the entire finishing process in one step. It doesn't get simpler than that!

Apply more pressure for excess removal and initial contouring; back it off for a naturally smooth surface. With no need to change instruments, you can stay focused on the procedure to achieve a better finish, faster.



Maximum gloss per system

After achieving a natural looking finish with Enhance®, create a high-end lustre with the Enhance® PoGo® diamond polishing system.



Ordering Information

Enhance® Finishing System Refills

- 40 Finishing Discs **624045**
- 40 Finishing Cups **624055**
- 40 Finishing Points **624065**

Enhance® mini Intro Kit **624385**

Enhance® mini 20 each: Discs, Cups, Points

Enhance® mini Finishing System Refills

- 40 Finishing Discs **624345**
- 40 Finishing Cups **624355**
- 40 Finishing Points **624365**

Enhance®/Enhance® PoGo® Complete Kit **662020**

Enhance® Finisher 10 each: Discs, Cups, Points
 Enhance® PoGo® Polisher 10 each: Discs, Cups, Points



Dentsply Sirona Pty Ltd

Australia Practice Consumables: 1300 552 929
 11-21 Gilby Road, Mount Waverley, VIC 3149
www.dentsplysirona.com.au

New Zealand Practice Consumables: 0800 336 877
 Unit 1, Building H, 4 Orbit Drive, Rosedale, Auckland 0632
www.dentsply.co.nz