

















Enhance® / Enhance® PoGo®

Finishing and Polishing Tips



<p>Site</p> <p>Steps</p>	<p>Occlusal Concave</p> 	<p>Proximal</p> 	<p>Incisal Edge Flat</p> 	<p>Convex</p> 
<p>Gross Reduction</p>	<p>Diamond or Carbide bur</p> 	<p>Diamond or Carbide bur</p> 	<p>Diamond or Carbide bur</p> 	<p>Diamond or Carbide bur</p> 
<p>Excess Removal + Contouring + Finishing</p>	<p>Finishing Point</p>  <p>Enhance®</p>	<p>Finishing Disc or Finishing Cup</p>  <p>Enhance®</p>	<p>Finishing Disc</p>  <p>Enhance®</p>	<p>Finishing Cup</p>  <p>Enhance®</p>
<p>Polishing High Gloss</p>	<p>Polishing Point</p>  <p>Enhance® PoGo®</p>	<p>Polishing Disc or Polishing Cup</p>  <p>Enhance® PoGo®</p>	<p>Polishing Disc</p>  <p>Enhance® PoGo®</p>	<p>Polishing Cup</p>  <p>Enhance® PoGo®</p>

Impressive Results

Ordering Information



Class II DO restoration with Ceram.x Spectra™ ST after finishing with Enhance®



Class II DO restoration with Ceram.x Spectra™ ST, final result after polishing with Enhance® PoGo®

Enhance® Finishing System Refills

40 Finishing Discs	624045
40 Finishing Cups	624055
40 Finishing Points	624065



Enhance® mini Intro Kit

Enhance® mini	624385
20 each: Discs, Cups, Points	



Enhance® mini Finishing System Refills

40 Finishing Discs	624345
40 Finishing Cups	624355
40 Finishing Points	624365

Enhance® PoGo® Polishing System Refills

40 Polishing Discs	662010
40 Polishing Cups	662023
40 Polishing Points	662022

