

Enhance® / Enhance® mini / Enhance® PoGo®

Procedure Comparison



	Excess Removal	Contouring	Finishing	Polishing High Gloss
Enhance® & Enhance® mini Finishing System Enhance® PoGo® Polishing System				
	Optional:			Optional:
Sof-Lex™, 1 Finishing and Polishing Systems (Discs) ²				
Sof-Lex™, 1 Finishing and Polishing Systems (Discs) ² Sof-Lex™, 1 Diamond Polishing System (Spiral)				<p>Pre-Polishing Spiral</p> <p>Diamond-Polishing Spiral</p>

¹ Not a registered trademark of Dentsply Sirona, Inc.

² "Skipping a grit size in the polishing sequence may reduce the quality of the restoration's polish." (Sof-Flex™ Instructions for Use)

Impressive Results Ordering Information



Class II DO restoration with Ceram.x Spectra™ ST after finishing with Enhance®



Class II DO restoration with Ceram.x Spectra™ ST, final result after polishing with Enhance® PoGo®

Enhance® Finishing System Refills

- 40 Finishing Discs **624045**
- 40 Finishing Cups **624055**
- 40 Finishing Points **624065**



Enhance® mini Intro Kit

- Enhance® mini **624385**
- 20 each: Discs, Cups, Points



Enhance® mini Finishing System Refills

- 40 Finishing Discs **624345**
- 40 Finishing Cups **624355**
- 40 Finishing Points **624365**

Enhance® PoGo® Polishing System Refills

- 40 Polishing Discs **662010**
- 40 Polishing Cups **662023**
- 40 Polishing Points **662022**

