

Kompatibilitätsübersicht inPost zum Scannen auf Multi-Unit Abutments von nt-trading

Implantat-hersteller	Implantatsystem	Plattform	Multi-Unit Abutments	nt-trading: 2-CONnect	SIRONA		
			Serie		inPost	REF.	
BioHorizons							
	Tapered plus, Tapered internal, Tapered internal tissue level, Internal dental implant, Single stage dental implants	■ 3.5 ■ 4.5	R-Serie	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483	
Biomet 3i							
	NanoTite Certain® PREVAIL, OSSEOTITE Certain PREVAIL, NanoTite Tapered Certain PREVAIL, NanoTite Parallel Walled Certain, OSSEOTITE Parallel Walled Certain, Full OSSEOTITE Parallel Walled Certain, OSSEOTITE XP Certain, Full OSSEOTITE XP Certain, NanoTite Tapered Certain, OSSEOTITE Tapered Certain, Full OSSEOTITE Tapered Certain	■ 3.4 ■ 4.1	H-Serie	H 800 S H 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483	
Dentsply Implants (AstraTech, Friadent)							
Astra Tech	OsseoSpeed TX	■ 3.5/4.0 ■ 4.5/5.0	S-Serie	S 800 S S 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639	
Friadent	Frialit / Xive	■ 3.4 ■ 3.8 ■ 4.5	T-Serie	T 800 S T 805 S T 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6559483 6551639 6551639	
MIS							
	Seven (internal hex.)	■ Standard ■ Wide	R-Serie	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483	
Nobel Biocare							
	NobelReplace® (Straight, Tapered, Platform Switching), Replace Select (Straight, Tapered, TC), NobelSpeedy Replace	■ RP ■ WP	E-Serie	E 810 S E 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639	
	Nobel Active, NobelReplace Conical Connection	■ NP ■ RP	F-Serie	F 800 S F 810 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639	
Osstem / Hiossen							
	Osstem TS (USA: Hiossen ET)	■ Standard USA:Regular	S-Serie	S 820 S	inPost for 2-CONnect KS61	6551639	
Straumann							
	Standard (Roxolid), Standard Plus, Tapered Effect	RN (4,8 mm) WN (6,5 mm)	N-Serie	N 810 S N 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639	
	Bone Level	NC (3,3 mm) RC (4,1/4,8 mm)	L-Serie	L 800 S L 810 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639	
Zimmer							
	Tapered Screw-Vent	■ 3.5 ■ 4.5	R-Serie	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483	

Technische Änderungen und Irrtümer vorbehalten, Best.-Nr. A91100-8945, Printed in Germany, Dispo-Nr. 04605, 1015.V0