

T2 S-Line/T3 S-Line - Overview of technical data

Gear		Ratio/Reduction	Speeds	Head diameter	Head height	Equipment	REF.
T2 S-Line AS 200 L	⚙️ ●	1:5	Up to approx. 200,000	approx. 8.7 mm	approx. 14.6 mm	Ø 1.6 FG	65 76 065
T2 S-Line AS 40 L	⚙️ ●	1:1	Up to approx. 40,000	approx. 8.7 mm	approx. 12.4 mm	Ø 2.35 WB	65 76 073
T2 S-Line AS 6 L	⚙️ ●	6:1	Up to approx. 6,000	approx. 8.7 mm	approx. 12.4 mm	Ø 2.35 WB	65 76 081

Contra-angle handpiece, red 1 : 5 standard with 4-nozzle spray, Contra-angle handpieces, blue 1 : 1 and green 6 : 1 with 1-nozzle spray.

Gear		Ratio/Reduction	Speeds	Head diameter	Head height	Equipment	REF.
T3 S-Line ES 200	●	1:5	Up to approx. 200,000	approx. 8.7 mm	approx. 14.6 mm	Ø 1.6 FG	65 76 123
T3 S-Line ES 40	●	1:1	Up to approx. 40,000	approx. 8.7 mm	approx. 12.4 mm	Ø 2.35 WB	65 76 131
T3 S-Line ES 6	●	6:1	Up to approx. 6,000	approx. 8.7 mm	approx. 12.4 mm	Ø 2.35 WB	65 76 149

Contra-angle handpiece, red 1 : 5 standard with 4-nozzle spray, Contra-angle handpieces, blue 1 : 1 and green 6 : 1 with 1-nozzle spray.

THE DENTAL
SOLUTIONS
COMPANY™



Subject to technical changes and errors. Order no. A91100-M45-C015-01-7600. Printed in Germany. Dispo-Nr. 04603-GI, OEW17 WS 0117V0





You don't have to be big to be great.

Outstanding ergonomics for the greatest precision. Meet the new line of contra-angle handpieces from Dentsply Sirona. The T2 S-Line and T3 S-Line instrument series meet these criteria, allowing you to experience outstanding ergonomics: With its short overall length, slim design, and reduced overall weight, each S-Line instrument lies perceptibly well-balanced in your hand. Its ideal balance allows you to work without strain or fatigue, even during longer treatments.

T2 S-Line / T3 S-Line - We make it shorter so you have it easier.

- Considerably shorter, thus improved ergonomics
- Perceptibly lighter thanks to compact size
- For fatigue-free work
- Smallest head diameter with \varnothing 8.7 mm
- With a universal ISO interface for motors with short coupling element
- Titanium sleeve for a comfortable grip
- Optimal view thanks to small head dimensions
- Optimal spray system
- Perfect hygiene - can be sterilized and thermally disinfected

Exclusively for the T2 S-Line

- 2-year warranty
- Adaptable for motors with light thanks to the ISOslider®
- Optimal illumination for an excellent view during treatment



Fatigue-free work thanks to the ergonomic design.