
















Overview Motors

Name	Image	REF (with tube)	Preis w/o VAT	REF (w/o tube)	Preis w/o VAT	Connection type	Tube REF. / coupling color (motor side)	Sanitation adapter	Sterilizabl* Brushless	Rotation per minute (min-1)	Torque	Gold Contact for ApexLocator function	Anti Retraction stop	Guarantee (in years)	Repexchange available	Suitable for	
Electric motor SL		Tube with LED: 6321231 M1/M1+/C1/C1+/C2/C2+/C4/C4+/C5+ 6321249 C3/C3+/C6 6320928 C8+	1,809 € 1,835 € 1,500 €	4176710	1,254 €	Suitable for T1 Classic system with Sirona-interface	With LED: 6319607 M1/M1+/C1/C1+/C2/C2+/C4/C4+/C5+ 6319623 C3/C3+/C6 6319631 C8 / 6319649 C8+	 4705732 6342377 (green)	- -	200 - 40,000	- 2.4	- -	- -	2 +	+	M1/M1+ / C/C+ Serie	
Electric motor BL with LED light		6320332	2,456 €	6319714	2,017 €	Suitable for T1 Classic system with Sirona-interface	6283118 (blue) 	4705732 6342377 (green)	+ +	90 - 40,000	- 2.4	+ +	+ +	2 +	+	Teneo / Sinius	
Electric moto BL E with LED		6557388	1,970 €	6546225	1,525 €	Suitable for T1 Classic system with Sirona-interface	6336700 (green) 	4705732 6342377 (green)	+ +	1,200 - 40,000 100 - 40,000 retrofit with eControl	- 2.4	+ +	+ +	2 +	+	Intego / Intego Pro	
Electric motor BL ISO C with LED light		6354356	2,611 €	6337518	2,102 €	Suitable for T1/T2/T3/T4 Line, T2-T3 S-Line system as well as other contra-angle handpieces with ISO 3964 interface or INTRAMatic* LUX	6336700 (green) 	5862920 6342385 (blue)	+ +	90 - 40,000	- 3.0	+ +	+ +	2 +	+	Teneo / Sinius	
Electric motor BL ISO E with LED light		Tube w/o ApexLocator function: 6496447 (SINIUS) 6445840 (incl. electronics) (INTEGO / INTEGO pro)	1,812 € 2,091 €	6445899	1,621 €	Suitable for T1/T2/T3/T4 Line, T2-T3 S-Line system as well as other contra-angle handpieces with ISO 3964 interface or INTRAMatic* LUX	Hose w/o ApexLocator function: 6336718 (green)  Hose with ApexLocator function: 6336700 (green)	5862920 6342385 (blue)	+ +	2,000 - 40,000 (Sinius) 1,200 - 40,000 (Intego / Intego Pro*) *100 - 40,000 retrofit with eControl	- 3.0	+ +	+ +	2 +	+	Sinuis / Intego / Intego Pro	
Electric motor BL Implant without light and spray, without outflow of air		6208198	2,011 €	6212208	1,631 €	Suitable for Implant 20:1 contra-angle (ISO 3964 interface)	6212182 (yellow) 	5862920 6342385 (blue)	+ +	90 - 40,000	- 5.0	+ -	- -	2 +	+	Teneo	
Air motor T4 BM/BO		n.a.	n.a.	BM 6439884 BO 6439892	270 €	Suitable for T4 Line system as well as other contra-angles and handpieces with ISO 3964 interface	Midwest 4-hole / Borden 2-hole	5862920 6342385 (blue)	+ +	90 - 18,000	-	- -	- -	- -	- -	-	Borden / Midwest connection


* Valid for motors

Overview Motors









Name	Image	REF (with hose)	Preis w/o VAT	REF (w/o tube)	Preis w/o VAT	Connection type	Tube REF. / coupling color (motor side)	Sanitation adapter	Sterilizable*	Brushless	Rotation per minute (min-1)	Adjustable torque	Gold Contact for ApexLocator function	Anti Retraction stop	Guarantee (in years)	Repexchange available	Suitable for
Adapter Basic Apex with light and without spray		n.a.	n.a.	5983072	365 €	ISO 3964 or INTRAMatic® LUX inter- face adapter to fit e.g. contra-angle Endo 6 L to motor SL or motor BL for ApexLocator function	n.a.	5862920 6342385 (blue)	+	n.a.	n.a.	n.a.	+	n.a.	2	n.a.	Motor BL / motor SL
Adapter ISO with light and spray		n.a.	n.a.	6000793	747 €	ISO 3964 or INTRAMatic® LUX interface adapter to fit e.g. T1 Line or contra- angle Endo 6 L to motor BL or motor SL (not suitable for ApexLocator function)	n.a.	5862920 6342385 (blue)	+	n.a.	n.a.	n.a.	n.a.	n.a.	2	n.a.	Motor BL / motor SL

* Valid for motors

Motors by treatment unit

for Teneo																		
Electric motor BL with LED light		6320332	2,456 €	6319714	2,017 €	Suitable for T1 Classic system with Sirona- interface	6283118 (blue)			4705732 6342377 (green)	+	+	90-40,000	- 2.4	+	+	2	+
Electric motor BL ISO C with LED light		6354356	2,611 €	6337518	2,102 €	Suitable for T1/T2/T3/ T4 Line, T2-T3 S-Line system as well as other contra-angles handpieces with ISO 3964 interface or INTRAMatic® LUX	6336700 (green)			5862920 6342385 (blue)	+	+	90-40,000	- 3.0	+	+	2	+
Electric motor BL Implant without light and spray, without outflow of air		6208198	2,011 €	6212208	1,631 €	Suitable for Implant 20:1 contra-angle (ISO 3964 interface)	6212182 (yellow)			5862920 6342385 (blue)	+	+	90-40,000	- 5.0	+	-	2	+
for Sinius																		
Electric motor BL with LED light		6320332	2,456 €	6319714	2,017 €	Suitable for T1 Classic system with Sirona- interface	6283118 (blue)			4705732 6342377 (green)	+	+	90-40,000	- 2.4	+	+	2	+
Electric motor BL ISO C with LED light		6354356	2,611 €	6337518	2,102 €	Suitable for T1/T2/T3/ T4 Line, T2-T3 S-Line system as well as other contra-angles handpieces with ISO 3964 interface or INTRAMatic® LUX	6336700 (green)			5862920 6342385 (blue)	+	+	90-40,000	- 3.0	+	+	2	+
Electric motor BL ISO E with LED light		Hose w/o ApexLocator function: 6496447	1,812 €	6445899	1,621 €	Suitable for T1/T2/T3/ T4 Line, T2-T3 S-Line system as well as other contra-angle handpieces with ISO 3964 interface or INTRAMatic® LUX	6336718 (green)			5862920 6342385 (blue)	+	+	2,000- 40,000	- 2.4	+	+	2	+

Overview Motors

Name	Image	REF (with tube)	Preis w/o VAT	REF (w/o tube)	Preis w/o VAT	Connection type	Tube REF. / coupling color (motor side)	Sanitation adapter	Sterilizable*	Brushless	Rotation per minute (min-1)	Torque	Gold Contact for ApexLocator function	Anti Retraction stop	Guarantee (in years)	Repexchange available
for Intego/Intego Pro																
Electric motor BL ISO E with LED light 	Tube w/o ApexLocator function 6445840 (incl. electronics)	2,091 €	6445899	1,621 €	Suitable for T1/T2/T3/ T4 Line, T2-T3 S-Line system as well as other contra-angles/handpieces with ISO 3964 interface or INTRAmatic® LUX	Hose w/o ApexLocator function 6336718 (green) 	5862920 6342385 (blue)	+	+	1,200- 40,000	- 3.0	+	+	2	+	
	Tube with ApexLocator function 6445857	2,120 €	Hose with ApexLocator function 6336700 (green) 	100 - 40,000 (Intego Pro: retrofit with eControl)												
Electric moto BL E with LED 	6557388	1,970€	6546225	1,525 €	Suitable for T1 Classic system with Sirona- interface	6336700 (green) 	4705732 6342377 (green)	+	+	1,200 - 40,000 100 - 40,000 retrofit with eControl	- 2.4	+	+	2	+	
for Midwest or Borden connection																
Air motor T4 BM/BO 	n.a.	n.a.	BM 6439884 BO 6439892	270 €	Suitable for T4 Line system as well as other contra angles hand- pieces with ISO 3964 interface	Midwest 4-Loch/ Borden 2-Loch	5862920 6342385 (blue)	+	+	90-18,000	-	-	-	-	-	
Additional																
Adapter Basic Apex with light and without spray 	n.a.	n.a.	5983072	365 €	ISO 3964 or INTRAmatic® LUX interface adapter to fit e.g. contra-angle Endo 6 L to motor SL or motor BL for ApexLocator function	n.a.	5862920 6342385 (blue)	+	n.a.	n.a.	n.a.	n.a.	+	n.a.	2	n.a.
Adapter ISO with light and spray 	n.a.	n.a.	6000793	746 €	ISO 3964 or INTRAmatic® LUX interface adapter to fit e.g. T1 Line or contra- angle Endo 6 L to motor BL or motor SL (not suitable for ApexLocator function)	n.a.	5862920 6342385 (blue)	+	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	2	n.a.

* Valid for motors