

inLab 3D V3.85 / V3.86 Model Milling Compatibility

To be able to utilize the model milling capability of the inLab 3D V3.85 & V3.86 software pack, a specially equipped inLab MC XL is required. Units with serial numbers 120000 and above are equipped with special milling motor assemblies and drain systems that are designed to mill models.

A conversion kit is required on inLab MC XL units below this serial number. As of June 1, 2011, the conversion kit availability and pricing are being established and will be released in the near future.

**IMPORTANT NOTE:**

Damage can occur to inLab MC XL units below serial number 120000 that do not have this conversion kit installed. As such, damages to inLab MC XL units in conjunction with milling models will not be covered under warranty by Sirona Dental Systems.

