

Hot Keys / ShortCuts for inLab SW4.02

ShortCuts		
Key(s)	Additional Key	Feature
Ctrl	N	Stratscreen
Ctrl	S	Save
Ctrl + Alt	S	Save As...
Ctrl	I	Import
Ctrl	E	Export
Ctrl	Q	Exit
Ctrl	F1	Adminphase
Ctrl	F2	Scanphase
Ctrl	F3	Modephase
Ctrl	F4	Designphase
Ctrl	Num 3	Buccal View
Ctrl	Num 5	Center
Ctrl	Num 2	Cervical View
Ctrl	Num 9	Distal View
Ctrl	Num 7	Lingual View
Ctrl	Num 1	Mesial View
Ctrl	Num 8	Occlusal View
Ctrl	O	Contact Occlusal & Approximal
Ctrl	M	Model Contacts
Ctrl	R	Colored restoration
Ctrl	B	Model Box
Ctrl	C	Slice
Ctrl	D	Cursor Details
Ctrl	U	Distance

ShortCuts		
Key(s)	Additional Key	Feature
Ctrl	G	Grid
Ctrl	P	Monocolor Models (Omnica)
Alt	V	ViewOptions
Alt	T	Toolmenu
Alt	D	Display Objects
Alt	A	Analyzing Tools
Alt	F	Formtool
Alt	M	Movetool
Alt	E	Edit Connector Lines
Alt	S	Shapetool
Alt	R	Reductetool
Alt	N	Recalculate
Alt	C	Contacttool
Alt	U	Unlock
Alt	H	Scaletool
Alt	B	Bioslider
Alt	X	Splittool
Alt	M	Connector Position
Alt	S	Connector Scale
Alt	F	Connector Shape

	Return	Ok, Yes, Continue, Start, Retry
	ESC	Cancel, No, Stop, Abort
	F1	Help

To toggle between commands use "Space-Bar".