

VITA Enamic Material Pack for CEREC SW \geq 4.0.3

Description

This software material pack updates existing CEREC SW \geq 4.0.3 installations and adds the VITA Enamic material as an additional material option. There are no other changes (features or fixes) in this material pack. This material pack is compatible with CEREC AC Bluecam (requires CEREC SW 4.0.3) and CEREC AC Omnicam (requires CEREC SW 4.0.5).

The material pack contains the following changes:

- **VITA Enamic material support**

The VITA Enamic material is added to the material section in the Administration phase of the software.