

SERVICE PACK inLab SW 4.0.2

Update Description

Service Pack inLab SW 4.0.2 updates inLab SW 4.0.0 and 4.0.1 installations and includes the following modifications:

■ **Compatibility with CEREC Omnicam data**

CEREC Omnicam files including the color information can be imported into the inLab SW 4.0.2:

- Files from CEREC SW 4.0.4 can be imported in *.dxd format.
- Files from CEREC Connect SW 4.1.1 can be imported via the Sirona Connect Portal.

■ **File size reduction in *.i.lab format**

Restoration file size saved in *.i.lab format has be reduced.

This Service Pack inLab SW 4.0.2 includes all changes from Service Pack inLab SW 4.0.1.

We recommend all inLab 4.0 users install this service pack.