

SERVICE PACK inLab 3D V3.88

# Update description

The Service Pack inLab 3D V3.88 contains the following functional enhancements:

■ **3M ESPE Lava™ Ultimate**

Under inLab 3D V3.88 it is possible to select the material 3M ESPE Lava™ Ultimate in the block sizes 12 and 14L. Please contact 3M ESPE directly for further details about this new material, as well as the scheduled launch dates in the various national markets.

The Service Pack inLab 3D V3.88 is available as a 64-bit version (inEos Blue) and a 32-bit version (inEos). Please select the relevant version. This new Service Pack requires inLab 3D V3.80 or higher. It embraces all previous service packs.

To receive regular information about our product updates you can subscribe to the inLab Newsletter at [www.sirona.com/inlab-newsletter](http://www.sirona.com/inlab-newsletter)