

# Your individual options

Here you will find the complete list of all individual options and you can see which equipment is already part of the basic units.

## Basic unit equipment options

VB CS TS

Options for the patient chair	VB	CS	TS
Double-articulated headrest	■	■	■
Motorized headrest	□	□	□
Premium upholstery	■	■	■
Premium thermo upholstery	□	□	□
Lounge upholstery, by request with colored stitches	□	□	□
Lumbar support function	□	□	□
Swiveling armrest, right or left	□	□	□
4-way foot control	■	■	■
Foot control, cable version	■	■	■
Wireless foot control	□	□	□
External device connection	□	□	□
Booster cushion for small patients (C seat cushion)	□	□	□
Installation plate to level out uneven floors	□	□	□

Options for the dentist element (5 handpiece positions, removable handle & instrument holder)	VB	CS	TS
Sprayvit M syringe with light	□	□	□
1st drive (for motor or turbine)	■	■	■
2nd drive (for motor or turbine)	■	■	■
3rd drive (for motor or turbine)	□	□	□
BL electric motor	□	□	□
BL ISO C electric motor	□	□	□
BL Implant electric motor	□	□	□
Implant 20:1 contra-angle handpiece	□	□	□
Turbine hose	□	□	□
SiroSonic TL ultrasonic handpiece	□	□	□
Hands-free operation with cursor control	■	■	■
Implantology function	□	□	□
Saline solution pump	□	□	□
Endodontic treatment function incl. reciprocal file systems	□	□	□
Contra-angle handpiece T1 Line Endo 6 L with light	□	□	□
ApexLocator function with display on the Sivision monitor	□	□	□
Swiveling tray holder for 1 or 2 standard size trays		□	□
Panoramic X-ray image viewer			□

Options for the assistant element (with 4 positions, removable holder, and silicone mat)	VB	CS	TS
Satelec mini LED curing light	□	□	□
Surgical suction device	□	□	□
2nd HVE kit	□	□	□
Sprayvit M syringe with light	□	□	□
Suction device	■	■	■
Saliva ejector	■	■	■

## Basic unit equipment options

VB CS TS

Options for the water unit	VB	CS	TS
Disinfection device with permanent disinfection and sanitization function	■	■	■
Suction hose cleaning with water	■	■	■
Suction hose cleaning with chemicals	□	□	□
Cuspidor valve for wet suction system	■	■	■
Automatic separator	□	□	□
Amalgam separator	□	□	□
Swiveling cuspidor	■	■	■
Without cuspidor	□	□	□
Tray with clip-on tray and handle	□		

Options for Sivision digital	VB	CS	TS
SiroCam AF	□	□	□
SiroCam AF+	□	□	□
22" monitor on the tray incl. tray	□		
22" monitor on the light support frame	□	□	□
Additional holder for camera	□		□
Camera power supply in the dentist element	□	□	□
Camera supply in the assistant element		□	□
Sivision ready installation kit with USB interface	□	□	□
Cable set for PC connection	□	□	□

Hugo working stool (optional long or short gas spring)	VB	CS	TS
Hugo Freehand (with foot release)	□	□	□
Hugo Manual (with manual release)	□	□	□
Hugo Manual Plus (with manual release and foot ring)	□	□	□

Carl working stool (optional long or short gas spring)	VB	CS	TS
Carl Manual (with manual release)	□	□	□
Carl Manual Plus (optional long or short gas spring)	□	□	□

Theo working stool (optional long or short gas spring)	VB	CS	TS
Saddle Stool Theo	□	□	□

Operating light	VB	CS	TS
LEDview Plus operating light (triple-jointed support arm, with sensor)	□	□	□
LEDview Plus ceiling model	□	□	□
Adaptation for Heliodent Plus ceiling combination	□	□	□
Adaptation for Heliodent Plus unit model	□	□	□

VB Sinius (sliding track version)  
 CS Sinius CS (whip arm version)  
 TS Sinius TS (OTP version)  
 ■ Included in the basic unit  
 □ Optional