


THE DENTAL
SOLUTIONS
COMPANY™



The best place to be
somewhere else.

Sinius.



When your patients are completely comfortable, they are sitting in a Sinius.

With its smart design and extraordinary comfort, Sinius is perfectly tailored to meet the individual needs of your patients. And thanks to the integrated functions such as endodontics, implantology, and patient communication, you can provide better, faster, and safer treatment. Discover Sinius – for your patients and yourself.

What patients think

"Knowing that I'm safe builds confidence."

Find out more about
Patient Communication
Comfort

"When my dentist feels good, I feel good as a patient."

Find out more about
User comfort

"I can see and feel top quality."

Find out more about
Hygiene Design

The ideal treatment concept for every practice



Sinius – the sliding track concept

- Unique ergonomic treatment concept
- Dentist element easy to move forward and backward manually
- Swivel-mounted dentist element ensures everything is within reach



Sinius TS – the swivel tray concept

- Ergonomic positioning with independently movable dentist element and tray system
- Functional tray system with detachable tray holder and multiple positioning options
- Ample instrument area on the dentist element and tray holder



Sinius CS – the whip arm concept

- Dentist element can be rotated by up to 240° for optimal positioning
- Maximum distance to instruments is 90 cm
- Ideal for the widest range of treatment situations for working alone or with a team

For just a moment,
I feel weightless.
Only the wind roars
past my ears –
until I pull my parachute.

The best place to be
somewhere else.

Sinius.




When your patients are completely comfortable, they are sitting in a Sinius.

Transparent communication and comfortable positioning at all times helps patients relax – while you focus on treatment.

Patient Communication

- With the SIDEXIS software, you can easily access patient data, X-ray images, intraoral images, and planning views.
- The SiroCam AF/AF+ autofocus camera combines superior image quality with easy operation.
- Windows Media Player and Microsoft PowerPoint allow you to inform your patients easily and intuitively directly via the EasyTouch user interface.



Comfort

1 Ergonomic patient positioning for a comfortable treatment position:

- With the combined movement of the backrest and the seat, the ErgoMotion function ensures ergonomic spine extension.
- The patient's vertical mouth point positioning is always the same. The motor-driven headrest automatically positions the patient's head – gently and without effort.

3 Armrests for practical comfort:

- The armrests are easy to swivel back and forth to allow comfortable entry and exit.
- The organic lines of the armrests fit elegantly into the high-quality design of the treatment center.

2 Maximum lowering for a comfortable entry position:

- The gentle lift drive for adjustable height ensures easy entry and exit, even for children and elderly patients.
- The highest position is 80 cm, the lowest 36 cm.

4 Lumbar support function for positioning adapted individually to patients:

- The backrest position can be adjusted individually to each patient's spine.
- Positioning that meets the needs of elderly or persons with physical impairments in particular.



User comfort

I make my way through
the jungle, reach the
secret cavern – and raise
the hidden treasure.

The best place to be
somewhere else.

Sinius.



denstplysirona.com/sinius



When your young patients are completely comfortable, they are sitting in a Sinius. You get them to trust you through a short, pleasant visit. Sinius supports your personal treatment concept with integrated treatment functions.

User comfort

Integrated implantology

The Sinius implantology function integrates the equipment needed for treatment and simplifies your practice workflow. It is very intuitive and easy to operate via the screen.

Set your usual work steps individually:

- Speed
- Torque
- Rotation direction
- Flow quantity of sterile cooling fluid


After every work step, you can directly access the next individually saved settings.



Integrated endodontics

The integrated endodontics function eliminates conversion times as there is no additional tabletop equipment. This allows better workflows – absolutely hygienic, safe, and efficient.

- **ApexLocator display** on the SIVISION monitor as a visual assistant
- **File libraries** with step-by-step instructions for preparing the root canal
- Better view with the **T1 Line Endo 6 L contra-angle handpiece** with light
- Integrated **reciprocal file systems**



The water's sparkling
in the morning sunlight.
I turn my board – and
ride the perfect wave.

The best place to be
somewhere else.

Sinius.



When your patients are completely comfortable, they are sitting in a Sinus. They feel at ease in its excellent design and its integrated hygiene system gives them the safety of the highest treatment standard.

Hygiene

Integrated disinfection system

Continuously adds Dentosept S disinfectant for water lines to the operating water. When used consistently, reliably prevents bacterial growth.

Built-in sanitization adapter

For regularly purging the water lines and monthly sanitization. At the press of a button, all water lines are automatically purged to be RKI compliant or sanitized with the Dentosept water disinfectant.

Automatic purging programs and dosage

Schnelles, einfaches und besonders sicheres Durchspülen der Wasserwege mit den Programmen Purge.

Efficient suction hose cleaning

The suction hoses can also be cleaned easily, quickly, and fully automatically.

Removable elements

All hygiene-critical elements can be removed and either thermally disinfected or sterilized.



The DGKH (German Society for Hospital Hygiene) has confirmed Dentsply Sirona's compliance with the highest hygiene safety standards in accordance with the recommendations of the Robert Koch Institute.

Design

Of course, root canal treatment is not a dream trip, but with its outstanding design, Sinius gives you the best possible conditions.



reddot design award
winner 2012



Best lying comfort

Comfortable material, soft filling, and a contoured surface for good support



Intelligent design

Innovative thermoforming method for special protection and excellent hygiene.



High-quality material

Pleasantly soft feel, easy to clean and disinfect.

Upholstery colors

The elegant world



Classic, timeless combinations for a truly sophisticated look.

The natural world



Soft, warm shades for a soothing treatment atmosphere.

The vital world



Vivid colors and bright tones for a modern, lively atmosphere.

* Available colors for Lounge upholstery.

Stitching colors



Standard Red Dark grey Light blue Rose

Paintwork colors



Rosé Platinum silver Matt gold Neutral white



Dentsply Sirona for you

We help dentists to provide better, safer, faster dental care. A person who performs at the top of their game every single day expects the same from their equipment. This is why we provide innovative dental solutions to improve oral health. For you, this means that with Dentsply Sirona, you are always working at the highest standard. You can rely on us - not only today, but in the future as well.

Dentsply Sirona

Sirona Dental Systems GmbH
Fabrikstraße 31, 64625 Bensheim, Germany
dentsplysirona.com

Procedural Solutions

Preventive
Restorative
Orthodontics
Endodontics
Implants
Prosthetics

Enabling Technologies

CAD/CAM
Imaging
Treatment Centers
Instruments