



Classic®

Mould Chart – Ultra Value Denture Teeth

For dental professionals and laboratories



Classic[®] Denture Teeth

Everyone deserves the best denture they can afford. So, even for patients on a tight budget, a great looking denture is still possible. Classic Ultra Value denture teeth offer natural-like aesthetics, a wide selection of moulds and the most popular shades.



Classic®















Anteriors Upper

Available in 24 shades and 1 bleach shade.

All dimensions in mm.
Images scaled to size.









Mould No.
(Lower
Articulation)

Mould No.
(Lower
Articulation)





Mould No.
(Lower
Articulation)



Classic®











Anteriors Lower

Available in 24 shades and 1 bleach shade.

All dimensions in mm.
Images scaled to size.

Mould No.
(Lower
Articulation)

Mould No.
(Lower
Articulation)










Classic® Posterior Teeth

Available in 24 shades and 1 bleach shade.

Classic® 0° Posteriors








Classic® 10° Posteriors



X = The average width of a 1x4 set of upper posteriors;
Y = the lower posteriors.
U = The average depth of the upper left first molar;
L = the lower left first molar.










ISO 22112: 2005. Processing instructions available upon request.

Classic® 20° Posteriors



All dimensions in mm.
Images scaled to size.

Classic® 33° Posteriors



LS indicates Long Buccal-Short Lingual moulds.

Denture Base

Lucitone® HIPA High Impact Pour Acrylic



Lucitone HIPA is a high impact pourable acrylic denture base. It exceeds the ISO standards for impact resistance. Lucitone HIPA is a self-curing denture base material designed for use in fabrication, repair, rebasing, or relining of full and partial dentures including implant overdentures or other dental appliances. It is available in shades: Original, Light, Light Reddish Pink, Dark Pink, and Original Opaque.

Lucitone® FRS® Flexible Dental Resin


Lucitone FRS Flexible Denture Resin is a temporary thermoplastic injectable resin that offers customizable flexibility, high impact resistance and transparency, delivering aesthetic removable partial dentures and clear or colored clasps.



- Semi-crystalline nylon composition for improved color stability and stain resistance. Provides aesthetics.
- Easy to finish and polish due to the microcrystalline structure. Saves time.
- Monomer - and metal-free. Non-irritating.
- Available in original, light pink, light reddish pink, dark pink and clear shades. Aesthetic shade options.
- Clear shade for fabrication of clear clasps, occlusal splints, and night guards. Offers greater versatility.
- Available in three (3) cartridge sizes: small, medium and large. Adapts to laboratory needs.

Classic® Anterior and Posterior Shades*:

CW2, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, 59C, 62C, 65C, 66C, 67C, 69C, 77C, 81C.
REF 905600



* The A1-D4 designations correspond to the VITA classical A1-D4* shade guide which is meant to be a guide, not a match. VITA classical A1-D4 is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.KG



Introducing Smart Cabinet Technology

Simplify your inventory management. Scan barcode and reorder Dentsply Sirona teeth through your preferred distributor on a mobile device. Download the app today on iOS and Android. You can now scan your cabinet and reorder Dentsply Sirona tooth cards with ease.



Scan Order

Quickly scan barcodes in your tooth cabinets to reorder Dentsply Sirona tooth stock with the push of a button.



Notification

Receive exclusive promotions through this portal as well as updates and news. Check back often for the latest notifications.



Order History

Review orders submitted through the Lab Smarter® app. See past orders down to the mould and shade. Past orders can be printed for your records.



How it works

Step-by-step instructions for using this app.



Contact my Dentsply Sirona Representative

A digital business card to call, text or email your local Dentsply Sirona sales representative.

Take control of your tooth stock in 3 easy steps:

1. Download

Go to iTunes or Google Play and download the Lab Smarter App on your iPhone or Android device.


2. Sign up

Complete the one-time registration and wait for the email confirmation. Sign into the app using your user name and password.

- a) If you are a Dentsply Sirona distributor, check the box at the bottom to gain access to the customer management mode.

3. Snap. Tap. Smile

Now scan your cabinet and reorder Dentsply Sirona teeth through your preferred dealer.



Google Play is a trademark of Google Inc. iPhone and iPad are registered trademarks of Apple Inc. App Store is a service mark of Apple Inc.



Notes

Combination Table

Select upper anterior mould form, and use table to determine recommended lower anterior and posterior options.

Classic® Anteriors		Classic® Posteriors			
Upper Mould No.	Articulates with Lower Mould	Combines with 33° Posterior Mould	Combines with 20° Posterior Mould	Combines with 10° Posterior Mould	Combines with 0° Posterior Mould
133	2C	30M	29M	F30	29M
134	3D	30M	31M	F33	31M
135	46	32M	33M	F33	33M
136	46	34M	33M	F33	33M
2D	2D	30L	29L	F32	29M
3D	3D	30M	31M	F33	31M
1N	46	30L	31M	F32	31M
2N	2N	30L	29L	F30	29M
2P	2P	30L	31L	F33	31M
A24	3N	30M	29M	F30	29M
A25	2E	30M	31M	F33	31M
A26	46	30M	31M	F33	31M
3M	3M	30M	29M	F30	29M
3N	3N	30L	31M	F32	31M
3P	3P	34M	33M	F33	33M
4M	3M	30M	29M	F30	29M
4N	2N	30L	31M	F30	31M
263	2N	30M	29L	F32	29M
264	2E	30M	31L	F33	31M
266	26	34L	33M	F33	33M
267	3R	34L	33M	F33	33M
268	3R	34L	31L	F33	33M

Dentsply Sirona Prosthetics
 570 West College Avenue
 York, PA 17401
 1-800-243-1942
prosthetics.dentsply.com | Ask your dealer for details.

Distributed By:
 Dentsply Canada
 161 Vinyl Court
 Woodbridge, ON L4L 4A3
 1-800-263-1437