



The Nupro® Polishing Portfolio

Put their
best smile forward.



“The majority of my patients believe providing that smooth and shiny feeling with the polishing IS their professional cleaning.”

Mandy Macachor, RDH
Wisconsin

The CDC recommends that dental handpieces be **heat sterilized** between patients.¹

57%
of patients suffer from dental hypersensitivity.²

Nearly **1/3**
of dental health professionals cite musculoskeletal disorders as the reason for early retirement.³

Delivering a pleasant polishing experience is an essential step in patient satisfaction. Our complete line of Nupro Polishing hygiene handpieces, disposable prophylaxis angles, and prophylaxis pastes are efficient, easy to use and effective so you can focus more on your patient.

Committed to providing the best polishing experience

At Dentsply Sirona, we are dedicated to providing innovative solutions that enhance infection prevention standards in your office, improve patient comfort during dental procedures, and deliver products designed for better ergonomics to protect the musculoskeletal health of clinicians.

Enhanced Infection Prevention

The Nupro Freedom® Cordless Prophy System with SmartMode® Technology features an autoclavable metal outer sheath for improved infection prevention.

Improved Patient Comfort

Our Nupro® Extra Care Prophylaxis Paste with NovaMin® is clinically proven to provide immediate relief from tooth sensitivity that lasts up to 28 days from just one application.

Better Ergonomics

Both the Nupro Freedom® and the Nupro RDH® Hygiene Handpieces are ergonomically designed and may help reduce the prevalence of musculoskeletal disorders among dental practitioners.

Experience convenience with CDC compliance



Nupro® Freedom® Cordless Prophy System

The best-selling Nupro® Freedom® Cordless Prophy System is designed to enhance convenience and comply with infection prevention standards. Autoclavable outer sheaths that replace the need to buy multiple handpieces for each operator.

The CDC sterilization guidelines state: Dental handpieces and associated attachments, including low-speed motors and reusable prophylaxis angles, should always be heat sterilized **between patients!**



Disposable polyethylene barrier protects the inner module from debris

Autoclavable, sealed outer sheath, protects the inner module from contamination


Nupro Freedom single-use disposable prophy angles guarantee performance and accurate cup speed

Additional features include:

- Operates with SmartMode® Technology
- 25% more speed and power⁴
- Lightweight, cordless design
- All-day battery life
- Rapid 15-minute recharge
- Quiet operation
- Enhanced control
- Reduced wrist fatigue
- Non-slip grip design
- Portable
- Enhanced convenience
- Optional wireless footpedal for choice of control

Choice of Control

Choose between SmartMode® Technology or the wireless foot pedal option. SmartMode technology responds to the clinician's touch with gentle pressure of the prophy cup against the tooth to increase speed. The wireless foot pedal provides and more traditional polishing experience.



ITEM #	DESCRIPTION	CONTENTS
90741L01	Nupro Freedom® Cordless Prophy System (Wireless Foot Pedal Sold Separately)	1 Inner Module 1 Metal Outer Sheath 1 Charging Base 1 Handpiece Tray Cradle 1 Power Supply 20 Nupro Freedom® Slim Contra Spiral Disposable Prophy Angles 1 Disposa-Shield® Trial Pack 1 Operation Manual
90741L03	Nupro Freedom® Cordless Prophy Package (Wireless Foot Pedal Sold Separately)	1 Nupro Freedom® Cordless Prophy System (listed above) PLUS: 2 Additional Metal Outer Sheaths 1 Box Disposa-Shield® Disposable Polyethylene Barriers (500 ct) 1 Box Nupro Freedom® Slim Contra Spiral Disposable Prophy Angles (200 ct)
90741L05	Nupro Freedom® Cordless Prophy Package with SmartMode® Technology and Wireless Foot Pedal	1 Nupro Freedom® Cordless Prophy Package (listed above) PLUS: 1 Wireless Foot Pedal

OPTIONAL ACCESSORIES

- 90796 Metal Outer Sheath
- 90799 Inner Module
- 90797 Charging Base⁵
- 90742 Wireless Foot Pedal⁵
- 90798 Power Supply
- A88019 Box of 500 Disposa-Shield® Disposable Polyethylene Barriers for Cordless Prophy System Inner Module
- 9070801 Handpiece Tray Cradle




A slim shape that makes polishing more pleasant for patients and clinicians.

Nupro Freedom® Slim Disposable Prophy Angles

Designed to help you deliver a better polishing experience for your patients, our Nupro Freedom® Slim Disposable Prophy Angles provide improved visibility and access with a 14% more slender angle.⁶ Additional features include:

- 16% shorter head⁶ for accessing hard-to-reach areas
- Reduced noise and vibration for a smooth, quiet operation
- Ergonomically designed for improved clinician comfort
- Available in a unique lavender pearl color

ITEM #	DESCRIPTION	PACKAGE (QTY)
965774	Slim Contra Spiral Cup	Box (200)
965772	Slim Soft Cup	Box (200)
965773	Slim Spiral Cup	Box (200)
965066	Slim Contra Spiral Cup	Box (1200)
965064	Slim Soft Cup	Box (1200)
965065	Slim Spiral Cup	Box (1200)



Nupro Freedom® Disposable Prophy Angles

Nupro Freedom® Disposable Prophy Angles are designed to operate with the Nupro Freedom® Cordless Prophy System with SmartMode® Technology for optimal efficiency.

LAVENDER ITEM #	PINK ITEM #	DESCRIPTION	PACKAGE (QTY)
965670	965673	Contra Spiral Cup	Box (200)
965671	965674	Contra Spiral Cup	Econo Box (1200)
965750	965752	Soft Cup	Box (200)
965751	965753	Soft Cup	Econo Box (1200)
965754	965756	Firm Cup	Box (200)
965755	965757	Firm Cup	Econo Box (1200)
965758	965760	Spiral Cup	Box (200)
965759	965761	Spiral Cup	Econo Box (1200)
965762	965764	Pedo Cup	Box (200)
965766	965767	Pointed Tip	Box (100)
965768	965769	Tapered Brush	Box (100)

Nupro Freedom Disposable Prophy Packs

An economical and convenient solution in one package- choose from three popular Nupro prophy paste flavors.

LAVENDER ITEM #	PINK ITEM #	ITEM DESCRIPTION	PACKAGE (QTY)
965685MM	965686MM	Contra Spiral Cup Mint Medium	Box (100)
965685MC	965686MC	Contra Spiral Cup Mint Coarse	Box (100)
965685BM	965686BM	Contra Spiral Cup BubbleExtreme Medium	Box (100)
965685BC	965686BC	Contra Spiral Cup BubbleExtreme Coarse	Box (100)
965716MM	965717MM	Soft Cup Mint Medium	Box (100)
965716MC	965717MC	Soft Cup Mint Coarse	Box (100)
965716RM	965717RM	Soft Cup Razzberry Medium	Box (100)
965716RC	965717RC	Soft Cup Razzberry Coarse	Box (100)
965716BM	965717BM	Soft Cup BubbleExtreme Medium	Box (100)
965716BC	965717BC	Soft Cup BubbleExtreme Coarse	Box (100)



Contra Spiral

Ergonomic angle for improved comfort during polishing
Unique cup designed to reduce splatter and improve interproximal cleaning



Soft

Light to moderate stain removal



Firm

Heavy stain removal



Spiral

Unique spiral cup is designed to reduce splatter with its internal web design
Features outer spiral ridges for improved interproximal cleaning



Pedo

Small tooth polishing



Pointed Tip

Cleaning around ortho appliances



Tapered Brush

Removal of pit and fissure stain before sealant placement and cleaning around orthodontic appliances

Exceptional ergonomics in the palm of your hand.

Nupro RDH® Hygienist Handpiece

Our best-selling Nupro RDH Hygienist Handpiece provides exceptional ergonomics in a lightweight design. The 360° cord swivels to improve maneuverability and access. Wide-back end design neutralizes hand position during procedures with specifically designed contours. And it protects both clinician and patient from splatter with a smooth, low-pressure start. Additional features include:

- Removable, personalized color bands
- Nonslip grip
- Reduces wrist strain
- Autoclavable parts control infection
- Reduces cord drag

Ergonomically designed dental instruments may help reduce the prevalence of musculoskeletal disorders among dental practitioners.⁷


ITEM #	ITEM DESCRIPTION	CONTENTS
740000	Nupro RDH® Hygienist Handpiece for Disposable Angles	1 Handpiece
740010	Nupro RDH® Hygienist Handpiece for Disposable Angles	3 Handpieces
740020	Nupro RDH® Hygienist Handpiece with Prophy Right Angle	1 Handpiece
740030	Nupro RDH® Hygienist Handpiece with Prophy Right Angle	3 Handpieces
740050	Nupro RDH® Featherweight Package - Gray	1 Handpiece 1 Featherweight Hose 1 Hose Installation Kit
740060	Nupro RDH® Featherweight Package - Beige	1 Handpiece 1 Featherweight Hose 1 Hose Installation Kit
740055	Nupro RDH® Featherweight Hose - Gray	1 Featherweight Hose 1 Hose Installation Kit
740065	Nupro RDH® Featherweight Hose - Beige	1 Featherweight Hose 1 Hose Installation Kit



Nupro® Revolv® Slim Disposable Prophylaxis Angles

Discover the Nupro® Revolv® Slim Disposable Prophylaxis Angles (DPAs), with a 25% shorter head and 38% more slender neck.⁵ Our slim angles provides improved visibility, access, and maneuverability—for a more pleasant polishing experience. Additional features include:

- 22% smaller head diameter than standard DPAs⁵
- Reduced noise and vibration
- Available in contra with spiral cup, or straight with choice of soft or spiral cup
- Comes in a unique lavender pearl color




ITEM #	ITEM DESCRIPTION	PACKAGE (QTY)
965262	Slim Contra Spiral (108° Angle)	Box (200)
965263	Slim Soft Cup	Box (200)
965264	Slim Spiral Cup	Box (200)
965069	Slim Contra Spiral (108° Angle)	Econo Box (1200)
965067	Slim Soft Cup	Econo Box (1200)
965068	Slim Spiral Cup	Econo Box (1200)

Nupro® Revolv® Disposable Prophylaxis Angles

Our Nupro Revolv Disposable Prophylaxis Angles are now available in a wide variety of multiple-cup options: soft, firm, pedo, tapered brush, pointed tip, and spiral-cup — for reduced splatter, interproximal cleaning, and improved access around crowns or orthodontics. Choose from assorted colors.

ITEM #	ITEM DESCRIPTION	PACKAGE (QTY)
965212	Contra Prophylaxis Angle (Spiral Cup)	Box (200)
965213	Contra Prophylaxis Angle (Spiral Cup)	Econo Box (1200)
965250	Soft Cup	Box (200)
965251	Soft Cup	Econo Box (1200)
965252	Firm Cup	Box (200)
965253	Firm Cup	Econo Box (1200)
965254	Spiral Cup	Box (200)
965255	Spiral Cup	Econo Box (1200)
965256	Pedo Cup	Box (200)
965258	Pointed Tip	Box (100)
965259	Tapered Brush	Box (100)

Disposable Prophylaxis Angles



Experience improved visibility and access.

9 out of **10** hygienists said they would recommend the Nupro® Revolv® Slim DPAs to their colleagues.⁵



Nupro Revolv Disposable Prophylaxis Packs

An economical and convenient solution in one package- choose from five popular Nupro prophylaxis paste flavors.

ITEM #	ITEM DESCRIPTION	PACKAGE (QTY)
96553MF	Soft Mint Fine	Box (100)
96553MM	Soft Mint Medium	Box (100)
96553MC	Soft Mint Coarse	Box (100)
96553MP	Soft Mint Plus	Box (100)
96553BM	Soft BubbleExtreme™ Medium	Box (100)
96553BC	Soft BubbleExtreme™ Coarse	Box (100)
96553RM	Soft Razzberry Medium	Box (100)
96553RC	Soft Razzberry Coarse	Box (100)
96553OM	Soft Orange Van. Swirl Medium	Box (100)
96553OC	Soft Orange Van. Swirl Coarse	Box (100)
96553FM	Soft Fruit Punch Medium	Box (100)
96553FC	Soft Fruit Punch Coarse	Box (100)
96522IMM	Contra Spiral Mint Medium	Box (100)
96522IMC	Contra Spiral Mint Coarse	Box (100)
96522IOM	Contra Spiral Orange Van. Swirl Medium	Box (100)
96522IOC	Contra Spiral Orange Van. Swirl Coarse	Box (100)
96522IFM	Contra Spiral Fruit Punch Medium	Box (100)
96522IFC	Contra Spiral Fruit Punch Coarse	Box (100)
96554MM	Soft Mint Medium	Econo Box (1000)
96554MC	Soft Mint Coarse	Econo Box (1000)

Designed to deliver consistent, splatter-free performance.

Nupro® Prophy Paste

With a wide variety of flavors to meet the needs of patient preferences, Nupro® Prophy Paste has a unique proprietary formulation designed to deliver consistent, splatter-free performance. We offer the largest selection that rinses easily and completely—featuring excellent stain removal and polishing, with color-coded grit options and 13 gluten-free flavors.

Grit Options:

Fine Grit (Gold) - Light stain removal

Medium Grit (Blue) - Normal cleaning and polishing; light plaque removal.

Coarse Grit (Red) - Heavy stain and plaque removal.

Plus Grit (Green) - Removal of extra heavy stain, including tobacco and heavy plaque deposits.


Nupro® Extra Care® Prophy Paste

Our dye-free Nupro® Extra Care® Prophy Paste with Fluoride and NovaMin® is clinically proven to offer immediate sensitivity relief that lasts up to 28 days, providing increased patient comfort with just one application. NovaMin® calcium phosphate technology releases the natural building blocks of teeth to create an enamel-like mineral layer that cleans and polishes in a single application.

Fluoride Cup Flavors 175 count box

FLAVOR	POLISH		STAIN REMOVAL	
	FINE - MEDIUM	MEDIUM - COARSE	MEDIUM - COARSE	MEDIUM - COARSE
Spearmint	801600		801601	
Peppermint	801602		801603	
Citrus Mint	801604		801605	
Orange	801606		801607	

Nupro Extra Care Prophy Paste



Non-Fluoride Cup Flavors 175 count box

FLAVOR	POLISH	
	FINE - MEDIUM	
Spearmint	801608	

Fluoride Cup Flavors 200 count box

FLAVOR	FINE	MEDIUM	COARSE	PLUS
Mint	801232	801222	801212	801299
Orange	801231	801221	801211	-
BubbleExtreme™	801302	801301	801300	-
Grape	801225	801226	801227	-
Razzberry	801311	801310	801309	-
Zinnamon	801314	801313	801312	-
Cookie Dough	801308	801307	801306	-
Orange Vanilla Swirl	801317	801318	801319	-
CherryBlast	-	801321	801322	-
Chocolate Bliss	801329	801330	801331	-
Strawberry Vanilla Swirl	801326	801327	801328	-
Pina Colada	-	801335	801336	-
Fruit Punch	801332	801333	801334	-

Non-Fluoride Cup Flavors 200 count box

FLAVOR	FINE	MEDIUM	COARSE
Mint	-	801223	801213
Orange	801235	-	-

Jars 340g (12oz)

FLAVOR	MEDIUM	COARSE
Mint	801122	801112
Grape	801126	801127
Orange	-	801111
CherryBlast	-	801171

Variety Packs 200 count box

FLAVOR	FLAVORS INCLUDED	ORDER #
Child - Medium Grit	Strawberry Vanilla Swirl, Mint, Razzberry, BubbleExtreme™	638027
Adult - Coarse Grit	CherryBlast, Mint, Razzberry, Zinnamon	638028

Nupro Prophy Paste



Prophy Grips

Box of 3 grips

801207

Nupro® Shimmer Restorative Polishing Paste

Restore the luster and shine of restorative materials (composites and amalgams) without causing damage to the surface with Nupro Shimmer Restorative Polishing Paste. Formulated with micro-fine aluminum oxide to protect restoration surface.



1 carton (2 syringes)

801351

Topex® Prep & Polish Paste

Convenient, fast, and easy way to clean and prepare tooth surfaces. For use prior to bonding or cementation of orthodontic bonding, placement of sealants and cementation of permanent bonded restorations. The Prep & Polish formulation does not contain additives, oils, waxes or flavors which can interfere with bonding surfaces.



100-2gm cups with prophy ring and 2 Upgrade prophy angles

AD30045

Creating a healthy practice takes a powerful partnership.

At Dentsply Sirona we strive to support your healthy practice by providing you with the most innovative, comprehensive portfolio of preventive products. Because the healthier the practice, the more healthy smiles you can provide.

To find out more about partnering with Dentsply Sirona Preventive to deliver a lifetime of oral health for your patients, visit dentsplysirona.com/preventive

Enhanced Patient Experiences

Efficient Workflows

Improved Predictability

Practice Growth



Soothe-Guard Air®
Lead-Free X-Ray Aprons

VoloWipes®
Disinfecting/Deodorizing/
Cleaning Wipes

Com-Fit Plush™
Natural Fit Masks

Disposashield®
Polyethelene Barriers

Flashtips®
Disposable Air/Water Syringe Tips

Rinn® XCP-DS® Fit
Universal Sensor Holders

Nupro Freedom®
Cordless Prophy System with
SmartMode® Technology

Nupro Freedom® Slim
Disposable Prophy Angles

Cavitron® 300 Series
Ultrasonic Scaling System

Captivate® by Nupro®
Tooth Whitening System

Purevac® HVE System
with Mirror Tips and Hose Adapter

Oraqix®
(lidocaine and prilocaine
periodontal gel) 2.5%/2.5%

Nupro®
White Varnish

Nupro® Extra Care
Prophylaxis Paste with NovaMin®

1. Centers for Disease Control and Prevention. "Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care." Atlanta, GA: US Department of Health and Human Services, Centers for Disease Control and Prevention, National Center for Chronic Disease Prevention and Health Promotion, Division of Oral Health; March 2016.
2. Saylor, C. D., & Overman, P. R. (2011). Hypersensitive Dentin Updates. Hypersensitive Dentin Updates, page 5.
3. OSHA Watch Updates July/August 2007 Volume 9, No. 4 page 2 - References The Quality America Dental OSHA Safety Program Manual Tab 4.
4. Compared to the Midwest® RDH® Freedom Cordless Prophy System
5. Data on file.
6. Compared to Nupro Freedom non-slim DPA
7. Dong H, Barr A, Loomer P, et al. The effects of periodontal instrument handle design on hand muscle load and pinch force. J Am Dent Assoc. 2006;137:1123-1130.

Learn more at dentsplysirona.com/preventive/polishing