

THE DENTAL
SOLUTIONS
COMPANY™



New: Ceiling model LEDview Plus
The better you see,
the more you detect

Natural lighting perfected.

[dentsplysirona.com](https://www.dentsplysirona.com)



New: Ceiling light LEDview Plus and ceiling combination with Heliodent Plus



A great combination

The ceiling model is available as a single ceiling light or as a ceiling combination with Heliodent Plus. The ceiling combination not only creates optimal viewing conditions but also enables you to take intraoral images directly on your treatment unit. This optimally supports the workflow, especially for endodontic and implant treatment. Available for ceiling heights of 2.5 m and 2.5 to 3.4 m.

Intuitive operation

The ceiling model features an integrated key pad where all functions can be operated independently of a treatment unit. On/off switch and composite mode can also be activated via an integrated no-touch sensor.

LEDview Plus

Optimal light for every treatment

LEDview Plus performs with the perfect interplay of light intensity, light color and light field size. It combines the latest LED technology with modern, ergonomically optimized design. Bring natural light into your practice so you can always see everything in perfect light.

Highest light quality

- Natural light from integrated multicolored LEDs
- Ideal working conditions thanks to individually adjustable brightness and color temperatures to match every kind of treatment
- Composite mode for processing filling materials without premature hardening

Optimal visibility

- Perfect and homogeneous illumination of the entire oral cavity with a 200x100 mm light field
- No shadows or glare through optimal arrangement of the light clusters

Easy handling and perfect hygiene

- Three movable joints for maximum flexibility in positioning
- Smooth surfaces and detachable handles for quick and reliable hygiene



More information at: LEDviewplus.com