



OSSIX™ Bone

Ossifying collagen  
sponge

# General guidelines for OSSIX™ Bone sizes and different indications

## Small (5x5x5mm)

- small socket preservation.
- Single unit sinus elevation.
- Immediate implant with small buccal/dehiscence.

## Medium (5x5x10mm)

- large socket preservation (multi-root).
- Multi-unit sinus elevation.
- Lateral ridge augmentation (1-2 units).

## Large (5x10x10mm)

- large defect after extraction of multiple teeth.
- Large defect after removal of cyst.
- Large span of lateral bone augmentation.

OSSIX™ Bone is available in 3 sizes



5x5x5mm (0.125 cc)



5x5x10mm (0.25 cc)



5x10x10mm (0.5 cc)

Order no.	3231 0009	3231 0010	3231 0011
-----------	-----------	-----------	-----------