

SureSmile®

Dentsply  
Sirona

The clear aligner experience  
that will make you and your  
patients smile



## Patient-inspired solutions

SureSmile® Clear Aligners offers solutions throughout the treatment journey to ensure patients are happy from their first aligner to their final smile. Now with SureSmile® VPro™ and SureSmile® Retainer, you can offer your patients more reasons to choose SureSmile® Clear Aligners to achieve their best smile.




## You determine the treatment goal — we provide the plan

With the user-friendly SureSmile® Treatment Planning Software, it's simple to upload digital impressions, X-rays, and patient photos from almost any source. You can rely on the experts at the SureSmile® Digital Lab to do the rest and send you an optimised treatment proposal. Or, if you prefer, you can customise to your own preferences.

**Our open platform fully integrates with all Dentsply Sirona scanners, as well as other commonly used scanners.**

## Review anytime, anywhere

When the SureSmile® Digital Lab has created the initial treatment proposal, you'll receive it for your review and approval. The initial setup is organised and easy. The SureSmile® Clear Aligner software is intuitive and offers a thorough review of each patient's case with a wide range of different tools. Our software is browser-based, so you can review your treatment plans whether you're at your practice, at home, or travelling.



Clear aligners with refinements throughout treatment time up to 3 years from aligner production with flat-rate pricing, now including SureSmile® VPro™\*.

Only pay for the number of aligners you need. You can order SureSmile® VPro™ for an additional fee. This is ideal for less complex cases.

## SureSmile technology has been used to treat over 1,000,000 patients\*\*


Through the course of treatment and beyond, you have the support from Dentsply Sirona with over 100 years in oral healthcare solutions. SureSmile® has more than 20 years of experience providing orthodontic treatment solutions. When you choose SureSmile® Clear Aligners, you can rest assured knowing you have the backing of decades of dental experience.

\* Based on doctor prescription. SureSmile VPro based on market availability.

\*\* SureSmile treated cases, data on file 2023. Not published.

## Clinical excellence

SureSmile® Clear Aligners are made by experts using the latest technology and proven materials. Unique features allow for personalisation for even the most challenging cases. Precise processing of the data and thorough quality control ensure outstanding fit, excellent durability, and aligners so transparent that they're practically undetectable.



### Comfort comes in a variable trimline


With SureSmile® Clear Aligners, you can customise the design of your trimline. An automated, high-precision, laser-cutting process is used to offer the ultimate in comfort. SureSmile® Clear Aligner edges are exceedingly smooth, whether you choose a straight or scalloped trimline.

The straight trimline follows the gingiva at a consistent distance, which you determine. A small overlap of the gingiva can have a positive effect on the force transmission and even reduce the number of necessary attachments. The scalloped trimline follows the course of the gingiva and is even less noticeable than the straight trimline.




## Predictable results

SureSmile Clear Aligners give you and your patients confidence in the treatment process. Each treatment proposal undergoes rigorous quality checks and is validated by a clinical expert from our Digital Lab. In a recent review, no refinements were required in 3 of 4 SureSmile Aligner cases.<sup>1</sup> So, SureSmile Clear Aligners help reduce chair time while delivering the results you and your patients expect.<sup>2</sup>




## SureSmile® Clear Aligners provides excellent force transmission

SureSmile Clear Aligners incorporate high performance materials to handle multiple attachments and are customisable.



## Enhanced treatment planning

3D tooth modification empowers the clinician to take control of the treatment plan. Clinicians have the flexibility to move teeth precisely where they want them. Modifying individual tooth movements enables more efficient treatment planning and approval.




## Smiles by design

When patient X-ray data is provided, a 3D model can be superimposed on photos to plan the final tooth position, in line with the natural facial features of your patient.



## Timelines and IPR scheduling

Our software shows a graphic animation of your patient's treatment plan. The timeline indicates when attachments have to be placed and at which appointments interproximal enamel reduction (IPR) has to be performed and on which teeth.



### Reference:

1. Based on an analysis of global SureSmile aligner cases conducted between Jan. 2019 – Dec. 2021. Data on File.
2. SureSmile Aligner Refinement Data, data on file 2020. Keim/Vogels/Vogels, 2020 JCO Study of Orthodontic Diagnosis and Treatment Procedures (JCO/October 2020 p. 581 et seq.).

# Seamless integration

We've taken our experience in digital planning for orthodontic treatment, combined with our experience as a dental solutions company, and provided you with a simple way to conduct initial setups and treatment plans that result in predictable patient outcomes and compelling results.

The SureSmile® Treatment Planning Workflow can be seamlessly integrated into virtually any practice. Our software is browser-based, which eliminates the need to download directly to your computer. We follow these steps:



## Collect records

Photos and a digital impression are taken at the first appointment. SureSmile® Treatment Planning Workflow accepts scans from Primescan® intraoral scanner and other commonly used scanners. Optional for clear aligner treatment: X-ray data for a more comprehensive diagnosis.



## Preview their new smile

Your patients will be able to see what their new smile could look like before they start treatment. Accessing the Outcome Simulator from DS Core\* facilitates a seamless start to your SureSmile® Clear Aligner workflow. Easily create 3D-visualizations using accurate intraoral scanner data while they're still in the chair.



## Submit the prescription

Rely on our Digital Lab Specialist or specify your detailed treatment preferences.



## Approval or modification of treatment plan

All treatment plans submitted to the SureSmile® Digital Lab go through a comprehensive planning process and a final quality inspection by a Digital Lab Specialist. This ensures excellent initial setups.



## Manufacture the aligners

Once you've approved the treatment plan, SureSmile® Clear Aligners are custom fabricated using an innovative process.



## Start treatment

When the Patient Kit, SureSmile® VPro™ and SureSmile® Clear Aligners arrive at your practice, it's time to start treatment.



## Finish treatment

When the treatment has ended, SureSmile® Retainers are essential to ensuring happy patients' new smiles stay straight.

\* Subscription to DS Core is required.

# Happier patients

SureSmile® VPro™ is a high-frequency vibration device that patients can use at home to help aligners fit properly. SureSmile VPro may reduce discomfort and may reduce treatment time when used five minutes a day with their aligners.<sup>1-3</sup>

## Helping your patients protect their investment

Don't forget to order the SureSmile® Retainer designed to prevent teeth relapse and protect patients' investment in their smile.

## SureSmile® VPro™ App

Patients can pair the SureSmile VPro app with their device in order to:

- Track daily SureSmile VPro usage
- Receive reminders to use device and change aligners as prescribed
- Monitor SureSmile VPro battery life
- Share results and progress photos



Partner with the clear aligner brand from Dentsply Sirona—the global leader in oral healthcare for over 100 years—and help spread the contagious power of smiling.

**References:** 1. Alansari S, Atique MI, Gomez JP, et al. The effects of brief daily vibration on clear aligner orthodontic treatment, J World Fed Ortho. 2018;7(4):134-140. doi:10.1016/j.ejwf.2018.10.002 2. Shipley T. Effects of high frequency acceleration device on aligner treatment—a pilot study. Dent J (Basel). 2018;6(3):32. doi:10.3390/dj6030032 3. Nicozisis J, Brigham G, Sparaga J, Shipley T. Effects of 120Hz high-frequency acceleration device on orthodontic discomfort. Preprints. Preprint posted online October 24, 2018. doi:10.20944/preprints201810.0557.v1

## Education, Training and Support every step of the way

As a SureSmile® provider, you will receive expert support tailored to your practice goals. Complete our Certification program including the SureSmileU® platform. Access ongoing support through DS Academy, an expansive library of educational resources and in-person training events. Receive in-office support through our team of experts across field sales, customer implementation specialists, and customer care associates who work together to help you achieve practice success. Promote the benefits of SureSmile® Clear Aligners with a variety of resources to support patient awareness and education.

### SureSmileU

Our SureSmileU Onboarding helps set you up for success from your very first case. With interactive content to support case submission, tips and tricks videos to lead patient interactions, and printable resources, your practice will be supported and guided during your onboarding. Deliver clear aligner treatment with confidence with SureSmileU.



### The SureSmile® Aligner Certification Program

Discover more about the SureSmile Certification Program, scan the QR code:



Scan the QR code to get in touch.

## SureSmile

[dentsplysirona.com/suresmile](https://dentsplysirona.com/suresmile)

