

Hemoban™ Gel 25% Aluiminum Choloride Hemostatic Solution

Hemoban™ Gel Complete Package 1 Hemoban™ Gel Syringe filled (30mL)

20 Empty Syringes (1.2mL)

20 Infusor Tips

Aquasil® Ultra+ Introductory Kit

1 digit® OLV/OXLV Regular Set Evaluator Kit 2 Cartridges Heavy Fast Set (50mL each)

1 Cartridge LV Fast Set (50mL) 1 Cartridge XLV Fast Set (50mL) 1 Bottle Caulk® Tray Adhesive (14mL) 12 Mixing Tips and 6 Intraoral Tips

Aquasil® Ultra+ Impression Material for digit® Targeted Delivery System

50 single-use cartridges, 50 intraoral mixing tips

	FAST SET		REGULAR SET		
	XLV	678645	Small Cartridge	678640	Small Cartridge
	XLV			678671	Large Cartridge
•	LV	678642	Small Cartridge	678640	Small Cartridge
	LV			678644	Large Cartridge

Primescan™ AC with Connect Software

For more information about Primescan™, consult your Dentsply Sirona Representative.

digit® Targeted Delivery System Dispenser

678007 digit® Syringe Dispenser package (2ct)

50mL Cartridge 4-pack Refills

4 cartridges (50mL ea), 12 mixing tips

	FAST SET	REGULAR SET	EXTENDED SET
XLV	678670	678668	678701
LV	678672	678671	
Medium	678674	678673	
Heavy	678676	678675	678703
Rigid	678678		

50mL Cartridge 24-pack Bulk Refills

24 cartridges (50mL ea)

	FAST SET	REGULAR SET
Medium	678601	678602
Heavy	678626	678625
Rigid	678627	

DECA™ Cartridge Standard Packs

2 Cartridges (380mL ea), Mixing Tips, 1 Bayonet Locking Ring

		FAST SET	REGULAR SET	EXTENDED SET
	Medium	678697	678682	
	Heavy	678684	678686	678704
	Rigid	678688		
•	Soft Putty		678622	

Hand Mix Putty

1 base Jar, 1 catalyst Jar, 2 measuring spoons

		FAST SET	REGULAR SET
	Soft Putty	678620 (900mL)	
		678708 (500mL)	
	Firm Putty		678621 (900mL)
			678709 (500mL)

For more information and to request a demo, contact your Dentsply Sirona Representative. *Some Items not available in all markets

Dentsply Sirona | 13320-B Ballantyne Corporate Place | Charlotte, NC 28277 | 1-800-532-2855 | www.dentsplysirona.com

Competitor brands are not registered trademarks of Dentsply Sirona. @2021 Dentsply Sirona Inc. All rights reserved ML661002B (R 3/26/21)





FINAL IMPRESSION SOLUTION

WORKFLOW SOLUTIONS

CAPTURE EVERY DETAIL QUICKLY AND RELIABLY



FINAL IMPRESSION SOLUTION

Every detail captured quickly and reliably

dentists have problems with crown fit due to impression errors.

A solution that works for your workflow.

Our Final Impression Solution creates a more efficient workflow by effectively stopping minor intraoral bleeding to create a clear field for both conventional and digital workflows. Our products help you capture the data you need so the lab can create the restoration you want.



Hemoban™ Gel 25% Aluminum Chloride Hemostatic Solution

Effectively stop intraoral bleeding—precisely and efficiently—using a transparent gel that stays where you place it and rinses off easily with no residue. Improve your productivity with:

- A clear gel that won't stain
- Performance that won't interfere with bonding
- Hemostasis without coagulation, reducing the risk of secondary bleeding when using cord
- Viscous yet spreadable non-drip gel for easy placement and removal



Cord soaked with liquid hemostatic agent

- Hemostatic agent absorbed fully into cord
- Textured surface



Cord coated with Hemoban™ Gel

- Hemostatic agent available on surface of cord for optimal delivery
- Smooth surface for easy insertion and removal without reoccurring bleeding
- Coated and soaked for easy packing

1. Data on file

CONVENTIONAL IMPRESSION WORKFLOW

Aquasil® Ultra+ Smart Wetting® Impression Material

Capture the right details every time

Only Aquasil® Ultra+ impression material has an optimized set of intraoral performance attributes—hydrophilicity, tear strength, and long work time—coupled with a precision delivery system for improved impression accuracy. Every time.

- Leading intraoral hydrophilicity for excellent wettability to minimize voids and bubbles while maximizing detail capture
- Leading intraoral tear strength for an added layer of confidence in even the most complex cases

	INTRAORAL HYDROPHILICITY (2" @ 80% RH)	CURED-FILM HYDROPHILICITY (5" @ 50% RH)	INTRAORAL TEAR STRENGTH (PSI @ 0.265mm thickness)	24 HOUR TEAR STRENGTH (PSI @ 0.265mm thickness)
AQUASIL® ULTRA+	15°	5°	607	678
IMPRINT™ 4	19°	6°	441	539
TAKE 1™ ADVANCED™	72°	88°	196	338
EXAFAST™ NDS	68°	66°	242	273
FLEXITIME®	55°	58°	291	399
IMPRINT™ 3	70°	43°	505	614
IMPREGUM™	49°	54°	165	259
V-POSIL	52°	42°	414	404



DIGITAL IMPRESSION WORKFLOW

Primescan™

The latest scanner from
Dentsply Sirona provides fast,
accurate, and easy-to-use digital
impressioning. With a large
field of view and sharpness up
to 20mm in depth, Primescan™
can easily capture full arch
impressions in 30 seconds or less.
Primescan™ enables restorative
workflows with your lab or
in-house manufacturing.

digit® Targeted Delivery System

Put precise placement at your fingertips.

Ergonomic Precision

The digit® Targeted Delivery System places your hand just 5.25 inches from the prepared tooth, up to 34% closer than a traditional 50mL cartridge dispenser. It's also significantly lighter, so you can control tip movement and place wash material accurately.

Optimum Work Time

Accurate material delivery and timing are critically important to the success of your final impression.

Our patented delivery system works with singleuse cartridges in multiple viscosities and set times, optimized for single- to multi-unit restoration cases, so you always have the right product for the right case.

