

# Your individual options.

Here is the full equipment list of the basic units as well as all optional features.

Patient chair	
Patient-specific chair positions*	■
Motorized headrest with patient-specific automatic headrest settings*	■
Ambient Light	□
Premium Thermo upholstery	□
Lounge upholstery, with colored stitching on request	□
Lumbar support	■
Massage function	□
Swiveling armrest, right	□
Armrest, left	□
4-way foot switch	■
Smart Control	■
Valve set for external devices	□
Booster cushion for small patients	□
ErgoMotion	■
Adjustable travel speed	■

Dentist element (with Smart Touch, 6 positions, removable handles and instrument holder)	
Smart Delivery	■
Smart Delivery Pro	□
Sprayvit with light	□
1st drive (for motor or turbine)	■
2nd drive (for motor or turbine)	■
3rd drive (for motor or turbine)	■
4th drive (for motor or turbine)	□
BL E electric motor	□
BL ISO F electric motor	□
Turbine hose	■
SiroSonic TL ultrasonic handpiece	□
Face function (fluorescence-aided caries excavation)	□
Implantology function	□
Integrated pump for saline solution	□
BL Implant E electric motor	□
Implant 20:1 contra-angle handpiece	□
Endodontic function incl. reciprocal file systems	□
ApexLocator with display on the Sivision monitor	□
T1 Line Endo 6 L contra-angle handpiece with light	□
USB interface	■
Hands-free operation with cursor control	■
Smart Device holder for tablet	□
Panoramic X-ray viewer	□

Assistant element (with touch display, 4 positions, removable holder)	
Sprayvit with light	□
Suction device	■
Saliva ejector	■
Surgical suction device	□
2nd HVE kit	□
Storage tray for SmartLite Pro curing light	□
Medical-grade 12 V connection	□

Water unit	
Disinfection device with permanent disinfection and sanitization function	■
Chemical cleaning of the suction hose and cuspidor	■
Cuspidor valve for wet suction system	■
Automatic separator	□
Amalgam separator	□
Motorized swiveling glass cuspidor bowl Removable	■
Hygiene system with removable cleaning modules and autopurge function	■
Aqua function (stand-alone water supply)	■

Patient communication	
SiroCam AF+	□
22" Sivision monitor, on tray incl. tray	□
22" Sivision monitor on light support column	□
Camera power supply in the dentist element	□
Cable set for PC connection	□

Working stools (optional long or short gas spring)	
Hugo Freehand (with foot control)	□
Hugo Manual (with manual control)	□
Hugo Manual Plus (with manual control and foot ring)	□
Carl Manual (with manual control)	□
Carl Manual Plus (with manual control and foot ring)	□
Theo with saddle seat	□

Additional options	
Sirona tray	□
LEDview Plus unit model	■
LEDview Plus ceiling model	■
Adapter for Heliident Plus / LEDview Plus ceiling combination	□
Adapter for Heliident Plus unit model	□
Patient mirror	□
Camera holder on the light	□

■ Included in the basic unit or package  
□ Optional

\* Available with Sidexis 4