


Intra-oral X-Ray Positioning Systems and Accessories For Digital Sensors



1. Xios AimRight Sensor Holder System (autoclavable)

Starter Kits

Aim Right Starter Kit for Size 2	B1074230
Aim Right Starter Kit for Size 1	B1074130

Individual Arms for Size 2 Sensors

Aim Right Universal Ring - Red	B1074051
Aim Right Horizontal Bitewing Holder Size 2 - Blue	B1074251
Aim Right Vertical Bitewing Holder Size 2 - Blue	B1074252
Aim Right Anterior Holder Size 2 - White	B1074253
Aim Right Periapical Holder Size 2 - Yellow	B1074254

Individual Arms for Size 1 Sensors

Aim Right Universal Ring - Red	B1074051
Aim Right Horizontal Bitewing Holder Size 1 - Blue	B1074151
Aim Right Vertical Bitewing Holder Size 1 - Blue	B1074152
Aim Right Anterior Holder Size 1 - White	B1074153
Aim Right Periapical Holder Size 1 - Yellow	B1074154



2. Xios Sensor Holder System (Adhesive)

Starter Kits


Xios XG Sensor Starter Kit for Size 2	6411297
Xios XG Sensor Starter Kit for Size 0 and 1	6411289

Individual Holders

Xios anterior sensor holder, blue (100 pcs)	6176510
Xios posterior sensor holder, yellow (100 pcs)	6176528
Xios bitewing sensor holder, red (100 pcs)	6176536
Xios endo universal sensor holder, green (50 pcs)	6176544
Xios endo sensor holder, grey (50 pcs)	6176551

Hygienic Covers for Xios XG Sensors

Hygienic covers, Xios XG, Size 0/1 (300 pcs)	6409960
Hygienic covers, Xios XG, Size 2 (300 pcs)	6409952




3. Rinn XCP-ORA-DS-Fit

XCP-ORA™ - One Ring & Arm Positioning System

The XCP-ORA - One Ring & Arm Positioning System reduces the number of components needed for accurately positioning radiographic media with RINN® XCP® or XCP-DS® bite-blocks for anterior, posterior and bitewings to just one arm, one ring.

XCP-ORA One Ring & Arm, 1-pack
XCP-ORA One Ring & Arm, 2-pack

550771
550772



XCP-DS FIT™

Reduce loading time during your digital radiograph procedures with the XCP-DS FIT Universal Autoclavable Sensor Holder.

- Universal design provides a custom fit to your sensor
- One holder works with both size 1 and size 2 sensors for less inventory
- Autoclavable for lower cost per use
- Rigid sensor support keeps sensor positioned accurately
- Silicone-cushioned edges help provide patient comfort and prevent sensor cover perforations
- Use with the XCP-ORA 3-in-1 Positioning System for easy paralleling



XCP-DS FIT Hygiene Kit


559900

XCP-ORA Arm and Ring, 8 XCP-DS FIT bite-blocks
 2 each: anterior, posterior, horizontal bitewing, vertical bitewing

XCP-DS FIT Complete Kit

559909

XCP-ORA Arm and Ring, Endo Ring and Arm,
 10 XCP-DS FIT bite-blocks
 2 each: anterior, posterior, horizontal bitewing, vertical bitewing, endo
 XCP-DS FIT Bite-block Refill Packs, 2 bite-blocks per pack



Anterior

559901

Posterior

559902

Horizontal Bitewing

559903

Vertical Bitewing

559904

Endo

559905

4. Rinn Snap-A-Ray DS

Snap-A-Ray® DS

Allows ability to quickly load and place digital sensor holders for accurate anterior and posterior periapical radiographs. An optional arm and ring is available for easy alignment.

Snap-A-Ray DS 1 pack	550321
Does not include the optional arm and ring	
Snap-A-Ray DS 3 pack	550323
Does not include the optional arm and ring	
Arm & Ring for Snap-a-Ray DS & DS Endo	550322



5. Rinn Uni-Grip

Uni-Grip® 360 Universal Digital Disposable Sensor Holder

Adjustable bite-block for anterior, posterior, horizontal bitewings, or vertical bitewings in the one holder, allowing the user to take a very fast full mouth series without the need to use multiple bite-blocks or holders. Universal fit for all digital sensors which locks in any size sensor cover for maximum stability - no need for tight fitting bags that make sensor removal hard.

Uni-Grip 360 Disposable Holder for Digital Sensors	550052
50 x Uni-Grips, 50 x Digital Sensor Barriers	
Uni-Grip 360 Combo Kit for Digital Sensors	550252
2 x 550052 kits (50 Uni-Grips and SD Digital Sensor Barrier), 1 x Uni-GripAR	

Uni-GripAR® 3-in-1 Arm & Ring Positioning System

Works with any Uni-Grip holder. Colour-coded for easy assembly & fully autoclavable.

Uni-GripAR Positioning Arm	550100
-----------------------------------	---------------




6. Rinn XCP – DS Fit Endo

XCP-DS FIT Endo Kit **559908**

1 x Endo Ring, 1 Endo Arm, 2 x XCP-DS FIT Endo Bite-blocks

XCP-DS FIT Bite-block Refill Packs, 2 bite-blocks per pack

Endo **559905**



7. Rinn Snap-A-Ray DS Endo

Snap-A-Ray® DS Endo

Similar to the Snap-A-Ray DS with a cut-out bite surface to accommodate rubber dam clamps.

Snap-A-Ray DS Endo 1 pack **550421**

Does not include the optional arm and ring

Snap-A-Ray DS Endo 3 pack **550423**

Does not include the optional arm and ring

Arm & Ring for Snap-a-Ray DS & DS Endo **550322**



Intra-oral X-Ray Positioning Systems and Accessories For Phosphor Plates



8. Rinn XCP

XCP® Kits for Film

Although the colour-coded XCP Arms and Rings have been superseded by the XCP-ORA™ (One Ring&Arm), they are still available to purchase along with the colour-coded XCP bite-blocks for film.


XCP Kit 2000 with Bitewing 542001

4 x Arms: 1 Anterior, 2 Posterior, 1 x Bitewing
4 x Aiming Rings: 2 Anterior/Bitewing, 2 x Posterior
17 x Bite Blocks: 4 Anterior, 8 Posterior, 5 x Bitewing
(2 x #2H, 1 x #0/1H, 1 x #2V, 1 x #1V)
1 x Instruction Poster



Arms sold individually

Anterior (Blue)	540857
Posterior (Yellow)	540858
Bitewing (Red)	540927
Endodontic (Green)	550597



Rings sold individually

Anterior (Blue)	540865
Posterior (Yellow)	540860
Bitewing (Red)	540934
Endodontic (Green)	550598



Bite-blocks

 Anterior (x25) (Blue)	540861
 Posterior (x25) (Yellow)	540862
 Ant/Post (x10/x10) (Blue/Yellow)	540863
 Endodontic (x3) (Green)	540935

Bitewing

#2 Horizontal (x6) (Red)	540929
#1 Horizontal (x6) (Red)	540933
#2 Vertical (x6) (Red)	540241
#1 Vertical (x6) (Red)	540922



Bitewing Bite-blocks for Kodak SURESOFT Film Packets

Dual-slot design accommodates both SURESOFT and standard #2 film packets.

XCP Horizontal Bitewing Bite-block for Kodak SURESOFT Film Packets (pk 2)	540629
XCP Vertical Bitewing Bite-block for Kodak SURESOFT Film Packets (pk 6)	540630



9. Rinn Snap-A-Ray Xtra

Snap-A-Ray® Xtra - Film & Phosphor Plate Holder

The Snap-A-Ray Xtra Film & Phosphor Plate Holder brings an extra measure of efficiency and comfort to anterior and posterior bisecting angle technique radiographs.

Snap-A-Ray Xtra Film & Phosphor Plate Holder, 1-pack 540501

Snap-A-Ray Xtra Film & Phosphor Plate Holder, 3-pack 540503



10. Rinn Uni-Grip

Uni-Grip® Disposable Radiography Holder for Film/PSP

Accommodates ANY #1 or #2 radiograph media; attaches in seconds. Reposition bitepieces in just a moment. Colour-coded channels ensure correct positioning. Compact size and rounded edges are easy on the patient.

Uni-Grip Disposable Holder for Film and Phosphor Plates 550055

50 x Uni-Grips, 50 x Digital Sensor Barriers

Uni-Grip Combo Kit for Film and Phosphor Plates 550255

2 x 550055 kits (50 Uni-Grips and SD Digital Sensor Barrier), 1 x Uni-GripAR

Uni-GripAR® 3-in-1 Arm & Ring Positioning System

Works with any Uni-Grip holder, colour-coded for easy assembly & fully autoclavable.

Uni-GripAR Positioning Arm 550100



11. Rinn XCP Endo Kit

XCP® Endo Kit 542045

Contains: XCP® Positioning Arm endo (1), XCP® Positioning Ring Endo (1), XCP® Biteblocks Endo (3)

3 x XCP® Biteblocks 540935

1 x XCP® Positioning Arms 550597

1x XCP® Positioning Rings 550598



12. Rinn Endo Ray II

EndoRay® II Film Holder

Prevent cone cutting with precise angulation and alignment of film to the x-ray cone. Fits over files, clamps and rubber dams.

EndoRay II, 2 pack 540303
Replacement Parts (1 per pack)

EndoRay II Arm 540310

EndoRay II Basket 540311

EndoRay II Ring 540312



Others – Digital Sensors

Universal Digital Sensor Covers

Protect your sensor with RINN's premium quality sensor covers

- Fits all #1 and #2 size sensors
- Rounded edges for comfort

Universal Digital Sensor Covers, (pk 500) 550500



RINN® XCP Kits For Digital Sensors

Although XCP-DS® has now been superseded by XCP-ORA (One Ring&Arm) and the universal, stretchy digital sensor holders, XCP-DS FIT™; the XCP-DS individual Arms & Rings and the dedicated bite-blocks for individual brands are still available.

The Arm, Ring and Bitepieces include:


Aiming Rings sold individually

Anterior (Blue)	540865
Posterior (Yellow)	540860
Bitewing (Red)	540934
Endo (Green)	550598






Arms sold individually

Anterior (Blue pins)	540857
Posterior (Yellow pins)	540858
Bitewing (Red pins)	540927
Endo (no pins)	550597



Bitepieces 3 per package

	Anterior (Blue)	550567
	Posterior (Yellow)	550568
	Bitewing (Red)	550570
	Endo (Green)	550562



Kits are available for the following brands of digital sensors and each kit includes arms, rings and colour-coded bite-blocks/sensor holders:

XCP-DS Kit

Schick CDR*	551900
Schick Elite*	552971
Sirona*	551150
Sirona XIOS*	551160
Dent X Eva +*	552101
Dexis*	552500
Gendex e2v*	551001
Kodak 5100/6100*	553001
Trophy Kodak RVG 5000/6000*	551070
Cygnus Black*	551010

*Each kit is suitable for a number of different digital sensors. Please refer to the table below.

XCP-DS Kit Name	Other systems the XCP-DS* Kit is suitable for
Schick CDR*	Schick CDR PlusWire, CDR Wireless, Belmont, BelSensor
Gendex e2v*	Gendex VisualiX/GX-S/HDI, Accent, Acteon Sopro, Cyber Medical XDR, Cygnus, Dentrax, Dürr, e2v, HDX, Mediadent, Progeny MPS/MPSe, and Visident systems
Trophy Kodak RVG 5000/6000*	RVG Access, RVG Ultimate
Trophy Kodak RVG 5*	RVG ui, RVG Reference, RVG Gold
Suni*	Suni metal-style only; not for Suni-Ray or Dr. Suni Plus Also suitable for CamSight, DMD, Escape, Lightyear, and Multimedia

*For full details and details of XCP-DS replacement baskets and bitepieces for specific digital sensors, please contact your local Dentsply Sirona Clinical Specialist.

Disposable Bite-Blocks, Bitewing Holders & Mouth Props

Disposable XCP® Bite-blocks, 100 pieces #0, #1 or #2 film	540870
EZ-Prop® Mouth Props, Styrofoam x 100 pieces (accommodates dental floss)	540900

