


Atlantis® abutments implant compatibility chart



| Manufacturer | Implant | Atlantis Healing Abutment | | Atlantis Abutment | | | Atlantis Crown ² | | Atlantis CustomBase solution - Atlantis Abutment & Atlantis Crown ² | | Atlantis Crown Abutment | | Atlantis Conus Abutment | |
|---|---|---------------------------|-------------------------|-------------------|-------------------------|-----------------------|-----------------------------|------------------------|--|----------|-------------------------|----------|-------------------------|--------------|
| | | | | | | | Cut-back or Full-contour | | Cut-back or Full-contour | | Cut-back | | Custom | Over-denture |
| | | Abutment material | | | | | | | | | | | | |
| | | Titanium Gold-shaded ti | Titanium Gold-shaded ti | Zirconia | Titanium Gold-shaded ti | Zirconia ¹ | Titanium Gold-shaded ti | Angulated screw access | Titanium | Zirconia | Titanium Gold-shaded ti | Titanium | | |
| Dentsply Sirona | INTERNAL | | | | | | | | | | | | | |
| | EV Connection | | | | | | | | | | | | | |
| | PrimeTaper EV 3.6 - (S) Purple | x | x | | x | | x | x | x | x | | x | x | |
| | PrimeTaper EV 4.2 - (M) Yellow | x | x | | x | | x | x | x | x | | x | x | |
| | PrimeTaper EV 4.8 - (L) Blue | x | x | | x | | x | x | x | x | | x | x | |
| | PrimeTaper EV 5.4 - (L) Blue | x | x | | x | | x | x | x | x | | x | x | |
| | OmniTaper EV 3.0 - (XS) Green | x | x | | x | | x | x | x | x | | | | |
| | OmniTaper EV 3.4 - (S) Purple | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | OmniTaper EV 3.8 - (M) Yellow | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | OmniTaper EV 4.5 - (L) Blue | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | OmniTaper EV 5.5 - (XL) Brown | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | Astra Tech Implant EV 3.0 - (XS) Green | x | x | | x | | x | x | x | | | | | |
| | Astra Tech Implant EV 3.6 - (S) Purple | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | Astra Tech Implant EV 4.2 - (M) Yellow | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | Astra Tech Implant EV 4.8 - (L) Blue | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | Astra Tech Implant EV 5.4 - (XL) Brown | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | Astra Tech Implant EV Profile | | | | | | | | | | | | | |
| | Astra Tech Implant EV Profile 4.2 Yellow | x | x | x | x | | x | x | x | | | | x | x |
| | Astra Tech Implant EV Profile 4.8 Blue | x | x | x | x | | x | x | x | | | | x | x |
| | Astra Tech Implant System TX | | | | | | | | | | | | | |
| | OsseoSpeed/OsseoSpeed TX 3.0 - Yellow | x | x | | x | | x | x | x | | | | | |
| | OsseoSpeed/OsseoSpeed TX 3.5, 4.0 - Aqua | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | OsseoSpeed/OsseoSpeed TX 4.5, 5.0 - Lilac | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | OsseoSpeed TX Profile 4.5, 5.0 | x | x | | x | | x | | x | | | | x | x |
| | OsseoSpeed 4.0 - Blue (Previous) | | x | | | | | | | | | | x | |
| | Xive / Frialit | | | | | | | | | | | | | |
| | Xive S 3.0 | x | x | | x | | x | x | x | | | | | |
| | Xive S 3.4, 3.8, 4.5, 5.5 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| Frialit 3.4, 3.8, 4.5, 5.5 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x | |
| Ankylos | | | | | | | | | | | | | | |
| Ankylos C/X - C/ 3.5, 4.5, 5.5, 7.0 (non-indexed) | | x | x | x | | | | | | | | | x | |
| Ankylos C/X - /X 3.5, 4.5, 5.5, 7.0 (indexed) | x | x | x | x | x ¹ | x | x | | | | | | x | x |
| MIS | MIS | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | C1 NP 3.3 | x | x | | x | | x | x | x | | | | x | x |
| | C1 SP 3.75, 4.2 | x | x | | x | | x | x | x | | | | x | x |
| | C1 WP 5.0 | x | x | | x | | x | x | x | | | | x | x |
| | M4 NP 3.3 | x | x | | x | | x | x | x | | | | x | x |
| | M4 SP 3.75, 4.2 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | M4 WP 5.0, 6.0 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | SEVEN NP 3.3 | x | x | | x | | x | x | x | | | | x | x |
| | SEVEN SP 3.75, 4.2 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | SEVEN WP 5.0, 6.0 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | x |
| | V3 NP 3.3 | x | x | | x | | x | x | x | | | | x | x |
| V3 SP 3.9, 4.3, 5.0 | x | x | | x | | x | x | x | | | | x | x | |

All trademarks and company names are the property of their respective owners.

Please turn over →

¹ Manufactured in Waltham, US. Turnaroundtime: 5-10 extra working days (excluding US Shipping).

² Includes zirconia crowns and Atlantis Temporary Crowns in PMMA.

| Manufacturer | Implant | Atlantis Healing Abutment | | Atlantis Abutment | | Atlantis Crown ² | | Atlantis CustomBase solution – Atlantis Abutment & Atlantis Crown ² | | Atlantis Crown Abutment | | Atlantis Conus Abutment | | |
|---|--|---------------------------|-------------------------|-------------------|-------------------------|-----------------------------|-------------------------|--|----------------|-------------------------|-------------------------|-------------------------|----------------|---|
| | | | | | | Cut-back or Full-contour | | Cut-back or Full-contour | | Cut-back | | Custom | Over-denture | |
| | | Abutment material | | | | | | | | | | | | |
| | | Titanium Gold-shaded ti | Titanium Gold-shaded ti | Zirconia | Titanium Gold-shaded ti | Zirconia ¹ | Titanium Gold-shaded ti | Angulated screw access | Titanium | Zirconia | Titanium Gold-shaded ti | Titanium | | |
| BioHorizons | BioHorizons | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | Internal 3.5, 4.5, 5.7 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | |
| | Tapered Internal 3.5, 4.5, 5.7 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | |
| Biomet 3i | Biomet 3i | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | Certain 4/3, 5/4 | x | x | x | x | x ¹ | x | | | x | x | x | x | |
| | Certain 3.25, 4.0, 5.0, 6.0 | x | x | x | x | x ¹ | x | | | x | x | x | x | |
| | Certain Prevail 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 | x | x | x | x | x ¹ | x | | | x | x | x | x | |
| | Certain XP 4/5, 5/6 | x | x | x | x | x ¹ | x | | | x | x | x | x | |
| | EXTERNAL | | | | | | | | | | | | | |
| | Miniplant 3.25 | | x | x | | | | | | x | x | x | x | |
| Osseotite 3.75, 4.0, 5.0, 6.0 | | x | x | | | | | | x | x | x | x | | |
| Osseotite XP 4/5, 5/6 | | x | x | | | | | | x | x | x | x | | |
| Camlog | Conelog | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | CONOLOG 3.3 | x | x | | x | | x | | | x | | | x | x |
| | CONOLOG 3.8, 4.3 | x | x | | x | | x | | | x | | | x | x |
| | CONOLOG 5.0 | x | x | | x | | x | | | x | | | x | x |
| | Screw-line Implant 3.3 | x | x | | x | | x | | | x | | | x | x |
| Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0 | x | x | x | x | x ¹ | x | | | x | x | | x | x | |
| Keystone Dental | Keystone Dental | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | PrimaConnex SD 3.3, 3.5, Genesis 3.8 | x ¹ | x ¹ | x ¹ | x ¹ | x ¹ | x ¹ | | | x ¹ | x ¹ | x ¹ | x ¹ | |
| | PrimaConnex RD 4.0, 4.1, Genesis 4.5 | x ¹ | x ¹ | x ¹ | x ¹ | x ¹ | x ¹ | | | x ¹ | x ¹ | x ¹ | x ¹ | |
| PrimaConnex WD 5.0, Genesis 5.5, 6.5 | x ¹ | x ¹ | x ¹ | x ¹ | x ¹ | x ¹ | | | x ¹ | x ¹ | x ¹ | x ¹ | | |
| Nobel Biocare | Nobel Biocare | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | NobelActive 3.0 | x | x | | x | | x | x | x | | | | | |
| | NobelActive/NobelReplace Conical Connection NP 3.5 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | |
| | NobelActive/NobelReplace Conical Connection RP 4.3, 5.0 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | |
| | NobelActive WP 5.5 | x | x | x | x | x ¹ | x | x | x | x | x | x | x | |
| | NobelReplace NP 3.5 (Tri-channel connection) | x | x | x | x | x ¹ | x | | x | x | x | x | x | |
| | NobelReplace RP 4.3 (Tri-channel connection) | x | x | x | x | x ¹ | x | | x | x | x | x | x | |
| | NobelReplace WP 5.0 (Tri-channel connection) | x | x | x | x | x ¹ | x | | x | x | x | x | x | |
| | NobelReplace 6.0 (Tri-channel connection) | x | x | x | x | x ¹ | x | | x | x | x | x | x | |
| | EXTERNAL | | | | | | | | | | | | | |
| Brånemark System NP 3.3 | | x | | x | | x | | | x | | | x | x | |
| Brånemark System RP 3.75, 4.0 | | x | x | x | | x | | | x | x | | x | x | |
| Brånemark System WP 5.0, 5.5 | | x | x | x | | x | | | x | x | | x | | |
| Osstem | Osstem | | | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | | | |
| | Osstem TS Mini 3.5 | x | x | x | x | | x | | | x | x | x | x | |
| Osstem TS Regular 4.0, 4.5, 5.0, 6.0, 7.0 | x | x | x | x | | x | | | x | x | x | x | | |

All trademarks and company names are the property of their respective owners.

Please turn over 

¹ Manufactured in Waltham, US. Turnaroundtime: 5-10 extra working days (excluding US Shipping).

² Includes zirconia crowns and Atlantis Temporary Crowns in PMMA.

| Manufacturer | Implant | Atlantis Healing Abutment | Atlantis Abutment | | Atlantis Crown ² | | Atlantis CustomBase solution – Atlantis Abutment & Atlantis Crown ² | | Atlantis Crown Abutment | | Atlantis Conus Abutment | |
|------------------------------------|---|---------------------------|-------------------------|----------|-----------------------------|-----------------------|--|------------------------|-------------------------|----------------|-------------------------|--------------|
| | | | | | Cut-back or Full-contour | | Cut-back or Full-contour | | Cut-back | | Custom | Over-denture |
| | | Abutment material | | | | | | | | | | |
| | | Titanium Gold-shaded ti | Titanium Gold-shaded ti | Zirconia | Titanium Gold-shaded ti | Zirconia ¹ | Titanium Gold-shaded ti | Angulated screw access | Titanium | Zirconia | Titanium Gold-shaded ti | Titanium |
| Straumann | Straumann | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | |
| | Bone Level 3.3 NC (Narrow CrossFit) | x | x | x | x | x ¹ | x | x | x | x | x | x |
| | Bone Level 4.1, 4.8 RC (Regular CrossFit) | x | x | x | x | x ¹ | x | x | x | x | x | x |
| | Standard/Standard Plus 4.8 RN (Regular Neck) | x | x | | x | | x | | x | | x | x |
| | Standard/Standard Plus 4.8 WN (Wide Neck) | x | x | | x | | x | | x | | x | x |
| | EXTERNAL | | | | | | | | | | | |
| Standard Plus 3.3 NN (Narrow Neck) | | x ¹ | | | | | | | | x ¹ | | |
| Zimmer Dental | Zimmer | | | | | | | | | | | |
| | INTERNAL | | | | | | | | | | | |
| | Tapered Screw-Vent 3.5, 4.5, 5.7 | x | x | x | x | x ¹ | x | x | x | x | x | x |
| | Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit | x | x | | x | | x | | x | | x | x |
| | Screw-Vent 3.3, 3.7, 4.5 | x | x | x | x | x ¹ | x | x | x | x | x | x |
| | Screw-Vent 3.3, 3.7, 4.5 friction fit | x | x | | x | | x | | x | | x | x |
| | SwissPlus 4.8 | x | x | | x | | x | | x | | x | x |

All trademarks and company names are the property of their respective owners.

¹ Manufactured in Waltham, US. Turnaroundtime: 5-10 extra working days (excluding US Shipping).

² Includes zirconia crowns and Atlantis Temporary Crowns in PMMA.