

Compatibility list Flex/Green/White Lid adapters for DAC Universal GUI/Text/D/S

Manufacturer	Current portfolio	Instrument	Ref. manufacturer	Type	Indication	Adapter Ref. Sirona	Adapter Sirona US Scaler tip	Adapter Sirona US Scaler handpiece (Sirosonic / L)	Adapter Sirona US Scaler handpiece (Sirosonic TL)	Adapter EMS US Scaler tip (Piezon®, Piezon® LED)	Adapter EMS US Scaler handpiece (Piezon®, Piezon® LED)	Adapter KaVo SONICflex Handpiece (2003/L, 2008/L, 2008S/LS)	Adapter KaVo SONICflex 2003 tip adapter (only for SONICflex 2003/L tips)	Adapter Satelec US Scaler tip	Adapter Sirona Sprayvit nozzle G	Adapter Sirona Sprayvit nozzle 4000	Adapter EMS AIR-FLOW-HANDY 3.0	Adapter Air-Flow Handpiece	Adapter Satelec Slim	Adapter Satelec Newton	Adapter Satelec Newton LED
EMS	yes	Instrument A	DS-001 A	Tip	Scaling	6610708				x											
EMS	yes	Instrument P	DS-011 A	Tip	Scaling	6610708				x											
EMS	yes	Instrument PS	DS-016 A	Tip	Scaling	6610708				x											
EMS	yes	Instrument B	DS-003 A	Tip	Scaling	6610708				x											
EMS	yes	Instrument C	DS-002 A	Tip	Scaling	6610708				x											
EMS	yes	Instrument PL1	DS-030 A	Tip	Perio	6610708				x											
EMS	yes	Instrument PL2	DS-031 A	Tip	Perio	6610708				x											
EMS	yes	Instrument PL3	DS-032 A	Tip	Perio	6610708				x											
EMS	yes	Instrument PL4	DS-033 A	Tip	Perio	6610708				x											
EMS	yes	Instrument PL5	DS-034 A	Tip	Perio	6610708				x											
EMS	yes	Instrument HPL3	DS-075 A	Tip	Perio	6610708				x											
EMS	yes	Instrument DPL3	DS-065 A	Tip	Perio	6610708				x											
EMS	yes	Piezon® LED Handpiece	EN-060	Handpiece	n.a.	6613538					x										
EMS	yes	Piezon® Handpiece	EN-061	Handpiece	n.a.	6613538					x										
EMS	yes	Piezon® Handpiece	EN-041	Handpiece	n. a.	6750090					x										
EMS	yes	AIR-FLOW Handpiece*	EL-308	Handpiece	n.a.	6623461													x		
EMS	yes	AIR-FLOW Handpiece*	EL-354	Handpiece	n.a.	6623461													x		
EMS	yes	AIR-FLOW Handy 3.0	EL-545	Handpiece nozzle	n.a.	6623453											x				
EMS	yes	AIR-FLOW Handy 3.0	EL-542	Handpiece nozzle	n.a.	6623453											x				
EMS	yes	AIR-FLOW Handy 3.0	EL-540	Handpiece nozzle	n.a.	6623453											x				
KaVo	yes	KaVo PiezoLED	1.007.3995	Handpiece	n.a.	6613538						x									
KaVo	yes	KaVo SONICflex 2003**	1.000.4246	Handpiece	n.a.	6732056						x									
KaVo	yes	KaVo SONICflex LUX 2003 L**	1.000.4243	Handpiece	n.a.	6732056						x									
KaVo	yes	KaVo SONICflex quick 2008**	1.005.9311	Handpiece	n.a.	6732056						x									
KaVo	yes	KaVo SONICflex quick 2008 L**	1.005.9310	Handpiece	n.a.	6732056						x									
KaVo	yes	KaVo SONICflex 2003 Spitzen**	several	Tip	check Dental-Information 09/2020	6735646							x								
Sirona	yes	Ultrasonic tip 1L	8924250	Tip	Scaling	6536127	x														
Sirona	yes	Ultrasonic tip 2L	8924268	Tip	Scaling	6536127	x														
Sirona	yes	Ultrasonic tip 3L	1894455	Tip	Scaling	6536127	x														
Sirona	yes	Ultrasonic tip 4L	1894463	Tip	Scaling	6536127	x														
Sirona	yes	Ultrasonic tip SiroRetro KL	5457887	Tip	Endodontics	6536127	x														
Sirona	yes	Ultrasonic tip SiroRetro KR	5457895	Tip	Endodontics	6536127	x														
Sirona	yes	Ultrasonic tip SiroRetro GL	5457903	Tip	Endodontics	6536127	x														
Sirona	yes	Ultrasonic tip SiroRetro GR	5457911	Tip	Endodontics	6536127	x														

