

THE DENTAL
SOLUTIONS
COMPANY™



Dentsply Sirona Treatment Centres

The focal point
of your practice

dentsplysirona.com/treatment-centres



Dentsply Sirona Treatment Centres

Treatment Centres are the focal point of every practice. Featuring a modern design, integrated functions and high-quality materials, our treatment centres transform the dentist and patient experience from start to finish. We make sure our units provide everything you need, so that your sole focus is on your patient and making them smile.

With a history of more than 130 years of innovation, we have a strong commitment to providing first-class quality for fundamental reliability. This is partnered with the continuous drive to deliver meaningful innovations that enhance comfort, workflow and experience.





1 Design

What makes Dentsply Sirona treatment centres unique is the harmony of its modern and ergonomic design concept in which form follows function. High-quality materials and smooth surfaces are designed to support the highest hygiene standards. The distinctive streamlined contours make our treatment centres instantly recognisable and leave a lasting impression.



2 Integration

You need a treatment centre with functions meticulously adapted to your workflow – and not the other way around. This is why our integrated solutions include Endo and Implant functions, an integrated infection prevention system and patient communication tools to support your treatment consultations. Seamless workflows and no additional set up time are just as important as our focus on ergonomics.



3 Experience

When your treatment centre fits your needs perfectly, you can focus fully on your patient. When everything comes together like this, work no longer feels like work, but like an experience. It is our goal to give you the tools to provide better dental care to your patients and ultimately make them smile. This is part of the reason why dentist and patient comfort is equally important to us. The Dentsply Sirona treatment centres provide numerous solutions for an ergonomic, painless and stress-free treatment, for the dentist, patient and team.

Teneo – Innovation in perfect harmony

Teneo takes an uncompromising approach: Only the best is good enough. Your treatment centre responds automatically and individually to all your needs. Carefully selected features assist you in every treatment situation. Comfortable features ensure that your patients feel more relaxed throughout the treatment, so you can focus entirely on your work. Why settle for less?

1 Motorised sliding track

To let you focus fully on your patients, Teneo and its motorised sliding track places everything you need at your fingertips – hands-free and effortless.

2 Dentist element

With a 7" EasyTouch Display and 6 + 1 instrument positions (incl. auxiliary storage tray position for the intraoral camera) you are ensured intuitive and self-explanatory operation.

3 Motorised

For optimum access any time, even with varying patient heights. Can be easily adjusted using the foot control for hands-free treatment, in addition to user interface operation and manual adjustment.

4 Integrated infection prevention


Maximum infection control in minimum time with minimum effort that saves your team valuable time.

5 Wireless foot control

Hands-free operation for all functions of the EasyTouch display, instruments and the motorised sliding track.

6 Motorised cuspidor

The motorised cuspidor automatically gets into the right position without the patient leaning over extensively.






Split concept for a better handling range

The dentist and assistant elements, as well as the tray, can be positioned independently allowing you to switch flexibly between team and solo treatment modes. You and your assistant can maintain an ergonomic working position at all times without the additional stress of worrying whether there is enough space on the tray and around the patient.



Holistic infection prevention concept

Easy-to-clean surfaces, removable components and built-in sanitisation adapters make hygiene processes efficient and safe. Instrument and suction hoses simply need to be placed on the integrated sanitisation adapters and the process starts at just the press of a button. The cleaning agent container for chemical flushing of the suction lines is integrated into the water unit and can be easily refilled. Additionally, the concentration of cleaning agent can be adjusted individually.



Integrated treatment functions


Integrated treatment functions for implantology and endodontics make setup times unnecessary, now that additional table equipment is no longer required. All functions can be operated with the wireless foot control. This allows a smooth transition between patient examination, education, and treatment.

Speed and torque are pre-set in an integrated library containing a variety of file and implant systems. Additional software functions can be added at a later time with plugins.

Sinius – More efficiency

Time is one of our most valuable assets, and the Sinius treatment centre delivers clear gains in terms of time and flexibility, as the unit can be configured to support any situation and need. Intuitive operation with the latest touch screen technology and numerous integrated functions provide further added value and benefit to your workflow.

Efficient workflows



Integrated endodontics and implantology

The integrated therapy functions make tabletop units a thing of the past. In addition to precise torque control, your endo function is optionally available with an ApexLocator for maximum clinical safety. Your implant treatments can be simplified with individual settings for speed, torque and flow rate of sterile cooling fluid for increased control at any time.

Holistic infection prevention concept

Easy-to-clean surfaces, removable components and built-in sanitisation adapters make hygiene efficient and safe. Instrument and suction hoses simply need to be placed on the integrated sanitisation adapters and the process starts at just the press of a button. The cleaning agent container for chemical flushing of the suction lines is integrated into the water unit and can be easily refilled. Additionally, the concentration of cleaning agent can be adjusted individually.

Treatment concepts



Sinius with manual sliding track



Sinius Traditional Style (TS) with over the patient arm design



Sinius Continental Style (CS) with whip arms

1 Dentist element

With a 4,3" EasyTouch Display and 5 + 1 instrument positions (incl. auxiliary storage tray position for the intraoral camera).

2 Motorised headrest


For optimum access any time, even with varying patient heights. Can be easily adjusted using the foot control for hands-free treatment, in addition to user interface operation and manual adjustment.

3 Integrated infection prevention

Maximum infection control in minimum time with minimum effort that saves your team valuable time.

4 Wireless foot control

With ultimate positioning freedom and no cable clutter, you can easily and freely navigate around the unit. Control all chair functions with the EasyTouch, completely hands-free.



Intego Pro – Performance you can rely on

The highlight of Intego Pro is undoubtedly its comfort water unit, bringing a comprehensive infection prevention concept to your practice. The integrated disinfection system reduces the workload supported by your team through automatic purging and sanitisation of the water lines, as well as a permanent disinfection with Dentosept S.

1 Dentist element

Available as TS or CS version. With the EasyPad or the 4,3" EasyTouch user interface, optional with torque control and ApexLocator, as well as 5+1 instrument positions.

2 Double-articulating or flat headrest

Double-articulating headrest for quick and flexible adaptations or flat headrest with a variety of positions facilitated by the magnetic cushion.

3 Comfort water unit and comfort assistant element

The comfort water unit and assistant element are based on a comprehensive infection prevention concept for safe and efficient hygiene processes.

4 Pneumatic or electronic foot control


Choose between the two foot controls for hands-free operation of your instruments.





Intego Ambidextrous

Optional Ambidextrous version with either Intego or Intego Pro to adjust the treatment centre for left-handed dentists. Fast, easy and elegant.



Integrated endodontic function

An ideal solution for increased safety, efficiency and simplicity, speed and torque can be preset for 6 different files. Additionally, certified WaveOne® Gold and RECIPROC® blue settings are integrated into the dentist element. The motor stops automatically when the set maximum torque is reached. An ApexLocator can also be integrated into the EasyTouch to allow for permanent control of the working length.



Integrated disinfection system

Integrated cleaning adapters make automatic purging of instrument hoses and sanitisation of the water lines easy. Furthermore, Intego Pro allows automatic sanitisation for regular intensive disinfection and permanent disinfection of all water lines with Dentosept S. It can also be used as a self-contained water supply.




Chemical suction hose cleaning

The cleaning agent container for chemical flushing of the suction lines is integrated into the water unit and can be easily refilled. Additionally, the concentration of cleaning agents can be individually adjusted. Moreover, the suction systems are compatible with all conventional central wet or dry suction systems, including an amalgam separator.

Intego - Performance you can rely on

If you are looking for a treatment centre to optimally support your personal workflows at a value, Intego is your choice. With simple operation, it provides all relevant functions when and where you need them. The ergonomic design offers high comfort for both you and your patients. The compact water unit with a fresh water bottle makes it easy to comply with infection prevention regulations.



Integrated endodontic function

As an ideal solution for increased safety, efficiency and simplicity, speed and torque can be preset for 6 different files. Additionally, certified WaveOne® Gold and RECIPROC® blue settings are integrated into the dentist element. The motor stops automatically when the set maximum torque is reached. An ApexLocator can also be integrated into the EasyTouch to allow for permanent control of the working length.



Fresh water bottle

The compact water unit can be adapted to fit many different needs and requirements through freely selectable suction systems and infection prevention options. A fresh water bottle can be integrated for controlled water quality. It can also be filled with Dentosept S solution for permanent disinfection and manual sanitisation.



Automatic suction hose cleaning

Cleaning adapter for automatic flushing of saliva ejector and HVE. The suction systems are compatible with all conventional wet and dry suction systems, including an amalgam separator.



Intego Ambidextrous

Optional Ambidextrous version with either Intego or Intego Pro to adjust the treatment centre for left-handed dentists. Fast, easy and elegant.

1 Dentist element

Available as TS or CS version. With the EasyPad or the 4,3" EasyTouch user interface, optional with torque control and ApexLocator, as well as 5+1 instrument positions to adapt to your needs.

2 Double-articulating or flat headrest

Double-articulating headrest for quick and flexible adaptations or the flat headrest with a variety of positions facilitated by the magnetic cushion.

3 Compact or comfort assistant element

Both assistant elements are compact, saving valuable space and ensuring maximum freedom of movement. With the comfort model, you can work comfortably even without an assistant.



Transforming the Patient Experience

No matter which treatment centre you choose – every Dentsply Sirona unit is tailored to the needs of your patients and as a result, can help you positively affect the way they remember the treatment. Patient experience is more than just comfortable positioning or friendly personnel. A dentist who knows which experiences their patients “collect” during their visit can improve their offering, fostering trust and loyalty.




Sidexis 4

A transparent treatment process contributes significantly to patient loyalty. With Sidexis 4, integrated and readily available at the unit, Teneo and Sinus provide clear on-site visualisation of the clinical situation and the ability to store individual patient positions. The software also allows easy access to patient and x-ray data for quick communication.



Patient communication

Establish a trust-based relationship with your patients through treatment education and clear communication. With the integrated monitor and the ability to display intraoral images captured with the intraoral camera, 3D x-ray images and presentations, patient consultation has never been simpler.



Intraoral camera

Optimise the visualisation and documentation of your treatment: The SiroCam F or AF+ combines excellent image quality with easy operation. The image is displayed on the 22" monitor in high-definition quality.



Ergonomics

Available for Teneo and Sinus, the motorised headrest is the key feature for optimum access to the treatment site. Adjustment is easy with the foot control, user interface or even manual adjustment, ensuring better positioning for any treatment.



Patient comfort

The ErgoMotion function provides combined movements of the backrest and seat for an ergonomic extension of the spine. Teneo even offers the OrthoMotion function, which is particularly relevant for treating elderly or disabled patients: The backrest and seat retain the same angle, protecting the lumbar vertebrae.



Your treatment centre. Your choice.

Working Stools



**Hugo –
the ergonomics expert**

Hugo promotes ergonomic and intuitive sitting and enables you to maintain a healthy posture throughout the treatment process. To suit the individual users' sizes and working habits, it is available in several variations.



**Carl –
the all-rounder**

The Carl round stool ensures maximum freedom of movement and meets the highest demands in comfort and quality. The round seat features comfortable upholstery and supports the highest hygiene standards.



**Theo –
the flexible one**

The saddle stool combines the healthy posture supported by standing with the advantages and stability that you get in a comfortable seated position. The elevated seat position allows optimal access to the patients.

Operating Lights



Intego /
Intego Pro
only



	LEDlight Plus
Illumination	approx. 5,000 – 30,000 lx in 5 levels
Color Temperature	approx. 5000 K not adjustable
CRI	>90
Light field	approx. 160 x 80 mm at a distance of 700 mm

	LEDview Plus
Illumination	approx. 5,000 – 40,000 lx in 7 levels
Color Temperature	Adjustable in 3 levels: approx. 4,600 – 6,200 K
CRI	Natural lighting: 95 at 5,400 K
Light field	Bigger light field: approx. 200 x 100 mm at a distance of 700 mm

Upholstery colors

The elegant world



Classic, timeless combinations for a truly sophisticated look.

The natural world




Soft, warm shades for a soothing treatment atmosphere.

The vital world



Vivid colors and bright tones for a modern, lively atmosphere.

Stitching colors**



Standard Red Dark grey Light blue Rose

Paintwork colors for Teneo and Sinius



Uncompromising quality



Premium upholstery

The high-quality, seamless upholstery ensures optimal comfort, even during long treatments.



Premium Thermo upholstery

Reduces the accumulated heat in the seat and back area and provides for patient relaxation as a result of this cooling effect.



Lounge upholstery

Its innovative material gives this upholstery option additional softness — for that bit of extra comfort.



Make your treatment centre one-of-a-kind. With a motorised headrest and lounge upholstery, featuring custom embroidery of two initials, you can create a truly unique and tailored feel in your practice environment.***


* Colors available for Lounge upholstery.

** With lounge upholstery

*** Optional. Available with a motorised headrest and Lounge upholstery.

You decide what fits your needs

Teneo



Feature	
Treatment concept	Motorised sliding track
Infection prevention	<ul style="list-style-type: none">• Holistic infection prevention concept with motorised cuspidor• Disinfection system for continuous disinfection• Integrated sanitisation adapters for automatic periodic cleaning of water lines• Integrated suction hose cleaning
Ergonomics	<ul style="list-style-type: none">• OrthoMotion• ErgoMotion• Vertical Mouth Point Positioning• Split Concept• Patient specific chair positions and automatic headrest adjustment• Lumbar support and massage function• Slow motion travel
Dentist element – Display	7" EasyTouch user interface
Dentist element – Positions	6+1 instruments positions
Assistant element	4 positions with operating keyboard
Headrest	Motorised
Foot control	Wireless
Left-handed version	No
Operating lights	LEDview Plus

Sinius



Feature	
Treatment concept	<ul style="list-style-type: none"> • Manual sliding track • Traditional Style (TS) with over the patient arm design • Continental Style (CS) with whip arm design
Infection prevention	<ul style="list-style-type: none"> • Holistic infection prevention concept with cuspidor • Disinfection system for continuous disinfection • Integrated sanitisation adapters for automatic periodic cleaning of water lines • Integrated suction hose cleaning
Ergonomics	<ul style="list-style-type: none"> • ErgoMotion • Lumbar support function • Split Concept • Vertical Mouth Point Positioning
Dentist element – Display	4,3" EasyTouch user interface
Dentist element – Positions	5+1 instruments positions
Assistant element	4 positions with operating keys
Headrest	<ul style="list-style-type: none"> • Motorised • Double-articulating
Foot control	Wireless or with cable
Left-handed version	No
Operating lights	LEDview Plus

Intego Pro



Feature	
Treatment concept	<ul style="list-style-type: none"> • Traditional Style (TS) with over the patient arm design • Ambidextrous Style (CS) with whip arm design
Infection prevention	<ul style="list-style-type: none"> • Holistic infection prevention with cuspidor • Disinfection system for continuous disinfection • Integrated sanitisation adapters for automatic periodic cleaning of water lines • Integrated suction hose cleaning
Ergonomics	ErgoMotion
Dentist element - Display	4.3" EasyPad or EasyTouch user interface
Dentist element - Positions	5 instruments positions (+1 with TS)
Assistant element	Comfort or Compact
Headrest	<ul style="list-style-type: none"> • Flat • Double-articulating
Foot control	Wireless or with cable
Left-handed version	Yes (Intego Pro Ambidextrous)
Operating lights	<ul style="list-style-type: none"> • LEDlight Plus • LEDview Plus

Intego



Feature	
Treatment concept	<ul style="list-style-type: none"> • Traditional Style (TS) with over the patient arm design • Ambidextrous Style (CS) with whip arm design
Infection prevention	<ul style="list-style-type: none"> • Fresh water bottle solution with cuspidor • Optional integrated suction hose cleaning system with water
Ergonomics	ErgoMotion
Dentist element - Display	4.3" EasyPad or EasyTouch user interface
Dentist element - Positions	5 instruments positions (+1 with TS)
Assistant element	Comfort or Compact
Headrest	<ul style="list-style-type: none"> • Flat • Double-articulating
Foot control	Wireless or with cable
Left-handed version	Yes (Intego Ambidextrous)
Operating lights	<ul style="list-style-type: none"> • LEDlight Plus • LEDview Plus

Australia ☎ 1300 552 929 🌐 www.dentsplysirona.com.au | **New Zealand** ☎ 0800 336 877 🌐 www.dentsplysirona.co.nz
ANZ ✉ clientservices@dentsplysirona.com 📍 11-21 Gilby Road, Mount Waverley, VIC 3149

