

Let us help you evolve your practice.

Contact your Dentsply Sirona representative
directly or call 844-848-0137.

dstreatmentcenters.com/teneo

Dentsply Sirona Treatment Centers
Better, Safer, Faster Dentistry

Enabling Technologies
CAD/CAM
Imaging
Treatment Centers
Instruments

Teneo
By Dentsply Sirona

The premium treatment center for intelligent dentistry

THE DENTAL
SOLUTIONS
COMPANY™

 **Dentsply
Sirona**

THE DENTAL
SOLUTIONS
COMPANY™

 **Dentsply
Sirona**

TENEO

Constantly evolving to meet your needs, Teneo is for clinicians who want the most modern and comfortable technology while maintaining the highest standards of quality and workflow. With a variety of equipment options for customization, Teneo is an investment in the future of your practice.

Integration

Teneo is the hub of your practice, enabling you to expand treatment options and keep more patients in-house.

Design

Featuring clean, flowing lines and premium materials, Teneo promotes safe, healthy, and comfortable treatment for patients and practitioners alike.

Experience

Teneo inspires comfort, ease, and confidence the moment patients settle into the chair, allowing doctor and staff to work with utmost efficiency and productivity.

Transition to implant or endo procedures with a **single touch**: Built-in implantology and endodontic functions eliminate bulky, obtrusive tabletop units and time-consuming setups while increasing productivity and profitability.

COMMUNICATION AND CONNECTIVITY




Information on-demand

Access to Windows Media Player, Sidexis 4, PowerPoint, intraoral camera images, and educational materials optimize workflows for increased productivity.




Digital connectivity

You can readily connect external USB devices to the treatment center, and also connect the treatment center to the office network for easy troubleshooting, simplified servicing, firmware/software updates and practice networking.



Sidexis 4

View your entire network through the HD monitor and control Sidexis 4 via the EasyTouch touchscreen for immediate access to patient records.




PROCEDURAL INTEGRATION



Integrated implantology

Integrated implant library streamlines workflow, while precise torque and speed control maximizes patient safety. Built-in peristaltic pump enables efficient delivery of sterile solution into surgical area.



Integrated endodontics

ApexLocator, precise torque and speed control, along with the WaveOne® reciprocal file system and integrated file library, simplifies workflow and helps ensure the highest degree of clinical efficacy.

TECHNOLOGY INTEGRATION



CEREC® AI

CEREC AI (Acquisition Integrated) gives you immediate access to all CEREC functions via a separate arm and PC.



Schick 33

Built-in USB port makes digital X-rays easy: Plug-and-play and you are ready to go.



Heliodent Plus

Smooth and steady tube positioning ensures clear, high-definition X-ray images.




Clean lines, smooth surfaces, and premium materials make for a simple yet sophisticated appearance that is stylish and practical.

Lounge upholstery


Soft-touch seating surfaces provide patients with a pleasant, cushioned, and supportive feel at any incline.




Basalt



Carbon




Dune




Orchid




Pacific



Rosewood



Carbon with red stitching



EasyTouch user interface

Fully customizable for the way you work, the easy-to-read touchscreen ensures intuitive and self-explanatory operation.



Motorized chair lift

Allows patients of all ages and sizes to comfortably get in and out of the chair.



Premium materials

Lightweight, high-strength construction offers long-term durability and simplified surface cleaning and disinfection.



6 + 1 instrument positions

Six instrument positions on the delivery allow you to adapt the configuration to your requirements with an optional intraoral camera external placement as the 7th position.

Ergonomically designed with integrated functions for a smooth, seamless workflow, Teneo ensures a healthy and comfortable procedural environment. Achieve optimal results from your first patient to your last.



INTUITIVE SEATING



Motorized track dentist element

Hands-free positioning of dentist element in the preferred ergonomic position for ideal workflow.



Wireless foot control and 4-way foot switch

For hands-free, multifunctional control of all EasyTouch display functions including endodontic and implant procedures.

COMFORTABLE POSITIONING



OrthoMotion & ErgoMotion

OrthoMotion permits the entire chair to be tilted without altering the angle between the seat and the backrest, while ErgoMotion reduces overstretching of the patient's head, and eliminates compression of the spine.



Patient position memory & massage function

Individual sitting positions can be saved for each patient, and the active lumbar support and massage functions help patients feel calm and relaxed.

OPTIMUM VISIBILITY



Motor-driven headrest


Allows optimal vertical mouth positioning for unencumbered access to treatment area and ensures a comfortable working position for the dentist and hygienist.



LEDview Plus light

3-axis movement for ideal positioning and neutral light output allows you to clearly see all areas of the mouth and teeth without color shift or distortion.

INTEGRATED WORKFLOW



Patient communication

SiroCam AF+ captures detailed digital images that can be shown and explained to patients via the 22" HD monitor for improved case acceptance.



Infection prevention

Built-in functions and adapters for automatic sanitization ensure maximum operatory safety at all times.

TENEQ EQUIPMENT OPTIONS

Options for the patient chair

Patient-specific chair positions*	<input checked="" type="checkbox"/>
Motor-driven headrest with automatic patient specific headrest adjustment*	<input checked="" type="checkbox"/>
Lounge upholstery	<input checked="" type="checkbox"/>
Up to 6 programmable patient treatment positions	<input checked="" type="checkbox"/>
Lumbar support / massage function	<input checked="" type="checkbox"/>
Armrest	<input checked="" type="checkbox"/>
4-way foot switch	<input checked="" type="checkbox"/>
Wireless electronic foot switch	<input checked="" type="checkbox"/>
Booster cushion for small patients	<input type="checkbox"/>

Water unit

Automatic sanitization system	<input checked="" type="checkbox"/>
Suction hose cleaning with chemicals	<input checked="" type="checkbox"/>
Cuspidor	<input type="checkbox"/>
Purge function with autopurge	<input checked="" type="checkbox"/>

■ Base Model

□ Optional

I = Implant Package
E = Endo Package
C = Communications Package

Dentist element

Sprayvit M syringe with light & water heater	<input checked="" type="checkbox"/>
BL ISO C electric motor	<input checked="" type="checkbox"/>
BL Implant electric motor	<input type="checkbox"/>
Additional BL ISO C electric motor	<input type="checkbox"/>
Turbine hose	<input checked="" type="checkbox"/>
SiroSonic TL ultrasonic handpiece	<input checked="" type="checkbox"/>
Implantology function	<input type="checkbox"/>
Endo function with WaveOne® reciprocal file system and integrated ApexLocator	<input type="checkbox"/>
ApexLocator with display on Sivation monitor*	<input type="checkbox"/>
Integrated pump for saline solution	<input checked="" type="checkbox"/>
USB interface	<input checked="" type="checkbox"/>
Hands-free operation with cursor control	<input checked="" type="checkbox"/>
Contra-angle handpiece T1 Line Endo 6 L with light	<input type="checkbox"/>
Implant 20:1 handpiece	<input type="checkbox"/>

Assistant element

Surgical suction device	<input checked="" type="checkbox"/>
Sprayvit M syringe with light	<input checked="" type="checkbox"/>
Suction device	<input checked="" type="checkbox"/>
Saliva ejector	<input checked="" type="checkbox"/>
HVE suction hose with slide	<input checked="" type="checkbox"/>

Digital communication

SiroCam AF+	<input type="checkbox"/>
22" monitor on the light support system	<input type="checkbox"/>
Cable set for PC connection	<input type="checkbox"/>
Sidexis 4	<input checked="" type="checkbox"/>
PC interface via EasyTouch	<input checked="" type="checkbox"/>
Remote diagnosis	<input checked="" type="checkbox"/>

Stools

1st stool, Hugo Manual**	<input type="checkbox"/>
2nd stool, Hugo Manual Plus**	<input type="checkbox"/>

Additional options


Tray with support arm	<input type="checkbox"/>
LEDview Plus operating light	<input type="checkbox"/>
Heliodent Plus mounted on light post***	<input type="checkbox"/>
CEREC AI mounted on light post***	<input type="checkbox"/>

*Available with Sidexis 4.

**Hugo stools can be substituted with the Carl model.

***Heliodent Plus and CEREC AI cannot be combined.

TENEQ OPTIONS



Hugo



Carl



Booster cushion



Cuspidor



“My fully integrated office allows me to provide specialty treatments at a high clinical level and provides the patient a better experience.”

— Dr. Sam Bollwinkel

When your dentistry is extraordinary, your treatment center is Teneo.

LEDview PLUS

With four different spotlights, seven intensity levels, three color temperature settings, and cure-safe mode, the LEDview Plus delivers accurate and precise illumination to maximize the view of the oral cavity from any position.



CONNECTIVITY

View your entire network through the monitor. Control Sidexis 4, Windows Media Player, and more via the EasyTouch touchscreen.



SANITIZATION

All instrument lines, including all vacuum, can be connected to the integrated sanitization system for periodic cleaning and monthly sanitization.



IMPLANTOLOGY

No need for a separate tabletop unit—with Teneo, implantology is part of the package—allowing you to seamlessly incorporate your Galileos-CEREC® integration workflow.




SIROCAM AF+ CAMERA

With its sophisticated autofocus technology, SiroCam AF+ delivers high resolution and excellent detail reproduction. At the touch of a button, the image is automatically focused in just a fraction of a second.



ENDODONTICS WITH APEXLOCATOR

Integrated endodontics makes tabletop units unnecessary and, along with precise torque control, also offers an ApexLocator for the highest degree of clinical safety.



“The Teneo story is pretty powerful. It reduces chairtime for the patient with a streamlined workflow.”

— Dr. Ryan O’Neill

