

GAC's In-Ovation® L

Leading the Way in
Self-Ligating Lingual Brackets



The Perfect Esthetic Solution

More and more, adults want to experience the benefits of a perfect smile, without the obvious appearance of metal or ceramic brackets. Now you have a completely undetectable solution to offer your patients with In-Ovation® L.

In-Ovation L Patient Benefits

- Absolutely undetectable
- Small, comfortable to wear
- Speak clearly and confidently
- Extended intervals between appointments
- Short treatment time
- Minimal chair time

Unprecedented Ease of Use

In-Ovation L leverages all the advantages of the In-Ovation self-ligating line of brackets to solve one of the biggest problems doctors have traditionally encountered with lingual systems — archwire ligation. The In-Ovation L clip opens and closes easily and quickly to make archwire ligation literally a “snap.”


Doctor and Staff Benefits:

- Interactive spring clip providing passive, expressive and active wire engagement for optimal control at each phase of treatment.
- Advanced manufacturing process leverages In-Ovation System success.
- Completely undetectable solution option for patients seeking invisible braces
- Low profile for patient comfort and avoiding gingival interference for better patient experience
- Beveled slot walls facilitates wire engagement and superior sliding mechanics
- True Twin Design offers optimal mesial, distal span for achieving superior rotation control



Experience the difference of In-Ovation L

GAC has translated the proven clinical benefits of In-Ovation into a lingual design. Only an interactive bracket provides full functionality throughout the course of treatment. This allows you to enjoy the early benefits of a free sliding, passive design and the precise control of an active one in the later stages of treatment. The In-Ovation L self-ligating bracket gives you a complete range of control options ... simply by changing the archwire.



Sentalloy Lingual Archwires


Size	Reference Number Small		Reference Number Medium		Reference Number Large		Reference Number X-Large	
	Maxillary	Mandibular	Maxillary	Mandibular	Maxillary	Mandibular	Maxillary	Mandibular
.012"	02-911-101	02-911-601	02-911-102	02-911-602	02-911-103	02-911-603	02-911-104	02-911-604
.014"	02-911-111	02-911-611	02-911-112	02-911-612	02-911-113	02-911-613	02-911-114	02-911-614

Neo Sentalloy Lingual Archwires

Size	Reference Number Small		Reference Number Medium		Reference Number Large		Reference Number X-Large	
	Maxillary	Mandibular	Maxillary	Mandibular	Maxillary	Mandibular	Maxillary	Mandibular
.017" x .025"	02-923-131	02-923-631	02-923-132	02-923-632	02-923-133	02-923-633	02-923-134	02-923-634

Resolve Beta Titanium Universal Archwires

Size	Reference Number
.016"	03-016-000
.016" x .022"	03-622-000



GAC's In-Ovation® L



Schedule Your Free Lunch and Learn with GAC!

We'll bring lunch to your office and provide a hands-on demonstration for you and your staff!¹

To sign up for your FREE GAC Lunch and Learn, visit:

dentsplysirona.com/orthodontics/ortholearn

1-800-645-5530

¹ Some geographical restrictions may apply.

The Dental Solutions Company

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies. We create industry-leading positions and platforms across consumables, equipment, technology, and specialty products. We are committed to introducing new, innovative products and complementary, end-to-end solutions to drive better, safer and faster dentistry.

To learn more, visit us online at **dentsplysirona.com**.



Dentsply GAC International
One CA Plaza, Suite 100
Islandia, NY 11749

©2016 Dentsply GAC Intl.
All Rights Reserved
0280-MM-120-090-01 Rev. 02