

## SELF-LIGATING BUCCAL TUBES

Maximize your In-Ovation treatment goals with In-Ovation self-ligating buccal tubes. They provide easy wire placement, quick snap closure and interactive technology.

Maxillary Tooth	T	A	Of	.018		.022	
				Right	Left	Right	Left
U6 Single	-10	0	14	88-161-05	88-261-05	88-162-05	88-262-05
U6 Single	-10	0	14	88-161-05RB	88-261-05RB	88-162-05RB	88-262-05RB
U6 Single	-10	0	14	87-161-05	87-261-05	87-162-05	87-262-05
U6 Single	-14	0	10	88-261-81	88-611-02	88-162-81	88-262-81
U6 Single	-14	0	10	88-261-81RB	88-611-02RB	88-162-81RB	88-262-81RB
U6 Single	-14	0	10	87-261-81	87-611-02	87-162-81	87-262-81
U6 Double	-20	0	20	88-161-62	88-261-62	88-162-62	88-262-62

Mandibular Tooth	T	A	Of				
				Right	Left	Right	Left
L6 Single	-25	0	4	88-361-09	88-461-09	88-362-09	88-462-09
6 Single	-25	0	4	88-361-09RB	88-461-09RB	88-362-09RB	88-462-09RB
6 Single	-25	0	4	87-361-09	87-461-09	87-362-09	87-462-09
L6 Single	-20	0	0	88-361-81	88-461-81	88-362-81	88-462-81
L6 Single	-20	0	0	88-361-81RB	88-461-81RB	88-362-81RB	88-462-81RB
L6 Single	-20	0	0	87-361-81	87-461-81	87-362-81	87-462-81
L6 Single	-14	0	4	88-361-62	88-461-62	88-362-62	88-462-62
L6 Double	-10	0	4	88-361-08	88-461-08	88-362-08	88-462-08



## In-Ovation® C by GAC Self-Ligating Ceramic Brackets

# Comprehensive Prescription Guide

## The Dental Solutions Company

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies. We create industry-leading positions and platforms across consumables, equipment, technology, and specialty products. We are committed to introducing new, innovative products and complementary, end-to-end solutions to drive better, safer and faster dentistry.

To learn more, visit us online at [dentsplysirona.com](http://dentsplysirona.com).



Dentsply GAC International  
One CA Plaza, Suite 100  
Islandia, NY 11749

©2017 Dentsply Sirona  
All Rights Reserved  
0280-MM-120-100-09 Rev. 05



## CCO

The CCO System is a comprehensive treatment philosophy designed to treat each and every case with goal directed treatment and pre-visualized results. The CCO Rx was developed to take full advantage of the bracket/archwire interaction when using an active clip and to achieve optimal tooth position at the end of treatment.

Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	12	5	0	H100-111-07	H100-211-07	H100-112-07	H100-212-07
U1/ltq	7	5	0	H100-111-08	H100-211-08	H100-112-08	H100-212-08
U2	10	9	0	H100-121-07	H100-221-07	H100-122-07	H100-222-07
U2/ltq	3	9	0	H100-121-08	H100-221-08	H100-122-08	H100-222-08
U3/hk	-7	10	2M	H100-131-17	H100-231-17	H100-132-17	H100-232-17
U3/hk ltq	0	10	2M	H100-131-18	H100-231-18	H100-132-18	H100-232-18
U4-5	-9	0	0	H100-541-07	H100-541-07	H100-542-007	H100-542-07
U4-5/hk	-9	0	0	H100-141-17	H100-241-17	H100-142-17	H100-242-17

Mandibular							
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	-6	0	0	H100-611-08	H100-611-08	H100-612-08	H100-612-08
L1-2/htq	-1	0	0	H100-611-07	H100-611-07	H100-612-07	H100-612-07
L3/hk	-8	3	0	H100-331-17	H100-431-17	H100-332-107	H100-432-17
L4†	-12	2	0	189-341-07	189-441-07	189-342-07	189-442-07
L4/hk†	-12	2	0	189-341-17	189-441-17	189-342-17	189-442-17
L5†	-17	-1	0	189-351-07	189-451-07	189-352-07	189-452-07
L5/hk†	-17	-1	0	189-351-17	189-451-17	189-352-17	189-452-1

## BASE Rx

Dentsply Sirona Orthodontics Base Rx has prescription values that are equivalent to the Roth prescription. ROTH is a registered trademark of Roth Licensing, LLC. All rights reserved.

Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	12	5	0	H100-111-00	H100-211-00	H100-112-00	H100-212-00
U2	8	9	0	H100-121-00	H100-221-00	H100-122-00	H100-222-00
U3	-2	13	4M	H100-131-00	H100-231-00	H100-132-00	H100-232-00
U3/hk	-2	13	4M	H100-131-10	H100-231-10	H100-132-10	H100-232-10
U4-5	-7	0	2D	H100-141-00	H100-241-00	H100-142-00	H100-242-00
U4-5/hk	-7	0	2D	H100-141-10	H100-241-10	H100-142-10	H100-242-10

Mandibular							
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	-1	2	0	H100-311-00	H100-411-00	H100-312-00	H100-412-00
L3	-11	7	2M	H100-331-00	H100-431-00	H100-332-00	H100-432-00
L3/hk	-11	7	2M	H100-331-10	H100-431-10	H100-332-10	H100-432-10
L4**	-17	-1	4D	H100-341-00	H100-441-00	H100-342-00	H100-442-00
L4/hk**	-17	-1	4D	H100-341-10	H100-441-10	H100-342-10	H100-442-10
L4/ob**	-17	-1	4D	H100-341-08	H100-441-08	H100-342-08	H100-442-08
L5**	-22	-1	4D	H100-351-00	H100-451-00	H100-352-00	H100-452-00
L5/hk	-22	-1	4D	H100-351-10	H100-451-10	H100-352-10	H100-452-10

\*\*In-Ovation Csr lower bicuspid brackets are made from zirconia. The high strength properties of Zirconia allow the bracket to be designed with low profile tie-wings that are swept to avoid occlusal interference.

## EURO

Euro Rx has prescription values that are equivalent to the \*MBT™ prescription.

Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	17	4	0	H100-111-80	H100-211-80	H100-112-80	H100-212-80
U2	10	8	0	H100-121-80	H100-221-80	H100-122-80	H100-222-80
U3/hk	-7	8	0	H100-131-81	H100-231-81	H100-132-81	H100-232-81
U4	-7	0	0	H100-541-80	H100-541-80	H100-542-80	H100-542-80
U4/hk	-7	0	0	H100-141-81	H100-241-81	H100-142-81	H100-242-81
U5	-7	0	0	H100-551-80	H100-551-80	H100-552-80	H100-552-80
U5/hk	-7	0	0	H100-151-81	H100-251-81	H100-152-81	H100-252-81

Mandibular							
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	-6	0	0	H100-611-80	H100-611-80	H100-612-80	H100-612-80
L3/hk	-6	3	0	H100-331-81	H100-431-81	H100-332-81	H100-432-81
L4†	-12	2	0	189-341-80	189-441-80	189-342-80	189-442-80
L4/hk†	-12	2	0	189-341-81	189-441-81	189-342-81	189-442-81
L5†	-17	2	0	189-351-80	189-451-80	189-352-80	189-452-80
L5/hk†	-17	2	0	189-351-81	189-451-81	189-352-81	189-452-81

\*MBT is a registered trademark of 3M Unitek. All rights reserved.

## RONCONE

In-Ovation C Roncone prescription features increased torque in the maxillary incisors and upper premolars, and reduced mandibular posterior torque. The prescription is based on the JSOP philosophy developed by Dr. Roncone.

Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	18	5	0	H100-111-60	H100-211-60	H100-112-60	H100-212-60
U2	10	8	0	H100-121-60	H100-221-60	H100-122-60	H100-222-60
U3	0	8	2M	H100-131-60	H100-231-60	H100-132-60	H100-232-60
U3/hk	0	8	2M	H100-131-61	H100-231-61	H100-132-61	H100-232-61
U4-5	-10	4	2D	H100-141-60	H100-241-60	H100-142-60	H100-242-60
U4-5/hk	-10	4	2D	H100-141-61	H100-241-61	H100-142-61	H100-242-61

Mandibular							
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	1	0	0	H100-611-60	H100-611-60	H100-612-60	H100-612-60
L3	-3	7	2M	H100-331-60	H100-431-60	H100-332-60	H100-432-60
L3/hk	-3	7	2M	H100-331-61	H100-431-61	H100-332-61	H100-432-61
L4†	-7	-1	2D	189-341-60	189-441-60	189-342-60	189-442-60
L4/hk†	-7	-1	2D	189-341-61	189-441-61	189-342-61	189-442-61
L5†	-7	-1	2D	189-351-60	189-451-60	189-352-60	189-452-60
L5/hk†	-7	-1	2D	189-351-61	189-451-61	189-352-61	189-452-61

\* Dentsply Sirona Orthodontics does not manufacture lower bicuspid ceramic brackets. We recommend utilizing In-Ovation R on the lower 4 and 5 bicuspids.

