

GAC's Ovation[®] C

Advanced Ceramic Esthetics



GAC's Ovation[®] C Ceramic Brackets

Precision Movement, Predictable Results

The Ovation C bracket is our most-advanced, traditional ceramic bracket. Engineered to deliver precision and control, Ovation C esthetically acclimates to its environment while offering optimal stain resistance. With a durability that defies ceramics, Ovation C brackets are strong enough to rebond as needed. All while delivering the predictable results you want, the patient comfort you expect and the esthetics you demand. It's everything you've come to expect from a GAC bracket.

Advanced Ceramic Technology

Esthetics

Translucent ceramic that blends with enamel while resisting stains.

Strength


Durability that defies traditional ceramics, strong enough to rebond.

Performance

100% mechanical base for predictable results and easy debonding.



Esthetics Without Compromise



Bracket ID System

for quick bracket identification and management before, during and after the procedure.




Debonding Indent

for easier debonding, indent provides positive seat for debonding instrument.



Smooth Swept Tie-Wings

for hassle-free ligation and enhanced patient comfort



Triple Chamfered Slot Walls

for easier engagement, reduced opportunity for archwire binding or crimping and increased inter-bracket span.



Rhomboid Shape

for better esthetics during treatment



100% Mechanical Base

for reliable bonding, debonding and ability to reposition.



Vertical Scribe Line

provides a visual reference for alignment with the long axis of the clinical crown

GAC's Ovation[®] C Ceramic Brackets



Schedule Your Free Lunch and Learn with GAC!

We'll bring lunch to your office and provide a hands-on demonstration for you and your staff!¹

To sign up for your FREE GAC Lunch and Learn, visit:

dentsplysirona.com/orthodontics/ortholearn

1-800-645-5530

¹ Some geographical restrictions may apply.

The Dental Solutions Company

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies. We create industry-leading positions and platforms across consumables, equipment, technology, and specialty products. We are committed to introducing new, innovative products and complementary, end-to-end solutions to drive better, safer and faster dentistry.

To learn more, visit us online at **dentsplysirona.com**.



Dentsply GAC International
One CA Plaza, Suite 100
Islandia, NY 11749

©2016 Dentsply GAC Intl.
All Rights Reserved
0280-MM-120-310-00 Rev. 04